Campylobacter Blood-Free Selective Agar (CCDA) (Dehydrated)*

REF CM0739B, CM0739R and CM0739T

EN

*These Instructions for Use (IFU) are intended to be read in conjunction with the IFU for Oxoid Limited's CCDA Selective Supplements (product codes: SR0155E, SR0155H).

Intended Use

Campylobacter Blood-Free Selective Agar (CCDA) (Dehydrated) (CM0739B, CM0739R or CM0739T) when supplemented with CCDA Selective Supplements (SR0155E or SR0155H), is a blood-free selective medium intended for the isolation of *Campylobacter* species from faecal samples.

Campylobacter Blood-Free Selective Agar (CCDA) (Dehydrated), with CCDA Selective Supplement added, is intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having campylobacteriosis.

The device is for professional use only, is not automated, nor is it a companion diagnostic.

Summary and Explanation

Campylobacter species are Gram-negative rods in appearance. Two of the most medically relevant *Campylobacter* species isolated from human samples are *Campylobacter coli* and *Campylobacter jejuni*.¹ The symptoms of *Campylobacter* spp. infections are characterised by vomiting, nausea, body aches, abdominal pain and in some cases bloody diarrhoea.² In acute stages, the symptoms of *Campylobacter* infections mimic ulcerative colitis and appendicitis.³ Long term complications of *Campylobacter* infections include arthritis, Guillain Barr Syndrome (GBS) and Irritable Bowel Syndrome (IBS).

The implications of *Campylobacter* infections can vary from person to person.⁴ In most cases, pregnant ladies, patients with AIDS, other blood disorders like thalassemia, patients undergoing chemotherapy and those with a weakened immune system are at high risk.⁵ Moreover, the resistance shown by the *Campylobacter* species to fluoroquinolones like ciprofloxacin complicates treatment.⁶ CCDA Selective Medium (PO0119A), and Campylobacter Blood-Free Selective Agar Base (Dehydrated) (CM0739B/R/T) with CCDA Selective Supplements (SR0155E/H) are therefore a key tool for the clinician to support the isolation of *Campylobacter* species from faecal samples.

Principle of Method

The formulation of Campylobacter Blood-Free Selective Agar (CCDA) (Dehydrated) (CM0739B, CM0739R and CM0739T) contains a nutrient broth comprised of peptone and Lab Lemco powder which with casein hydrolysate are sources of nitrogen and other nutrients. Ferrous sulphate, sodium pyruvate and charcoal replace the blood used in other formulations of campylobacter culture media. Charcoal acts as an adsorbent for toxic compounds. Sodium desoxycholate is added to suppress the growth of Gram-positive organisms. CCDA Selective Supplements (SR0155E and SR0155H) added to DCM contains cefoperazone, a cephalosporin that acts as a broad- spectrum antimicrobial with activity against a wide range of microorganism including Enterobacteriaceae and *Pseudomonas* species. Amphotericin B is added to inhibit fungi, this is required when medium is incubated at 37°C to prevent overgrowth of yeasts.

Typical Formula

	grams per litre
Nutrient Broth No. 2	25.0
Activated carbon	4.0
Casein hydrolysate	3.0
Sodium desoxycholate	1.0
Iron (II) sulphate	0.25
Sodium pyruvate	0.25
Agar	12.0

Materials Provided

- CM0739B: 500g pot
- CM0739R: 2.5kg pot
- CM0739T: 5kg pot

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Selective supplements (SR0155E or SR0155H)
- Petri Dishes

Storage

Store dehydrated media in its original packaging between 10°C and 30°C. Keep out of direct sunlight. Keep away from moisture.

Once reconstituted, store media between 2°C and 12°C.

Warnings and Precautions

- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Materials of Animal Origin

Campylobacter Blood-Free Selective Agar (CCDA) (Dehydrated) (CM0739B, CM0739R, CM0739T) contains peptone manufactured from porcine, bovine and microbial raw materials.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local the recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) Q 5.

Procedure

- Suspend 22.75g in 500ml of distilled water.
- Bring to the boil to dissolve completely.
- Sterilize by autoclaving at 121°C for 15 minutes.
- Cool to 50°C.
- Aseptically add the contents of 1 vial of CCDA Selective Supplement (SR0155E or SR0155H) reconstituted as directed.
- Mix well and dispense into sterile Petri dishes.

ATCC® 25922™	
<i>Escherichia coli</i> ATCC® 8739™	No growth
<i>Staphylococcus aureus</i> ATCC® 25923™	No growth

Limitations

Identifications are presumptive and should be confirmed.

Plates of C.C.D.A should not be overdried. Non target microorganisms may grow on this medium if they are resistant to the antimicrobials present. This includes cefaperazone-resistant Enterobacteriaceae and Pseudomonas species. *Campylobacter* species vary in their sensitivity to antimicrobial agents and therefore may not grow on this medium. Incubation at 42°C may prevent the growth of some *Campylobacter* species such as *C. jejuni* subspecies doylei and *C. fetus* sub species intestinalis.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of *Campylobacter* species is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the device.

The precision of Campylobacter Blood-Free Selective Agar (CCDA) (Dehydrated) (CM0739B, CM0739R, CM0739T) was demonstrated by an overall pass rate of 100% over six months of testing (July 2021 - January 2022). This shows that the performance is reproducible.

Interpretation

Campylobacter spp. will show as 0.5 - 2mm grey colonies. Identification is presumptive and should be confirmed.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 48 h @ 37° ± 2°C microaerophilic

Positive Controls	
Colony count is ≥ 50% of the control medium count	
Inoculum level 10 – 100 cfu	
<i>Campylobacter jejuni</i> ATCC® 33560™	0.5-2.0 mm grey colonies
A satisfactory result is represented by growth in accordance with the specification.	
Inoculum level 10 ⁴ – 10 ⁶ cfu	
<i>Campylobacter lari</i> ATCC® 35221™	0.5-2.0 mm grey colonies
Negative Controls	
<i>Escherichia coli</i>	No growth


Campylobacter Blood-Free Selective Agar (CCDA) (Dehydrated) (CM0739) device is tested in-house as part of the QC process since the products were launched in 1996. For target organisms, when using 10-100 cfu inoculum of *Campylobacter jejuni* ATCC® 33560 and *and 10⁴-10⁶ cfu Campylobacter lari* ATCC® 35221 and incubating the device at 37 ± 2°C for 48 hours, the user can recover organisms with colony size and morphology as listed in this document. For non-target organisms, when using 10⁴-10⁶ cfu inoculum of *Escherichia coli* ATCC® 25922, *Escherichia coli* ATCC® 8739 and *Staphylococcus aureus* ATCC® 25923 incubating the device at 37 ± 2°C for 48 hours, the user can recover organisms with colony size and morphology as listed in this document.

Bibliography

1. Public Health England. 2014 'Identification of Campylobacter species ID_23i3.1 I. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
2. Skirrow M.B. (1977) *Brit. Med. Jour.* 2, 9-11.
3. Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. (2018). The possible relationship between Campylobacter spp./Arcobacter spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
4. Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A., & Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with Campylobacter jejuni infection. *Neurology- Neuroimmunology Neuroinflammation*, 8(2).
5. Centers for Disease Control and Prevention. 2019. 'Campylobacter (Campylobacteriosis). <https://www.cdc.gov/campylobacter/technical.html>
6. Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, D. W., Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, S. D. (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant Campylobacter jejuni infections. *Frontiers in Public Health*, 9, 820.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limitation (storage temp.)
	Use by (expiration date) YYYY-MM
	Keep away from sunlight
	Do not re-use
	Consult Instructions for Use

	Contains sufficient for <n> tests
	Do not use if packaging damaged
	Manufacturer
	Authorized representative in the European Community
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of issue and modifications introduced
2.0	2023-12-15

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
- Public Health England. 2014 'Identification of Campylobacter species ID_23i3.1 I. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
- Skirrow M.B. (1977) *Brit. Med. Jour.* 2, 9-11.
- Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. (2018). The possible relationship between Campylobacter spp./Arcobacter spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
- Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A., & Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with Campylobacter jejuni infection. *Neurology- Neuroimmunology Neuroinflammation*, 8(2).
- Centers for Disease Control and Prevention. 2019. 'Campylobacter (Campylobacteriosis). <https://www.cdc.gov/campylobacter/technical.html>
- Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, D. W., Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, S. D. (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant Campylobacter jejuni infections. *Frontiers in Public Health*, 9, 820.

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
Litteratur

- Public Health England. 2014 'Identification of Campylobacter species ID_23i3.1 I. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
- Skirrow M.B. (1977) *Brit. Med. Jour.* 2, 9-11.
- Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. (2018). The possible relationship between Campylobacter spp./Arcobacter spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531-538.
- Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A., & Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with Campylobacter jejuni infection. *Neurology- Neuroimmunology Neuroinflammation*, 8(2).
- Centers for Disease Control and Prevention. 2019. 'Campylobacter (Campylobacteriosis). <https://www.cdc.gov/campylobacter/technical.html>
- Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, DW, Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, SD (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant Campylobacter jejuni infections. *Frontiers in Public Health*, 9, 820.

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Literaturverzeichnis

- Public Health England. 2014 „Identification of *Campylobacter* species“ ID_23i3.1 l. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
- Skirrow M.B. (1977) *Brit. Med. Jour.* 2, 9-11.
- Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. 2018(2018). The possible relationship between *Campylobacter* spp./*Arcobacter* spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
- Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Diones, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A. & Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with *Campylobacter jejuni* infection. *Neurology- Neuroimmunology Neuroinflammation*, 8(2).
- Centers for Disease Control and Prevention. 2019. „*Campylobacter* (Campylobacteriosis). <https://www.cdc.gov/campylobacter/technical.html>
- Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, D. W., Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, S. D. (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant *Campylobacter jejuni* infections. *Frontiers in Public Health*, 9, 820.

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- Public Health England. 2014 'Identification of *Campylobacter* species ID_23i3.1 I. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
- Skirrow M.B. (1977) *Brit. Med. Jour.* 2, 9-11.
- Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ates, I., & Deniz, K. (2018). The possible relationship between *Campylobacter* spp./*Arcobacter* spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
- Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A., & Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with *Campylobacter jejuni* infection. *Neurology- Neuroimmunology Neuroinflammation*, 8(2).
- Centers for Disease Control and Prevention. 2019. '*Campylobacter*(*Campylobacteriosis*). <https://www.cdc.gov/campylobacter/technical.html>
- Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, DW, Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, S. D. (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant *Campylobacter jejuni* infections. *Frontiers in Public Health*, 9, 820.

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- Public Health England. 2014 'Identification of Campylobacter species ID_23i3.1 I. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
- Skirrow M.B. (1977) *Brit. Med. Jour.* 2, 9-11.
- Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. (2018). The possible relationship between *Campylobacter* spp./*Arcobacter* spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
- Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A., & Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with *Campylobacter jejuni* infection. *Neurology-Neuroimmunology Neuroinflammation*, 8(2).
- Centers for Disease Control and Prevention. 2019. '*Campylobacter* (Campylobacteriosis). <https://www.cdc.gov/campylobacter/technical.html>
- Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, D. W., Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, S. D. (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant *Campylobacter jejuni* infections. *Frontiers in Public Health*, 9, 820.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Reino Unido

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Bibliographie

- Public Health England. 2014 "Identification of *Campylobacter* species ID_23i3.1 I". <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
- Skirrow M. B. (1977) *Brit. Med. Jour.* 2, 9-11.
- Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. (2018). The possible relationship between *Campylobacter* spp./*Arcobacter* spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
- Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A., & Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with *Campylobacter jejuni* infection. *Neurology- Neuroimmunology Neuroinflammation*, 8(2).
- Centers for Disease Control and Prevention. 2019. "*Campylobacter* (Campylobacteriosis)". <https://www.cdc.gov/campylobacter/technical.html>
- Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, D. W., Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, S. D. (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant *Campylobacter jejuni* infections. *Frontiers in Public Health*, 9, 820.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Royaume-Uni

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
Bibliografija

- Public Health England. 2014 'Identification of *Campylobacter* species ID_23i3.1 I. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
- Skirrow M.B. (1977) Brit. Med. Jour. 2, 9-11.
- Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. (2018). The possible relationship between *Campylobacter* spp./*Arcobacter* spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
- Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A., & Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with *Campylobacter jejuni* infection. *Neurology- Neuroimmunology Neuroinflammation*, 8(2).
- Centers for Disease Control and Prevention. 2019. '*Campylobacter* (Campylobacteriosis). <https://www.cdc.gov/campylobacter/technical.html>
- Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, D. W., Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, S. D. (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant *Campylobacter jejuni* infections. *Frontiers in Public Health*, 9, 820.

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- Public Health England. 2014 'Identification of *Campylobacter* species ID_23i3.1 l. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
- Skirrow M.B. (1977) *Brit. Med. Jour.* 2, 9-11.
- Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. (2018). The possible relationship between *Campylobacter* spp./*Arcobacter* spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
- Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A., & Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with *Campylobacter jejuni* infection. *Neurology- Neuroimmunology Neuroinflammation*, 8(2).
- Centers for Disease Control and Prevention. 2019. '*Campylobacter* (Campylobacteriosis). <https://www.cdc.gov/campylobacter/technical.html>
- Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, D.W., Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, S. D. (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant *Campylobacter jejuni* infections. *Frontiers in Public Health*, 9, 820.

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
Bibliografia

- Public Health England. 2014 'Identification of *Campylobacter* species ID_23i3.1 I. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
- Skirrow M.B. (1977) *Brit. Med. Jour.* 2, 9-11.
- Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. (2018). The possible relationship between *Campylobacter* spp./*Arcobacter* spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531-538.
- Ramos, A. P., Leonhard, SE, Halstead, SK., Cuba, m. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A. e Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with *Campylobacter jejuni* infection. *Neurology-Neuroimmunology Neuroinflammation*, 8(2).
- Centers for Disease Control and Prevention. 2019. '*Campylobacter* (Campylobacteriosis). <https://www.cdc.gov/campylobacter/technical.html>
- Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, DW, Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T. e Manning, SD (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant *Campylobacter jejuni* infections. *Frontiers in Public Health*, 9, 820.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Regno Unito

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- Public Health England. 2014 'Identification of Campylobacter species ID_23i3.1 I. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
- Skirrow M.B. (1977) Brit. Med. Jour. 2, 9-11.
- Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. (2018). The possible relationship between Campylobacter spp./Arcobacter spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
- Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A., & Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with Campylobacter jejuni infection. *Neurology- Neuroimmunology Neuroinflammation*, 8(2).
- Centers for Disease Control and Prevention. 2019. 'Campylobacter (Campylobacteriosis). <https://www.cdc.gov/campylobacter/technical.html>
- Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, D. W., Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, S. D. (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant Campylobacter jejuni infections. *Frontiers in Public Health*, 9, 820.

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Bibliografi

- Public Health England. 2014 'Identification of Campylobacter species ID_23i3.1 l.
<https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
- Skirrow M.B. (1977) *Brit. Med. Jour.* 2, 9-11.
- Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. (2018). The possible relationship between *Campylobacter* spp./*Arcobacter* spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
- Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A., & Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with *Campylobacter jejuni* infection. *Neurology- Neuroimmunology Neuroinflammation*, 8(2).
- Centers for Disease Control and Prevention. 2019. 'Campylobacter (Campylobacteriosis).
<https://www.cdc.gov/campylobacter/technical.html>
- Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, D. W., Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, S. D. (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant *Campylobacter jejuni* infections. *Frontiers in Public Health*, 9, 820.

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- Public Health England. 2014 'Identification of Campylobacter species ID_23i3.1 I. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
- Skirrow M.B. (1977) *Brit. Med. Jour.* 2, 9-11.
- Akar, M., Aydın, F., Yurci, M. A., Abay, S., Ateş, İ. i Deniz, K. (2018). The possible relationship between *Campylobacter* spp./*Arcobacter* spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
- Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A. i Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with *Campylobacter jejuni* infection. *Neurology-Neuroimmunology Neuroinflammation*, 8(2).
- Centers for Disease Control and Prevention. 2019. 'Campylobacter (Campylobacteriosis). <https://www.cdc.gov/campylobacter/technical.html>
- Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, D. W., Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T. i Manning, SD (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant *Campylobacter jejuni* infections. *Frontiers in Public Health*, 9, 820.

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Bibliografia

- Public Health England. 2014 'Identification of *Campylobacter* species ID_23i3.1 l. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
- Skirrow M.B. (1977) *Brit. Med. Jour.* 2, 9-11.
- Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. (2018). The possible relationship between *Campylobacter* spp./*Arcobacter* spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
- Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A., & Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with *Campylobacter jejuni* infection. *Neurology-Neuroimmunology Neuroinflammation*, 8(2).
- Centers for Disease Control and Prevention. 2019. '*Campylobacter* (Campylobacteriosis). <https://www.cdc.gov/campylobacter/technical.html>
- Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, D. W., Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, S. D. (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant *Campylobacter jejuni* infections. *Frontiers in Public Health*, 9, 820.

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Bibliografie

- Public Health England. 2014 'Identification of *Campylobacter* species ID_23i3.1 I. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
- Skirrow M.B. (1977) *Brit. Med. Jour.* 2, 9-11.
- Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. (2018). The possible relationship between *Campylobacter* spp./*Arcobacter* spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
- Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A., & Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with *Campylobacter jejuni* infection. *Neurology- Neuroimmunology Neuroinflammation*, 8(2).
- Centers for Disease Control and Prevention. 2019. 'Campylobacter (Campylobacteriosis). <https://www.cdc.gov/campylobacter/technical.html>
- Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, D. W., Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, S. D. (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant *Campylobacter jejuni* infections. *Frontiers in Public Health*, 9, 820.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Regatul Unit

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Bibliografia

- Public Health England. 2014 'Identification of *Campylobacter* species ID_23i3.1 I. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
- Skirrow M.B. (1977) *Brit. Med. Jour.* 2, 9-11.
- Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. (2018). The possible relationship between *Campylobacter* spp./*Arcobacter* spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
- Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A. & Gourlay, D. (2021). Prevalence of Guillain-Barré syndrome in Peru in 2019 associated with *Campylobacter jejuni*. *Neurology- Neuroimmunology Neuroinflammation*, 8(2).
- Centers for Disease Control and Prevention. 2019. '*Campylobacter* (Campylobacteriosis). <https://www.cdc.gov/campylobacter/technical.html>
- Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, D. W., Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, S. D. (2021). Epidemiological associations between *Campylobacter jejuni* resistant to tetracycline and fluoroquinolone. *Frontiers in Public Health*, 9, 820.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Storbritannien