www.thermofisher.com

Thermo
SCIENTIFIC

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers.
- Incubators.
- Quality control organisms.
- Petri dishes.
- Supplements (SR0158E).
- Antimicrobial Susceptibility Testing Discs.

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Haemophilus Test Medium

EN

REF CM0898B

Intended Use

Haemophilus Test Medium (HTM) Base (CM0898B) is a medium formulated for the susceptibility testing of *Haemophilus* species isolated from clinical samples.

The device is used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having *Haemophilus* infections.

The devices are for professional use only, are not automated, nor are they companion diagnostics.

Summary and Explanation

The genus *Haemophilus* is comprised of a wide variety of clinically relevant species capable of causing infection in humans¹. These organisms are Gram-negative, film-forming, coccobacilli bacteria with marked pleomorphism and an ideal growth temperature of 35 to 37°C¹.

Members of *Haemophilus* are typically cultured on chocolate as these media allow for easy access to X (haemin) and V (nicotinamide adenine dinucleotide) blood factors required for *Haemophilus* growth². This was first documented in 1893 by Richard Pfeiffer, who identified haemoglobin as the most important constituent of blood to culture influenza bacillus³. *H. influenzae* requires both X and V factors² however, these complex growth conditions have complicated the routine susceptibility testing due to antagonism between some essential nutrients and certain antimicrobial agents. Studies conducted by Jorgensen *et al.* (1987) led to the development of Haemophilus Test Medium (HTM), which consists of Mueller-Hinton agar supplemented with both X and V factors along with yeast extract. The major advantage of this medium is its transparency, which allows for zones of inhibition to be easily read through the bottom of the plate. Low levels of antimicrobial antagonists also permit testing of trimethoprim/sulphamethoxazole with greater confidence.

Principle of Method

The Haemophilus Test Medium (HTM) Base (CM0898B) and HTM Supplement (SR0158E) devices facilitate the susceptibility testing of *Haemophilus* species isolated from clinical samples. The medium consists of Mueller-Hinton Agar with added yeast extract and supplemented with X (haemin) and V (nicotinamide adenine dinucleotide; NAD) factors required for the growth of *Haemophilus* species. Yeast extract provides an additional source of growth factors. Beef extract and casein hydrolysate supply amino acids, nitrogenous substances, vitamins, and minerals necessary for growth. Agar is the solidifying agent.

Typical Formula

	grams per litre
Mueller Hinton Agar	38.0
Yeast Extract	5.0

Materials Provided

CM0898B: 500g of Haemophilus Test Medium.

500g of dehydrated Haemophilus Test Medium yields approximately 11.6L after reconstitution.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 12 and Q 5.

Procedure

Suspend 21.5g of Haemophilus Test Medium Base in 500ml of distilled water. Bring to the boil to dissolve. Sterilise by autoclaving at 121° C for 15 minutes. Cool to 50°C and aseptically add the contents of 1 vial of Haemophilus Test Medium Supplement SR0158E, as directed. Mix well and pour into sterile Petri dishes.

- Prepare the inoculum in Mueller-Hinton Broth CM0405 or 0.9% saline, to match the turbidity of 0.5 McFarland standard.
- Using a swab saturated with the above inoculum suspension, inoculate the surface of a Haemophilus Test Medium Agar plate to give confluent growth.
- Apply the antimicrobial discs on to the surface of the Haemophilus Test Medium plate.

Incubate the plates at 35°C in 5-7% carbon dioxide for 16-18 hours and measure the zones of inhibition.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Controls	
Inoculum level: 10 – 100 cfu (colony forming units). Colony count is ≥ 70% of the control medium count. Incubation Conditions: 5% CO ₂ at 35°C, 16 – 18 hours incubation	
<i>Haemophilus influenzae</i> ATCC® 9334™	0.5 – 2 mm Straw colonies.
<i>Haemophilus influenzae</i> ATCC® 35056™	0.5 – 2 mm Straw colonies.
<i>Haemophilus influenzae</i> ATCC® 10211™	0.5 – 2 mm Straw colonies.
Incubation Conditions: 5% CO ₂ at 35°C, 20 - 24 hours incubation	
<i>Haemophilus influenzae</i> ATCC® 49247™	AMP10 – Ampicillin, 13 – 21 mm
	CTX30 – Cefotaxime, 31 – 39 mm
	C30 – Chloramphenicol, 31 – 40 mm
	TE30 – Tetracycline, 14 – 22 mm
	STX25 – Trimethoprim/Sulphamethox, 24 – 32 mm
<i>Haemophilus influenzae</i> ATCC® 49766™	CFM5 – Cefixime, 25 – 33 mm
	CEC30 – Cefaclor, 25 – 31 mm
<i>Haemophilus influenzae</i> ATCC® 49766™	CXM30 – Cefuroxime, 28 – 36 mm

Limitations

To obtain accurate and reproducible results it is necessary to ensure that the method is strictly followed (refer to the CLSI documents for further information). Cultures of *Haemophilus* spp. must be pure and identity confirmed before carrying out susceptibility testing.

Performance Characteristics

Accuracy has been demonstrated through review of internal QC data. Testing of well-characterised isolates in the QC processes is performed as part of the manufacture of each batch of the devices.

Accuracy has been demonstrated through review of the QC data. Correct detection of *Haemophilus* strains is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the devices. The precision of Haemophilus Test Medium (HTM) Base (CM0898B) was demonstrated by an overall pass rate of 100% obtained for the product over 1 year and 2 months of testing (21.05.2021 – 24.07.2022; 10 batches). This shows that the performance is reproducible. Haemophilus Test Medium (HTM) Base (CM0898B) has been tested in-house as part of the QC process since the product was launched in 1996. For target organisms, when using 10-100 cfu inoculum of *Haemophilus influenzae* (ATCC® 9334™, ATCC® 35056™, and ATCC® 10211™) and incubating the devices at 35 ± 2°C for 16-18 hours, the user can recover organisms with colony morphology as listed in this document. For *Haemophilus influenzae* (ATCC® 49247™ and ATCC® 49766™), when incubating the devices with antimicrobial susceptibility discs at 35 ± 2°C for 16-18 hours, the user can recover organisms with zone diameters as listed in this document

Bibliography

1. Baron, S. 1996. Medical Microbiology. 4th edition. Galveston, Texas, USA.
2. Public Health England. 2021. 'Identification of Haemophilus species and the HACEK group of organisms.' ID 12 (4). <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>.
3. Pfeiffer, R. 1893. Die aetologie der influenza. Zeitschrift Für Hygiene Und Infektionskrankheiten, 13(1), 357–386.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Consult instructions for use or consult electronic instructions for use

	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
CLSI is a trademark of the Clinical Laboratory and Standards Institute.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke,
Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of modifications introduced
3.0	2023-06-07

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