Campylobacter Karmali Agar Base

EN

REF CM0935B

Intended Use

Campylobacter Karmali Agar Base (CM0935B), when supplemented with SR0167E and SR0205E devices, is a selective medium for the isolation of *Campylobacter* species from faecal samples. CM0935B, with SR0167E and SR0205E devices added, is intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having *Campylobacter* infections. Can also be used for testing food samples.

The device is for professional use only, is not automated and nor is it a companion diagnostic.

Summary and Explanation

Campylobacter infections are among the most common bacterial infections, producing both gastrointestinal and systemic illnesses. *Campylobacter jejuni* is the most commonly reported bacterial cause of foodborne infection in the United States. *Campylobacter* are small, non-spore forming, helical, curved gram-negative rods that require a microaerobic environment and cannot tolerate drying¹.

Principle of Method

Isolation of *Campylobacter* species is achieved through the inclusion of a combination of selective and non-selective agents. Vancomycin suppresses the growth of gram-positive microorganisms, whilst cefoperazone inhibits the growth of a wide range of microorganisms including *Enterobacteriaceae* and *Pseudomonas* species. Karmali Campylobacter Selective Supplement (SR0167E) uses cycloheximide as an antifungal agent, whilst Modified Karmali Selective Supplement (SR0205E) uses amphotericin B. Sodium pyruvate is added to improve the recovery of stressed cells, whilst charcoal and haemin allow the adsorption of toxic compounds produced by growing microorganisms.

Typical Formula

	<u>grams per litre</u>
Columbia Agar Base	39.0
Activated charcoal	4.0
Haematin	0.032

Materials Provided

500g of dehydrated Karmali Agar Base powder that yields approximately 23.25L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Selective supplements (SR0167E)
- Petri dish

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 8°C.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.

- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 1, ID 23 and S 7.

Procedure

Suspend 21.5 g in 500ml of distilled water and bring to the boil to dissolve. Sterilize by autoclaving at 121°C for 15 minutes. Cool to 50°C. Aseptically add 1 vial of *Campylobacter* Selective Supplement (Karmali) (SR0167E) reconstituted as directed. Mix well and pour into sterile Petri dishes.

Interpretation

Once the medium is reconstituted:

Grey colonies indicate *Campylobacter jejuni* and *Campylobacter coli*.

Feather edged grey colonies indicate *Candida albicans*

Negative strains are inhibited.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 48h @ 42°C microaerophilic

Positive Controls	
Medium is challenged with 10-100 colony forming units Colony count is $\geq$ 50% of the control medium count	
<i>Campylobacter jejuni</i> ATCC® 29428	0.5-2mm grey colonies
<i>Campylobacter jejuni</i> ATCC® 33560	0.5-2mm grey colonies
<i>Campylobacter jejuni</i> ATCC® 33291	0.5-2mm grey colonies
<i>Campylobacter coli</i> ATCC® 43478	0.5-2mm grey colonies
For <i>Candida albicans</i> ATCC® 10231, a satisfactory result is represented by recovery equal to or greater than 40% of the control medium.	
<i>Candida albicans</i> ATCC® 10231	Pinpoint-0.5mm feather-edged, grey colonies
Medium is challenged with 1E+04 to 1E+06 colony forming units	
<i>Campylobacter lari</i> ATCC® 35221	0.5-2mm grey colonies
Negative Controls	

<i>Pseudomonas aeruginosa</i> ATCC® 27853	No growth
<i>Escherichia coli</i> ATCC® 8739	No growth
<i>Staphylococcus aureus</i> ATCC® 25923	No growth
Positive Control Testing performed in accordance with current CLSI M22 A Tested with the addition of Campylobacter Selective Supplement (Karmali) SR0167 Medium is challenged with 10-100 colony forming units	
<i>Campylobacter jejuni</i> ATCC® 33291	0.5-2mm grey/brown colonies
Negative Control	
<i>Escherichia coli</i> ATCC® 25922	No growth

Limitations

Identifications are presumptive and colonies should be confirmed using appropriate methods.

Performance Characteristics


Accuracy has been demonstrated through review of the QC data. Correct detection of *Campylobacter* species is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the devices. The precision of Karmali Agar Base (CM0935B) device was demonstrated by an overall pass rate of 97.5% obtained for the product over several years of testing (21.FEB.2017 – 01.AUG.2022; 30 batches). This shows that the performance is reproducible.


Karmali Agar Base (CM0935B) device is tested in-house as part of the QC process since at least 1995. For target organisms, when using 10-100 cfu inoculum of *Campylobacter jejuni* ATCC® 29428™, *Campylobacter jejuni* ATCC® 33560™, *Campylobacter jejuni* ATCC® 33291™, *Campylobacter coli* ATCC® 43478™, or *Candida albicans* ATCC® 10231™ and incubating the device at 42°C for 48 hours in microaerophilic conditions, or when using 10⁴ – 10⁶ cfu inoculum of *Campylobacter lari* ATCC® 352211™ the user can recover organisms with colony size and morphology as listed in this document. For non-target organisms, When using 10⁴ – 10⁶ cfu inoculum of *Pseudomonas aeruginosa* ATCC® 27853™ the user can expect no growth.

Bibliography

1. Karmali M.A., Simor A.E., Roscoe M., Fleming P.C, Smith S.S. and Lane J. (1986) J.Clin.Micro. 23. 456-459

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use

	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of issue and modifications introduced
2.0	2023-07-10 Minor formatting changes

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©2022 Thermo Fisher Scientific Inc. Všetky práva vyhradené.
Všetky ostatné ochranné známky sú vlastníctvom Thermo Fisher Scientific Inc. a jej dcérskych spoločností.
Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Spojené kráľovstvo