	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Aztreonam Discs 30 µg (ATM30)

REF CT0264B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Aztreonam Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent aztreonam. The discs are labelled on both sides with details of the antimicrobial (ATM) and amount present (µg): ATM30 (30µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Aztreonam AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which aztreonam has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling, and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100,^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- N/A

Gram-Negative

- *Enterobacterales*
- *Pseudomonas aeruginosa*
- *Haemophilus influenzae* and *Haemophilus parainfluenzae*

Species with published breakpoints according to current CLSIM100 literature:

Gram-Positive

- N/A

Gram-Negative

- *Enterobacterales*
- *Pseudomonas aeruginosa*
- *Haemophilus influenzae* and *Haemophilus parainfluenzae*

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- N/A

Gram-Negative

- *Enterobacterales*
- *Pseudomonas* spp.
- *Aeromonas* spp.

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Aztreonam AST Discs may be used in the semi-quantitative agar diffusion test method for *in vitro* susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology please refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Aztreonam AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Aztreonam AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Aztreonam AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Component Description	Material Description
Cartridge with spring, cap and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Pale to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Reagent	Description of Function
Aztreonam	White crystalline powder. Aztreonam exerts its mechanism of action through the binding of the penicillin binding protein 3, subsequently leading to the inhibition of the cell wall synthesis.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Aztreonam AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Aztreonam AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Results Summary for Tested QC Isolates in Accordance with CLSI Methodology.

Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3164330	<i>Escherichia coli</i>	25922	MHA	28	36	32	32	32	32	0 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	23	29	26	26	26	26	0 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	30	38	34	31	31	31	3 ± 0(CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV) (%)
							1	2	3	
3179506	<i>Escherichia coli</i>	25922	MHA	28	36	32	31	31	31	1 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	23	29	26	25	25	25	1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	30	38	34	30	30	30	4 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV) (%)
3195892	<i>Escherichia coli</i>	25922	MHA	28	36	32	28	28	28	4 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	23	29	26	24	24	24	2 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	30	38	34	30	30	31	3.66 ± 0.577 (CV=0.15)

Table 4. Results Summary for Tested QC Isolates in Accordance with EUCAST Methodology.

Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3164330	<i>Escherichia coli</i>	25922	MHA	28	36	32	32	32	32	0 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	23	29	26	26	26	26	0 ± 0 (CV=0)
	<i>Klebsiella pneumoniae</i>	700603	MHA	9	17	13	11	11	11	2 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV) (%)
							1	2	3	
3179506	<i>Escherichia coli</i>	25922	MHA	28	36	32	31	31	31	1 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	23	29	26	25	25	25	1 ± 0 (CV=0)
	<i>Klebsiella pneumoniae</i>	700603	MHA	9	17	13	9	9	9	4 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV) (%)
							1	2	3	
3195892	<i>Escherichia coli</i>	25922	MHA	28	36	32	28	28	28	4 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	23	29	26	24	24	24	2 ± 0 (CV=0)
	<i>Klebsiella pneumoniae</i>	700603	MHA	9	17	13	9	9	9	4 ± 0 (CV=0)

Clinical Performance Characteristics

A study was performed to identify the antimicrobial susceptibility patterns of *K. pneumoniae* and *E. coli* isolates. The study demonstrated that 64% of the positive controls (ESBL producers) were susceptible to aztreonam.⁴ Antibiotics including cefotetan, cefoxitin, ceftriaxone and cefepime demonstrated 100% susceptibility. In contrast, the susceptibility rate of aztreonam among the ESBL producing isolates was higher than cefsulodin, ceftazidime, and cefpodoxime. Oxoid Aztreonam (30µg) ATM30 Antimicrobial Susceptibility Test Discs (CT0264B) were found, through the Jarlier test, to show 100% specificity to all isolates.

Another study collected clinical isolates from 11 hospitals in Shanghai between April 2000 and March 2001 resulting in a total of 3,905 non-fermentative isolates.⁵ Of these, 1,790 were *P. aeruginosa* and 169 were other *Pseudomonas* spp. isolates. The isolates were identified through routine microbiological methods and susceptibility testing was performed through the disc diffusion method for a variety of antibiotics, including the Oxoid Aztreonam AST Discs. The study did not report any inconsistent or out of specification results when using the discs.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References


1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing

- c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. F. M. MacKenzie, C. A. Miller and I. M. Gould. 2002. Comparison of screening methods for TEM- and SHV-derived extended-spectrum beta-lactamase detection. *Clinical Microbiology and Infection* 8(11): 715-724
5. W. Fu, Z. Demei, W. Shi, H. Fupin and Z. Yingyuan. 2003. The susceptibility of non-fermentative Gram-negative bacilli to cefoperazone and sulbactam compared with other antibacterial agents. *International Journal of Antimicrobial Agents* 22(4): 444-448

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark


The UKCA logo, consisting of the letters "UK" stacked above "CA" in a bold, black, sans-serif font.	UK Conformity Mark
---	--------------------

Revision Information

Version	Date of issue and modifications introduced
1.0	2023-01-13

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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Litteratur

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).

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Viited

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (-)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing

3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. F. M. MacKenzie, C. A. Miller and I. M. Gould. 2002. Comparison of screening methods for TEM- and SHV-derived extended-spectrum beta-lactamase detection. *Clinical Microbiology and Infection* 8(11): 715-724
5. W. Fu, Z. Demei, W. Shi, H. Fupin and Z. Yingyuan. 2003. The susceptibility of non-fermentative Gram-negative bacilli to cefoperazone and sulbactam compared with other antibacterial agents. *International Journal of Antimicrobial Agents* 22(4): 444-448

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Viitteet

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (-)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).

4. F. M. MacKenzie, C. A. Miller and I. M. Gould. 2002. Comparison of screening methods for TEM- and SHV-derived extended-spectrum beta-lactamase detection. *Clinical Microbiology and Infection* 8(11): 715-724
5. W. Fu, Z. Demei, W. Shi, H. Fupin and Z. Yingyuan. 2003. The susceptibility of non-fermentative Gram-negative bacilli to cefoperazone and sulbactam compared with other antibacterial agents. *International Journal of Antimicrobial Agents* 22(4): 444-448

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Références

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (-)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.

2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. F. M. MacKenzie, C. A. Miller and I. M. Gould. 2002. Comparison of screening methods for TEM- and SHV-derived extended-spectrum beta-lactamase detection. *Clinical Microbiology and Infection* 8(11): 715-724
5. W. Fu, Z. Demei, W. Shi, H. Fupin and Z. Yingyuan. 2003. The susceptibility of non-fermentative Gram-negative bacilli to cefoperazone and sulbactam compared with other antibacterial agents. *International Journal of Antimicrobial Agents* 22(4): 444-448

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Referenzen

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (-)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method

- b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. F. M. MacKenzie, C. A. Miller and I. M. Gould. 2002. Comparison of screening methods for TEM- and SHV-derived extended-spectrum beta-lactamase detection. *Clinical Microbiology and Infection* 8(11): 715-724
5. W. Fu, Z. Demei, W. Shi, H. Fupin and Z. Yingyuan. 2003. The susceptibility of non-fermentative Gram-negative bacilli to cefoperazone and sulbactam compared with other antibacterial agents. *International Journal of Antimicrobial Agents* 22(4): 444-448

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Nuorodos

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- F. M. MacKenzie, C. A. Miller and I. M. Gould. 2002. Comparison of screening methods for TEM- and SHV-derived extended-spectrum beta-lactamase detection. *Clinical Microbiology and Infection* 8(11): 715-724
- W. Fu, Z. Demei, W. Shi, H. Fupin and Z. Yingyuan. 2003. The susceptibility of non-fermentative Gram-negative bacilli to cefoperazone and sulbactam compared with other antibacterial agents. *International Journal of Antimicrobial Agents* 22(4): 444-448

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1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (-)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. F. M. MacKenzie, C. A. Miller and I. M. Gould. 2002. Comparison of screening methods for TEM- and SHV-derived extended-spectrum beta-lactamase detection. *Clinical Microbiology and Infection* 8(11): 715-724
5. W. Fu, Z. Demei, W. Shi, H. Fupin and Z. Yingyuan. 2003. The susceptibility of non-fermentative Gram-negative bacilli to cefoperazone and sulbactam compared with other antibacterial agents. *International Journal of Antimicrobial Agents* 22(4): 444-448

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