	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Azithromycin Discs 15 µg (AZM15)

REF CT0906B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Azithromycin Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent (Product name). The discs are labelled on both sides with details of the antimicrobial (AZM) and amount present (µg): AZM (15µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Azithromycin AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Azithromycin AST Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which azithromycin has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling, and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100,^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- *Staphylococcus* spp.
- *Streptococcus pneumoniae*

Gram-Negative

- *Haemophilus influenzae*
- *Haemophilus parainfluenza*

Species with published breakpoints according to current CLSI M100 literature:

Gram-Positive

- *Streptococcus* spp. β-Hemolytic Group
- *Streptococcus* spp. Viridans Group
- *Streptococcus pneumoniae*
- *Staphylococcus* spp.

Gram-Negative

- *Enterobacterales*
- *Acinetobacter* spp.
- *Haemophilus influenzae* and *Haemophilus parainfluenzae*
- *Neisseria gonorrhoeae*
- *Neisseria meningitidis*

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- N/A

Gram-Negative

- *Salmonella Typhi*
- *Shigella* spp.

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Azithromycin AST Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology please refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Azithromycin AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Azithromycin AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Azithromycin AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with CT0906B	
Component Description	Material Description
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Pale beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge, 5 cartridges per pack.

Table 2. Description of Azithromycin AST Disc Reagents	
Reagent	Description of Function
Azithromycin	White-to-yellow powder that is water soluble. Azithromycin is a broad-spectrum macrolide antibiotic which inhibits bacterial growth through the inhibition of protein synthesis and translation. Azithromycin achieves this through binding to the 23S rRNA of the bacterial 50S ribosomal subunit which prevents the assembly of the 50S ribosomal subunit.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Azithromycin AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Azithromycin AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Raw data of zone sizes, in accordance with CLSI methodology, taken from the 3 most recent batches of Azithromycin (15 µg) AZM15 Antimicrobial Susceptibility Test Discs (CT0906B).

CLSI Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
330089	<i>Staphylococcus aureus</i>	25923	MHA	21	26	23.5	22	22	22	-1.5 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	13	21	17	21	21	21	4 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA+5 %SB	19	25	22	23	23	23	1 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3259907	<i>Staphylococcus aureus</i>	25923	MHA	21	26	23.5	24	24	24	0.5 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	13	21	17	21	21	21	4 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA+5 %SB	19	25	22	25	25	25	3 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3164397	<i>Staphylococcus aureus</i>	25923	MHA	21	26	23.5	25	25	25	1.5 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	13	21	17	20	20	20	3 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA+5 %SB	19	25	22	22	22	22	0 ± 0 (CV=0)

Clinical Performance Characteristics

In 2020, a study was conducted to determine the reliability and performance of disc diffusion against agar dilution for susceptibility testing of *N. gonorrhoeae* isolates to azithromycin.⁴ The study used 189 *N. gonorrhoeae* clinical isolates. The disc diffusion method was performed (CLSI MO2-A13) with Oxoid Azithromycin AZM15 AST Discs and a competitor's discs. Reproducibility was tested through use of *N. gonorrhoeae* ATCC® 49226 and nine World Health Organisation reference strains. Each strain was tested three times each day for three days.

The disc diffusion method was compared to the agar dilution method using a linear regression plot. This illustrated an agreement and kappa index of 96.8% and 0.933 ($P < 0.001$) for Oxoid Azithromycin AZM15 AST Discs compared to a kappa index of 95.2% and 0.899 ($P < 0.001$) for the competitor's discs. Very major discrepancies (VMDs) and major discrepancies (MDs) were identified for both manufacturer's discs. The VMD and MD values for Oxoid, 4.0% and 2.6% respectively, and for the competitor's discs, 9.4% and 1.7%, demonstrate the ability of Oxoid discs to provide more reliable results in accordance with the CLSI.

Mal et al. conducted a study to assess the breakpoints for antimicrobial agents as determined by the BSAC and calibrated dichotomous sensitivity (CDS) for isolates of *N. gonorrhoeae* strains to determine resistance.⁵ The disc diffusion method was performed with Oxoid Azithromycin (15 µg) AZM15 AST Discs. Breakpoints were performed to the guidelines of the CLSI, BSAC and CDS. Etests demonstrated that the *N. gonorrhoeae* isolates showed a 99% susceptibility to azithromycin. The BSAC and CDS disc diffusion method for Oxoid Azithromycin AST Discs demonstrated a 99% susceptibility rate and a 1% resistance, agreeing with the Etest results completely. This demonstrated the Oxoid Azithromycin AST Discs performed to the ideal specification.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Gianecini, Ricardo, Lucia Irazu, Marcelo Rodríguez, Paula Cristaldo, Claudia Oviedo, Marisa Turco, Veronica Rodrigo, et al. 2020. "Testing of In Vitro Susceptibility of Neisseria Gonorrhoeae to Azithromycin: Comparison of Disk Diffusion and Reference Agar Dilution Methods." *Journal of Clinical Microbiology* 58 (11). <https://doi.org/10.1128/JCM.01398-20>.
5. Mal, Pushpa Bhawan, Kauser Jabeen, Joveria Farooqi, Magnus Unemo, and Erum Khan. 2016. "Antimicrobial Susceptibility Testing of Neisseria Gonorrhoeae Isolates in Pakistan by Etest Compared to Calibrated Dichotomous Sensitivity and Clinical Laboratory Standards Institute Disc Diffusion Techniques." *BMC Microbiology* 16 (1). <https://doi.org/10.1186/s12866-016-0707-6>.

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date

	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision Information

Version	Date of issue and modifications introduced
1.0	2023-01-12.

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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Referencije

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Gianecini, Ricardo, Lucia Irazu, Marcelo Rodríguez, Paula Cristaldo, Claudia Oviedo, Marisa Turco, Veronica Rodrigo, et al. 2020. "Testing of In Vitro Susceptibility of *Neisseria Gonorrhoeae* to Azithromycin: Comparison of Disk Diffusion and Reference Agar Dilution Methods." *Journal of Clinical Microbiology* 58 (11). <https://doi.org/10.1128/JCM.01398-20>.
- Mal, Pushpa Bhawan, Kauser Jabeen, Joveria Farooqi, Magnus Unemo, and Erum Khan. 2016. "Antimicrobial Susceptibility Testing of *Neisseria Gonorrhoeae* Isolates in Pakistan by Etest Compared to Calibrated Dichotomous Sensitivity and Clinical Laboratory Standards Institute Disc Diffusion Techniques." *BMC Microbiology* 16 (1). <https://doi.org/10.1186/s12866-016-0707-6>.

© 2023. Thermo Fisher Scientific Inc. Sva prava pridržana. ATCC® je zaštitni znak društva ATCC. Svi ostali zaštitni znakovi vlasništvo su društva Thermo Fisher Scientific Inc. i njegovih podružnica. Ove informacije nisu namijenjene za poticanje na uporabu tih proizvoda na način kojim se mogu povrijediti prava intelektualnog vlasništva drugih društava.

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Reference

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Gianecini, Ricardo, Lucia Irazu, Marcelo Rodríguez, Paula Cristaldo, Claudia Oviedo, Marisa Turco, Veronica Rodrigo, et al. 2020. "Testing of In Vitro Susceptibility of Neisseria Gonorrhoeae to Azithromycin: Comparison of Disk Diffusion and Reference Agar Dilution Methods." *Journal of Clinical Microbiology* 58 (11). <https://doi.org/10.1128/JCM.01398-20>.
- Mal, Pushpa Bhawan, Kauser Jabeen, Joveria Farooqi, Magnus Unemo, and Erum Khan. 2016. "Antimicrobial Susceptibility Testing of Neisseria Gonorrhoeae Isolates in Pakistan by Etest Compared to Calibrated Dichotomous Sensitivity and Clinical Laboratory Standards Institute Disc Diffusion Techniques." *BMC Microbiology* 16 (1). <https://doi.org/10.1186/s12866-016-0707-6>.

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1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Gianecini, Ricardo, Lucia Irazu, Marcelo Rodríguez, Paula Cristaldo, Claudia Oviedo, Marisa Turco, Veronica Rodrigo, et al. 2020. "Testing of In Vitro Susceptibility of *Neisseria Gonorrhoeae* to Azithromycin: Comparison of Disk Diffusion and Reference Agar Dilution Methods." *Journal of Clinical Microbiology* 58 (11). <https://doi.org/10.1128/JCM.01398-20>.
5. Mal, Pushpa Bhawan, Kauser Jabeen, Joveria Farooqi, Magnus Unemo, and Erum Khan. 2016. "Antimicrobial Susceptibility Testing of *Neisseria Gonorrhoeae* Isolates in Pakistan by Etest Compared to Calibrated Dichotomous Sensitivity and Clinical Laboratory Standards Institute Disc Diffusion Techniques." *BMC Microbiology* 16 (1). <https://doi.org/10.1186/s12866-016-0707-6>.

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Referências

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Gianecini, Ricardo, Lucia Irazu, Marcelo Rodríguez, Paula Cristaldo, Claudia Oviedo, Marisa Turco, Veronica Rodrigo, et al. 2020. "Testing of In Vitro Susceptibility of *Neisseria Gonorrhoeae* to Azithromycin: Comparison of Disk Diffusion and Reference Agar Dilution Methods." *Journal of Clinical Microbiology* 58 (11). <https://doi.org/10.1128/JCM.01398-20>.
5. Mal, Pushpa Bhawan, Kausar Jabeen, Joveria Farooqi, Magnus Unemo, and Erum Khan. 2016. "Antimicrobial Susceptibility Testing of *Neisseria Gonorrhoeae* Isolates in Pakistan by Etest Compared to Calibrated Dichotomous Sensitivity and Clinical Laboratory Standards Institute Disc Diffusion Techniques." *BMC Microbiology* 16 (1). <https://doi.org/10.1186/s12866-016-0707-6>.

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