EN

Oxoid™ Ertapenem Discs 10 µg (ETP10)

REFCT1761B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the Oxoid AST generic instructions for use supplied with the product and available online.

Ertapenem Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent ertapenem. The discs are labelled on both sides with details of the antimicrobial (ETP) and amount present (µg): ETP10 (10µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Ertapenem AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which ertapenem has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling, and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100,^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Negative

- Enterobacteriaceae
- *Haemophilus influenzae*
- *Haemophilus parainfluenza*

Gram-Positive

- *Staphylococcus* spp.

Species with published breakpoints according to current CLSI M100 literature:

Gram-Negative

- Enterobacterales
- *Haemophilus influenzae*
- *Haemophilus parainfluenza*

Gram-Positive

- N/A

Species with published breakpoints according to current EUCAST breakpoint tables:

Gram-Negative

- Enterobacterales
- *Haemophilus influenzae*

- *Moraxella catarrhalis*

Gram-Positive

- N/A

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Ertapenem AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Ertapenem AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Ertapenem AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Ertapenem AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with CT1761B	
Component Description	Material Description
Cartridge with spring, cap and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Pale beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Ertapenem AST Disc Reagent	
Reagent	Description of Function
Ertapenem	Antimicrobial agent. White/off white powder. Inhibition of cell wall synthesis mediated through binding to penicillin binding proteins.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Ertapenem AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Ertapenem AST Discs have a shelf-life of 12 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Raw data of zone sizes, in accordance with CLSI methodology, taken from 3 batches of Ertapenem (10 µg) ETP10 Antimicrobial Susceptibility Test Discs (CT1761B).

CLSI Methodology										
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3245377	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	24	31	27.5	27	27	27	-0.5 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC® 25922	MHA	29	36	32.5	33	33	33	0.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	13	21	17	17	17	17	0 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49247	HTM	20	28	24	28	28	28	4 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49766	HTM	27	33	30	28	28	28	-2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA+ 5%SB	28	35	31.5	35	35	35	3.5 ± 0 (CV=0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3245376	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	24	31	27.5	30	30	30	2.5 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC® 25922	MHA	29	36	32.5	33	33	33	0.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	13	21	17	20	20	20	3 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49247	HTM	20	28	24	24	24	24	0 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49766	HTM	27	33	30	31	31	31	1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA+ 5%SB	28	35	31.5	32	32	32	0.5 ± 0 (CV=0)

Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3235301	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	24	31	27.5	27	27	27	-0.5 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC® 25922	MHA	29	36	32.5	33	33	33	0.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	13	21	17	20	20	20	3 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49247	HTM	20	28	24	25	25	25	1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49766	HTM	27	33	30	32	32	32	2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA+ 5%SB	28	35	31.5	33	33	33	1.5 ± 0 (CV=0)

Table 4. Raw data of zone sizes, in accordance with EUCAST methodology, taken from 3 batches of Ertapenem (10 µg) ETP10 Antimicrobial Susceptibility Test Discs (CT1761B).

EUCAST Methodology										
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3245377	<i>Escherichia coli</i>	ATCC® 25922	MHA	29	36	32.5	33	33	33	0.5 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49766	MHA-F	27	33	30	27	27	27	-3 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA-F	28	34	31	34	34	34	3 ± 0 (CV=0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3245376	<i>Escherichia coli</i>	ATCC® 25922	MHA	29	36	32.5	33	33	33	0.5 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49766	MHA-F	27	33	30	28	28	28	-2 ± 0 (CV=0)

	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA-F	28	34	31	31	31	31	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3235301	<i>Escherichia coli</i>	ATCC® 25922	MHA	29	36	32.5	33	33	33	0.5 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49766	MHA-F	27	33	30	32	32	32	2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA-F	28	34	31	32	32	32	1 ± 0 (CV=0)

Clinical Performance Characteristics

A study by Almeida et al. 2012, aimed to evaluate the positive predictive value and the specificity of the ertapenem disc to predict the presence of KPC-producing Gram-negative bacilli (GNB). During the study period, January 2010 to June 2011, a total of 21839 cultures from different biological materials were collected and processed at Hospital Municipal Dr Moisés Deautch, in M'Boi Mirim, SP, Brazil. Of these 21839, 3010 had positive culture, with GNB isolated in 708 (23.52%) of the positive cultures. These GNB were tested for sensitivity to three antibiotic discs including Ertapenem (10 µg) ETP10 Antimicrobial Susceptibility Test Discs (CT1761B) among the other discs (Oxoid), following the Technical Note published by ANVISA in 2010 guidelines. Results showed 111 of the 708 isolates (15.7%) were considered resistant or intermediate to ertapenem (Almeida et al. 2012). In conclusion, this article shows that Ertapenem Oxoid discs are effective in the clinical analysis of Gram-negative clinical isolates, given the low level of resistance and no out of specification or inconsistent results were reported.

A study by Vading et al. 2011, aimed to compare CLSI and EUCAST MIC and disc diffusion carbapenem breakpoints for the detection of carbapenemase-producing *Klebsiella pneumoniae*. During the study period a total of 51 clinical isolates of carbapenem-non-susceptible *K. pneumoniae* with known KPC production (n= 31) or VIM production (n= 20) were collected from microbiological laboratories in Sweden, Greece, the USA and Norway; isolates had earlier been genotypically characterised. Antimicrobial susceptibility testing was carried using the disc diffusion method in accordance with CLSI and EUCAST methodologies. Isolates were tested for sensitivity to three antibiotic discs including Ertapenem (10 µg) ETP10 Antimicrobial Susceptibility Test Discs (CT1761B) among the other discs (Oxoid). Results using the CLSI breakpoints showed three KPC-producing isolates (n=31) were susceptible to ertapenem (9.7%). For the EUCAST breakpoints, no isolates were categorised as susceptible to ertapenem. The study concluded that when using the EUCAST disk diffusion breakpoints, all carbapenemase producers were discovered with the ertapenem discs, whereas the CLSI breakpoints saw three susceptible isolates to ertapenem. Additionally, it was reported that ertapenem and meropenem performed better in the separation of KPC/VIM producers and wild-type strains than imipenem. Furthermore, it was concluded that disk diffusion breakpoints identified a much higher proportion of the carbapenemase producers than the MIC breakpoints (Vading et al. 2011). In conclusion, this article demonstrates that Oxoid Ertapenem discs are effective in the clinical analysis of *K. pneumoniae*, given that it has detected resistance and has been used to identify specific resistance genotypes (KPC/VIM production) and no out of specification or inconsistent results were reported.

A study by Martinez, Garzón, and Mattar 2012, aimed to describe the presence of CTX-M-1 phylogenetic subgroup of extended-spectrum β-lactamases (ESBL), associated with TEM and SHV genes, and the gene encoding cephalosporinase, CMY-2 in *Escherichia coli* and *Klebsiella pneumoniae* isolates from community-acquired urinary tract infections (UTIs). During the study period of two months (February to March 2011), 123 Enterobacteriaceae isolates were collected from patients with culture proven UTI in four clinical microbiology laboratories in Valledupar, Colombia. Of these 123, 102 were *E. coli* and 21 were *K. pneumoniae*, isolates were identified using conventional methods. Antibiotic susceptibility testing was performed using the disk diffusion method and was interpreted in accordance with CLSI guidelines. Isolates were tested for sensitivity against a panel of antibiotic discs (Oxoid) including Ertapenem (10 µg) ETP10 Antimicrobial Susceptibility Test Discs (CT1761B). Quality control strains of *E. coli* ATCC 25922 and *K. pneumoniae* ATCC 700603 were used. Results showed all *E. coli* (100%) and all *K. pneumoniae* (100%) isolates were susceptible to ertapenem.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing

3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Almeida, Lais Pinto de, Fabiana Puerro de Carvalho, Alexandre Gímenes Marques, Andrea dos Santos Pereira, Renata Puzzo Bortoleto, and Marinês Dalla Valle Martino. 2012. "Ertapenem Disk Performance to Predict Klebsiella Pneumoniae Carbapenemase Produced by Gram-Negative Bacilli Isolated in a São Paulo City Public Hospital." *Einstein (São Paulo, Brazil)* 10 (4): 439–41. <https://doi.org/10.1590/s1679-45082012000400008>.
5. Vading, M., Samuelsen, B. Haldorsen, A. S. Sundsfjord, and C. G. Giske. 2011. "Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing Klebsiella Pneumoniae with the EUCAST and CLSI Breakpoint Systems." *Clinical Microbiology and Infection* 17 (5): 668–74. <https://doi.org/10.1111/j.1469-0691.2010.03299.x>.
6. Martinez, Pedro, Denisse Garzón, and Salim Mattar. 2012. "CTX-M-Producing Escherichia Coli and Klebsiella Pneumoniae Isolated from Community-Acquired Urinary Tract Infections in Valledupar, Colombia." *Brazilian Journal of Infectious Diseases* 16 (5): 420–25. <https://doi.org/10.1016/j.bjid.2012.05.001>.

© 2025 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of Symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limitation (storage temp.)
	Lot / Batch Number
	Catalog Number
	Do not re-use
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Expiry Date
	Period-after-opening – the number of days the product can be used after opening
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Do not use if package is damaged
	Authorized representative in the European Community/ European Union
	Unique device identifier
	European Conformity Mark

	UK Conformity Mark
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

Revision Information

Version	Date of issue and modifications introduced
2.0	2025-01-31 Typographical error corrected Minor formatting changes Shelf life corrected

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- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:

- a) EUCAST Disk Diffusion Method
- b) Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Almeida, Lais Pinto de, Fabiana Puerro de Carvalho, Alexandre Gimenes Marques, Andrea dos Santos Pereira, Renata Puzzo Bortoleto, and Marinês Dalla Valle Martino. 2012. "Ertapenem Disk Performance to Predict Klebsiella Pneumoniae Carbapenemase Produced by Gram-Negative Bacilli Isolated in a São Paulo City Public Hospital." *Einstein (São Paulo, Brazil)* 10 (4): 439–41. <https://doi.org/10.1590/s1679-45082012000400008>.
5. Vading, M., Samuelsen, B. Haldorsen, A. S. Sundsfjord, and C. G. Giske. 2011. "Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing Klebsiella Pneumoniae with the EUCAST and CLSI Breakpoint Systems." *Clinical Microbiology and Infection* 17 (5): 668–74. <https://doi.org/10.1111/j.1469-0691.2010.03299.x>.
6. Martinez, Pedro, Denisse Garzón, and Salim Mattar. 2012. "CTX-M-Producing Escherichia Coli and Klebsiella Pneumoniae Isolated from Community-Acquired Urinary Tract Infections in Valledupar, Colombia." *Brazilian Journal of Infectious Diseases* 16 (5): 420–25. <https://doi.org/10.1016/j.bjid.2012.05.001>.

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Reference

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Almeida, Lais Pinto de, Fabiana Puerro de Carvalho, Alexandre Gimenes Marques, Andrea dos Santos Pereira, Renata

- Puzzo Bortoletto, and Marinês Dalla Valle Martino. 2012. "Ertapenem Disk Performance to Predict Klebsiella Pneumoniae Carbapenemase Produced by Gram-Negative Bacilli Isolated in a São Paulo City Public Hospital." *Einstein (São Paulo, Brazil)* 10 (4): 439–41. <https://doi.org/10.1590/s1679-45082012000400008>.
5. Vading, M., Samuelsen, B. Haldorsen, A. S. Sundsfjord, and C. G. Giske. 2011. "Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing Klebsiella Pneumoniae with the EUCAST and CLSI Breakpoint Systems." *Clinical Microbiology and Infection* 17 (5): 668–74. <https://doi.org/10.1111/j.1469-0691.2010.03299.x>.
 6. Martinez, Pedro, Denisse Garzón, and Salim Mattar. 2012. "CTX-M-Producing Escherichia Coli and Klebsiella Pneumoniae Isolated from Community-Acquired Urinary Tract Infections in Valledupar, Colombia." *Brazilian Journal of Infectious Diseases* 16 (5): 420–25. <https://doi.org/10.1016/j.bjid.2012.05.001>.

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3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Almeida, Lais Pinto de, Fabiana Puerro de Carvalho, Alexandre Gimenes Marques, Andrea dos Santos Pereira, Renata Puzzo Bortoleto, and Marinês Dalla Valle Martino. 2012. "Ertapenem Disk Performance to Predict Klebsiella Pneumoniae Carbapenemase Produced by Gram-Negative Bacilli Isolated in a São Paulo City Public Hospital." *Einstein (São Paulo, Brazil)* 10 (4): 439–41. <https://doi.org/10.1590/s1679-45082012000400008>.
5. Vading, M., Samuelsen, B. Haldorsen, A. S. Sundsfjord, and C. G. Giske. 2011. "Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing Klebsiella Pneumoniae with the EUCAST and CLSI Breakpoint Systems." *Clinical Microbiology and Infection* 17 (5): 668–74. <https://doi.org/10.1111/j.1469-0691.2010.03299.x>.
6. Martinez, Pedro, Denisse Garzón, and Salim Mattar. 2012. "CTX-M-Producing Escherichia Coli and Klebsiella Pneumoniae Isolated from Community-Acquired Urinary Tract Infections in Valledupar, Colombia." *Brazilian Journal of Infectious Diseases* 16 (5): 420–25. <https://doi.org/10.1016/j.bjid.2012.05.001>.

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Referenzen

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Almeida, Lais Pinto de, Fabiana Puerro de Carvalho, Alexandre Gímenes Marques, Andrea dos Santos Pereira, Renata Puzzo Bortoleto, and Marinês Dalla Valle Martino. 2012. "Ertapenem Disk Performance to Predict Klebsiella Pneumoniae Carbapenemase Produced by Gram-Negative Bacilli Isolated in a São Paulo City Public Hospital." *Einstein (São Paulo, Brazil)* 10 (4): 439–41. <https://doi.org/10.1590/s1679-45082012000400008>.
- Vading, M., Samuelsen, B. Haldorsen, A. S. Sundsfjord, and C. G. Giske. 2011. "Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing Klebsiella Pneumoniae with the EUCAST and CLSI Breakpoint Systems." *Clinical Microbiology and Infection* 17 (5): 668–74. <https://doi.org/10.1111/j.1469-0691.2010.03299.x>.
- Martínez, Pedro, Denisse Garzón, and Salim Mattar. 2012. "CTX-M-Producing Escherichia Coli and Klebsiella Pneumoniae Isolated from Community-Acquired Urinary Tract Infections in Valledupar, Colombia." *Brazilian Journal of Infectious Diseases* 16 (5): 420–25. <https://doi.org/10.1016/j.bjid.2012.05.001>.

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- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Almeida, Lais Pinto de, Fabiana Puerro de Carvalho, Alexandre Gímenes Marques, Andrea dos Santos Pereira, Renata Pozzo Bortoleto, and Marinês Dalla Valle Martino. 2012. "Ertapenem Disk Performance to Predict *Klebsiella Pneumoniae* Carbapenemase Produced by Gram-Negative Bacilli Isolated in a São Paulo City Public Hospital." *Einstein (São Paulo, Brazil)* 10 (4): 439–41. <https://doi.org/10.1590/s1679-45082012000400008>.
- Vading, M., Samuelsen, B. Haldorsen, A. S. Sundsfjord, and C. G. Giske. 2011. "Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing *Klebsiella Pneumoniae* with the EUCAST and CLSI Breakpoint Systems." *Clinical Microbiology and Infection* 17 (5): 668–74. <https://doi.org/10.1111/j.1469-0691.2010.03299.x>.
- Martínez, Pedro, Denisse Garzón, and Salim Mattar. 2012. "CTX-M-Producing *Escherichia Coli* and *Klebsiella Pneumoniae* Isolated from Community-Acquired Urinary Tract Infections in Valledupar, Colombia." *Brazilian Journal of Infectious Diseases* 16 (5): 420–25. <https://doi.org/10.1016/j.bjid.2012.05.001>.

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Viited

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).

4. Almeida, Lais Pinto de, Fabiana Puerro de Carvalho, Alexandre Gimenes Marques, Andrea dos Santos Pereira, Renata Pozzo Bortoleto, and Marinês Dalla Valle Martino. 2012. "Ertapenem Disk Performance to Predict Klebsiella Pneumoniae Carbapenemase Produced by Gram-Negative Bacilli Isolated in a São Paulo City Public Hospital." *Einstein (São Paulo, Brazil)* 10 (4): 439–41. <https://doi.org/10.1590/s1679-45082012000400008>.
5. Vading, M., Samuelsen, B. Haldorsen, A. S. Sundsfjord, and C. G. Giske. 2011. "Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing Klebsiella Pneumoniae with the EUCAST and CLSI Breakpoint Systems." *Clinical Microbiology and Infection* 17 (5): 668–74. <https://doi.org/10.1111/j.1469-0691.2010.03299.x>.
6. Martinez, Pedro, Denisse Garzón, and Salim Mattar. 2012. "CTX-M-Producing Escherichia Coli and Klebsiella Pneumoniae Isolated from Community-Acquired Urinary Tract Infections in Valledupar, Colombia." *Brazilian Journal of Infectious Diseases* 16 (5): 420–25. <https://doi.org/10.1016/j.bjid.2012.05.001>.

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Références

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Almeida, Lais Pinto de, Fabiana Puerro de Carvalho, Alexandre Gímenes Marques, Andrea dos Santos Pereira, Renata Pozzo Bortoleto, and Marinês Dalla Valle Martino. 2012. "Ertapenem Disk Performance to Predict Klebsiella Pneumoniae Carbapenemase Produced by Gram-Negative Bacilli Isolated in a São Paulo City Public Hospital." *Einstein (São Paulo, Brazil)* 10 (4): 439–41. <https://doi.org/10.1590/s1679-45082012000400008>.
- Vading, M., Samuelsen, B. Haldorsen, A. S. Sundsfjord, and C. G. Giske. 2011. "Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing Klebsiella Pneumoniae with the EUCAST and CLSI Breakpoint Systems." *Clinical Microbiology and Infection* 17 (5): 668–74. <https://doi.org/10.1111/j.1469-0691.2010.03299.x>.
- Martínez, Pedro, Denisse Garzón, and Salim Mattar. 2012. "CTX-M-Producing Escherichia Coli and Klebsiella Pneumoniae Isolated from Community-Acquired Urinary Tract Infections in Valledupar, Colombia." *Brazilian Journal of Infectious Diseases* 16 (5): 420–25. <https://doi.org/10.1016/j.bjid.2012.05.001>.

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Referencije

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a) M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b) M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c) M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.


2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a) EUCAST Disk Diffusion Method
 - b) Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Almeida, Lais Pinto de, Fabiana Puerro de Carvalho, Alexandre Gímenes Marques, Andrea dos Santos Pereira, Renata Puzzo Bortoleto, and Marinês Dalla Valle Martino. 2012. "Ertapenem Disk Performance to Predict *Klebsiella Pneumoniae* Carbapenemase Produced by Gram-Negative Bacilli Isolated in a São Paulo City Public Hospital." *Einstein (São Paulo, Brazil)* 10 (4): 439–41. <https://doi.org/10.1590/s1679-45082012000400008>.
5. Vading, M., Samuelsen, B. Haldorsen, A. S. Sundsfjord, and C. G. Giske. 2011. "Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing *Klebsiella Pneumoniae* with the EUCAST and CLSI Breakpoint Systems." *Clinical Microbiology and Infection* 17 (5): 668–74. <https://doi.org/10.1111/j.1469-0691.2010.03299.x>.
6. Martinez, Pedro, Denisse Garzón, and Salim Mattar. 2012. "CTX-M-Producing *Escherichia Coli* and *Klebsiella Pneumoniae* Isolated from Community-Acquired Urinary Tract Infections in Valledupar, Colombia." *Brazilian Journal of Infectious Diseases* 16 (5): 420–25. <https://doi.org/10.1016/j.bjid.2012.05.001>.

© 2025. Thermo Fisher Scientific Inc. Sva prava pridržana. ATCC® je zaštitni znak društva ATCC. Svi ostali zaštitni znakovi vlasništvo su društva Thermo Fisher Scientific Inc. i njegovih podružnica. Ove informacije nisu namijenjene za poticanje na uporabu tih proizvoda na način kojim se mogu povrijediti prava intelektualnog vlasništva drugih društava.

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
- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Almeida, Lais Pinto de, Fabiana Puerro de Carvalho, Alexandre Gimenes Marques, Andrea dos Santos Pereira, Renata Puzzo Bortoleto, and Marinês Dalla Valle Martino. 2012. "Ertapenem Disk Performance to Predict Klebsiella Pneumoniae Carbapenemase Produced by Gram-Negative Bacilli Isolated in a São Paulo City Public Hospital." *Einstein (São Paulo, Brazil)* 10 (4): 439–41. <https://doi.org/10.1590/s1679-45082012000400008>.
- Vading, M., Samuelsen, B. Haldorsen, A. S. Sundsfjord, and C. G. Giske. 2011. "Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing Klebsiella Pneumoniae with the EUCAST and CLSI Breakpoint Systems." *Clinical Microbiology and Infection* 17 (5): 668–74. <https://doi.org/10.1111/j.1469-0691.2010.03299.x>.
- Martinez, Pedro, Denisse Garzón, and Salim Mattar. 2012. "CTX-M-Producing Escherichia Coli and Klebsiella Pneumoniae Isolated from Community-Acquired Urinary Tract Infections in Valledupar, Colombia." *Brazilian Journal of Infectious Diseases* 16 (5): 420–25. <https://doi.org/10.1016/j.bjid.2012.05.001>.

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Referințe

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method

- b) Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
 - Almeida, Lais Pinto de, Fabiana Puerro de Carvalho, Alexandre Gimenes Marques, Andrea dos Santos Pereira, Renata Puzzo Bortoleto, and Marinês Dalla Valle Martino. 2012. "Ertapenem Disk Performance to Predict Klebsiella Pneumoniae Carbapenemase Produced by Gram-Negative Bacilli Isolated in a São Paulo City Public Hospital." *Einstein (São Paulo, Brazil)* 10 (4): 439–41. <https://doi.org/10.1590/s1679-45082012000400008>.
 - Vading, M., Samuelsen, B. Haldorsen, A. S. Sundsfjord, and C. G. Giske. 2011. "Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing Klebsiella Pneumoniae with the EUCAST and CLSI Breakpoint Systems." *Clinical Microbiology and Infection* 17 (5): 668–74. <https://doi.org/10.1111/j.1469-0691.2010.03299.x>.
 - Martinez, Pedro, Denisse Garzón, and Salim Mattar. 2012. "CTX-M-Producing Escherichia Coli and Klebsiella Pneumoniae Isolated from Community-Acquired Urinary Tract Infections in Valledupar, Colombia." *Brazilian Journal of Infectious Diseases* 16 (5): 420–25. <https://doi.org/10.1016/j.bjid.2012.05.001>.

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Referencie

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Almeida, Lais Pinto de, Fabiana Puerro de Carvalho, Alexandre Gimenes Marques, Andrea dos Santos Pereira, Renata Puzzo Bortoleto, and Marinês Dalla Valle Martino. 2012. "Ertapenem Disk Performance to Predict *Klebsiella pneumoniae* Carbapenemase Produced by Gram-Negative Bacilli Isolated in a São Paulo City Public Hospital." *Einstein (São Paulo, Brazil)* 10 (4): 439–41. <https://doi.org/10.1590/s1679-45082012000400008>.
- Vading, M., Samuelsen, B. Haldorsen, A. S. Sundsfjord, and C. G. Giske. 2011. "Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing *Klebsiella pneumoniae* with the EUCAST and CLSI Breakpoint Systems." *Clinical Microbiology and Infection* 17 (5): 668–74. <https://doi.org/10.1111/j.1469-0691.2010.03299.x>.
- Martinez, Pedro, Denisse Garzón, and Salim Mattar. 2012. "CTX-M-Producing *Escherichia coli* and *Klebsiella pneumoniae* Isolated from Community-Acquired Urinary Tract Infections in Valledupar, Colombia." *Brazilian Journal of Infectious Diseases* 16 (5): 420–25. <https://doi.org/10.1016/j.bjid.2012.05.001>.

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Referenser

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a) M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b) M100, Performance Standards for Antimicrobial Susceptibility Testing

- c) M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a) EUCAST Disk Diffusion Method
 - b) Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Almeida, Lais Pinto de, Fabiana Puerro de Carvalho, Alexandre Gimenes Marques, Andrea dos Santos Pereira, Renata Puzzo Bortoleto, and Marinês Dalla Valle Martino. 2012. "Ertapenem Disk Performance to Predict Klebsiella Pneumoniae Carbapenemase Produced by Gram-Negative Bacilli Isolated in a São Paulo City Public Hospital." *Einstein (São Paulo, Brazil)* 10 (4): 439–41. <https://doi.org/10.1590/s1679-45082012000400008>.
5. Vading, M., Samuelsen, B. Haldorsen, A. S. Sundsfjord, and C. G. Giske. 2011. "Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing Klebsiella Pneumoniae with the EUCAST and CLSI Breakpoint Systems." *Clinical Microbiology and Infection* 17 (5): 668–74. <https://doi.org/10.1111/j.1469-0691.2010.03299.x>.
6. Martinez, Pedro, Denisse Garzón, and Salim Mattar. 2012. "CTX-M-Producing Escherichia Coli and Klebsiella Pneumoniae Isolated from Community-Acquired Urinary Tract Infections in Valledupar, Colombia." *Brazilian Journal of Infectious Diseases* 16 (5): 420–25. <https://doi.org/10.1016/j.bjid.2012.05.001>.

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