
	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Tigecycline Discs 15µg (TGC15)

REF CT1841B

Antimicrobial Susceptibility Test Discs

EN

NOTE: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Tigecycline Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs that contain specific amounts of the antimicrobial agent tigecycline. The discs are labelled on both sides with details of the antimicrobial (TGC) and drug content (15µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Tigecycline AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which tigecycline has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling, and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- *Staphylococcus aureus* (including Methicillin-resistant isolates)
- *Streptococcus* spp. (other than *Streptococcus pneumoniae*)
- *Streptococcus pneumoniae*
- *Enterococcus faecalis* (vancomycin-susceptible isolates)

Gram-Negative

- Enterobacteriaceae
- *Haemophilus influenzae*

Species with published breakpoints according to current CLSI M100 literature:

Gram-Positive

- N/A

Gram-Negative

- N/A

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- *Staphylococcus* spp.
- *Enterococcus* spp.
- *Streptococcus* groups A, B, C, and G

Gram-Negative

- Enterobacterales (*E. coli* and *C. koseri*)

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principles of Method

Tigecycline AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology, refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Tigecycline AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Tigecycline AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Tigecycline AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to table 2.

Table 1. Materials Provided with CT1841B	
Component Description	Material Description
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Tigecycline AST Disc Reagents	
Reagent	Description of Function
Tigecycline	Orange Crystalline Powder. Antimicrobial agent. The drug mechanism of action involves the inhibition of protein translation by binding to the 30S ribosomal subunit that blocks the elongation of peptide chains.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Tigecycline AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Tigecycline AST discs have a shelf-life of 12 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2°C to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this document.

Analytical Performance Characteristics

Table 3. Raw data of zone sizes, in accordance with CLSI methodology, taken from the 3 most recent batches of Tigecycline (15 µg) TGC15 Antimicrobial Susceptibility Test Discs (CT1841B).

CLSI Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3304645	<i>Escherichia coli</i>	ATCC 25922	MHA	20	27	23.5	20	20	20	3.5 ± 0 (CV = 0)
	<i>Staphylococcus aureus</i>	ATCC 25923	MHA	20	25	22.5	23	22	22	0.17 ± 0.58 (CV = 3.46)
	<i>Pseudomonas aeruginosa</i>	ATCC 27853	MHA	9	13	11	11	11	11	0 ± 0 (CV = 0)
	<i>Haemophilus influenzae</i>	ATCC 49247	HTM	23	31	27	25	26	26	1.33 ± 0.58 (CV = 0.43)
	<i>Neisseria gonorrhoeae</i>	ATCC 49226	GCV	30	40	35	36	36	36	1 ± 0 (CV = 0)
	<i>Streptococcus pneumoniae</i>	ATCC 49619	MHA + 5% SB	23	29	26	28	28	28	2 ± 0 (CV = 0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3300906	<i>Escherichia coli</i>	ATCC 25922	MHA	20	27	23.5	21	21	21	2.5 ± 0 (CV = 0)
	<i>Staphylococcus aureus</i>	ATCC 25923	MHA	20	25	22.5	24	24	24	1.5 ± 0 (CV = 0)
	<i>Pseudomonas aeruginosa</i>	ATCC 27853	MHA	9	13	11	10	10	10	1 ± 0 (CV = 0)
	<i>Haemophilus influenzae</i>	ATCC 49247	HTM	23	31	27	25	25	25	2 ± 0 (CV = 0)
	<i>Neisseria gonorrhoeae</i>	ATCC 49226	GCV	30	40	35	33	33	33	2 ± 0 (CV = 0)
	<i>Streptococcus pneumoniae</i>	ATCC 49619	MHA + 5% SB	23	29	26	26	26	26	0 ± 0 (CV = 0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3281500	<i>Escherichia coli</i>	ATCC 25922	MHA	20	27	23.5	24	24	23	0.17 ± 0.57 (CV = 3.46)
	<i>Staphylococcus aureus</i>	ATCC 25923	MHA	20	25	22.5	20	20	20	2.5 ± 0 (CV = 0)
	<i>Pseudomonas aeruginosa</i>	ATCC 27853	MHA	9	13	11	9	9	9	2 ± 0 (CV = 0)
	<i>Haemophilus influenzae</i>	ATCC 49247	HTM	23	31	27	24	24	24	3 ± 0 (CV = 0)

	<i>Neisseria gonorrhoeae</i>	ATCC 49226	GCV	30	40	35	34	35	35	0.33 ± 0.58 (CV = 1.73)
	<i>Streptococcus pneumoniae</i>	ATCC 49619	MHA + 5% SB	23	29	26	27	27	27	1 ± 0 (CV = 0)

Table 4. Raw data of zone sizes, in accordance with EUCAST methodology, taken from the 3 most recent batches of Tigecycline (15 µg) TGC15 Antimicrobial Susceptibility Test Discs (CT1841B).

EUCAST Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3304645	<i>Escherichia coli</i>	ATCC 25922	MHA	20	27	23.5	20	20	20	3.5 ± 0 (CV = 0)
	<i>Staphylococcus aureus</i>	ATCC 29213	MHA	19	25	22	20	20	21	1.67 ± 0.58 (CV = 0.35)
	<i>Enterococcus faecalis</i>	ATCC 29212	MHA	20	26	23	22	22	22	1 ± 0 (CV = 0)
	<i>Streptococcus pneumoniae</i>	ATCC 49619	MHA-F	24	30	27	25	25	25	2 ± 0 (CV = 0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3300906	<i>Escherichia coli</i>	ATCC 25922	MHA	20	27	23.5	21	21	21	2.5 ± 0 (CV = 0)
	<i>Staphylococcus aureus</i>	ATCC 29213	MHA	19	25	22	21	21	21	1 ± 0 (CV = 0)
	<i>Enterococcus faecalis</i>	ATCC 29212	MHA	20	26	23	22	22	22	1 ± 0 (CV = 0)
	<i>Streptococcus pneumoniae</i>	ATCC 49619	MHA-F	24	30	27	24	24	24	3 ± 0 (CV = 0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3281500	<i>Escherichia coli</i>	ATCC 25922	MHA	20	27	23.5	24	24	23	0.17 ± 0.58 (CV = 3.46)
	<i>Staphylococcus aureus</i>	ATCC 29213	MHA	19	25	22	21	21	21	1 ± 0 (CV = 0)
	<i>Enterococcus faecalis</i>	ATCC 29212	MHA	20	26	23	21	21	21	2 ± 0 (CV = 0)
	<i>Streptococcus pneumoniae</i>	ATCC 49619	MHA-F	24	30	27	24	24	24	3 ± 0 (CV = 0)

Clinical Performance Characteristics

Begum et al., in 2013, aimed to determine the prevalence of multi-drug resistant *A. baumannii* isolated from patients of the Surgical Intensive Care Units (SICU) of the Pakistan Institute of Medical Sciences.⁴ A total of 91 bacterial isolates were collected and characterised as *A. baumannii* as they appeared as non-lactose fermenters and non-pigmented on MacConkey and blood agar respectively. All 91 *A. baumannii* strains isolated demonstrated zero resistance to tigecycline (as determined by Oxoid Tigecycline (15 µg) TGC15 AST Discs utilising CLSI guidelines). This study showed that Oxoid Tigecycline AST Discs are effective in the clinical analysis of *A. baumannii* as there were no out-of-specification results and high susceptibility was displayed.

Hafeez et al., 2016 aimed to determine the frequency and antimicrobial susceptibility of 526 samples isolated from the intensive care units (ICU) from a Military hospital in Rawalpindi.⁵ Of the isolates, *A. baumannii* made up the highest proportion (31.08%), then *K. pneumoniae* (19.18%), *P. aeruginosa* (16.75%), *E. coli* (8.37%). Susceptibility among the pathogens was highest to tigecycline with 89.65%, 85.91% and 80.64% of *A. baumannii*, *K. pneumoniae* and *E. coli* isolates being susceptible. This study demonstrated that Oxoid Tigecycline (15 µg) TGC15 AST Discs are effective in the in vitro analysis of clinical isolates and demonstrated high susceptibility among Gram-negative pathogens.

In 2007, Hope et al., determined zone breakpoints to be used with Oxoid Tigecycline (15 µg) TGC15 AST Discs (in accordance with BSAC methods) utilising bacterial isolates from the 2004 BSAC bacteraemia surveillance programme.⁶ The collected isolates consisted of *Streptococci*, *Staphylococci*, *Citrobacter* spp., *Enterobacter* spp., *Klebsiella* spp., among others. The study concluded that the Oxoid Tigecycline (15 µg) TGC15 AST Discs accurately differentiated between susceptible, intermediate, and resistant isolates for most species, excluding *Klebsiella* spp. for which differentiation of the susceptible/intermediate border was challenging.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References


1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests.
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing.
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Begum, Shahzeera, Fariha Hasan, Shagufta Hussain, and Aamer Ali Shah. 2013. "Prevalence of Multi Drug Resistant Acinetobacter Baumannii in the Clinical Samples from Tertiary Care Hospital in Islamabad, Pakistan." *Pakistan Journal of Medical Sciences* 29 (5): 1253–58.
5. Hafeez, Amira, Tehmina Munir, Sara Najeeb, Sabahat Rehman, Mehreen Gilani, Maliha Ansari, and Nadia Saad. 2016. "ICU Pathogens: A Continuous Challenge." *Journal of the College of Physicians and Surgeons Pakistan* 26 (7): 577–80. <https://doi.org/2374>.
6. Hope, Russell, T. Parsons, S. Mushtaq, D. James, and D. M. Livermore. 2007. "Determination of Disc Breakpoints and Evaluation of Etests for Tigecycline Susceptibility Testing by the BSAC Method." *Journal of Antimicrobial Chemotherapy* 60 (4): 770–74. <https://doi.org/10.1093/jac/dkm297>.

© 2025 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number

	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Period-after-opening – the number of days the product can be used after opening
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

Revision Information

Version	Date of issue and modifications introduced
2.0	2025-02-12 Shelf life corrected

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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Reference

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests.
 - M100, Performance Standards for Antimicrobial Susceptibility Testing.
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Begum, Shahzeera, Fariha Hasan, Shagufta Hussain, and Amer Ali Shah. 2013. "Prevalence of Multi Drug Resistant Acinetobacter Baumannii in the Clinical Samples from Tertiary Care Hospital in Islamabad, Pakistan." *Pakistan Journal of Medical Sciences* 29 (5): 1253–58.
- Hafeez, Amira, Tehmina Munir, Sara Najeeb, Sabahat Rehman, Mehreen Gilani, Maliha Ansari, and Nadia Saad. 2016. "ICU Pathogens: A Continuous Challenge." *Journal of the College of Physicians and Surgeons Pakistan* 26 (7): 577–80. <https://doi.org/2374>.
- Hope, Russell, T. Parsons, S. Mushtaq, D. James, and D. M. Livermore. 2007. "Determination of Disc Breakpoints and Evaluation of Etests for Tigecycline Susceptibility Testing by the BSAC Method." *Journal of Antimicrobial Chemotherapy* 60 (4): 770–74. <https://doi.org/10.1093/jac/dkm297>.

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Litteratur

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests.
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing.
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Begum, Shahzeera, Fariha Hasan, Shagufta Hussain, and Aamer Ali Shah. 2013. "Prevalence of Multi Drug Resistant Acinetobacter Baumannii in the Clinical Samples from Tertiary Care Hospital in Islamabad, Pakistan." *Pakistan Journal of Medical Sciences* 29 (5): 1253–58.
5. Hafeez, Amira, Tehmina Munir, Sara Najeeb, Sabahat Rehman, Mehreen Gilani, Maliha Ansari, and Nadia Saad. 2016. "ICU Pathogens: A Continuous Challenge." *Journal of the College of Physicians and Surgeons Pakistan* 26 (7): 577–80. <https://doi.org/2374>.
6. Hope, Russell, T. Parsons, S. Mushtaq, D. James, and D. M. Livermore. 2007. "Determination of Disc Breakpoints and Evaluation of Etests for Tigecyclin Susceptibility Testing by the BSAC Method." *Journal of Antimicrobial Chemotherapy* 60 (4): 770–74. <https://doi.org/10.1093/jac/dkm297>.

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1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests.
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing.
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Begum, Shahzeera, Fariha Hasan, Shagufta Hussain, and Aamer Ali Shah. 2013. "Prevalence of Multi Drug Resistant Acinetobacter Baumannii in the Clinical Samples from Tertiary Care Hospital in Islamabad, Pakistan." Pakistan Journal of Medical Sciences 29 (5): 1253–58.
5. Hafeez, Amira, Tehmina Munir, Sara Najeeb, Sabahat Rehman, Mehreen Gilani, Maliha Ansari, and Nadia Saad. 2016. "ICU Pathogens: A Continuous Challenge." Journal of the College of Physicians and Surgeons Pakistan 26 (7): 577–80. <https://doi.org/2374>.
6. Hope, Russell, T. Parsons, S. Mushtaq, D. James, and D. M. Livermore. 2007. "Determination of Disc Breakpoints and Evaluation of Etests for Tigecycline Susceptibility Testing by the BSAC Method." Journal of Antimicrobial Chemotherapy 60 (4): 770–74. <https://doi.org/10.1093/jac/dkm297>.

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Atsauces

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests.
 - M100, Performance Standards for Antimicrobial Susceptibility Testing.
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Begum, Shahzeera, Fariha Hasan, Shagufta Hussain, and Aamer Ali Shah. 2013. "Prevalence of Multi Drug Resistant Acinetobacter Baumannii in the Clinical Samples from Tertiary Care Hospital in Islamabad, Pakistan." *Pakistan Journal of Medical Sciences* 29 (5): 1253–58.
- Hafeez, Amira, Tehmina Munir, Sara Najeeb, Sabahat Rehman, Mehreen Gilani, Maliha Ansari, and Nadia Saad. 2016. "ICU Pathogens: A Continuous Challenge." *Journal of the College of Physicians and Surgeons Pakistan* 26 (7): 577–80. <https://doi.org/2374>.
- Hope, Russell, T. Parsons, S. Mushtaq, D. James, and D. M. Livermore. 2007. "Determination of Disc Breakpoints and Evaluation of Etests for Tigecycline Susceptibility Testing by the BSAC Method." *Journal of Antimicrobial Chemotherapy* 60 (4): 770–74. <https://doi.org/10.1093/jac/dkm297>.

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Referencie

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests.
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing.
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Begum, Shahzeera, Fariha Hasan, Shagufta Hussain, and Amer Ali Shah. 2013. "Prevalence of Multi Drug Resistant Acinetobacter Baumannii in the Clinical Samples from Tertiary Care Hospital in Islamabad, Pakistan." *Pakistan Journal of Medical Sciences* 29 (5): 1253–58.
5. Hafeez, Amira, Tehmina Munir, Sara Najeeb, Sabahat Rehman, Mehreen Gilani, Maliha Ansari, and Nadia Saad. 2016. "ICU Pathogens: A Continuous Challenge." *Journal of the College of Physicians and Surgeons Pakistan* 26 (7): 577–80. <https://doi.org/2374>.
6. Hope, Russell, T. Parsons, S. Mushtaq, D. James, and D. M. Livermore. 2007. "Determination of Disc Breakpoints and Evaluation of Etests for Tigecycline Susceptibility Testing by the BSAC Method." *Journal of Antimicrobial Chemotherapy* 60 (4): 770–74. <https://doi.org/10.1093/jac/dkm297>.

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