Sabouraud Dextrose Agar with Chloramphenicol

EN

REF PO0161A, PO0358A (Deep Fill), PO1001A (Extra Deep Fill) and PB1219E

* This instructions for use (IFU) document is intended to be read in conjunction with the IFU for Columbia CAP Selective Agar with Horse Blood. Product code: PB1161A See additional IFUs available at www.thermofisher.com.

Intended Use

Sabouraud Dextrose Agar with Chloramphenicol (PO0161A, PO0358A, PO1001A and PB1219E) is an acidic pH selective medium for the isolation of dermatophytes, other fungi, and yeasts from skin, hair, nails, genital, respiratory and urine samples. Sabouraud Dextrose Agar with Chloramphenicol devices (PO0161A, PO0358A, PO1001A and PB1219E) are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having fungal infections.

The devices are for professional use only, are not automated nor are they companion diagnostics.

Summary and Explanation

Candidiasis can have a range of manifestations ranging from superficial to life-threatening.¹⁻¹² Less invasive manifestations of *Candida* spp. infections include superficial dermatological infections such as athlete's foot or onychomycosis, to more complicated issues such as bloodstream infections or endocarditis.^{4,10} Eye infections can typically exhibit themselves as mild, but quickly become infectious or sight-threatening if treatment is ineffective.⁸ Keratitis is an infection of the cornea that requires rapid and thorough investigation but is caused by a range of fungi including *Aspergillus* and *Candida* species.^{13,14} Similarly, infectious endophthalmitis is a potentially sight-hindering infection of the intraocular tissue and surrounding fluids that can be caused by yeast and mould infections. Due to the range of consequences caused, it is important to distinguish between *Candida* spp. in clinical samples as susceptibility to antifungal treatment can differ between species and ineffective treatment can be detrimental to recovery.

Principle of Method

These devices contain mycological peptone and dextrose to support the growth of fungi. The agar is slightly acidic to help inhibit bacterial growth in favour of fungi growth. Chloramphenicol is added to the media because it has bacteriostatic activity against Gram-negative and Gram-positive bacteria, further favouring fungal growth and increasing selectivity of the medium.

Typical Formulas*

PO0161A, PO0358A, PO1001A, PB1219E

	<u>grams per litre</u>
Mycological peptone	10.0
Glucose	40.0
Agar	15.0

Additions

Chloramphenicol	100 mg
-----------------	--------

*Adjusted as required to meet performance standards.

Physical Appearance

PO0161A

Colour	Straw
Clarity	Opaque
pH	5.6 ± 0.2
Fill weight	19.0 ± 2.0

PO0358A

Colour	Straw
Clarity	Opaque
pH	5.6 ± 0.2
Fill weight	22.75 ± 0.75

PO1001A

Colour	Straw 2
Clarity	Opaque
pH	6.0 ± 0.2
Fill weight	25.5 ± 1

PB1219E

Colour	Straw 1, Straw 1-2 or Straw 2
Clarity	Clear
pH	5.6 ± 0.2
Total fill weight	19.0 ± 2.0 g

Materials Provided

PO0161A: 10 x 90mm Sabouraud Dextrose Agar with Chloramphenicol plates

PO0358A: 10 x 90mm Sabouraud Dextrose Agar with Chloramphenicol (deep fill) plates

PO1001A: 10 x 90mm Sabouraud Dextrose Agar with Chloramphenicol (deep fill) plates

PB1219E: 10 x 90mm Staph Strep/ Sabouraud Dextrose Agar + Chloramphenicol biplates

Each plate should only be used once.

Materials Required but Not Supplied

- Inoculating loops
- Swabs
- Collection containers
- Incubators
- Quality control organisms

Storage

- Store product in its original packaging at 2–10°C until used.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

Warnings and Precautions

- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or plates.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Material Safety Data Sheet (MSDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local the recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 1, B 2, B 5, B 9, B11, B 28, B 39, B 41, and B 57.

Procedure

The sample preparation and inoculation procedure may vary with the specimen; an example is given in the steps below.

- Inoculate with the sample (skin dandruff, skin scrapings, hair, nail pieces, etc.) directly or after fragmentation and dilution.
- Incubate plates aerobically at 23° to 27°C for 5 days depending on the type of sample and type of fungus sought.
- Examine the plates every 2 to 3 days.
- Identify isolates by performing morphological (microscopic and macroscopic), biochemical, and/or serological tests as required.

Interpretation

Candida albicans is indicated by the presence of cream coloured, domed colonies.

Aspergillus brasiliensis is indicated by the presence of white mycelium, black spores.

For PO1001A, *Trichophyton rubrum* is indicated by the presence of white mycelia, buff spores.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

PO0161A, PO0358A, PO1001A, PB1219E:

Positive Controls Inoculum 10 -100 colony forming units Colony count shall be equal to or greater than 50% of the control medium	
Results after incubation at 23-27°C for 5 days	
<i>Candida albicans</i> ATCC® 10231™	Cream, domed colonies
<i>Aspergillus brasiliensis</i> ATCC® 16404™	White mycelia, black spores
Results after incubation at 28-32°C for 7 days	
<i>Trichophyton rubrum</i> ATCC® 28191™*	White mycelia, buff spores
Negative Controls Inoculum 10,000 – 100,000 colony forming units Negative strains are restricted in growth or are inhibited.	
<i>Escherichia Coli</i> ATCC® 8739™	No growth

* for PO1001A only

Limitations

Devices PO0161A and PO0358A may not be suitable for the isolation of all fungi. Some microorganisms, for example *Nocardia spp.* (filamentous bacteria) may be inhibited by chloramphenicol. Due to variation in nutritional requirements some strains may grow poorly or fail to grow such as *Malassezia furfur* which requires a lipid supplement. Bacteria may grow if they are resistant to chloramphenicol. Medium choice and incubation protocol will be dependent upon the specimen type. All identifications are presumptive and should be confirmed using appropriate methods. A non-selective medium should be inoculated alongside the selective medium.

Identifications of fungi, dermatophytes and yeasts are presumptive and should be confirmed.

Performance

Accuracy has been demonstrated through review of the QC data. Correct detection of fungi, yeast and dermatophyte strains is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the devices

The precision of the devices was demonstrated by an overall pass rate of 100% for all products; ten batches per product:

- PO0358A: between 26.10.2021 – 29.12.2021
- PO1001A: between 04.07.2021 – 27.12.2021
- PO0161A: between 06.12.2021 – 10.01.2022
- PB1219E: between 15.10.2021 – 28.12.2021

The devices are tested in-house as part of the QC process since the products were launched in 1997 (PO0358A, PO0161A), 2002 (PO1001A) and 2011 (PB1219E). For target organisms, when using 10-100 cfu of *Candida albicans*, *Aspergillus brasiliensis* and *Trichophyton rubrum* and incubating the device at 23-27°C for 5 days (*Candida albicans*, *Aspergillus brasiliensis*) or 28-32°C for 7 days (PO1001A) (*Trichophyton rubrum* ATCC® 28191™), the user can recover organisms with colony size and morphology as listed in this document. For non-target organisms, when using >10⁴ cfu for *Escherichia* and incubating the device at 32°C for 5 days (PO0161A) or 23-27°C for 5 days (PB1219E, PO0358A, PO1001A), (28-32°C for 7 days PO1001A), the users can expect no growth.


Bibliography

1. Public Health England. 2013. 'Sexually Transmitted Infections'. UK Standards for Microbiology Investigations S 6 (1.2). <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>.
2. Public Health England. 2015a. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
3. Public Health England. 2015b. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
4. Public Health England. 2016. 'Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
5. Public Health England. 2017a. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
6. Public Health England. 2017b. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
7. Public Health England. 2017c. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
8. Public Health England. 2017d. 'Investigation of Bacterial Eye Infections'. UK Standards for Microbiology Investigations B 2 (6.1). <https://www.gov.uk/government/publications/smi-b-2-investigation-of-bacterial-eye-infections>.
9. Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical>.

wound-swabs.

10. Public Health England. 2019a. 'Investigation of Blood Cultures (for Organisms Other than Mycobacterium Species)'. UK Standards for Microbiology Investigations B 37 (8.2). <https://www.gov.uk/government/publications/smi-b-37-investigation-of-blood-cultures-for-organisms-other-than-mycobacterium-species>.
11. Public Health England. 2019b. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
12. Public Health England. 2019c. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
13. Pflugfelder SC, Flynn HW, Jr. Infectious endophthalmitis. *Infectious Disease Clinics of North America*, 6, 859-73
14. Brod R.D., Flynn H.W.Jr. (1993). Endophthalmitis: current approaches to diagnosis and therapy. *Current Opinion in Infectious Diseases*, 6, 628-37

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

CLSI is a trademark of the Clinical Laboratory and Standards Institute. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, England

For technical assistance please contact your local distributor.

Revision information

Version	Date of issue and modifications introduced
2.0	2024-06-12

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1. Public Health England. 2013. 'Sexually Transmitted Infections'. UK Standards for Microbiology Investigations S 6 (1.2). <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>.
2. Public Health England. 2015a. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
3. Public Health England. 2015b. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
4. Public Health England. 2016. 'Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
5. Public Health England. 2017a. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
6. Public Health England. 2017b. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
7. Public Health England. 2017c. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.

8. Public Health England. 2017d. 'Investigation of Bacterial Eye Infections'. UK Standards for Microbiology Investigations B 2 (6.1). <https://www.gov.uk/government/publications/smi-b-2-investigation-of-bacterial-eye-infections>.
9. Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
10. Public Health England. 2019a. 'Investigation of Blood Cultures (for Organisms Other than Mycobacterium Species)'. UK Standards for Microbiology Investigations B 37 (8.2). <https://www.gov.uk/government/publications/smi-b-37-investigation-of-blood-cultures-for-organisms-other-than-mycobacterium-species>.
11. Public Health England. 2019b. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
12. Public Health England. 2019c. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
13. Pflugfelder SC, Flynn HW, Jr. Infectious endophthalmitis. *Infectious Disease Clinics of North America*, 6, 859-73
14. Brod R.D., Flynn H.W.Jr. (1993). Endophthalmitis: current approaches to diagnosis and therapy. *Current Opinion in Infectious Diseases*, 6, 628-37

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Sabouraud Dextrose Agar With Chloramphenicol

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
Referenser

1. Public Health England. 2013. 'Sexually Transmitted Infections'. UK Standards for Microbiology Investigations S 6 (1.2). <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>.
2. Public Health England. 2015a. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
3. Public Health England. 2015b. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
4. Public Health England. 2016. 'Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
5. Public Health England. 2017a. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
6. Public Health England. 2017b. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
7. Public Health England. 2017c. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
8. Public Health England. 2017d. 'Investigation of Bacterial Eye Infections'. UK Standards for Microbiology Investigations B 2 (6.1). <https://www.gov.uk/government/publications/smi-b-2-investigation-of-bacterial-eye-infections>.
9. Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non->

- surgical- wound-swabs.
10. Public Health England. 2019a. 'Investigation of Blood Cultures (for Organisms Other than Mycobacterium Species)'. UK Standards for Microbiology Investigations B 37 (8.2). <https://www.gov.uk/government/publications/smi-b-37-investigation-of-blood-cultures-for-organisms-other-than-mycobacterium-species>.
 11. Public Health England. 2019b. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
 12. Public Health England. 2019c. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
 13. Pflugfelder SC, Flynn HW, Jr. Infectious endophthalmitis. *Infectious Disease Clinics of North America*, 6, 859-73
 14. Brod R.D., Flynn H.W.Jr. (1993). Endophthalmitis: current approaches to diagnosis and therapy. *Current Opinion in Infectious Diseases*, 6, 628-37

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