REF R8311006.....
Σ 20

1. INTENDED USE

RapID ONE System is a qualitative micromethod using enzyme reactions to identify clinical isolates grown on agar of oxidase negative, Gram-negative Enterobacteriaceae. Used in a diagnostic workflow to aid clinicians in treatment options for patients suspected of having bacterial infections. The device is not automated, is for professional use only and is not a companion diagnostic.

A complete listing of the organisms addressed by the RapID ONE system is provided in the RapID ONE Differential Chart.

2. SUMMARY AND EXPLANATION

RapID ONE System is comprised of (1) RapID ONE Panels and (2) RapID ONE Reagent. RapID ONE Panels are disposable plastic trays with 18 reaction cavities, which contain dehydrated reactants. The panel allows for the simultaneous inoculation of each cavity with a predetermined amount of inoculum. A suspension of the test organism in RapID Inoculation Fluid is used as the inoculum which rehydrates and initiates test reactions. After incubation of the panel, each test cavity is examined for reactivity by noting the development of a color. In some cases, reagents must be added to the test cavities to provide a color change. The resulting pattern of positive and negative test scores is used as the basis for identification of the test isolate by comparison of test results with the probability values in the Differential Chart (Table 4), or by use of the RapID ERIC™ software.

3. PRINCIPLE

The tests used in the RapID ONE System are based on microbial degradation of specific substrates detected by various indicator systems. The reactions employed are a combination of conventional tests and single-substrate chromogenic tests, described in Table 1.

4. REAGENTS

RapID ONE Reagent (provided with kit)	(15 ml/Btl)
Reactive ingredient per liter:	
p-Dimethylaminocinnamaldehyde	0.06 g
RapID Inoculation Fluid (R8325106, supplied separately)	(2 ml/tube)
KCl	6.0 g
CaCl ₂	0.5 g
Demineralized Water	1000.0 ml
RapID Spot Indole Reagent (R8309002, supplied separately)	(15 ml/Btl)
p-Dimethylaminocinnamaldehyde	10.0 g
Hydrochloric Acid	100.0 ml
Demineralized Water	900.0 ml

5. PRECAUTIONS

This product is for *in vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, media, and test panels after use. Directions should be read and followed carefully.

Do not use reagents beyond the printed expiration dates.

Do not use if there is evidence of contamination or other signs of deterioration.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction do not use device.

Caution!

1. RapID ONE Reagent is toxic and may cause harm to the environment. Harmful by inhalation, contact with skin or eyes, or if swallowed. May impair fertility or cause harm to unborn child.
2. RapID Spot Indole Reagent may cause irritation to skin, eyes, and respiratory system.
3. Refer to Safety Data Sheet for detailed information on reagent chemicals.

Composition / information on ingredients

2-Methoxyethanol 109-86-4
Acetic acid 64-19-7
Hydrochloric acid 7647-01-0

WARNING! This product contains a chemical known in the State of California to cause birth defects or other reproductive harm.

Emergency Telephone Number

INFOTRAC - 24 Hour Number: 1-800-535-5053

Outside of the United States, call 24 Hour Number: 001-352-323-3500 (Call Collect)

DANGER

US ONLY

US & EU

6. STORAGE

2°C

8°C

RapID ONE System and RapID Spot Indole Reagent should be stored in their original containers at 2-8°C until use. Allow products to equilibrate to room temperature before use. DO NOT interchange reagents among different RapID systems. Remove only the number of panels necessary for testing. Reseal the plastic pouch and promptly return to 2-8°C. Panels must be used the same day they are removed from storage. RapID Inoculation Fluid should be stored in its original container at room temperature (20-25°C) until used.

7. PRODUCT DETERIORATION

This product should not be used if (1) the expiration date has passed, (2) the plastic tray is broken or the lid is compromised, or (3) there are other signs of deterioration.

8. SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected and handled following recommended guidelines.^{8,9}

9. MATERIALS SUPPLIED

- 20 RapID ONE Panels
- 20 report forms
- RapID ONE Reagent (one plastic dropper-bottle containing reagent sufficient for 20 panels),
- 2 chipboard incubation trays
- 1 color guide
- Instructions for use (IFU).

10. CONTENTS SYMBOLS

ONE Panels	ONE Panels
Report Forms	RapID Report Forms
ONE Reagent	ONE Reagent
Incubation Trays	Incubation Trays

11. MATERIALS REQUIRED BUT NOT SUPPLIED

- Loop sterilization device
- Inoculating loop, swab, collection containers,
- Incubators, alternative environmental systems,
- Supplemental media
- Quality control organisms
- Gram stain reagents
- Microscope slides
- Oxidase reagent
- Cotton swabs
- RapID Inoculation Fluid - 2 ml (R8325106)
- McFarland #2 turbidity standard or equivalent (R20412)
- Pipettes
- RapID Spot Indole Reagent (R8309002)
- ERIC (R8323600) (optional).

12. PROCEDURE

Inoculum Preparation:

1. Test organisms must be grown in pure culture and examined by Gram stain and oxidase prior to use in the system.

Notes:

- Only oxidase-negative, Gram-negative bacilli should be tested using the RapID ONE System. Oxidase-positive bacilli should be tested using the RapID NF Plus System (R8311005).
- The oxidase test should be interpreted with caution when using bacterial growth from differential agars that contain dyes which may interfere with interpretation.

2. Test organisms may be removed from a variety of selective and nonselective agar growth media. The following types of media are recommended:

Nonselective Media: Tryptic Soy Agar (TSA) with blood.

Differential or Selective Media: Hektoen Enteric (HE) Agar; MacConkey Agar; Eosin Methylene Blue (EMB) Agar; Desoxycholate Agar; Salmonella-Shigella (SS) Agar.

Notes:

- Plates used for inoculum preparation should preferably be 18-24 hours old. Slow-growing isolates may be tested using 48-hour plates.
- The use of media other than those recommended may compromise test performance.

3. Using a cotton swab or inoculating loop, suspend sufficient growth from the agar plate culture in RapID Inoculation Fluid (2 ml) to achieve a visual turbidity equal to a #2 McFarland turbidity standard or equivalent.

Notes:

- Suspensions significantly less turbid than a #2 McFarland standard will result in aberrant reactions.
- Bacterial suspensions that are slightly more turbid than a #2 McFarland standard will not affect test performance and are recommended for stock cultures and quality control strains. However, suspensions prepared with a turbidity far greater than a #2 McFarland standard will compromise test performance.
- Suspensions should be mixed thoroughly and vortexed if required.
- Suspensions should be used within 15 minutes of preparation.
- 4. An agar plate may be inoculated for purity and any additional testing that may be required using a loopful of the test suspension from the inoculation fluid tube. Incubate the plate for at least 18-24 hours at 35-37°C.

Inoculation of RapID ONE Panels:

1. Peel back the lid of the panel over the inoculation port by pulling the tab marked "Peel to Inoculate" up and to the left.
2. Using a pipette, gently transfer the entire contents of the Inoculation Fluid tube into the upper right-hand corner of the panel. Reseal the inoculation port of the panel by pressing the peel-back tab back in place.
3. After adding the test suspension, and while keeping the panel on a level surface, tilt the panel back away from the reaction cavities at approximately 45° (see below).

4. While tilted back, gently rock the panel from side to side to evenly distribute the inoculum along the rear baffles as illustrated below.


Table 1. Principles and Components of the RapID ONE System

Cavity #	Test Code	Reactive Ingredient	Quantity	Principle	Bibliography #
1	URE	Urea	0.25%	Hydrolysis of urea produces basic products which raise the pH and change the indicator	1, 3
2	ADH	Arginine	1.0%	Hydrolysis of the amino acid substrate produces basic products which raise the pH and change the indicator	1, 3
3	ODC	Ornithine	1.0%		
4	LDC	Lysine	1.0%		
5	TET	Aliphatic thiol	0.2%	Utilization of the thiol compound produces acidic products which lower the pH and change the indicator	4
6	LIP	Fatty acid ester	1.0%	Hydrolysis of the fatty acid ester releases acidic products which lower the pH and change the indicator	3, 4
7	KSF	Sugar aldehyde	1.0%		
8	SBL	Sorbitol	1.0%		
9	GUR	<i>p</i> -Nitrophenyl-β, D-glucuronide	0.1%		
10	ONPG	<i>o</i> -Nitrophenyl-β, D-galactoside	0.1%		
11	βGLU	<i>p</i> -Nitrophenyl-β,D-glucoside	0.1%		
12	βXYL	<i>p</i> -Nitrophenyl-β,D-xyloside	0.1%		
13	NAG	<i>p</i> -Nitrophenyl-N-acetyl-β, D-glucosaminide	0.1%		
14	MAL	Malonate	0.5%	Utilization of malonate produces basic products which raise the pH and change the indicator	1, 3
15	PRO	Proline-β-naphthylamide	0.1%		
16	GGT	γ-Glutamyl- β-naphthylamide	0.1%		
17	PYR	Pyrrolidonyl- β-naphthylamide	0.1%		
18	ADON	Adonitol	1.0%	Utilization of the carbohydrate substrate produces acidic products which lower the pH and change the indicator	1, 3
19	IND	Tryptophane	0.4%	Utilization of tryptophane results in the formation of indole which is detected with RapID Spot Indole Reagent.	1-3

Scoring of RapID ONE Panels:

RapID ONE panels contain 18 reaction cavities that provide 19 test scores. Test cavity 18 is bifunctional, containing two separate tests in the same cavity. Bifunctional tests are first scored before the addition of reagent providing the first test result, and then the same cavity is scored again after the addition of reagent to provide the second test result. Test cavities 15 through 17 require RapID ONE Reagent and are designated with a box drawn around them. Bifunctional test 18, which uses RapID Spot Indole Reagent, is designated with a box drawn around the reagent-requiring test.

1. While firmly holding the RapID ONE panel on the bench top, peel off the label lid over the reaction cavities by pulling the lower right hand tab up and to the left.

2. Add 2 drops of RapID ONE Reagent to cavities 15 (PRO) through 17 (PYR).

3. Read and score cavities 1 (URE) through 18 (ADON) from left to right using the interpretation guide presented in Table 2 and the color guide. Panels should be read by looking down through the reaction wells against a white background. Record test scores in the appropriate boxes on the report form using the test code above the box for the bifunctional test.

4. Add 2 drops of RapID Spot Indole Reagent to cavity 18 (ADON/IND).

Note: Only RapID Spot Indole Reagent (R8309002) should be used. Kovacs' or Ehrlich's Indole reagent will not provide satisfactory results.

5. Allow at least 10 seconds but no more than 2 minutes for color development.
6. Read and score test cavity 18 (IND). Record the score in the appropriate box of the report form.
7. Reference the microcode obtained on the report form in ERIC for the identification.

13. RESULTS AND RANGE OF EXPECTED VALUES

The RapID ONE Differential Chart (Table 4) illustrates the expected results for the RapID ONE System. Differential Chart results are expressed as a series of positive percentages for each system test. This information statistically supports the use of each test and provides the basis, through numerical coding of digital test results, for a probabilistic approach to the identification of the test isolate.

Identifications are made using individual test scores from RapID ONE panels in conjunction with other laboratory information (e.g., Gram stain and oxidase) to produce a pattern that statistically resembles known reactivity for taxa recorded in the RapID System database. These patterns are compared through the use of the RapID ONE Differential Chart, or by derivation of a microcode and the use of ERIC.

RapID ONE Panel Test Location

Cavity #	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18
Test Code	URE	ADH	ODC	LDC	TET	LIP	KSF	SBL	GUR	ONPG	βGLU	βXYL	NAG	MAL	PRO	GGT	PYR	ADON

</tbl_r

Table 3. Quality Control Chart for RapID ONE Panels

Organism	URE	ADH	ODC	LDC	TET	LIP	KSF	SBL	GUR	ONPG	βGLU	βXYL	NAG	MAL	PRO	GGT	PYR	ADON	IND	OXI
<i>Proteus vulgaris</i> ATCC® 6380	+	V	-	-	+	-	-	-	-	V	-	-	-	-	+	-	-	+	-	-
<i>Escherichia coli</i> ATCC® 25922	-	-	+	(+)	-	-	-	(+)	+	+	-	-	-	-	-	-	-	+	-	-
<i>Pseudomonas aeruginosa</i> ATCC® 27853	V	+	V	V	+	-	-	-	-	-	-	-	V	+	+	V	-	-	+	-
<i>Enterobacter aerogenes</i> ATCC® 13048	-	-	+	+	-	-	+	+	-	+	+	+	+	+	V	V	+	+	-	-

+, positive; -, negative; V, variable; (+), usually positive

^a Key indicator strains demonstrate acceptable performance of the most labile substrate in the system and reactivity in a significant number of wells, according to Clinical and Laboratory Standards Institute recommendations for streamlined quality control.¹⁴

14. QUALITY CONTROL

All lot numbers of the RapID ONE System have been tested using the following quality control organisms (Table 3) and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported. Table 3 lists expected results for the selected battery of test organisms.

Notes:

- RapID reagent quality control is accomplished by obtaining the expected reactions for tests requiring the addition of the reagents (cavities 15-18).
- Organisms which have been repeatedly transferred on agar media for prolonged periods may provide aberrant results.
- Quality control strains should be stored frozen or lyophilized. Prior to use, quality control strains should be transferred 2-3 times from storage on an agar medium that is recommended for use with the RapID ONE System.
- Formulations, additives, and ingredients of culture media vary from manufacturer to manufacturer and may vary from batch to batch. As a result, culture media may influence constitutive enzymatic activity of designated quality control strains. If quality control strain results differ from the patterns indicated, a subculture onto medium from a different batch or from another manufacturer will often resolve quality control discrepancies.

15. LIMITATIONS

- The use of RapID ONE System and the interpretation of results requires the knowledge of a competent microbiologist, familiar with laboratory procedures, who is trained in general microbiological methods and who judiciously makes use of training, experience, specimen information, and other pertinent procedures before reporting the identification obtained using this system.
- Specimen source, oxidase reaction, Gram stain characteristics, and growth on selective agars should be considered when using the RapID ONE System.

- The RapID ONE System must be used with pure cultures of test organisms. The use of mixed microbial populations or direct testing of clinical material without culture will result in aberrant results.
- The RapID ONE System is designed for use with the taxa listed in the RapID ONE Differential Chart. The use of organisms not specifically listed may lead to misidentifications.
- Expected values listed for RapID ONE System tests may differ from conventional test results or previously reported information.
- The accuracy of the RapID ONE System is based upon the statistical use of a multiplicity of specially designed tests and an exclusive, proprietary database. The use of any single test found in the RapID ONE System to establish the identification of a test isolate is subject to the error inherent in that test alone.

16. PERFORMANCE CHARACTERISTICS

The performance characteristics of RapID ONE System have been established by laboratory testing of reference and stock cultures at Remel, and by clinical evaluations using fresh clinical and stock isolates.¹⁰⁻¹³

17. BIBLIOGRAPHY

- Ewing, W.H. 1986. Edwards and Ewing's Identification of Enterobacteriaceae. 4th ed. Elsevier Science Publishing Company, Inc., New York, NY.
- Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy. 1991. Manual of Clinical Microbiology. 5th ed. ASM, Washington, D.C.
- MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Lippincott Williams & Wilkins, Philadelphia, PA.
- Holt J.G., N.R. Krieg, P.H.A. Sneath, J.T. Staley, and S.T. Williams. 2000. Bergey's Manual of Determinative Bacteriology. 9th ed. Lippincott Williams and Wilkins. Philadelphia, PA.
- Eriquez, L.A., N.E. Hodinka, and K.A. Hornsby. 1993. Abstract C-294. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.

18. PACKAGING

REF R8311006 RapID ONE System..... 20 Tests/Kit

19. SYMBOL LEGEND

REF	Catalogue Number
IVD	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Contains sufficient for <N> tests
	Do not use if package is damaged
	Do not re-use
LOT	Batch Code (Lot Number)
	Use By (Expiration Date)
	Importer
UDI	Unique Device Identifier
EC REP	Authorized representative in the European Community
UK CA	UK Conformity Assessed
CE	European Conformity Assessment
	Manufacturer

Rapid™ and ERIC™ are trademarks of Thermo Fisher Scientific and its subsidiaries.

ATCC® is a registered trademark of American Type Culture Collection.

Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, USA
www.thermofisher.com/microbiology
Tel: (800) 255-6730 • International: (913) 888-0939
www.oxoid.com/IFU
Europe +800 135 79 135 • US 1 855 2360 190
CA 1 855 805 8539 • ROW +31 20 794 7071

Version	Date of modifications introduced
IFU8311006	August 2023 Updated to meet IVDR requirements

Printed in the UK

Table 4 - RapID ONE Differential Chart

Organism	URE	ADH	ODC	LDC	TET	LIP	KSF	SBL	GUR	ONPG	βGLU	βXYL	NAG	MAL	PRO	GGT	PYR	ADON	IND	OXI
<i>Acinetobacter calcoaceticus</i>	5	29	11	7	2	94	0	1	0	0	0	0	2	21	3	5	2	0	0	0
<i>Burkholderia cepacia</i> ^a	14	12	31	84	3	95	2	0	0	14	9	0	19	16	71	95	0	0	0	62
<i>Cedecea davisae</i> (EG-15)	2	27	51	0	0	71	88	0	0	51	98	96	99	92	1	82	0	0	0	0
<i>Cedecea lapageli</i>	0	38	0	7	0	90	95	0	0	90	92	0	99	98	1	98	0	0	0	0
<i>Cedecea neteri</i>	2	46	2	5	0	93	98	95	0	98	98	0	99	98	0	99	0	0	0	0
<i>Cedecea</i> sp. 3	0	90	0	0	0	95	90	0	0	91	97	0	98	0	0	97	0	0	0	0
<i>Cedecea</i> sp. 5	0	91	35	2	0	81	90	98	0	90	98	0	98	0	0	95	0	0	0	0
<i>Citrobacter amalonaticus</i>	11	11	96	3	98	0	95	96	0	95	83	1	9	2	2	92	98	2	98	0
<i>Citrobacter freundii</i>	11	18	24	0	88	0	94	92	0	92	5	1	0	5	3	88	96	1	1	0
<i>Citrobacter koseri</i> ^b	5	20	95	2	0	95	95	95	0	95	92	0	1	90	4	95	96	80	96	0
<i>Cronobacter sakazakii</i> ^c	0	89	92	3	0	0	98	1	0	96	84	96	98	19	0	95	38	0	20	0
<i>Edwardsiella hoshinae</i>	0	0	79	98	11	0	8	0	0	0	0	2	90	69	99	5	0	0	78	0
<i>Edwardsiella tarda</i>	0	0	83	89	90	0	2	0	0	3	0	2	96	0	98	4	0	0	92	0
<i>Enterobacter asburiae</i> (EG-17)	7	6	89	0	0	0	98	94	0	98	94	30	97	2	0	82	9	0	0	0
<i>Enterobacter cancerogenus</i> ^d	0	91	95	1	0	0	98	0	0	98	84	2	98	98	0	98	96	0	0	0
<i>Enterobacter cloacae</i>	9	92	92	1	0	2	95	96	0	98	82	92	94	76	1	94	21	17	0	0
<i>Escherichia coli</i>	1	15	75	84	2	0	1	90	92	95	1	0	0	1	3	44	1	2	98	0
<i>Escherichia fergusonii</i>	0	5	88	96	0	0	96	0	0	86	2	0	0	14	2	90	98	82	94	0
<i>Escherichia hermannii</i> (EG-11)	0	2	87	2	98	0	98	3	0	85	67	0	1	0	4	96	92	0	96	0
<i>Ewingella americana</i> ^e	0	0	0	0	0	0	56	0	0	92	44	0	94	0	0	21	99	0	0	0
<i>Hafnia alvei</i>	0	13	92	96	5	1	87	0	0	31	7	0	77	20	97	96	11	0	2	0
<i>Klebsiella aerogenes</i> ^f	4	7	97	97	1	1	96	97	0	98	98	97	92	1	93	98	88	0	0	0
<i>Klebsiella oxytoca</i>	93	12	9	97	9	2	98	96	0	95	98	84	0	95	0	96	82			

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- Ewing, W.H. 1986. Edwards and Ewing's Identification of Enterobacteriaceae. 4th ed. Elsevier Science Publishing Company, Inc., New York, NY.
- Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadowy. 1991. Manual of Clinical Microbiology. 5th ed. ASM, Washington, D.C.
- MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Lippincott Williams & Wilkins, Philadelphia, PA.
- Holt J.G., N.R. Krieg, P.H.A. Sneath, J.T. Staley, and S.T. Williams. 2000. Bergey's Manual of Determinative Bacteriology. 9th ed. Lippincott Williams and Wilkins. Philadelphia, PA.
- Eriquez, L.A., N.E. Hodinka, and K.A. Hornsby. 1993. Abstract C-294. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Kilian, M. and P. Bulow. 1976. Acta Pathol. Microbiol. Scand. B. 84:245-251.
- Humble, W.M., A. King, and I. Phillips. 1977. J. Clin. Pathol. 30:275-277.
- K.C. Carroll, M.A. Pfaller, M.L. Landry, A.J. McAdam, R. Patel, S.S. Richter, D.W. Warnock. 2019. Manual of Clinical Microbiology. 12th Edition. ASM Press.
- Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott's Diagnostic Microbiology. 12th ed. Mosby Elsevier, St. Louis, MO.
- Jerris, R., A. Kabant, S. Peroulas, T. Lee, K.A. Hornsby, and D. Lockhart. 1993. Abstract C-304. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Kitch, T., M.R. Jacobs, and P.C. Appelbaum. 1993. Abstract C-303. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Schreckenberger, P., M. Montero, and N. Heldt. 1993. Abstract C-309. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Eriquez, L.A., A.P. Jones, and N.E. Hodinka. 1992. Abstract C-5. Abstracts of the 92nd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Clinical and Laboratory Standards Institute (CLSI). 2008. Quality Control for Commercial Microbial Identification Systems; Approved Guideline. M50-A. CLSI, Wayne, PA.

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- Kitch, T., M.R. Jacobs, and P.C. Appelbaum. 1993. Abstract C-303. *Abstracts of the 93rd General Meeting of the American Society for Microbiology*. ASM, Washington, D.C.

- Schreckenberger, P., M. Montero, and N. Heldt. 1993. Abstract C-309. *Abstracts of the 93rd General Meeting of the American Society for Microbiology*. ASM, Washington, D.C.

- Eriquez, L.A., A.P. Jones, and N.E. Hodinka. 1992. Abstract C-5. *Abstracts of the 92nd General Meeting of the American Society for Microbiology*. ASM, Washington, D.C.

- Clinical and Laboratory Standards Institute (CLSI). 2008. *Quality Control for Commercial Microbial Identification Systems*; Approved Guideline. M50-A. CLSI, Wayne, PA.

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ATCC™ je registrovaná ochranná známka sbírky American Type Culture Collection.

Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, USA

www.thermofisher.com/microbiology

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- Jerris, R. A. Kabant, S. Peroulas, T. Lee, K.A. Hornsby, and D. Lockhart. 1993. Abstract C-304. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Kitch, T., M.R. Jacobs, and P.C. Appelbaum. 1993. Abstract C-303. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Schreckenberger, P., M. Montero, and N. Heldt. 1993. Abstract C-309. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Eriuez, L.A., A.P. Jones, and N.E. Hodinka. 1992. Abstract C-5. Abstracts of the 92nd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Clinical and Laboratory Standards Institute (CLSI). 2008. Quality Control for Commercial Microbial Identification Systems; Approved Guideline. M50-A. CLSI, Wayne, PA.

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- Ewing, W.H. 1986. Edwards and Ewing's Identification of Enterobacteriaceae. 4th ed. Elsevier Science Publishing Company, Inc., New York, NY.
- Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy. 1991. Manual of Clinical Microbiology. 5th ed. ASM, Washington, D.C.
- MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Lippincott Williams & Wilkins, Philadelphia, PA.
- Holt J.G., N.R. Krieg, P.H.A. Sneath, J.T. Staley, and S.T. Williams. 2000. Bergey's Manual of Determinative Bacteriology. 9th ed. Lippincott Williams and Wilkins. Philadelphia, PA.
- Eriksen, L.A., N.E. Hodinka, and K.A. Hornsby. 1993. Abstract C-294. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Kilian, M. and P. Bulow. 1976. Acta Pathol. Microbiol. Scand. B. 84:245–251.
- Humble, W.M., A. King, and I. Phillips. 1977. J. Clin. Pathol. 30:275–277.
- K.C. Carroll, M.A. Pfaller, M.L. Landry, A.J. McAdam, R. Patel, S.S. Richter, D.W. Warnock. 2019. Manual of Clinical Microbiology. 12th Edition. ASM Press.

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- Jerris, R., A. Kabant, S. Peroulas, T. Lee, K.A. Hornsby, and D. Lockhart. 1993. Abstract C-304. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.

11. Kitch, T., M.R. Jacobs, and P.C. Appelbaum. 1993. Abstract C-303. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.

- Schreckenberger, P., M. Montero, and N. Heldt. 1993. Abstract C-309. Abstracts of the 93rd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.

13. Eriquez, L.A., A.P. Jones, and N.E. Hodinka. 1992. Abstract C-5. Abstracts of the 92nd General Meeting of the American Society for Microbiology. ASM, Washington, D.C.

- Clinical and Laboratory Standards Institute (CLSI, Klinikai és Laboratóriumi Minősítő Intézet). 2008. Quality Control for Commercial Microbial Identification Systems; Approved Guideline. M50-A. CLSI, Wayne, PA.

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17. BIBLIOGRAFIA

- Ewing, W.H. 1986. *Edwards and Ewing's Identification of Enterobacteriaceae*. 4th ed. Elsevier Science Publishing Company, Inc., New York, NY.
- Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy. 1991. *Manual of Clinical Microbiology*. 5th ed. ASM, Washington, D.C.
- MacFaddin, J.F. 2000. *Biochemical Tests for Identification of Medical Bacteria*. 3rd ed. Lippincott Williams & Wilkins, Philadelphia, PA.

4. Holt J.G., N.R. Krieg, P.H.A. Sneath, J.T. Staley, and S.T. Williams. 2000. *Bergey's Manual of Determinative Bacteriology*. 9th ed. Lippincott Williams and Wilkins. Philadelphia, PA.

- Eriquez, L.A., N.E. Hodinka, and K.A. Hornsby. 1993. Abstract C-294. *Abstracts of the 93rd General Meeting of the American Society for Microbiology*. ASM, Washington, D.C.
- Kilian, M. and P. Bulow. 1976. *Acta Pathol. Microbiol. Scand. B*. 84:245-251.
- Humble, W.M., A. King, and I. Phillips. 1977. *J. Clin. Pathol.* 30:275-277.
- K.C. Carroll, M.A. Pfaller, M.L. Landry, A.J. McAdam, R. Patel, S.S. Richter, D.W. Warnock. 2019. *Manual of Clinical Microbiology*. 12th Edition. ASM Press.
- Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. *Bailey and Scott's Diagnostic Microbiology*. 12th ed. Mosby Elsevier, St. Louis, MO.

10. Jerris, R., A. Kabant, S. Peroulas, T. Lee, K.A. Hornsby, and D. Lockhart. 1993. Abstract C-304. *Abstracts of the 93rd General Meeting of the American Society for Microbiology*. ASM, Washington, D.C.

- Kitch, T., M.R. Jacobs, and P.C. Appelbaum. 1993. Abstract C-303. *Abstracts of the 93rd General Meeting of the American Society for Microbiology*. ASM, Washington, D.C.
- Schreckenberger, P., M. Montero, and N. Heldt. 1993. Abstract C-309. *Abstracts of the 93rd General Meeting of the American Society for Microbiology*. ASM, Washington, D.C.

13. Eriquez, L.A., A.P. Jones, and N.E. Hodinka. 1992. Abstract C-5. *Abstracts of the 92nd General Meeting of the American Society for Microbiology*. ASM, Washington, D.C.

- Clinical and Laboratory Standards Institute (CLSI). 2008. *Quality Control for Commercial Microbial Identification Systems; Approved Guideline*. M50-A. CLSI, Wayne, PA.

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5. Eriequez, L.A., N.E. Hodinka, and K.A. Hornsby. 1993. Abstract C-294. *Abstracts of the 93rd General Meeting of the American Society for Microbiology*. ASM, Washington, D.C.

- Kilian, M. and P. Bulow. 1976. *Acta Pathol. Microbiol. Scand. B*. 84:245–251.
- Humble, W.M., A. King, and I. Phillips. 1977. *J. Clin. Pathol.* 30:275–277.
- K.C. Carroll, M.A. Pfaller, M.L. Landry, A.J. McAdam, R. Patel, S.S. Richter, D.W. Warnock. 2019. *Manual of Clinical Microbiology*. 12th Edition. ASM Press.
- Forbes, B.A., D.F. Sahn, and A.S. Weissfeld. 2007. *Bailey and Scott's Diagnostic Microbiology*. 12th ed. Mosby Elsevier, St. Louis, MO.
- Jerris, R., A. Kabant, S. Peroulas, T. Lee, K.A. Hornsby, and D. Lockhart. 1993. Abstract C-304. *Abstracts of the 93rd General Meeting of the American Society for Microbiology*. ASM, Washington, D.C.
- Kitch, T., M.R. Jacobs, and P.C. Appelbaum. 1993. Abstract C-303. *Abstracts of the 93rd General Meeting of the American Society for Microbiology*. ASM, Washington, D.C.
- Schreckenberger, P., M. Montero, and N. Heldt. 1993. Abstract C-309. *Abstracts of the 93rd General Meeting of the American Society for Microbiology*. ASM, Washington, D.C.
- Eriequez, L.A., A.P. Jones, and N.E. Hodinka. 1992. Abstract C-5. *Abstracts of the 92nd General Meeting of the American Society for Microbiology*. ASM, Washington, D.C.
- Clinical and Laboratory Standards Institute (CLSI). 2008. *Quality Control for Commercial Microbial Identification Systems; Approved Guideline*. M50-A. CLSI, Wayne, PA.

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