Urea 40% Solution

EN

REF SR0020K

* This instructions for use (IFU) document is intended to be read in conjunction with the IFU for Thermo Scientific™ Oxoid™ Urea Agar Base (Dehydrated) (product code: CM0053B) and Thermo Scientific™ Oxoid™ Urea Broth Base (Dehydrated) (product code: CM0071B)

Intended Use

Thermo Scientific™ Urea 40% Solution (SR0020K) is intended to be added to Thermo Scientific™ Oxoid™ Urea Agar Base (Dehydrated) (CM0053B) for the preparation of Christensen's medium. Christensen's medium is intended to be used to detect rapid urease activity of the *Proteae* and non-rapid urease activity of some *Enterobacteriaceae* in faecal and urine samples. Urea 40% Solution (SR0020K) is also intended to be used in the preparation of Thermo Scientific™ Oxoid™ Urea Broth Base (Dehydrated) (CM0071B), a liquid version of Christensen's medium intended for the differentiation of urease-producing *Enterobacteriaceae* from members of the *Salmonella* and *Shigella* groups. Urea 40% Solution (SR0020K) used in the preparation of Christensen's media, is intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having *Enterobacteriaceae* infections. The device is for professional use only, is not automated and nor is it a companion diagnostic.

Summary and Explanation

The urease test is used to differentiate urease-positive organisms (e.g. *Proteus*) from other organisms and to differentiate between *Candida albicans* and *Cryptococcus neoformans*¹. The urease test is used to determine the ability of an organism to split urea through the production of the enzyme urease¹. Quality control organisms to be used as a positive control are *Proteus mirabilis* National Collection of Type Cultures (NCTC) NCTC 10975 and *Cryptococcus neoformans* National Collection of Pathogenic Fungi (NCPF) NCPF 3168, organisms to be used as negative control are *Escherichia coli* NCTC 10418 or NCTC 12241 and *Candida albicans* NCPF 3281¹. Due to rapid urease production the test is used for the presumptive identification of *Cryptococcus neoformans* and *Cryptococcus gattii* species¹. Some strains of *Enterobacter* and *Klebsiella* species and other species of *Cryptococcus* can give a positive result for the urease test¹. *Candida albicans* does not produce urease, therefore giving a negative urease result¹.

Principle of Method

Urea 40% Solution (SR0020K) is a solution intended to be used in the preparation of Urea Agar Base (Dehydrated) (CM0053B) and Urea Broth Base (Dehydrated) (CM0071B). When urea is hydrolysed by urease positive organisms, ammonia is produced making the reaction alkaline, Urea Agar Base (Dehydrated) (CM0053B) turns pink/red colour and Urea Broth Base (Dehydrated) (CM0071B) turns pink/red colour.

Typical Formula

	Milligrams per 5 ml vial
Urea Purified	2000 mg

Materials Provided

SR0020K: 10x vials each for 95 ml of medium

Materials Required but Not Supplied

1. Urea Agar Base (Dehydrated) product code CM0053B
2. Urea Broth Base (Dehydrated) product code CM0071B
3. Inoculating loops
4. Swabs
5. Collection containers
6. Incubators
7. Quality control organisms

Storage

- Store product in its original packaging between 2°C and 8°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow product to equilibrate to room temperature before use.

Warnings and Precautions

- Each vial is for single use. Do not re-use.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or vials.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

- Supplements may be dispensed on the bench, wear disposable gloves and safety glasses.
- Disposable gloves must be worn for rehydrating, dispensing and disposing of ALL Supplements (i.e. wear disposable gloves whilst handling supplements).
- All waste must be emptied into the HAZARDOUS WASTE container situated in the Microbiology laboratory.
- The vial and rubber bung must be rinsed into the HAZARDOUS WASTE container.
- The metal lid and rinsed rubber bung should be placed in the waste bin.
- The rinsed glass vial should then be placed in the glass bin.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) Q5.

Procedure


1. Allow supplement to equilibrate to room temperature before use.
2. Suspend 2.4 g of Urea Agar Base (Dehydrated) (CM0053B) or 0.9 g of Urea Broth Base (Dehydrated) (CM0071B) in 95 ml of distilled water. Bring to the boil to dissolve completely.
3. Sterilise by autoclaving at 115°C for 20 minutes.
4. Cool to 50°C for CM0053B and 55°C for CM0071B and aseptically add one vial of sterile Urea 40% Solution (SR0020K).
5. Mix well. Aseptically dispense 10 ml amounts into sterile containers and allow to set in the slope position for CM0053B.

Bibliography


1. Public Health England. 2019. 'Urease Test.' UK Standards for Microbiology Investigations. UK Standards for Microbiology Investigations TP 36 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/794218/TP_36i4.pdf

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union

	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

For technical assistance please contact your local distributor.

Revision information

Version	Date of issue and modifications introduced
2.0	2024-05-20

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Storbritannien