www.thermofisher.com

C-N Supplement

REF SR0102E

EN

This instructions for use (IFU) document is intended to be read in conjunction with the IFU for Pseudomonas Agar Base (Dehydrated). Product codes CM0559B/R/T.

Intended Use **IVD**

C-N Supplement (SR0102E) is a selective supplement that is intended to be used in the preparation of Pseudomonas Agar Base (Dehydrated) (CM0559B, CM0559R and CM0559T) for the isolation of *Pseudomonas* spp. from a wide range of clinical samples including urine, wound swabs, and sputum. The devices can also be used for testing food, water, and environmental samples. C-N Supplement (SR0102E) is recommended for the selective isolation of *Pseudomonas aeruginosa*.

C-N Supplement (SR0102E) is used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having Pseudomonas infections.

The device is for professional use only, is not automated nor is it a companion diagnostic.

Summary and Explanation

Pseudomonas species are motile, gram negative bacteria that are straight or curved in appearance¹. *Pseudomonas* species can cause a variety of infections in humans, ranging from lungs (pneumonia), wound infections, septicaemia, urinary tract and gastrointestinal infections².

Pseudomonas aeruginosa is a non-fermentive, strictly aerobic, catalase and oxidase positive gram-negative rod. It is widely distributed in soil and water and may cause disease in animals, plants and insects. *P. aeruginosa* is known as the key respiratory pathogen in patients with cystic fibrosis³. In addition to it, people with, chronic obstructive pulmonary disease (COPD), pregnant ladies and people with a weakened immune system are most affected by the *Pseudomonas* infection². *P. aeruginosa* are known to cause a condition known as pyelonephritis particularly in patients with catheters⁴. In general, infections caused by *Pseudomonas aeruginosa* are treated with antibiotics. However, with an increase in the number of multi drug resistant strains of *Pseudomonas aeruginosa*, presumptive identification of the strains is highly important. For the identification of best antibiotics to treat the infection clinical specimens are sent to a laboratory and are grown against a set of antibiotics like fluoroquinolones².

Principle of Method

C-N Supplement (SR0102E) when used with the Pseudomonas Agar Base (Dehydrated) (CM0559B, CM0559R and CM0559T), is used for the isolation of *P. aeruginosa*, which is achieved through the inclusion of cetrimide and sodium nalidixate to inhibit the growth of competing flora.

Typical Formula

	Milligrams per 500ml
Cetrimide	100.0mg
Sodium nalidixate	7.5mg

Materials Provided

SR0102E: 10x freeze dried vials each for 500ml of medium.

Materials Required but Not Supplied

- Pseudomonas Agar Base (dehydrated) product code CM0559B/R/T
- Inoculating loops
- Swabs
- Collection containers
- Incubators
- Quality control organisms
- Petri dishes

Storage

- Store product in its original packaging between 2°C and 8°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Warnings and Precautions

Signal: Word Danger

Hazard Statements

- H302 - Harmful if swallowed
- H315 - Causes skin irritation
- H317 - May cause an allergic skin reaction
- H318 - Causes serious eye damage
- H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled
- H335 - May cause respiratory irritation
- H373 - May cause damage to organs through prolonged or repeated exposure
- H400 - Very toxic to aquatic life

Precautionary Statements

- P302 + P352 - IF ON SKIN: Wash with plenty of soap and water
- P362 - Take off contaminated clothing and wash before reuse
- P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
- P310 - Immediately call a POISON CENTER or doctor/physician
- P304 + P340 - IF INHALED: Remove person to fresh air and keep comfortable for breathing
- P280 - Wear protective gloves/protective clothing/eye protection/face protection

Each vial is for single use. Do not re-use.

For in vitro diagnostic use only.

For professional use only.

Inspect the product packaging before first use.

Do not use the product if there is any visible damage to the packaging or vials.

Do not use the product beyond the stated expiry date.

Do not use the device if signs of contamination are present.

Do not use the device if the colour has changed or there are other signs of deterioration.

It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local the recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 17, B 11, B 41 and B 57.

Procedure

- Allow supplement to equilibrate to room temperature before use.
- Aseptically add 2ml of 1:1 ethanol:sterile distilled water to 1 vial and mix gently to dissolve.
- Aseptically add the vial contents to 500ml of sterile Pseudomonas Basal Medium, Prepared from Pseudomonas Agar Base (CM0559B/R/T) plus 5ml glycerol and cooled to 50°C.
- Mix well and aseptically dispense into suitable containers (i.e. petri dishes).

Bibliography

- 1) Public Health England 'Identification of Pseudomonas species and other NonGlucose Fermenters.' ID17/3 24.04.2015.
- 2) Centers for Disease control and Prevention (CDC). 2019. 'Pseudomonas aeruginosa in Healthcare Settings.
- 3) Langan, K. M., Kotsimbos, T., & Peleg, A. Y. (2015). Managing Pseudomonas aeruginosa respiratory infections in cystic fibrosis. *Current Opinion in Infectious Diseases*, 28(6), 547–556.
- 4) Newman, J. W., Floyd, R. v. & Fothergill, J. L. (2017). The contribution of Pseudomonas aeruginosa virulence factors and host factors in the establishment of urinary tract infections. *FEMS Microbiology Letters*, 364(15).

Symbol Legend

Symbol	Definition
REF	Catalogue number

	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2023 Thermo Fisher Scientific Inc. All rights reserved. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Modifications introduced
3.0	2024-05-24

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©2023 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.