www.thermofisher.com

Modified Karmali Selective Supplement

[REF] SR0205E

[EN]

This instructions for use (IFU) document is intended to be read in conjunction with the IFU for *Campylobacter* Agar Base, product code: CM0935B.

Intended Use

Modified Karmali Selective Supplement (SR0205E) is used in the preparation of *Campylobacter* (Karmali) medium for the isolation of *Campylobacter* species from faecal samples. These devices are intended for in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having *Campylobacter* infections. Can also be used for testing food samples.

The device is for professional use only, is not automated and nor is it a companion diagnostic.

Summary and Explanation

Campylobacter infections are among the most common bacterial infections, producing both gastrointestinal and systemic illnesses. *Campylobacter jejuni* is the most commonly reported bacterial cause of foodborne infection in the United States. *Campylobacter* are small, non-spore forming, helical, curved gram-negative rods that require a microaerobic environment and cannot tolerate drying. *Campylobacter* Karmali Selective Medium is based on the formulation described by Karmali et al¹.

Principle of Method

The original *Campylobacter* Blood Free medium in the Oxoid product range contains sodium pyruvate in the agar base. *Campylobacter* medium (Karmali) incorporates this ingredient into the selective supplement. The original medium also contains sodium deoxycholate for the inhibition of Gram positive organisms, whereas, with *Campylobacter* Medium (Karmali) suppression of Gram positives is achieved by the inclusion of Vancomycin. Modified Karmali Selective Supplement (SR0205E) offers a cycloheximide free alternative by using amphotericin B as an antifungal agent.

Typical Formula

	mg/500ml
Sodium pyruvate	50.0mg
Cefoperazone	16.0mg
Vancomycin	10.0mg
Amphotericin B	5.0mg

Materials Provided

SR0205E: 10x freeze dried vials, each for supplementing 500ml of medium.

Materials Required but Not Supplied

- CM0935B: 500g of Karmali Agar Base.
- Inoculating loops, swabs, collection containers.

- Petri dishes.

Storage

- Store product in its original packaging between 2°C and 8°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Warnings and Precautions

Signal Word: Danger

Hazard Statements

H317 - May cause an allergic skin reaction
H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled

Precautionary Statements

P280 - Wear protective gloves/protective clothing
P302 + P352 - IF ON SKIN: Wash with plenty of soap and water
P333 + P313 - If skin irritation or rash occurs: Get medical advice/attention
P285 - In case of inadequate ventilation wear respiratory protection
P342 + P311 - If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician
P261 - Avoid breathing dust
P304 + P340 - IF INHALED: Remove person to fresh air and keep comfortable for breathing

Each vial is for single use. Do not re-use.

For in vitro diagnostic use only.

For professional use only.

Inspect the product packaging before first use.

Do not use the product if there is any visible damage to the packaging or vials.

Do not use the product beyond the stated expiry date.

Do not use the device if signs of contamination are present.

Do not use the device if the colour has changed or there are other signs of deterioration.

It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

There are no specimen collection and/or transport materials provided with this device. Specimens should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI).

Procedure

Aseptically add 2ml of 1:1 ethanol:sterile distilled water to 1 vial. Mix gently to dissolve. Aseptically add the vial contents to 500ml of sterile Karmali Agar Base (CM0935B) prepared as directed and cooled to 50°C. Mix well and pour into sterile Petri dishes. SR0205E is only suitable for use with Karmali Agar Base (CM0935B).

Bibliography

1. Karmali M.A., Simor A.E., Roscoe M., Fleming P.C, Smith S.S. and Lane J. (1986) J.Clin.Micro. 23. 456-459.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment

	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Modifications introduced
2.0	2024-01-12

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Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, Regatul Unit