	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Cefixime Discs 5µg (CFM5)

REF CT0653B

Antimicrobial Susceptibility Test Discs

Note: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Cefixime Antimicrobial Susceptibility Test (AST) Discs are 6mm paper discs which contain specific amounts of the antimicrobial agent cefixime. The discs are labelled on both sides with details of the antimicrobial (CFM) and amount present (µg): CFM5 (5µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Cefixime AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which cefixime has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate, or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling, and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100^{1ac}, or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- N/A

Gram-Negative

- Enterobacteriaceae
- *Haemophilus influenzae* and *parainfluenzae*
- *Neisseria gonorrhoeae*

Species with published breakpoints according to current CLSI literature:

Gram-Positive

- N/A

Gram-Negative

- Enterobacteriaceae
- *Haemophilus influenzae* and *parainfluenzae*
- *Neisseria gonorrhoeae*

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- N/A

Gram-Negative

- Enterobacteriaceae (Uncomplicated UTI only)
- *Haemophilus influenzae*
- *Moraxella catarrhalis*

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Cefixime AST Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology, refer to the

relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Cefixime AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Cefixime AST Discs consist of 6mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Cefixime AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with CT0653B		Table 2. Description of Cefixime AST Disc Reagent	
Component Description	Material Description	Reagent	Description of Function
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.	Cefixime	Off white to pale yellow powder. Cefixime is a semisynthetic antibiotic that belongs to the third generation of cephalosporins. It acts by binding specific penicillin-binding proteins located on the bacterial cell wall, thus, inhibiting the mucopeptide synthesis in the third and last stage of the bacterial cell wall's production. Cefixime is also highly stable in presence of beta-lactamase enzymes.
Desiccant tablet (x5)	Pale beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.		
Foil	Foil individually sealing each cartridge with its desiccant.		
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.		

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Cefixime AST Discs can be dispensed using a Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Cefixime AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Table 3. Raw data of zone sizes, in accordance with CLSI methodology, taken from recent batches of Cefixime (5 µg) CFM5 Antimicrobial Susceptibility Test Discs (CT0653B).

CLSI Methodology																	
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)										Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	4	5	6	7	8	9	10	
3326512	<i>Escherichia coli</i>	25922	MHA	20	26	23	20	20	20	N/A	N/A	N/A	N/A	N/A	N/A	N/A	-3 ± 0 (CV=0)
	<i>Neisseria gonorrhoeae</i>	49226	GCV	37	45	41	45	44	44	N/A	N/A	N/A	N/A	N/A	N/A	N/A	3.3 ± 0.6 (CV=0.2)
	<i>Haemophilus influenzae</i>	49247	HTM	25	33	29	25	26	26	26	25	25	26	26	26	26	25.7 ± 0.5 (CV=0.02)
	<i>Streptococcus pneumoniae</i>	49619	MHA+5 %SB	16	23	19.5	23	23	23	N/A	N/A	N/A	N/A	N/A	N/A	N/A	3.5 ± 0 (CV=0)
Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)									Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)	
							1			2			3				
3315578	<i>Escherichia coli</i>	25922	MHA	20	26	23	20			20			20			-3 ± 0 (CV=0)	
	<i>Neisseria gonorrhoeae</i>	49226	GCV	37	45	41	40			40			41			-0.7 ± 0.6 (CV= -0.9)	
	<i>Haemophilus influenzae</i>	49247	HTM	25	33	29	29			29			28			-0.3 ± 0.6 (CV= -1.7)	
	<i>Streptococcus pneumoniae</i>	49619	MHA+5 %SB	16	23	19.5	21			21			21			1.5 ± 0 (CV=0)	
Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)									Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)	
							1			2			3				
3232843	<i>Escherichia coli</i>	25922	MHA	20	26	23	21			21			21			-2 ± 0 (CV=0)	
	<i>Neisseria gonorrhoeae</i>	49226	GCV	37	45	41	37			38			39			-3 ± 1 (CV= -0.3)	
	<i>Haemophilus influenzae</i>	49247	HTM	25	33	29	28			28			28			-1 ± 0 (CV=0)	
	<i>Streptococcus pneumoniae</i>	49619	MHA+5 %SB	16	23	19.5	19			19			19			-0.5 ± 0 (CV=0)	

Table 4. Raw data of zone sizes, in accordance with EUCAST methodology, taken from recent batches of Cefixime (5 µg) CFM5 Antimicrobial Susceptibility Test Discs (CT0653B).

EUCAST Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3326512	<i>Escherichia coli</i>	25922	MHA	20	26	23	20	20	20	-3 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49766	MHA-F	29	35	32	30	30	30	-2 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3315578	<i>Escherichia coli</i>	25922	MHA	20	26	23	20	20	20	-3 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49766	MHA-F	29	35	32	29	29	29	-3 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3232843	<i>Escherichia coli</i>	25922	MHA	20	26	23	21	21	21	-2 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49766	MHA-F	29	35	32	33	33	34	1.3 ± 0.6 (CV=0.4)

Clinical Performance Characteristics

A study was conducted to compare antimicrobial susceptibility patterns of *Typhoidal Salmonella* isolates using antibiotic drugs prescribed for conventional treatment, newer drugs such as fluoroquinolones, and third generation cephalosporins⁴. 200 samples were collected, then isolated and tested for *Typhoidal Salmonella*. Once positivity was confirmed (as either *Salmonella typhi* or *Salmonella paratyphi A*), the isolates were assessed for antimicrobial susceptibility testing using disc diffusion through the Kirby Bauer method and according to CLSI guidelines. Although a significant percentage of the *S. paratyphi A* isolates were sensitive to drugs conventionally used for treatment against typhoidal fever such as ampicillin, chloramphenicol, and co-trimoxazole, the sensitivity results for *S. typhi* isolates against the same antibiotics were not as encouraging. For the fluoroquinolones, the susceptibility of the isolates was quite low for both species, while all 200 isolates tested showed 100% susceptibility against the third generation of cephalosporins including cefixime, thus demonstrating the effectiveness of the Oxoid Cefixime (5 µg) AST Discs as an excellent tool for *in vitro* analysis of *S. typhi*.

A study held in 2020 assessed the risk factors and the prevalence of STIs with a specific focus on *Neisseria gonorrhoeae*, *Trichomonas vaginalis* and *Treponema pallidum*.⁵ The antimicrobial susceptibility patterns for *N. gonorrhoeae* were also

determined. Three samples were collected from 189 students. The first sample from each student was sent to a laboratory to isolate and identify *N. gonorrhoeae*. Once confirmed, antimicrobial susceptibility testing was performed through disc diffusion against six antibiotic discs including Cefixime (5 µg) CFM5; the results were interpreted against CLSI guidelines. The second sample collected was smeared for Gram stain identification, while the third sample was processed for the presence of *T. vaginalis* and/or *T. pallidum* using saline wet mount preparation and a commercially available screening test kit. Out of the 189 subjects participating in the study, 27 were confirmed positive for one or more of the infections under testing, corresponding to 14% of the total. For *N. gonorrhoeae* isolates, the antimicrobial susceptibility patterns demonstrated that 85.7% of the strains were sensitive to ceftriaxone and cefixime, while 78.6% were susceptible to cefoxitin. The susceptibility for ciprofloxacin was lower (50%), while all *N. gonorrhoeae* were resistant to penicillin and tetracycline. The discs proved to be a valuable tool showing that the resistance of *N. gonorrhoeae* against the third generation of cephalosporins, such as cefixime, is slowly increasing, highlighting the importance of continuous surveillance for antimicrobial resistance with this type of infection.

In 2020, a study was also conducted to examine the effectiveness of conventional antibiotics in *Salmonella typhi*.⁶ A total of 100 blood samples were isolated from patients, then tested for the presence of *S. Typhi*. Once *S. typhi* was confirmed, antimicrobial susceptibility testing was performed through disc diffusion and results interpreted according to CLSI guidelines. Three groups of antibiotics were tested – first choice treatments, second choice drugs (fluoroquinolones), and third generation cephalosporins including cefixime. Cefixime (5µg) CFM5 AST Discs clearly demonstrated the effectiveness of the third generation cephalosporins for the treatment of non-resistant typhoidal fever cases. The results of this study proved that the antibiotics, which are part of the third group of cephalosporins, are a valid alternative for the cure of the disease, especially in children whom constitute the vast majority of the patients.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established. In the event of malfunction do not use device.

References

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Shujat, Umer, Aamer Ikram, Inam Qadir Javaid Hashmi, Shahid Ahmed Abbasi, Amna Afzal, and Muhammad Ayyub. 2016. Current Antimicrobial Sensitivity Pattern of Typhoidal Salmonellae in a Referral Diagnostic Centre. *Microbiologia Medica* 31 (1). <https://doi.org/10.4081/MM.2016.4669>.
5. Tamrat, Rahel, Tesfaye Kasa, Zewdineh Sahilemariam, and Mulatu Gashaw. 2020. Prevalence and Factors Associated with Sexually Transmitted Infections among Jimma University Students, Southwest Ethiopia. *International Journal of Microbiology* 2020. <https://doi.org/10.1155/2020/8859468>.
6. Khan, Aftab Ahmad, Rizwan Uppal, Gul E. Rehan, Farhat Khurshid, Khurshid Ahmad, and Hina Zaib. 2020. "Effectiveness of Conventional Antibiotics in *Salmonella Typhi* Positive Blood Cultures." *Journal of Medical Sciences (Peshawar)* 28 (4): 318–22. <https://doi.org/10.52764/JMS.20.28.4.2>.

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of Symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision Information

Version	Date of issue and modifications introduced
1.0	2023-01-13

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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Referencije

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Shujat, Umer, Aamer Ikram, Inam Qadir Javaid Hashmi, Shahid Ahmed Abbasi, Amna Afzal, and Muhammad Ayyub. 2016. Current Antimicrobial Sensitivity Pattern of Typhoidal *Salmonellae* in a Referral Diagnostic Centre. *Microbiologia Medica* 31 (1). <https://doi.org/10.4081/MM.2016.4669>.
5. Tamrat, Rahel, Tesfaye Kasa, Zewdineh Sahilemariam, and Mulatu Gashaw. 2020. Prevalence and Factors Associated with Sexually Transmitted Infections among Jimma University Students, Southwest Ethiopia. *International Journal of Microbiology* 2020. <https://doi.org/10.1155/2020/8859468>.
6. Khan, Aftab Ahmad, Rizwan Uppal, Gul E. Rehan, Farhat Khurshid, Khurshid Ahmad, and Hina Zaib. 2020. "Effectiveness of Conventional Antibiotics in *Salmonella Typhi* Positive Blood Cultures." *Journal of Medical Sciences (Peshawar)* 28 (4): 318–22. <https://doi.org/10.52764/JMS.20.28.4.2>.

© 2023. Thermo Fisher Scientific Inc. Sva prava pridržana. ATCC® je zaštitni znak društva ATCC. Svi ostali zaštitni znakovi vlasništvo su društva Thermo Fisher Scientific Inc. i njegovih podružnica. Ove informacije nisu namijenjene za poticanje na uporabu tih proizvoda na način kojim se mogu povrijediti prava intelektualnog vlasništva drugih društava.

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Reference

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Shujat, Umer, Aamer Ikram, Inam Qadir Javaid Hashmi, Shahid Ahmed Abbasi, Amna Afzal, and Muhammad Ayyub. 2016. Current Antimicrobial Sensitivity Pattern of Typhoidal Salmonellae in a Referral Diagnostic Centre. *Microbiologia Medica* 31 (1). <https://doi.org/10.4081/MM.2016.4669>.

5. Tamrat, Rahel, Tesfaye Kasa, Zewdineh Sahilemariam, and Mulatu Gashaw. 2020. Prevalence and Factors Associated with Sexually Transmitted Infections among Jimma University Students, Southwest Ethiopia. *International Journal of Microbiology* 2020. <https://doi.org/10.1155/2020/8859468>.
6. Khan, Aftab Ahmad, Rizwan Uppal, Gul E. Rehan, Farhat Khurshid, Khurshid Ahmad, and Hina Zaib. 2020. "Effectiveness of Conventional Antibiotics in *Salmonella Typhi* Positive Blood Cultures." *Journal of Medical Sciences (Peshawar)* 28 (4): 318–22. <https://doi.org/10.52764/JMS.20.28.4.2>.

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- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Shujat, Umer, Aamer Ikram, Inam Qadir Javaid Hashmi, Shahid Ahmed Abbasi, Amna Afzal, and Muhammad Ayyub. 2016. Current Antimicrobial Sensitivity Pattern of Typhoidal Salmonellae in a Referral Diagnostic Centre. *Microbiologia Medica* 31 (1). <https://doi.org/10.4081/MM.2016.4669>.
- Tamrat, Rahel, Tesfaye Kasa, Zewdineh Sahilemariam, and Mulatu Gashaw. 2020. Prevalence and Factors Associated with Sexually Transmitted Infections among Jimma University Students, Southwest Ethiopia. *International Journal of Microbiology* 2020. <https://doi.org/10.1155/2020/8859468>.
- Khan, Aftab Ahmad, Rizwan Uppal, Gul E. Rehan, Farhat Khurshid, Khurshid Ahmad, and Hina Zaib. 2020. "Effectiveness of Conventional Antibiotics in *Salmonella Typhi* Positive Blood Cultures." *Journal of Medical Sciences (Peshawar)* 28 (4): 318–22. <https://doi.org/10.52764/JMS.20.28.4.2>.

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Referencias

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Shujat, Umer, Aamer Ikram, Inam Qadir Javaid Hashmi, Shahid Ahmed Abbasi, Amna Afzal, and Muhammad Ayyub. 2016. Current Antimicrobial Sensitivity Pattern of Typhoidal Salmonellae in a Referral Diagnostic Centre. *Microbiology Medica* 31 (1). <https://doi.org/10.4081/MM.2016.4669>.
5. Tamrat, Rahel, Tesfaye Kasa, Zewdineh Sahilemariam, and Mulatu Gashaw. 2020. Prevalence and Factors Associated with Sexually Transmitted Infections among Jimma University Students, Southwest Ethiopia. *International Journal of Microbiology* 2020. <https://doi.org/10.1155/2020/8859468>.
6. Khan, Aftab Ahmad, Rizwan Uppal, Gul E. Rehan, Farhat Khurshid, Khurshid Ahmad, and Hina Zaib. 2020. "Effectiveness of Conventional Antibiotics in *Salmonella Typhi* Positive Blood Cultures." *Journal of Medical Sciences (Peshawar)* 28 (4): 318–22. <https://doi.org/10.52764/JMS.20.28.4.2>.

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