EN

Oxoid™ Ceftriaxone Discs 30 µg (CRO30)

REFCT0417B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the Oxoid AST generic instructions for use supplied with the product and available online.

Ceftriaxone Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent ceftriaxone. The discs are labelled on both sides with details of the antimicrobial (CRO) and amount present (µg): CRO30 (30µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Ceftriaxone AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which ceftriaxone has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling, and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA, CLSI M100,^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- *Streptococcus* spp β- Hemolytic Group
- *Streptococcus* spp Viridans Group

Gram-Negative

- *Enterobacteriaceae*
- *Haemophilus influenzae* and *parainfluenzae*
- *Neisseria gonorrhoeae*
- *Neisseria meningitidis*

Species with published breakpoints according to current CLSI M100 literature:

Gram-Positive

- *Streptococcus* spp. β-Hemolytic Group
- *Streptococcus* spp. Viridans Group

Gram-Negative

- *Enterobacterales*
- *Acinetobacter* spp.
- *Haemophilus influenzae* and *Haemophilus parainfluenzae*
- *Neisseria gonorrhoeae*
- *Neisseria meningitidis*

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- Viridans group *Streptococci*

Gram-Negative

- *Enterobacterales*
- *Haemophilus influenzae*
- *Moraxella catarrhalis*
- *Kingella kingae*

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Ceftriaxone AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Ceftriaxone AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Ceftriaxone AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Ceftriaxone AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to table 2.

Table 1. Materials Provided with CT0417B	
Component Description	Material Description
Cartridge with spring, cap and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Pale beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Ceftriaxone AST Disc Reagent	
Reagent	Description of Function
Ceftriaxone	White to yellow powder that is soluble in water. Ceftriaxone is a β -lactam antimicrobial agent that causes bacterial cell lysis through cell wall inhibition. Ceftriaxone achieves this through binding to, and subsequently inactivating, penicillin-binding proteins found on the cell membrane.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Ceftriaxone AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Ceftriaxone AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant, or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Results Summary for Tested QC Isolates in Accordance with CLSI Methodology.

CLSI Methodology										
Product Batch	Organism	ATCC® number™	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3203453	<i>Staphylococcus aureus</i>	ATCC® 25923™	MHA	22	28	25	26	26	26	1 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC® 25922	MHA	29	35	32	29	29	29	-3 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853™	MHA	17	23	20	23	23	23	3 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49247™	HTM	31	39	35	33	34	34	-1.34 ± 0.58 (CV=0.43)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619™	MHA+ 5%SB	30	35	32.5	35	35	35	2.5 ± 0 (CV=0)
	<i>Neisseria gonorrhoeae</i>	ATCC® 49226	GCV	39	51	45	44	44	44	-1 ± 0 (CV=0)
Product Batch	Organism	ATCC® number™	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3179509	<i>Staphylococcus aureus</i>	ATCC® 25923™	MHA	22	28	25	27	27	27	2 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC® 25922	MHA	29	35	32	30	30	30	-2 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853™	MHA	17	23	20	21	21	21	1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49247™	HTM	31	39	35	31	31	31	-4 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619™	MHA+ 5%SB	30	35	32.5	32	32	32	-0.5 ± 0 (CV=0)
	<i>Neisseria gonorrhoeae</i>	ATCC® 49226	GCV	39	51	45	41	41	41	-4 ± 0 (CV=0)
Product Batch	Organism	ATCC® number™	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3164331	<i>Staphylococcus aureus</i>	ATCC® 25923™	MHA	22	28	25	25	25	25	0 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC® 25922	MHA	29	35	32	34	34	34	2 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853™	MHA	17	23	20	21	21	21	1 ± 0 (CV=0)

	<i>Haemophilus influenzae</i>	ATCC® 49247™	HTM	31	39	35	34	34	34	-1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619™	MHA+ 5%SB	30	35	32.5	33	33	33	0.5 ± 0 (CV=0)
	<i>Neisseria gonorrhoeae</i>	ATCC® 49226	GCV	39	51	45	43	43	43	-2 ± 0 (CV=0)

Table 4. Raw Data of Zone Sizes in Accordance with EUCAST Methodology.

EUCAST Methodology										
Product Batch	Organism	ATCC® number™	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3203453	<i>Escherichia coli</i>	ATCC® 25922™	MHA	29	35	32	29	29	29	-3 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49766™	MHA-F	33	41	37	38	38	38	1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619™	MHA-F	32	38	35	37	37	37	2 ± 0 (CV=0)
	<i>Klebsiella pneumoniae</i>	ATCC® 700603™	MHA	16	22	19	17	17	17	-2 ± 0 (CV=0)
Product Batch	Organism	ATCC® number™	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3179509	<i>Escherichia coli</i>	ATCC® 25922™	MHA	29	35	32	30	30	30	-2 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49766™	MHA-F	33	41	37	40	40	40	3 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619™	MHA-F	32	38	35	35	35	35	0 ± 0 (CV=0)
	<i>Klebsiella pneumoniae</i>	ATCC® 700603™	MHA	16	22	19	19	19	19	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® number™	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3164331	<i>Escherichia coli</i>	ATCC® 25922™	MHA	29	35	32	34	34	34	2 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49766™	MHA-F	33	41	37	40	40	40	3 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619™	MHA-F	32	38	35	35	35	35	0 ± 0 (CV=0)
	<i>Klebsiella pneumoniae</i>	ATCC® 700603™	MHA	16	22	19	16	16	16	-3 ± 0 (CV=0)

Clinical Performance Characteristics

A study was designed to determine the prevalence and antimicrobial susceptibility of the *Acinetobacter baumannii* infection using the disc diffusion and E-test methods in intensive care patients.⁴

Of the 521 clinical samples collected from the patients, 95 (18.2%) were present for *A. baumannii*. Of the 95 clinical isolates, 76 (80%) demonstrated sensitivity to ceftriaxone, 12 (12.6%) exhibited intermediate sensitivity, and 7 (7.4%) were resistant. The clinical isolates demonstrated 100% susceptibility (sensitive and intermediate) to antimicrobial agents such as colistin, cefepime, imipenem and tazocin. In contrast, antimicrobial resistance was found at a high level to amikacin and minocycline with 75.7% and 82.1% resistance, respectively.

The low level of ceftriaxone resistance in *A. baumannii* shown in the study is in keeping with previous studies. The study proves the importance of comprehending prevalence and the antimicrobial susceptibility pattern of *A. baumannii* amongst hospital patients to avoid widespread drug resistance and to improve treatment. The study did not report any consistent or out of specification results for the Oxoid Ceftriaxone AST Discs and demonstrates that the discs are an effective method for the in vitro analysis of *A. baumannii*.

A study was conducted to determine the bacterial profile of non-chlamydia conjunctival bacteria in patients who then go on to have trachomatous trichiasis surgery.⁵ The antimicrobial susceptibility of the determined bacteria to various routinely used antibiotics was then found using the disc diffusion method.

Of the total number of enrolled patients, 438 (31%) presented with non-chlamydia pathogenic bacteria. There were many different bacterial profiles determined including 120 *S. aureus*, 110 coagulase-negative *Staphylococci* (CONS), 40 *Streptococcus* spp, 25 *S. pneumoniae*, 20 *E. coli*, 20 *Klebsiella* spp. and 20 *Enterobacter*. All *E. coli*, *Citrobacter* spp. and *Proteus* spp. isolates showed susceptibility to ceftriaxone. A majority of the *S. aureus*, CONS, *Streptococcus* spp., *S. pneumoniae*, *Klebsiella* spp. and *Enterobacter* spp. isolates demonstrated resistance to ceftriaxone. Of all Gram-positive isolates, 87 (29.5%) demonstrated ceftriaxone resistance, and 14 (10%) of Gram-negative isolates showed resistance.

This study concluded that there are multiple clinical isolates that cause non-chlamydia conjunctival bacteria, many of which demonstrated unacceptable resistance to chloramphenicol and azithromycin, two routinely administered antibiotics for prevention of said disease. The study found many antimicrobial agents that exhibit reduced resistance including ceftriaxone, demonstrating the need for further surveillance of resistance among patients undergoing trachomatous trichiasis surgery to improve treatment and reduce resistance. Oxoid Ceftriaxone AST Discs are an effective antimicrobial device to clinically analyse a variety of bacteria in non-chlamydia conjunctival bacterial infection.

A 2017 study investigated clinical isolates of *Staphylococcus aureus* to determine its nasal colonisation and subsequent antimicrobial susceptibility using the disc diffusion method of various antibiotics.⁶

Of the 400 collected nasal swabs, 52 (13%) clinical isolates of *Staphylococcus aureus* were found prevalent among the children. Susceptibility testing demonstrated that 86.6% of isolates were susceptible to ceftriaxone. All isolates were found to show some sensitivity to all the antimicrobial agents, with 100% sensitivity to drugs including oxacillin, cefoxitin, and gentamycin. Isolates demonstrated high rates of resistance to penicillin and ampicillin.

The study found that the frequency of *Staphylococcus aureus* in the children was relatively low compared to previous studies, at just 13% of the 400 subjects being prevalent for the bacteria. It was also determined that the clinical isolates were highly sensitive to certain antimicrobial drugs including ceftriaxone but demonstrated resistance to the β -lactam antimicrobial drugs, ampicillin and penicillin. This study, therefore, demonstrates that Oxoid Ceftriaxone AST Discs are an effective method for the in vitro analysis of *S. aureus*.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.


References

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Memon, Memon, and Memon (2017). Drug Resistance Patterns of *Acinetobacter baumannii* Infection in Intensive Care Unit of a Tertiary Care Hospital in Sindh, *Medical Forum Monthly* 28(4): 100-103. ISSN: 1029385X.
5. B. Abera and M. Kibret (2014.) Azithromycin, Fluoroquinolone and Chloramphenicol Resistance of Non-Chlamydia Conjunctival Bacteria in Rural Community of Ethiopia, *Indian Journal of Ophthalmology*, 62 (2): 236-239. DOI: 10.4103/0301-4738.99974.
6. Reta, M. Wubie and G. Mekuria (2017). Nasal Colonization and Antimicrobial Susceptibility Pattern of *Staphylococcus aureus* Among Pre-school Children in Ethiopia, *Journal: BMC Research Notes*, 10 (1). DOI: 10.1186/s13104-017-3079-6.


© 2025 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of Symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limitation (storage temp.)
	Lot / Batch Number
	Catalog Number
	Do not re-use
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Expiry Date
	Period-after-opening – the number of days the product can be used after opening
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Do not use if package is damaged
	Authorized representative in the European Community/ European Union
	Unique device identifier
	European Conformity Mark
	UK Conformity Mark
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

Revision Information

Version	Date of issue and modifications introduced
2.0	2025-03-14 Minor formatting changes Test media in table 4 updated

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Referencias

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Memon, Memon, and Memon (2017). Drug Resistance Patterns of *Acinetobacter baumannii* Infection in Intensive Care Unit of a Tertiary Care Hospital in Sindh, *Medical Forum Monthly* 28(4): 100-103. ISSN: 1029385X.
- B. Abera and M. Kibret (2014.) Azithromycin, Fluoroquinolone and Chloramphenicol Resistance of Non-Chlamydia Conjunctival Bacteria in Rural Community of Ethiopia, *Indian Journal of Ophthalmology*, 62 (2): 236-239. DOI: 10.4103/0301-4738.99974.
- Reta, M. Wubie and G. Mekuria (2017). Nasal Colonization and Antimicrobial Susceptibility Pattern of *Staphylococcus aureus* Among Pre-school Children in Ethiopia, *Journal: BMC Research Notes*, 10 (1). DOI: 10.1186/s13104-017-3079-6.

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Referanser

1. Clinical and Laboratory Standards Institute (CLSI). Se siste utgave:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Se siste versjon:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Memon, Memon, and Memon (2017). Drug Resistance Patterns of *Acinetobacter baumannii* Infection in Intensive Care Unit of a Tertiary Care Hospital in Sindh, *Medical Forum Monthly* 28(4): 100-103. ISSN: 1029385X.
5. B. Abera and M. Kibret (2014.) Azithromycin, Fluoroquinolone and Chloramphenicol Resistance of Non-Chlamydia Conjunctival Bacteria in Rural Community of Ethiopia, *Indian Journal of Ophthalmology*, 62 (2): 236-239. DOI: 10.4103/0301-4738.99974.
6. Reta, M. Wubie and G. Mekuria (2017). Nasal Colonization and Antimicrobial Susceptibility Pattern of *Staphylococcus aureus* Among Pre-school Children in Ethiopia, *Journal: BMC Research Notes*, 10 (1). DOI: 10.1186/s13104-017-3079-6.

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