
EN

Thermo Scientific™ Oxoid™ X and V Factor Diagnostic Discs **REF (DD0003B^Δ, DD0003T^Δ, DD0004B*, DD0004T*, DD0005B* and DD0005T*)**

Intended Use

The X and V Factor Diagnostic Discs are impregnated with X and/or V factor and are used in a diagnostic workflow to aid microbiology professionals in identifying potential pathogenic bacteria suspected in patients having a *Haemophilus* spp. infection. These discs enable the growth of *Haemophilus* spp. and provide a qualitative result.

The device is not automated, is for professional use only and is not a companion diagnostic.

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of the Test

X and V Factor Diagnostic Discs are used in the qualitative agar diffusion test method for *in vitro* susceptibility testing. A suitable therapeutic agent for *in vivo* use can be determined using filter paper discs impregnated with specified concentrations of growth factor agents placed on the surface of a suitable test medium. Pure cultures of clinical isolates are inoculated onto the test medium and the Diagnostic Disc (DD) placed on the surface. The growth factors within the disc diffuses through the agar to form a gradient. After incubation, the devices perform to the current performance criteria of the presence or absence of growth around the discs¹.

Components

X and V Factor Diagnostic Discs Diagnostic Discs consist of 6mm diameter paper discs impregnated concentrations of haemin and/or nicotinamide adenine dinucleotide (NAD). The discs are marked on both sides with the corresponding code for the antimicrobial agent (X DD, V DD or XV DD).

X and V Factor Diagnostic Discs Diagnostic Discs Diagnostic Discs are supplied in cartridges of 50 discs.

There is one cartridge per pack for DD0003T/DD0004T/DD0005T.

There are five cartridges per pack for DD0003B/DD0004B/DD0005B.

Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

X Factor Diagnostic Discs

DD0003T.....	Σ 50
DD0003B.....	Σ 250

V Factor Diagnostic Discs

DD0004T.....	Σ 50
DD0004B.....	Σ 250

X and V Factor Diagnostic Discs

DD0005T.....	Σ 50
DD0005B.....	Σ 250

Materials Provided

DD0003T/0004T/0005T

- 1 x Spring cap and plunger
- 1 x Desiccant tablet per cartridge
- 50 x discs in each cartridge
- 1 x Cartridge per pack
- Foil sealing each cartridge

DD0003B/0004B/0005B

- 1 x Spring cap and plunger per cartridge
- 1 x Desiccant tablet in each cartridge
- 50 x discs in each cartridge
- 5 x Cartridges per pack
- Foil sealing each cartridge

X Factor Diagnostic Discs contain the active ingredient Haemin an iron-containing porphyrin, typically a dark brown powder when in dry form

V Factor Diagnostic Discs contain the active ingredient Nicotinamide-adenine dinucleotide (NAD) a naturally occurring pyridine nucleotide, typically a white powder when in dry form.

X and V Factor Diagnostic Discs contain the active ingredient Haemin and NAD, typically a white powder when in dry form.

Materials Provided

X and V Factor Diagnostic Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. X and V Factor Diagnostic Discs are supplied in cartridges of 50 discs. There is one cartridge per pack for DD0003T/DD0004T/DD0005T.

There are five cartridges per pack for DD0003B/DD0004B/DD0005B. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with DD0003B/T, DD0004B/T, DD0005B/T	
Component Description	Material Description
Cartridge with spring, cap, and plunger (x1) DD0003T/DD0004T/DD0005T (x5) DD0003B/DD0004B/DD0005B.	Assembly components and plastic cartridge containing 50x Diagnostic discs.
Desiccant tablet (x1) DD0003T/DD0004T/DD0005T (x5) DD0003B/DD0004B/DD0005B.	Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
DD0003T/DD0004T/DD0005T. Diagnostic Discs (x50)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 1 cartridge per pack.
DD0003B/DD0004B/DD0005B. Diagnostic Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of X and V Factor Diagnostic Disc Reagents	
Reagent	Description of Function
Nicotinamide-adenine dinucleotide (Supplier: Ubichem)	Nicotinamide-adenine dinucleotide (NAD) is a naturally occurring pyridine nucleotide found in living cells. NAD acts as coenzymes and substrates involved in multiple cellular metabolic processes, including energy production. NAD has shown clinical importance as drug targets, supplementations and in the development of antibiotics.
Haemin (Supplier: Merck Life Sciences Ltd)	Haemin, is also known as protoporphyrin IX, is an iron-containing porphyrin that contains chlorine and is naturally made by red blood cells. Haemin is often administered as an injection in acute porphyria treatment.

The concentration of the active ingredient on the Diagnostic disc is analysed for every batch and is controlled using internal specifications. The actual concentration is detailed on the Certificate of Analysis.

X and V Factor Diagnostic Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Materials Required but Not Supplied

- Agar plates with appropriate media, inoculum suspension medium, sterile loops and swabs, sterile forceps, McFarland turbidity standards, incubator, modified atmosphere environments.
- Antibiotic disc dispensers.
- Quality control (QC) strains, apparatus to measure zone sizes, and interpretative criteria for local standard methods.

Warnings and Precautions

- Not intended for susceptibility testing.
- This product is for *in vitro* diagnostic use only.
- Follow directions for use.
- Observe aseptic techniques and established precautions against all microbiological hazards throughout all procedures.
- Cultures, containers and other contaminated materials must be sterilized after use in accordance with guidelines for the handling and disposal of biohazardous waste.
- Refer to the Material Safety Data Sheet (MSDS) for safe handling and disposal of the product (found on www.thermofisher.com).
- Inspect the product packaging before first use. Do not use product if there is any visible damage to the foil seal.
- Do not use beyond the stated expiry date.
- Once the cartridge is open, ensure it is stored in an opaque desiccated environment to prevent degradation.
- If the discs do not produce the expected reaction with recommended control microorganisms, review the entire procedure.
- In the event of malfunction do not use the device.
- The device does not contain any carcinogenic, mutagenic or toxic to reproduction, endocrine disrupting substances or materials likely to cause sensitisation or an allergic reaction under normal use.

Shelf-life and Storage Conditions

Unopened cartridges of X Factor (DD0003B, DD0003T) Diagnostic discs have a shelf-life of 12 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Unopened cartridges of V Factor (DD0004B, DD0004T) Diagnostic discs have a shelf-life of 12 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to -10°C until required.

Unopened cartridges of X and V Factor (DD0005B, DD0005T) Diagnostic discs have a shelf-life of 12 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to -10°C until required.

Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Method

1. Specimens should be collected and handled following recommended guidelines.
2. Prepare Blood Agar Base plates (CM0055) as directed.
3. Dip a sterile cotton wool swab into the suspension, rotate the swab several times, and press firmly on the inside wall of the tube to remove excess fluid. Spread the organisms evenly over the entire surface of the plate three times, in three directions by rotating the plate approximately 60° at a time. As a final step swab the rim of the agar.
4. Ensure there is no visible moisture on the surface of the plates before applying the discs. Discs should be applied to the surface of the agar within 15 minutes of inoculation.
5. Using sterile forceps, remove the disc from the cartridge and apply to the surface of the agar plate in the following positions around the periphery of the plate (approximately 1 or 2 cm in from the edge of the medium). X Factor disc at 12 o'clock, V Factor discs at 4 o'clock, and X + V Factor discs at 8 o'clock ensuring they are evenly spaced.
6. Lightly touch the disc with the sterile forceps to ensure complete contact between the disc and agar surface.
7. V factor diffuses more readily than X factor, do not place discs too close together as the V factor may diffuse towards the X factor discs, leading to growth due to the X factor rather than V factor.
8. Invert and incubate plates at 35-37°C for 18-24 hours.
9. After incubation inspect for the presence or absence of growth around the discs.

Summary:

Haemophilus influenzae requires both the X and the V factor to grow

Haemophilus parainfluenzae only requires the V factor to grow

Haemophilus haemoglobinophilus only requires the X factor to grow

Identification is presumptive and should be confirmed.

Interpretation

See Table 3. below for an Interpretation of the results.

Table 3. Interpretation of results

Product	Organism ATCC® number™/NCTC	Expected Result
X Factor Discs DD0003B DD0003T	<i>Haemophilus influenzae</i> ATCC® 19418™	No growth
	<i>Haemophilus haemoglobinophilus</i> NCTC 8540	Growth
	<i>Haemophilus influenzae</i> ATCC® 49247™	No growth
	<i>Haemophilus parainfluenzae</i> ATCC® 33392™	No growth
V Factor Discs DD0004B DD0004T	<i>Haemophilus influenzae</i> ATCC® 19418™	No growth
	<i>Haemophilus haemoglobinophilus</i> NCTC 8540	No growth
	<i>Haemophilus influenzae</i> ATCC® 49247™	No growth
	<i>Haemophilus parainfluenzae</i> ATCC® 33392™	Growth
X and V Factors DD0005B DD0005T	<i>Haemophilus influenzae</i> ATCC® 19418™	Growth
	<i>Haemophilus haemoglobinophilus</i> NCTC 8540	Growth
	<i>Haemophilus influenzae</i> ATCC® 49247™	Growth
	<i>Haemophilus parainfluenzae</i> ATCC® 33392™	Growth

User Quality Control Procedures

See the strains as listed in Table 4. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

Analytical Performance Characteristics

Repeatability and reproducibility were measured internally by testing the organisms in Table 4. Refer to Table 5 for the number of batches tested for each device and respective pass rates.

Table 4. QC Testing Panel

Product Code	Organism	ATCC® number™ /NCTC
DD0003B DD0003T	<i>Haemophilus influenzae</i>	ATCC® 19418™
	<i>Haemophilus haemoglobinophilus</i>	NCTC 8540
	<i>Haemophilus influenzae</i>	ATCC® 49247™
DD0004B DD0004T	<i>Haemophilus parainfluenzae</i>	ATCC® 33392™
	<i>Haemophilus influenzae</i>	ATCC® 19418™
	<i>Haemophilus haemoglobinophilus</i>	NCTC 8540
	<i>Haemophilus influenzae</i>	ATCC® 49247™
DD0005B DD0005T	<i>Haemophilus parainfluenzae</i>	ATCC® 33392™
	<i>Haemophilus influenzae</i>	ATCC® 19418™
	<i>Haemophilus influenzae</i>	ATCC® 9334™
	<i>Haemophilus influenzae</i>	ATCC® 49247™

Table 5. Repeatability and Reproducibility of Tested Batches

Product Code	Batches Tested	Time Period	Pass Rate
DD0003B	176	2001 – 2021	99.43%
	7	17/MAR/2023 – 22/APR/2024*	100%
DD0003T	148	1998 – 2021	95.27%
	3	24/AUG/2023 – 16/APR/2024*	100%
DD0004B	163	2001 – 2021	95.09%
	7	21/JUL/2023 – 16/APR/2024*	100%
DD0004T	118	1998 – 2021	98.31%
	2	21/JUL/2023 – 20/OCT/2023*	100%
DD0005B	113	2001 – 2021	95.58%
	5	29/JUN/2023 – 02/OCT/2023*	100%
DD0005T	142	1998 – 2021	97.06%
	3	13/JUN/2023 – 16/MAY/2024*	100%

*Quality Control test organisms were updated to include *Haemophilus haemoglobinophilus* NCTC 8540 and remove *Haemophilus influenzae* ATCC® 19418™

Clinical Performance Characteristics

A study by Su et al² in 2020 used the Thermo Scientific™ Oxoid™ X and V Factor Diagnostic Discs to explore the epidemiology of over 1,400 *H. influenzae* strains in Taiwan. Their results found antimicrobial resistant *H. influenzae* at a rate that remained constant over the years. The isolates were frequently resistant to the common drug of choice, ampicillin. This signifies the need for correct identification of the isolates for effective treatment.

Thermo Scientific™ Oxoid™ X and V Factor Diagnostic Discs were used by Unger et al³ in 2015 to determine the impact of IPTp with SPAZ on the nasopharyngeal carriage of *H. influenzae* among pregnant women in Papua New Guinea. IPTp with SPAZ was compared with the control of IPTp with SPCQ and the Thermo Scientific™ Oxoid™ X and V Factor Diagnostic Discs successfully identified the incidence of *H. influenzae*. Their results demonstrated that IPTp with SPAZ significantly reduced the nasopharyngeal carriage of *H. influenzae* compared with IPTp with SPCQ (Unger et al. 2015).

In 2005, Poyrazoglu et al⁴ collected oropharyngeal samples from over 800 young children in Turkey and the Thermo Scientific™ Oxoid™ X and V Factor Diagnostic Discs were used to accurately identify *H. influenzae*. These isolates were subsequently confirmed as *H. influenzae* type B and carriage rate was highest in the older children with antimicrobial resistance emerging (Poyrazoglu et al. 2005).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

1. Public Health England. 2019. "UK Standards for Microbiology Investigations-Standards-for-Microbiology-Investigations-Smi-Quality-and-Consistency-in-Clinical-Laboratories PHE Publications Gateway Number: 2018381 UK Standards for Microbiology Investigations Are Produced in Association With." <https://www.rcpath.org/>.
2. Su, P, A Huang, C Lai, H Lin, T Lin, and C Ho. 2020. "Extensively Drug-Resistant Haemophilus Influenzae - Emergence, Epidemiology, Risk Factors, and Regimen." *BMC Microbiology* 20 (1). <https://doi.org/10.1186/S12866-020-01785-9>
3. Unger, Holger W., Celestine Aho, Maria Ome-Kaius, Regina A. Wangnapi, Alexandra J. Umbers, Wanda Jack, Alice Lafana, et al. 2015. "Impact of Intermittent Preventive Treatment in Pregnancy with Azithromycin-Containing Regimens on Maternal Nasopharyngeal Carriage and Antibiotic Sensitivity of Streptococcus Pneumoniae, Haemophilus Influenzae, and Staphylococcus Aureus: A Cross-Sectional Survey at Delivery." *Journal of Clinical Microbiology* 53 (4): 1317–23. <https://doi.org/10.1128/JCM.03570-14>.
4. Poyrazoglu, S, S Komec, G Gokcay, and B Ongen. 2005. "Haemophilus Influenzae Type b Carriage among 3- to 24-Month-Old Turkish Children." *Epidemiology and Infection* 133 (6): 1113–17. <https://doi.org/10.1017/S0950268805004>

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
IVD	In Vitro Diagnostic Medical Device
	Temperature limitation (storage temp.)
LOT	Lot / Batch Number
REF	Catalog Number
	Do not re-use
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Expiry Date
	Period-after-opening – the number of days the product can be used after opening
R_x only	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Do not use if package is damaged
EC REP	Authorized representative in the European Community/ European Union
UDI	Unique device identifier
CE	European Conformity Mark
UK CA	UK Conformity Mark
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

Revision Information

Version	Date of issue
3.0	2024-08-15. Updated Principle of test Method and QC procedures to align with current methodology. Updated Clinical Performance Characteristics to align with IVDR file

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1. Public Health England. 2019. "UK Standards for Microbiology Investigations-Standards-for-Microbiology-Investigations-Smi-Quality-and-Consistency-in-Clinical-Laboratories PHE Publications Gateway Number: 2018381 UK Standards for Microbiology Investigations Are Produced in Association With." <https://www.rcpath.org/>.
2. Su, P, A Huang, C Lai, H Lin, T Lin, and C Ho. 2020. "Extensively Drug-Resistant Haemophilus Influenzae - Emergence, Epidemiology, Risk Factors, and Regimen." *BMC Microbiology* 20 (1). <https://doi.org/10.1186/S12866-020-01785-9>
3. Unger, Holger W., Celestine Aho, Maria Ome-Kaius, Regina A. Wangnapi, Alexandra J. Umbers, Wanda Jack, Alice Lafana, et al. 2015. "Impact of Intermittent Preventive Treatment in Pregnancy with Azithromycin-Containing Regimens on Maternal Nasopharyngeal Carriage and Antibiotic Sensitivity of Streptococcus Pneumoniae, Haemophilus Influenzae, and Staphylococcus Aureus: A Cross-Sectional Survey at Delivery." *Journal of Clinical Microbiology* 53 (4): 1317–23. <https://doi.org/10.1128/JCM.03570-14>.
4. Poyrazoglu, S, S Komec, G Gokcay, and B Ongen. 2005. "Haemophilus Influenzae Type b Carriage among 3- to 24-Month-Old Turkish Children." *Epidemiology and Infection* 133 (6): 1113–17. <https://doi.org/10.1017/S09502688050004>

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Literatura

1. Public Health England. 2019. "UK Standards for Microbiology Investigations-Standards-for-Microbiology-Investigations-Smi-Quality-and-Consistency-in-Clinical-Laboratories PHE Publications Gateway Number: 2018381 UK Standards for Microbiology Investigations Are Produced in Association With." <https://www.rcpath.org/>.
2. Su, P, A Huang, C Lai, H Lin, T Lin, and C Ho. 2020. "Extensively Drug-Resistant *Haemophilus Influenzae* - Emergence, Epidemiology, Risk Factors, and Regimen." *BMC Microbiology* 20 (1). <https://doi.org/10.1186/S12866-020-01785-9>
3. Unger, Holger W., Celestine Aho, Maria Ome-Kaius, Regina A. Wangnapi, Alexandra J. Umbers, Wanda Jack, Alice Lafana, et al. 2015. "Impact of Intermittent Preventive Treatment in Pregnancy with Azithromycin-Containing Regimens on Maternal Nasopharyngeal Carriage and Antibiotic Sensitivity of *Streptococcus Pneumoniae*, *Haemophilus Influenzae*, and *Staphylococcus Aureus*: A Cross-Sectional Survey at Delivery." *Journal of Clinical Microbiology* 53 (4): 1317–23. <https://doi.org/10.1128/JCM.03570-14>.
4. Poyrazoglu, S, S Komec, G Gokcay, and B Ongen. 2005. "Haemophilus Influenzae Type b Carriage among 3- to 24-Month-Old Turkish Children." *Epidemiology and Infection* 133 (6): 1113–17. <https://doi.org/10.1017/S0950268805004>

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Referencer

1. Public Health England. 2019. "UK Standards for Microbiology Investigations-Standards-for-Microbiology-Investigations-Smi-Quality-and-Consistency-in-Clinical-Laboratories PHE Publications Gateway Number: 2018381 UK Standards for Microbiology Investigations Are Produced in Association With." <https://www.rcpath.org/>.
2. Su, P, A Huang, C Lai, H Lin, T Lin, and C Ho. 2020. "Extensively Drug-Resistant *Haemophilus Influenzae* - Emergence, Epidemiology, Risk Factors, and Regimen." *BMC Microbiology* 20 (1). <https://doi.org/10.1186/S12866-020-01785-9>
3. Unger, Holger W., Celestine Aho, Maria Ome-Kaius, Regina A. Wangnapi, Alexandra J. Umbers, Wanda Jack, Alice Lafana, et al. 2015. "Impact of Intermittent Preventive Treatment in Pregnancy with Azithromycin-Containing Regimens on Maternal Nasopharyngeal Carriage and Antibiotic Sensitivity of *Streptococcus Pneumoniae*, *Haemophilus Influenzae*, and *Staphylococcus Aureus*: A Cross-Sectional Survey at Delivery." *Journal of Clinical Microbiology* 53 (4): 1317–23. <https://doi.org/10.1128/JCM.03570-14>.

4. Poyrazoglu, S, S Komec, G Gokcay, and B Ongen. 2005. "Haemophilus Influenzae Type b Carriage among 3- to 24-Month-Old Turkish Children." *Epidemiology and Infection* 133 (6): 1113–17. <https://doi.org/10.1017/S0950268805004>

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Literatura

- Public Health England. 2019. "UK Standards for Microbiology Investigations-Standards-for-Microbiology-Investigations-Smi-Quality-and-Consistency-in-Clinical-Laboratories PHE Publications Gateway Number: 2018381 UK Standards for Microbiology Investigations Are Produced in Association With." <https://www.rcpath.org/>.
- Su, P, A Huang, C Lai, H Lin, T Lin, and C Ho. 2020. "Extensively Drug-Resistant *Haemophilus Influenzae* - Emergence, Epidemiology, Risk Factors, and Regimen." *BMC Microbiology* 20 (1). <https://doi.org/10.1186/S12866-020-01785-9>
- Unger, Holger W., Celestine Aho, Maria Ome-Kaius, Regina A. Wangnapi, Alexandra J. Umbers, Wanda Jack, Alice Lafana, et al. 2015. "Impact of Intermittent Preventive Treatment in Pregnancy with Azithromycin-Containing Regimens on Maternal Nasopharyngeal Carriage and Antibiotic Sensitivity of *Streptococcus Pneumoniae*, *Haemophilus Influenzae*, and *Staphylococcus Aureus*: A Cross-Sectional Survey at Delivery." *Journal of Clinical Microbiology* 53 (4): 1317–23. <https://doi.org/10.1128/JCM.03570-14>.
- Poyrazoglu, S, S Komec, G Gokcay, and B Ongen. 2005. "Haemophilus Influenzae Type b Carriage among 3- to 24-Month-Old Turkish Children." *Epidemiology and Infection* 133 (6): 1113–17. <https://doi.org/10.1017/S09502688050004>

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
Bibliografia

1. Public Health England. 2019. "UK Standards for Microbiology Investigations-Standards-for-Microbiology-Investigations-Smi-Quality-and-Consistency-in-Clinical-Laboratories PHE Publications Gateway Number: 2018381 UK Standards for Microbiology Investigations Are Produced in Association With." <https://www.rcpath.org/>.
2. Su, P, A Huang, C Lai, H Lin, T Lin, and C Ho. 2020. "Extensively Drug-Resistant *Haemophilus Influenzae* - Emergence, Epidemiology, Risk Factors, and Regimen." *BMC Microbiology* 20 (1). <https://doi.org/10.1186/S12866-020-01785-9>

3. Unger, Holger W., Celestine Aho, Maria Ome-Kaius, Regina A. Wangnapi, Alexandra J. Umbers, Wanda Jack, Alice Lafana, et al. 2015. "Impact of Intermittent Preventive Treatment in Pregnancy with Azithromycin-Containing Regimens on Maternal Nasopharyngeal Carriage and Antibiotic Sensitivity of Streptococcus Pneumoniae, Haemophilus Influenzae, and Staphylococcus Aureus: A Cross-Sectional Survey at Delivery." *Journal of Clinical Microbiology* 53 (4): 1317–23. <https://doi.org/10.1128/JCM.03570-14>.
4. Poyrazoglu, S, S Komec, G Gokcay, and B Ongen. 2005. "Haemophilus Influenzae Type b Carriage among 3- to 24-Month-Old Turkish Children." *Epidemiology and Infection* 133 (6): 1113–17. <https://doi.org/10.1017/S0950268805004>

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- Public Health England. 2019. „UK Standards for Microbiology Investigations-Standards-for-Microbiology-Investigations-Smi-Quality-and-Consistency-in-Clinical-Laboratories PHE Publications Gateway Number: 2018381 UK Standards for Microbiology Investigations Are Produced in Association With.“ <https://www.rcpath.org/>.
- Su, P, A Huang, C Lai, H Lin, T Lin, and C Ho. 2020. „Extensively Drug-Resistant *Haemophilus Influenzae* - Emergence, Epidemiology, Risk Factors, and Regimen.“ *BMC Microbiology* 20 (1). <https://doi.org/10.1186/S12866-020-01785-9>
- Unger, Holger W., Celestine Aho, Maria Ome-Kaius, Regina A. Wangnapi, Alexandra J. Umbers, Wanda Jack, Alice Lafana, et al. 2015. „Impact of Intermittent Preventive Treatment in Pregnancy with Azithromycin-Containing Regimens on Maternal Nasopharyngeal Carriage and Antibiotic Sensitivity of *Streptococcus Pneumoniae*, *Haemophilus Influenzae*, and *Staphylococcus Aureus*: A Cross-Sectional Survey at Delivery.“ *Journal of Clinical Microbiology* 53 (4): 1317–23. <https://doi.org/10.1128/JCM.03570-14>.
- Poyrazoglu, S, S Komec, G Gokcay, and B Ongen. 2005. „*Haemophilus Influenzae* Type b Carriage among 3- to 24-Month-Old Turkish Children.“ *Epidemiology and Infection* 133 (6): 1113–17. <https://doi.org/10.1017/S09502688050004>

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
Atsauces

- Public Health England. 2019. "UK Standards for Microbiology Investigations-Standards-for-Microbiology-Investigations-Smi-Quality-and-Consistency-in-Clinical-Laboratories PHE Publications Gateway Number: 2018381 UK Standards for Microbiology Investigations Are Produced in Association With." <https://www.rcpath.org/>.
- Su, P, A Huang, C Lai, H Lin, T Lin, and C Ho. 2020. "Extensively Drug-Resistant *Haemophilus Influenzae* - Emergence, Epidemiology, Risk Factors, and Regimen." *BMC Microbiology* 20 (1). <https://doi.org/10.1186/S12866-020-01785-9>
- Unger, Holger W., Celestine Aho, Maria Ome-Kaius, Regina A. Wangnapi, Alexandra J. Umbers, Wanda Jack, Alice Lafana, et al. 2015. "Impact of Intermittent Preventive Treatment in Pregnancy with Azithromycin-Containing Regimens on Maternal Nasopharyngeal Carriage and Antibiotic Sensitivity of *Streptococcus Pneumoniae*, *Haemophilus Influenzae*, and *Staphylococcus Aureus*: A Cross-Sectional Survey at Delivery." *Journal of Clinical Microbiology* 53 (4): 1317–23. <https://doi.org/10.1128/JCM.03570-14>.
- Poyrazoglu, S, S Komec, G Gokcay, and B Ongen. 2005. "*Haemophilus Influenzae* Type b Carriage among 3- to 24-Month-Old Turkish Children." *Epidemiology and Infection* 133 (6): 1113–17. <https://doi.org/10.1017/S09502688050004>

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NO

Thermo Scientific™ Oxoid™ X and V Factor Diagnostic Discs

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Referanser

1. Public Health England. 2019. "UK Standards for Microbiology Investigations-Standards-for-Microbiology-Investigations-Smi-Quality-and-Consistency-in-Clinical-Laboratories PHE Publications Gateway Number: 2018381 UK Standards for Microbiology Investigations Are Produced in Association With." <https://www.rcpath.org/>.
2. Su, P, A Huang, C Lai, H Lin, T Lin, and C Ho. 2020. "Extensively Drug-Resistant *Haemophilus Influenzae* - Emergence, Epidemiology, Risk Factors, and Regimen." *BMC Microbiology* 20 (1). <https://doi.org/10.1186/S12866-020-01785-9>
3. Unger, Holger W., Celestine Aho, Maria Ome-Kaius, Regina A. Wangnapi, Alexandra J. Umbers, Wanda Jack, Alice Lafana, et al. 2015. "Impact of Intermittent Preventive Treatment in Pregnancy with Azithromycin-Containing Regimens on Maternal Nasopharyngeal Carriage and Antibiotic Sensitivity of *Streptococcus Pneumoniae*, *Haemophilus Influenzae*, and *Staphylococcus Aureus*: A Cross-Sectional Survey at Delivery." *Journal of Clinical Microbiology* 53 (4): 1317–23. <https://doi.org/10.1128/JCM.03570-14>.
4. Poyrazoglu, S, S Komec, G Gokcay, and B Ongun. 2005. "Haemophilus Influenzae Type b Carriage among 3- to 24-Month-Old Turkish Children." *Epidemiology and Infection* 133 (6): 1113–17. <https://doi.org/10.1017/S09502688050004>

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Referințe

1. Public Health England. 2019. "UK Standards for Microbiology Investigations-Standards-for-Microbiology-Investigations-Smi-Quality-and-Consistency-in-Clinical-Laboratories PHE Publications Gateway Number: 2018381 UK Standards for Microbiology Investigations Are Produced in Association With." <https://www.rcpath.org/>.
2. Su, P., A Huang, C Lai, H Lin, T Lin, and C Ho. 2020. "Extensively Drug-Resistant *Haemophilus Influenzae* - Emergence, Epidemiology, Risk Factors, and Regimen." *BMC Microbiology* 20 (1). <https://doi.org/10.1186/S12866-020-01785-9>
3. Unger, Holger W., Celestine Aho, Maria Ome-Kaius, Regina A. Wangnapi, Alexandra J. Umbers, Wanda Jack, Alice Lafana, et al. 2015. "Impact of Intermittent Preventive Treatment in Pregnancy with Azithromycin-Containing Regimens on Maternal Nasopharyngeal Carriage and Antibiotic Sensitivity of *Streptococcus Pneumoniae*, *Haemophilus Influenzae*, and *Staphylococcus Aureus*: A Cross-Sectional Survey at Delivery." *Journal of Clinical Microbiology* 53 (4): 1317–23. <https://doi.org/10.1128/JCM.03570-14>.

4. Poyrazoglu, S, S Komec, G Gokcay, and B Ongen. 2005. "Haemophilus Influenzae Type b Carriage among 3- to 24-Month-Old Turkish Children." *Epidemiology and Infection* 133 (6): 1113–17. <https://doi.org/10.1017/S0950268805004>

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
Referenser

1. Public Health England. 2019. "UK Standards for Microbiology Investigations-Standards-for-Microbiology-Investigations-Smi-Quality-and-Consistency-in-Clinical-Laboratories PHE Publications Gateway Number: 2018381 UK Standards for Microbiology Investigations Are Produced in Association With." <https://www.rcpath.org/>.
2. Su, P, A Huang, C Lai, H Lin, T Lin, and C Ho. 2020. "Extensively Drug-Resistant *Haemophilus Influenzae* - Emergence, Epidemiology, Risk Factors, and Regimen." *BMC Microbiology* 20 (1). <https://doi.org/10.1186/S12866-020-01785-9>
3. Unger, Holger W., Celestine Aho, Maria Ome-Kaius, Regina A. Wangnapi, Alexandra J. Umbers, Wanda Jack, Alice Lafana, et al. 2015. "Impact of Intermittent Preventive Treatment in Pregnancy with Azithromycin-Containing Regimens on Maternal Nasopharyngeal Carriage and Antibiotic Sensitivity of *Streptococcus Pneumoniae*, *Haemophilus Influenzae*, and *Staphylococcus Aureus*: A Cross-Sectional Survey at Delivery." *Journal of Clinical Microbiology* 53 (4): 1317–23. <https://doi.org/10.1128/JCM.03570-14>.
4. Poyrazoglu, S, S Komec, G Gokcay, and B Ongen. 2005. "*Haemophilus Influenzae* Type b Carriage among 3- to 24-Month-Old Turkish Children." *Epidemiology and Infection* 133 (6): 1113–17. <https://doi.org/10.1017/S09502688050004>

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