<https://www.thermofisher.com>

EN

Oxoid™ Fosfomycin Trometamol 200 µg (FOT200)

REF CT0758B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the Oxoid AST generic instructions for use supplied with the product and available online.

Fosfomycin Trometamol Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent fosfomycin trometamol. The discs are labelled on both sides with details of the antimicrobial (FOT) and amount present (µg): FOT200 (200µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Fosfomycin Trometamol AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which Fosfomycin and Trometamol has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- *Enterococcus faecalis*

Gram-Negative

- Enterobacterales

Species with published breakpoints according to current CLSI M100 literature:

Gram-Positive

- *Enterococcus* spp.

Gram-Negative

- Enterobacterales

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- N/A

Gram-Negative

- Enterobacterales (Fosfomycin iv, *E. coli* only)
- Enterobacterales (Fosfomycin oral, uncomplicated UTI only, *E. coli*)

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Oxoid Fosfomycin Trometamol AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of

inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Oxoid Fosfomycin Trometamol AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Oxoid Fosfomycin Trometamol AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Fosfomycin Trometamol AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with CT0758B	
Component Description	Material Description
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Pale beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Fosfomycin Trometamol AST Disc Reagents	
Reagent	Description of Function
Fosfomycin	White to almost white powder. Fosfomycin is a broad-spectrum antimicrobial agent belonging to the phosphonic class of antibiotics. Fosfomycin exhibits its mechanism of action through the inhibition of peptidoglycan synthesis, subsequently inhibiting cell wall synthesis and ultimately leading to cell lysis.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Oxoid Fosfomycin Trometamol AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Fosfomycin Trometamol AST discs have a shelf-life of 24 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Raw data of zone sizes, in accordance with CLSI/EUCAST methodology, taken from the 3 most recent batches of Fosfomycin Trometamol (200 µg) FOT200 Antimicrobial Susceptibility Test Discs (CT0758B).

CLSI Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3330349	<i>Staphylococcus aureus</i>	25923	MHA	25	33	29	25	25	25	-4 ± 0 (CV=0)
3312219	<i>Staphylococcus aureus</i>	25923	MHA	25	33	29	25	25	25	-4 ± 0 (CV=0)
3312220	<i>Staphylococcus aureus</i>	25923	MHA	25	33	29	27	27	27	-2 ± 0 (CV=0)
EUCAST Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3330349	<i>Escherichia coli</i>	25922	MHA	26	34	30	27	27	27	-3 ± 0 (CV=0)
3312219	<i>Escherichia coli</i>	25922	MHA	26	34	30	28	28	28	-2 ± 0 (CV=0)
3312220	<i>Escherichia coli</i>	25922	MHA	26	34	30	28	28	28	-2 ± 0 (CV=0)

Clinical Performance Characteristics

A study by Garcia-Fulguerias et al. (2021) collected urine samples from children aged six to 14 years of age with UTIs in a hospital in Uruguay.⁴ The isolates were identified using an automated technique; susceptibility testing was performed via disc diffusion method for a variety of antibiotics, including Oxoid Fosfomycin Trometamol AST Discs (CT0758B). Susceptibility testing was in accordance with CLSI guidelines. Polymerase chain reaction was performed to identify fosfomycin resistant associated genes in the fosfomycin resistant strains.

Significant bacteriuria was identified in 286 urine samples; *E. coli* was the most frequently isolated pathogen in 255 samples. Susceptibility testing found 97.9% of isolates to be susceptible to fosfomycin. The rate of susceptibility was found to be higher than all other studied antibiotics except cefuroxime and ceftriaxone, which were susceptible in three more isolates (99%). Lower susceptibility among the isolates were found for ampicillin, ciprofloxacin, nitrofurantoin, and gentamicin. More specifically, fosfomycin was susceptible in 100% of *P. mirabilis* and *K. pneumoniae* isolates and in 99.6% of *Enterobacterales* isolates. Five of the six fosfomycin resistant isolates were found to have natural resistance to fosfomycin and one of the six had acquired resistance.

This study found relatively high susceptibility for multiple orally administered antibiotics, including fosfomycin, and concluded that fosfomycin exhibited an advantageous therapy option for treatment of paediatric UTIs due to its low side effects, high activity, and once daily dosing. No inconsistent or out of specification results for this device were reported. The study showed that Oxoid Fosfomycin Trometamol AST Disc are an effective method for the in vitro analysis of *Enterobacterales*.

Choudhury, Yeng, and Krishnan (2015) collected carbapenemase-producing *Enterobacteriaceae* (CPE) isolates from clinical samples (urine, blood, skin and soft tissue infection) from September 2012 to March 2015 in a hospital in Singapore.⁵ The isolates were identified, and CPE associated genes were confirmed. Antimicrobial susceptibility testing was performed for fosfomycin using Oxoid Fosfomycin Trometamol AST Discs (CT0758B) and tigecycline. Susceptibility testing was undertaken according to CLSI guidelines.

Thirty-six CPE isolates were collected. Most frequently identified strains were *Klebsiella pneumoniae* and *E. coli*. Genotyping found four distinct genotypes. Susceptibility testing found 88.8% of isolates to be susceptible to fosfomycin; all urinary isolated pathogens were susceptible. Tigecycline susceptibility was found in 66.6% of all isolates with 71.0% of skin and soft tissue infection derived isolates being susceptible.

The study demonstrated higher susceptibility to fosfomycin. The study shows potential treatment benefits for both antibiotics for site specific infections, but more specifically the use of fosfomycin for urinary tract infections caused by CPE isolates. No inconsistent or out of specification results for this device were reported. This study demonstrated that Oxoid Fosfomycin Trometamol AST Discs are an effective method for the in vitro analysis of *Enterobacteriaceae* pathogens.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Test
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. V. Garcia-Fulgueiras, L. Caiata, I. Bado, G. Giachetto and L. Robino (2021). Antibiotic susceptibility and fosfomycin resistance characterization in a cohort of children older than 6 years of age with urinary tract infection. *Revista Argentina de Microbiologia*. DOI: 10.1016/J.RAM.2021.04.002
5. S. Choudhury, J. Yeng and P. Krishnan (2015). In vitro susceptibilities of clinical isolates of carbapenemase-producing *Enterobacteriaceae* to fosfomycin and tigecycline. *Clinical Microbiology and Infection* 21(10): 75-76. DOI: 10.1016/J.CMI.2015.06.005

© 2025 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limitation (storage temp.)
	Lot / Batch Number
	Catalog Number
	Do not re-use
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Expiry Date
	Period of use after opening 7 days
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Do not use if package is damaged
	Authorized representative in the European Community/ European Union
	Unique device identifier
	European Conformity Mark
	UK Conformity Mark
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union

Made in the United Kingdom	Made in the United Kingdom
-------------------------------	----------------------------

Revision Information

Version	Date of issue and modifications introduced
3.0	2025-02-04 QC results amended for E. coli

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- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- V. Garcia-Fulgueiras, L. Caiata, I. Bado, G. Giachetto and L. Robino (2021). Antibiotic susceptibility and fosfomycin resistance characterization in a cohort of children older than 6 years of age with urinary tract infection. *Revista Argentina de Microbiologia*. DOI: 10.1016/J.RAM.2021.04.002
- S. Choudhury, J. Yeng and P. Krishnan (2015). In vitro susceptibilities of clinical isolates of carbapenemase-producing Enterobacteriaceae to fosfomycin and tigecycline. *Clinical Microbiology and Infection* 21(10): 75-76. DOI: 10.1016/J.CMI.2015.06.005

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Reference

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. V. Garcia-Fulgueiras, L. Caiata, I. Bado, G. Giachetto and L. Robino (2021). Antibiotic susceptibility and fosfomycin resistance characterization in a cohort of children older than 6 years of age with urinary tract infection. *Revista Argentina de Microbiologia*. DOI: 10.1016/J.RAM.2021.04.002
5. S. Choudhury, J. Yeng and P. Krishnan (2015). In vitro susceptibilities of clinical isolates of carbapenemase-producing *Enterobacteriaceae* to fosfomycin and tigecycline. *Clinical Microbiology and Infection* 21(10): 75-76. DOI: 10.1016/J.CMI.2015.06.005

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Referenzen

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. V. Garcia-Fulgueiras, L. Caiata, I. Bado, G. Giachetto and L. Robino (2021). Antibiotic susceptibility and fosfomycin resistance characterization in a cohort of children older than 6 years of age with urinary tract infection. *Revista Argentina de Microbiología*. DOI: 10.1016/J.RAM.2021.04.002
5. S. Choudhury, J. Yeng and P. Krishnan (2015). In vitro susceptibilities of clinical isolates of carbapenemase-producing *Enterobacteriaceae* to fosfomycin and tigecycline. *Clinical Microbiology and Infection* 21(10): 75-76. DOI: 10.1016/J.CMI.2015.06.005

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Referencias

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. V. García-Fulgueiras, L. Caiata, I. Bado, G. Giachetto and L. Robino (2021). Antibiotic susceptibility and fosfomycin resistance characterization in a cohort of children older than 6 years of age with urinary tract infection. *Revista Argentina de Microbiología*. DOI: 10.1016/J.RAM.2021.04.002
5. S. Choudhury, J. Yeng and P. Krishnan (2015). In vitro susceptibilities of clinical isolates of carbapenemase-producing *Enterobacteriaceae* to fosfomycin and tigecycline. *Clinical Microbiology and Infection* 21(10): 75-76. DOI: 10.1016/J.CMI.2015.06.005

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Références

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. V. Garcia-Fulgueiras, L. Caiata, I. Bado, G. Giachetto and L. Robino (2021). Antibiotic susceptibility and fosfomycin resistance characterization in a cohort of children older than 6 years of age with urinary tract infection. *Revista Argentina de Microbiología*. DOI: 10.1016/J.RAM.2021.04.002
5. S. Choudhury, J. Yeng and P. Krishnan (2015). In vitro susceptibilities of clinical isolates of carbapenemase-producing Enterobacteriaceae to fosfomycin and tigecycline. *Clinical Microbiology and Infection* 21(10): 75-76. DOI: 10.1016/J.CMI.2015.06.005

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1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. V. Garcia-Fulguerias, L. Caiata, I. Bado, G. Giachetto and L. Robino (2021). Antibiotic susceptibility and fosfomicin resistance characterization in a cohort of children older than 6 years of age with urinary tract infection. *Revista Argentina de Microbiologia*. DOI: 10.1016/J.RAM.2021.04.002
5. S. Choudhury, J. Yeng and P. Krishnan (2015). In vitro susceptibilities of clinical isolates of carbapenemase-producing *Enterobacteriaceae* to fosfomicin and tigecycline. *Clinical Microbiology and Infection* 21(10): 75-76. DOI: 10.1016/J.CMI.2015.06.005

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Bibliografia

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. V. Garcia-Fulguerias, L. Caiata, I. Bado, G. Giachetto and L. Robino (2021). Antibiotic susceptibility and fosfomycin resistance characterization in a cohort of children older than 6 years of age with urinary tract infection. *Revista Argentina de Microbiologia*. DOI: 10.1016/J.RAM.2021.04.002
5. S. Choudhury, J. Yeng and P. Krishnan (2015). In vitro susceptibilities of clinical isolates of carbapenemase-producing *Enterobacteriaceae* to fosfomycin and tigecycline. *Clinical Microbiology and Infection* 21(10): 75-76. DOI: 10.1016/J.CMI.2015.06.005

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Atsauces

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. V. Garcia-Fulguerias, L. Caiata, I. Bado, G. Giachetto and L. Robino (2021). Antibiotic susceptibility and fosfomycin resistance characterization in a cohort of children older than 6 years of age with urinary tract infection. *Revista Argentina de Microbiologia*. DOI: 10.1016/J.RAM.2021.04.002
5. S. Choudhury, J. Yeng and P. Krishnan (2015). In vitro susceptibilities of clinical isolates of carbapenemase-producing Enterobacteriaceae to fosfomycin and tigecycline. *Clinical Microbiology and Infection* 21(10): 75-76. DOI: 10.1016/J.CMI.2015.06.005

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