	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Nitrofurantoin Discs 100 µg/300 µg (F100/F300)

REF CT0034B & CT0036B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the Oxoid AST generic instructions for use supplied with the product and available online.

Thermo Scientific™ Oxoid™ Nitrofurantoin Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent Nitrofurantoin. The discs are labelled on both sides with details of antimicrobial and amount present (µg): F100 (100µg) or F300 (300µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Oxoid Nitrofurantoin AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disk should only be used once.

Intended Use

Oxoid Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which Nitrofurantoin has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Species with published breakpoints according to current FDA literature (all associated with 300 µg):

Gram-Positive:

- Staphylococcus spp.
- Enterococcus spp.

Gram-Negative:

- Enterobacteriaceae

Species with published breakpoints according to current CLSI M100 (all associated with 300 µg):

Gram-Positive:

- Staphylococcus spp.
- Enterococcus spp.

Gram-Negative:

- Enterobacterales

Species with published breakpoints according to current EUCAST breakpoint tables (all associated with 100 µg):

Gram-Positive:

- Staphylococcus spp. (Uncomplicated UTI only) (S. saprophyticus)
- Enterococcus spp. (Uncomplicated UTI only) (E. faecalis)
- Streptococcus spp. (Uncomplicated UTI only) (S. agalactiae, group B streptococci)
- Aerococcus sanguinicola and urinae (Uncomplicated UTI only)

Gram-Negative:

- Enterobacterales (Uncomplicated UTI only) (E. coli)

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principles of Method

Oxoid Nitrofurantoin AST Discs may be used in the semi-quantitative agar diffusion test method for *in vitro* susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology please refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological traceability of calibrator and control material values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Oxoid Nitrofurantoin AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Oxoid Nitrofurantoin AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Nitrofurantoin AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, please refer to table 2.

Table 1. Materials provided with CT0034B & CT0036B	
Component description	Material description
Cartridge with spring, cap and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs
Desiccant tablet (x5)	Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge, 5 cartridges per pack.
Cartridge with spring, cap and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs

Table2. Description of Nitrofurantoin AST Disc reagents	
Reagent	Description of Function
Nitrofurantoin	Yellow crystalline powder. Antimicrobial agent. Mechanism of action is inactivation of or alteration of bacterial ribosomal proteins and other related macromolecules.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Oxoid Nitrofurantoin AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life & Storage Conditions

Unopened cartridges of Nitrofurantoin AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided (with an unsaturated (orange) desiccant), or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2 - 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Results summary for tested QC isolates Nitrofurantoin F100 (100µg) Antimicrobial Susceptibility Disks (CT0034B), in accordance with EUCAST methodology

Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3116889	<i>Escherichia coli</i>	ATCC® 25923	MHA	17	23	20	21	21	21	1 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC® 25922	MHA	17	23	20	21	21	21	1 ± 0 (CV=0)
	<i>Enterococcus faecalis</i>	ATCC® 27853	MHA+5 %SB	18	24	21	21	21	21	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA-F	25	31	28	28	28	28	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3111881	<i>Escherichia coli</i>	ATCC® 25923	MHA	17	23	20	20	20	20	0 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC® 25922	MHA	17	23	20	22	22	22	0 ± 0 (CV=0)
	<i>Enterococcus faecalis</i>	ATCC® 27853	MHA+5 %SB	18	24	21	22	22	22	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA-F	25	31	28	28	28	28	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3108139	<i>Escherichia coli</i>	ATCC® 25923	MHA	17	23	20	20	20	20	0 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC® 25922	MHA	17	23	20	22	22	22	0 ± 0 (CV=0)
	<i>Enterococcus faecalis</i>	ATCC® 27853	MHA+5 %SB	18	24	21	22	22	22	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA-F	25	31	28	29	29	29	0 ± 0 (CV=0)

Table 4. Results summary for tested QC isolates Nitrofurantoin F300 (300µg) Antimicrobial Susceptibility Disks (CT0036B), in accordance with CLSI methodology

Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3108140	<i>Escherichia coli</i>	ATCC® 25922	MHA	20	25	22.5	23	23	23	0.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	18	22	20	21	21	21	1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA+5%SB	23	29	26	29	29	29	3 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3104594	<i>Escherichia coli</i>	ATCC® 25922	MHA	20	25	22.5	23	23	23	0.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	18	22	20	22	22	22	2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA+5%SB	23	29	26	26	26	26	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3001685	<i>Escherichia coli</i>	ATCC® 25922	MHA	20	25	22.5	21	21	21	-1.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	18	22	20	22	22	22	2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA+5%SB	23	29	26	26	26	26	0 ± 0 (CV=0)

Clinical performance characteristics

A study by O'Halloran et al. 2018, aimed to determine the comparability of CLSI, EUCAST and Stokes antimicrobial susceptibility testing (AST) methods for determining susceptibility of uropathogenic *Escherichia coli* to a variety of antimicrobial agents including, Nitrofurantoin⁴. During the study period (January-February 2017), urine samples were collected, from samples submitted to Cork University Hospital laboratory, from patients attending GP surgeries. Overall, 100 *Escherichia coli* isolates were obtained. The CLSI, EUCAST and Stokes methods were performed using Oxoid AST discs including, Nitrofurantoin F100 [100µg] (CT0034B), and F300 [300µg] (CT0036B) Antimicrobial Susceptibility Disks. *E. coli* ATCC®25922 was used as a quality control strain. It was concluded that nitrofurantoin showed 97% level of accordance across the CLSI and EUCAST methods.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- O'Halloran, C., N. Walsh, M. C. O'Grady, L. Barry, C. Hooton, G. D. Corcoran, and B. Lucey. 2018. "Assessment of the Comparability of CLSI, EUCAST and Stokes Antimicrobial Susceptibility Profiles for *Escherichia Coli* Uropathogenic Isolates." *British Journal of Biomedical Science* 75 (1): 24–29. <https://doi.org/10.1080/09674845.2017.1392736>.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier

The symbol consists of a large 'R' with a diagonal line through it, and the word 'only' in a smaller font to the right.	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
The CE mark is a stylized 'C' followed by a 'E'.	European Conformity Mark
The UKCA mark consists of the letters 'UK' stacked above 'CA'.	UK Conformity Mark

Revision information

Version	Date of issue and modifications introduced
1.0	2023-01-17

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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Viited

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. O'Halloran, C., N. Walsh, M. C. O'Grady, L. Barry, C. Hooton, G. D. Corcoran, and B. Lucey. 2018. "Assessment of the Comparability of CLSI, EUCAST and Stokes Antimicrobial Susceptibility Profiles for *Escherichia Coli* Uropathogenic Isolates." *British Journal of Biomedical Science* 75 (1): 24–29. <https://doi.org/10.1080/09674845.2017.1392736>.

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Referenzen

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- O'Halloran, C., N. Walsh, M. C. O'Grady, L. Barry, C. Hooton, G. D. Corcoran, and B. Lucey. 2018. "Assessment of the Comparability of CLSI, EUCAST and Stokes Antimicrobial Susceptibility Profiles for *Escherichia coli* Uropathogenic Isolates." *British Journal of Biomedical Science* 75 (1): 24–29. <https://doi.org/10.1080/09674845.2017.1392736>.

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Referencie

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- O'Halloran, C., N. Walsh, M. C. O'Grady, L. Barry, C. Hooton, G. D. Corcoran, and B. Lucey. 2018. "Assessment of the Comparability of CLSI, EUCAST and Stokes Antimicrobial Susceptibility Profiles for Escherichia Coli Uropathogenic Isolates." British Journal of Biomedical Science 75 (1): 24–29. <https://doi.org/10.1080/09674845.2017.1392736>.

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© 2023 Thermo Fisher Scientific Inc. Всички права запазени. ATCC® е търговска марка на ATCC. Всички други търговски марки са собственост на Thermo Fisher Scientific Inc. и неговите дъщерни дружества. Тази информация няма за цел да насърчава използването на тези продукти по начин, който може да наруши правата на интелектуална собственост на други.