	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Oxacillin Discs 1 µg (OX1)

REF CT0159B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Oxacillin Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent oxacillin. The discs are labelled on both sides with details of the antimicrobial (OX) and amount present (µg): OX1 (1µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Oxacillin AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which oxacillin has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling, and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100,^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- *Staphylococcus* spp.

Gram-Negative

- N/A

Species with published breakpoints according to current CLSI M100 literature:

Gram-Positive

- *Staphylococcus epidermidis* (30 µg cefoxitin surrogate test for oxacillin)
- *Staphylococcus pseudintermedius*
- *Staphylococcus schleiferi*
- *Staphylococcus lugdunensis* (30 µg cefoxitin surrogate test for oxacillin)
- *Staphylococcus aureus* (30 µg cefoxitin surrogate test for oxacillin).
- Other *Staphylococcus* spp., excluding *S. aureus*, *S. lugdunensis*, *S. epidermidis*, *S. pseudintermedius*, and *S. schleiferi* (30 µg cefoxitin surrogate test for oxacillin).
- *Streptococcus pneumoniae*

Gram-Negative

- N/A

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- *Staphylococcus pseudintermedius* (for screening for methicillin resistance only)
- *Staphylococcus schleiferi* (for screening for methicillin resistance only)
- *Streptococcus pneumoniae* (for screening for beta-lactam resistance only)

Gram-Negative

- N/A

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Oxacillin AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology please refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Oxacillin AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Oxacillin AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Oxacillin AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with CT0159B	
Component Description	Material Description
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Oxacillin AST Disc Reagent	
Reagent	Description of Function
Oxacillin	White crystalline powder. Oxacillin is a water-soluble penicillinase resistant β -lactam antimicrobial agent utilized for the treatment of staphylococci infections. The agent exerts its effects by binding to penicillin-binding proteins inside the bacterial cell wall, subsequently inhibiting the synthesis of the cell wall, ultimately resulting in cell lysis.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Oxacillin AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Raw data of zone sizes, in accordance with CLSI methodology, taken from the 3 most recent batches of Oxacillin (1 µg) OX1 Antimicrobial Susceptibility Test Discs (CT0159B).

CLSI Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3100812	<i>Staphylococcus aureus</i>	25923	MHA	18	24	21	19	19	19	2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA+ 5%SB	0	12	6	11	11	11	5 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV) (%)
							1	2	3	
2981452	<i>Staphylococcus aureus</i>	25923	MHA	18	24	21	24	24	24	3 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA+ 5%SB	0	12	6	12	12	12	6 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV) (%)
							1	2	3	
2964633	<i>Staphylococcus aureus</i>	25923	MHA	18	24	21	21	21	21	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA+ 5%SB	0	12	6	10	10	10	4 ± 0 (CV=0)

Table 4. Raw data of zone sizes, in accordance with EUCAST methodology, taken from the 3 most recent batches of Oxacillin (1 µg) OX1 Antimicrobial Susceptibility Test Discs (CT0159B).

EUCAST Methodology										
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3100812	<i>Streptococcus pneumoniae</i>	49619	MHA+5 %SB	8	14	11	13	13	13	2 ± 0 (CV=0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV) (%)
							1	2	3	
2981452	<i>Streptococcus pneumoniae</i>	49619	MHA+5 %SB	8	14	11	13	13	13	2 ± 0 (CV=0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV) (%)
							1	2	3	
2964633	<i>Streptococcus pneumoniae</i>	49619	MHA+5 %SB	8	14	11	11	11	11	0 ± 0 (CV=0)

Clinical Performance Characteristics

The purpose of this study was to identify and evaluate the performance of antimicrobial susceptibility test discs from nine manufacturers in 2014 and 2017. The study demonstrated that more than 99% of all the Oxoid AST Discs were within the QC range in 2014 and 100% in 2017. More specifically, 100% of the Oxoid Oxacillin (1 µg) OX1 AST Discs (CT0159B) from both years were within ± 1 mm of the QC target for all organisms. The Oxoid Antimicrobial Susceptibility Test Discs performed excellently and consistently over the three years. This study demonstrated that the Oxoid Oxacillin (1 µg) OX1 AST Discs are an accurate and reliable method for the in vitro analysis of bacteria and are more competitively attractive than other manufactures.⁴

Another study compared methods for the identification of oxacillin-resistant *Staphylococcus aureus* (ORSA) in clinical samples.⁵ A total of 500 *S. aureus* isolates were collected. The results of the oxacillin susceptibility testing demonstrated that 103 (20.6%) isolates were resistant to oxacillin according to CLSI breakpoints. Of these 103, the PCR identification of the *mecA* gene confirmed 88 of isolates to be ORSA. When testing the current diagnostic method routinely used in the hospitals, all of the isolates were oxacillin susceptible. This study demonstrated that the use of Oxoid Oxacillin AST Discs provides more accurate detection of ORSA than current methods used in the hospitals, highlighting the need to review the current routinely used ORSA identification methods. The study also concluded that the identification of the *mecA* gene can exhibit varied results as ORSA isolates have been shown to be *mecA* negative. This study, therefore, demonstrated the Oxoid Oxacillin AST Discs can be used in the identification of oxacillin-resistant *S. aureus* in clinical samples.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References


- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing

3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. J. Ahman, E. Matuschek and G. Kahlmeter (2019), The Quality of Antimicrobial Discs from Nine Manufacturers – EUCAST Evaluations in 2014 and 2017, *Clinical Microbiology and Infection* 25(3): 346-352. DOI: 10.1016/j.cmi.2018.05.021.
5. A M. Momenah et al., (2005) Accuracy of Current Oxacillin Sensitivity Tests Routinely Used in Hospitals in Western Saudi Arabia, *Saudi Medical Journal* 26 (4): 584-587. ISSN: 03795284.

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision Information

Oxoid Limited, Wade Road, Basingstoke,
Hampshire, RG24 8PW, UK

Version	Date of issue and modifications introduced
1.0	2023-01-17

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).

4. J. Ahman, E. Matuschek and G. Kahlmeter (2019), The Quality of Antimicrobial Discs from Nine Manufacturers – EUCAST Evaluations in 2014 and 2017, *Clinical Microbiology and Infection* 25(3): 346-352. DOI: 10.1016/j.cmi.2018.05.021.
5. A M. Momenah et al., (2005) Accuracy of Current Oxacillin Sensitivity Tests Routinely Used in Hospitals in Western Saudi Arabia, *Saudi Medical Journal* 26 (4): 584-587. ISSN: 03795284.

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Reference

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:

- a. EUCAST Disk Diffusion Method
- b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. J. Ahman, E. Matuschek and G. Kahlmeter (2019), The Quality of Antimicrobial Discs from Nine Manufacturers – EUCAST Evaluations in 2014 and 2017, *Clinical Microbiology and Infection* 25(3): 346-352. DOI: 10.1016/j.cmi.2018.05.021.
5. A M. Momenah et al., (2005) Accuracy of Current Oxacillin Sensitivity Tests Routinely Used in Hospitals in Western Saudi Arabia, *Saudi Medical Journal* 26 (4): 584-587. ISSN: 03795284.

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Litteratur

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- J. Ahman, E. Matuschek and G. Kahlmeter (2019), The Quality of Antimicrobial Discs from Nine Manufacturers – EUCAST Evaluations in 2014 and 2017, *Clinical Microbiology and Infection* 25(3): 346-352. DOI: 10.1016/j.cmi.2018.05.021.
- A M. Momenah et al., (2005) Accuracy of Current Oxacillin Sensitivity Tests Routinely Used in Hospitals in Western Saudi Arabia, *Saudi Medical Journal* 26 (4): 584-587. ISSN: 03795284.

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Viited

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. J. Ahman, E. Matuschek and G. Kahlmeter (2019), The Quality of Antimicrobial Discs from Nine Manufacturers – EUCAST Evaluations in 2014 and 2017, *Clinical Microbiology and Infection* 25(3): 346-352. DOI: 10.1016/j.cmi.2018.05.021.

5. A M. Momenah et al., (2005) Accuracy of Current Oxacillin Sensitivity Tests Routinely Used in Hospitals in Western Saudi Arabia, *Saudi Medical Journal* 26 (4): 584-587. ISSN: 03795284.

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Referenzen

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing

- c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. J. Ahman, E. Matuschek and G. Kahlmeter (2019), The Quality of Antimicrobial Discs from Nine Manufacturers – EUCAST Evaluations in 2014 and 2017, *Clinical Microbiology and Infection* 25(3): 346-352. DOI: 10.1016/j.cmi.2018.05.021.
5. A M. Momenah et al., (2005) Accuracy of Current Oxacillin Sensitivity Tests Routinely Used in Hospitals in Western Saudi Arabia, *Saudi Medical Journal* 26 (4): 584-587. ISSN: 03795284.

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Referințe

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.

2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. J. Ahman, E. Matuschek and G. Kahlmeter (2019), The Quality of Antimicrobial Discs from Nine Manufacturers – EUCAST Evaluations in 2014 and 2017, *Clinical Microbiology and Infection* 25(3): 346-352. DOI: 10.1016/j.cmi.2018.05.021.
5. A M. Momenah et al., (2005) Accuracy of Current Oxacillin Sensitivity Tests Routinely Used in Hospitals in Western Saudi Arabia, *Saudi Medical Journal* 26 (4): 584-587. ISSN: 03795284.

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Referenser

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- J. Ahman, E. Matuschek and G. Kahlmeter (2019), The Quality of Antimicrobial Discs from Nine Manufacturers – EUCAST Evaluations in 2014 and 2017, *Clinical Microbiology and Infection* 25(3): 346-352. DOI: 10.1016/j.cmi.2018.05.021.

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