	https://www.thermofisher.com	
		
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Rifampicin Discs 5 µg (RD5)

REF CT0207B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Rifampicin Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent rifampicin. The discs are labelled on both sides with details of the antimicrobial (RD) and amount present (µg): RD5 (5µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Rifampicin AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Rifampicin AST Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which rifampicin has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture, Oxoid Rifampicin Discs are for professional use only.

The test provides information to categorise organisms as either resistant, intermediate, or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling, and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100,^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- N/A

Gram-Negative

- *Neisseria meningitidis*

Species with published breakpoints according to current CLSI M100 literature:

Gram-Positive

- *Staphylococcus* spp.
- *Enterococcus* spp.
- *Streptococcus pneumoniae*

Gram-Negative

- *Haemophilus influenzae*
- *Haemophilus parainfluenzae*
- *Neisseria meningitidis*

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- *Staphylococcus* spp.
- *Streptococcus* groups A, B, C and G
- *Streptococcus pneumoniae*
- *Corynebacterium* spp.
- *Aerococcus sanguinicola* and *urinae*

Gram-Negative

- *Haemophilus influenzae*
- *Kingella kingae*

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Rifampicin AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology, refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Rifampicin AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Rifampicin AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Rifampicin AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with CT0207B	
Component Description	Material Description
Cartridge with spring, cap and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge, 5 cartridges per pack.

Table 2. Description of Rifampicin AST Disc Reagent	
Reagent	Description of Function
Rifampicin	Rifampicin is an antibiotic with a broad antibacterial spectrum: it acts suppressing the initiation of RNA synthesis by forming a stable complex with the DNA-dependent RNA polymerase enzyme. It has a bactericidal action on a both intracellular and extracellular level.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Rifampicin AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Rifampicin AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Raw data of zone sizes, in accordance with CLSI methodology, taken from the three most recent batches of Rifampicin (5 µg) RD5 Antimicrobial Susceptibility Test Discs (CT0207B).

CLSI Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3315571	<i>Staphylococcus aureus</i>	25923	MHA	26	34	30	32	32	32	2 ± 0 (CV=0)
	<i>Escherichia coli</i>	25922	MHA	8	10	9	8	8	8	-1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	22	30	26	23	23	23	-3 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA+5 %SB	25	30	27.5	28	28	28	0.5 ± 0 (CV=0)
	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
3274008	<i>Staphylococcus aureus</i>	25923	MHA	26	34	30	29	29	29	-1 ± 0 (CV=0)
	<i>Escherichia coli</i>	25922	MHA	8	10	9	8	8	8	-1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	22	30	26	25	25	25	-1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA+5 %SB	25	30	27.5	27	27	28	-0.2 ± 0.6 (CV=-3.5)
	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
3245239	<i>Staphylococcus aureus</i>	25923	MHA	26	34	30	31	31	31	1 ± 0 (CV=0)
	<i>Escherichia coli</i>	25922	MHA	8	10	9	8	8	8	-1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	22	30	26	22	22	22	-4 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA+5 %SB	25	30	27.5	26	26	26	-1.5 ± 0 (CV=0)
	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)

Table 4. Raw data of zone sizes, in accordance with EUCAST methodology, taken from the three most recent batches of Rifampicin (5 µg) RD5 Antimicrobial Susceptibility Test Discs (CT0207B).

EUCAST Methodology										
Product Batch	Organism	ATCC® umber	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid- point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3315571	<i>Staphylococcus aureus</i>	29213	MHA	30	36	33	32	32	32	-1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49766	MHA-F	21	27	24	24	24	24	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA-F	26	32	29	28	28	28	-1 ± 0 (CV=0)
	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
3274008	<i>Staphylococcus aureus</i>	29213	MHA	30	36	33	33	33	33	0 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49766	MHA-F	21	27	24	25	25	25	1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA-F	26	32	29	29	29	29	0 ± 0 (CV=0)
	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
3245239	<i>Staphylococcus aureus</i>	29213	MHA	30	36	33	35	35	35	2 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49766	MHA-F	21	27	24	23	23	23	-1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA-F	26	32	29	28	28	28	-1 ± 0 (CV=0)
	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)

Clinical Performance Characteristics

A study investigated antimicrobial susceptibility patterns for methicillin-resistant *Staphylococcus aureus* (MRSA) isolates taken from various body parts. 58 samples were confirmed positive for *S. aureus*, then classified according to their clinical specimens, and tested for resistance to methicillin.⁴ Antimicrobial testing was performed through disc diffusion using 10 antibiotics including Oxoid Rifampicin RD5 AST Discs. The study confirmed the majority of MRSA isolates came from open wounds, and that isolates that showed resistance against other antibiotics, present in greater numbers in the MRSA group rather than the methicillin-susceptible staph (MSSA) group. Oxoid Rifampicin AST Discs proved to be a valuable tool for in vitro analysis of MRSA susceptibility and monitoring of antimicrobial resistance in MRSA isolates.

Another study analysed antimicrobial susceptibility patterns of methicillin-resistant *Staphylococcus aureus* nasal isolates.⁵ Nasal swab samples were collected from 182 medical students. 32 samples (control group) were from students not exposed to the hospital environment. Isolates were selected for *S. aureus*. It was found that the control group students were carrying only MSSA isolates, which confirmed the hypothesis of 'acquired' MRSA isolates following exposure in a contaminated environment.

Antimicrobial susceptibility testing using the disc diffusion method proved that the *S. aureus* strains were susceptible to four antibiotics, including the Oxoid Rifampicin RD5 AST Discs.

An investigation was conducted to study antibiotic resistance on *Helicobacter pylori* isolates collected from antral biopsy specimens.⁶ Antimicrobial susceptibility disc diffusion testing was performed on the isolates positive for *H. pylori* against amoxicillin, tetracycline, clarithromycin, metronidazole, and rifampicin. For rifampicin, Oxoid Rifampicin RD5 AST Discs were used. Minimum Inhibitory Concentrations (MIC) for each antibiotic were determined through Etests. The testing demonstrated increasing resistance to metronidazole, however other antibiotics, including rifampicin, amoxicillin, or tetracycline were shown to be a good choice of treatment. The Oxoid Rifampicin AST Discs were proven to be a valuable tool for in vitro analysis and monitoring of *H. pylori* susceptibility and for surveillance of antimicrobial resistance.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority on which the user and/or the patient is established.

References

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. A. A. Agboke and A. A. Attama (2015). Antibiotics Susceptibility Pattern of Methicillin-Resistant *Staphylococcus aureus* (MRSA) In Enugu State, South-East Region of Nigeria. *BMRJ*, 10(1): 1-7, 2015; DOI: [10.9734/BMRJ/2015/18219](https://doi.org/10.9734/BMRJ/2015/18219).
5. S. A. Zakai (2015). Prevalence of methicillin-resistant *Staphylococcus aureus* nasal colonization among medical students in Jeddah, Saudi Arabia. *Saudi Med J* 2015; Vol. 36 (7): 807-812. DOI: [10.15537/smj.2015.7.11609](https://doi.org/10.15537/smj.2015.7.11609).
6. A. I. Al-Omar, M. Arjunan and A. A. Al-Salamah (2013). Prevalence of Antibiotic Resistance and Minimum Inhibitory Concentrations of *Helicobacter pylori* Isolates from Clinical Gastritis Patients, Saudi Arabia. *Global Journal of Gastroenterology and Hepatology*, 2013; 1, 101-1062013. E-ISSN: 2308-6483/13.

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date

	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision Information

Version	Date of issue and modifications introduced
1.0	2023-01-18

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. A. A. Agboke and A. A. Attama (2015). Antibiotics Susceptibility Pattern of Methicillin-Resistant *Staphylococcus aureus* (MRSA) In Enugu State, South-East Region of Nigeria. *BMRJ*, 10(1): 1-7, 2015; DOI: [10.9734/BMRJ/2015/18219](https://doi.org/10.9734/BMRJ/2015/18219).
5. S. A. Zakai (2015). Prevalence of methicillin-resistant *Staphylococcus aureus* nasal colonization among medical students in Jeddah, Saudi Arabia. *Saudi Med J* 2015; Vol. 36 (7): 807-812. DOI: 10.15537/smj.2015.7.11609.
6. A. I. Al-Omar, M. Arjunan and A. A. Al-Salamah (2013). Prevalence of Antibiotic Resistance and Minimum Inhibitory Concentrations of *Helicobacter pylori* Isolates from Clinical Gastritis Patients, Saudi Arabia. *Global Journal of Gastroenterology and Hepatology*, 2013; 1, 101-106/2013. E-ISSN: 2308-6483/13.

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Referencije

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. A. A. Agboke and A. A. Attama (2015). Antibiotics Susceptibility Pattern of Methicillin-Resistant *Staphylococcus aureus* (MRSA) In Enugu State, South-East Region of Nigeria. *BMRJ*, 10(1): 1-7, 2015; DOI: [10.9734/BMRJ/2015/18219](https://doi.org/10.9734/BMRJ/2015/18219).
5. S. A. Zakai (2015). Prevalence of methicillin-resistant *Staphylococcus aureus* nasal colonization among medical students in Jeddah, Saudi Arabia. *Saudi Med J* 2015; Vol. 36 (7): 807-812. DOI: 10.15537/smj.2015.7.11609.
6. A. I. Al-Omar, M. Arjunan and A. A. Al-Salamah (2013). Prevalence of Antibiotic Resistance and Minimum Inhibitory Concentrations of *Helicobacter pylori* Isolates from Clinical Gastritis Patients, Saudi Arabia. *Global Journal of Gastroenterology and Hepatology*, 2013; 1, 101-106. E-ISSN: 2308-6483/13.

© 2023. Thermo Fisher Scientific Inc. Sva prava pridržana. ATCC® je zaštitni znak društva ATCC. Svi ostali zaštitni znakovi vlasništvo su društva Thermo Fisher Scientific Inc. i njegovih podružnica. Ove informacije nisu namijenjene za poticanje na uporabu tih proizvoda na način kojim se mogu povrijediti prava intelektualnog vlasništva drugih društava.

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Reference

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. A. A. Agboke and A. A. Attama (2015). Antibiotics Susceptibility Pattern of Methicillin-Resistant *Staphylococcus aureus* (MRSA) In Enugu State, South-East Region of Nigeria. *BMRJ*, 10(1): 1-7, 2015; DOI: [10.9734/BMRJ/2015/18219](https://doi.org/10.9734/BMRJ/2015/18219).
5. S. A. Zakai (2015). Prevalence of methicillin-resistant *Staphylococcus aureus* nasal colonization among medical students in Jeddah, Saudi Arabia. *Saudi Med J* 2015; Vol. 36 (7): 807-812. DOI: 10.15537/smj.2015.7.11609.
6. A. I. Al-Omar, M. Arjunan and A. A. Al-Salamah (2013). Prevalence of Antibiotic Resistance and Minimum Inhibitory Concentrations of *Helicobacter pylori* Isolates from Clinical Gastritis Patients, Saudi Arabia. *Global Journal of Gastroenterology and Hepatology*, 2013; 1, 101-1062013. E-ISSN: 2308-6483/13.

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Atsauces

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. A. A. Agboke and A. A. Attama (2015). Antibiotics Susceptibility Pattern of Methicillin-Resistant *Staphylococcus aureus* (MRSA) In Enugu State, South-East Region of Nigeria. *BMRJ*, 10(1): 1-7, 2015; DOI: [10.9734/BMRJ/2015/18219](https://doi.org/10.9734/BMRJ/2015/18219).
5. S. A. Zakai (2015). Prevalence of methicillin-resistant *Staphylococcus aureus* nasal colonization among medical students in Jeddah, Saudi Arabia. *Saudi Med J* 2015; Vol. 36 (7): 807-812. DOI: [10.15537/smj.2015.7.11609](https://doi.org/10.15537/smj.2015.7.11609).
6. A. I. Al-Omar, M. Arjunan and A. A. Al-Salamah (2013). Prevalence of Antibiotic Resistance and Minimum Inhibitory Concentrations of *Helicobacter pylori* Isolates from Clinical Gastritis Patients, Saudi Arabia. *Global Journal of Gastroenterology and Hepatology*, 2013; 1, 101-106. E-ISSN: 2308-6483/13.

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1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. A. A. Agboke and A. A. Attama (2015). Antibiotics Susceptibility Pattern of Methicillin-Resistant *Staphylococcus aureus* (MRSA) In Enugu State, South-East Region of Nigeria. *BMRJ*, 10(1): 1-7, 2015; DOI: [10.9734/BMRJ/2015/18219](https://doi.org/10.9734/BMRJ/2015/18219).
5. S. A. Zakai (2015). Prevalence of methicillin-resistant *Staphylococcus aureus* nasal colonization among medical students in Jeddah, Saudi Arabia. *Saudi Med J* 2015; Vol. 36 (7): 807-812. DOI: 10.15537/smj.2015.7.11609.
6. A. I. Al-Omar, M. Arjunan and A. A. Al-Salamah (2013). Prevalence of Antibiotic Resistance and Minimum Inhibitory Concentrations of *Helicobacter pylori* Isolates from Clinical Gastritis Patients, Saudi Arabia. *Global Journal of Gastroenterology and Hepatology*, 2013; 1, 101-106 2013. E-ISSN: 2308-6483/13.

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Referencie

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- A. A. Agboke and A. A. Attama (2015). Antibiotics Susceptibility Pattern of Methicillin-Resistant *Staphylococcus aureus* (MRSA) In Enugu State, South-East Region of Nigeria. *BMRJ*, 10(1): 1-7, 2015; DOI: [10.9734/BMRJ/2015/18219](https://doi.org/10.9734/BMRJ/2015/18219).
- S. A. Zakai (2015). Prevalence of methicillin-resistant *Staphylococcus aureus* nasal colonization among medical students in Jeddah, Saudi Arabia. *Saudi Med J* 2015; Vol. 36 (7): 807-812. DOI: 10.15537/smj.2015.7.11609.
- A. I. Al-Omar, M. Arjunan and A. A. Al-Salamah (2013). Prevalence of Antibiotic Resistance and Minimum Inhibitory Concentrations of *Helicobacter pylori* Isolates from Clinical Gastritis Patients, Saudi Arabia. *Global Journal of Gastroenterology and Hepatology*, 2013; 1, 101-1062013. E-ISSN: 2308-6483/13.

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Referenser

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. A. A. Agboke and A. A. Attama (2015). Antibiotics Susceptibility Pattern of Methicillin-Resistant *Staphylococcus aureus* (MRSA) In Enugu State, South-East Region of Nigeria. BMRJ, 10(1): 1-7, 2015; DOI: 10.9734/BMRJ/2015/18219.
5. S. A. Zakai (2015). Prevalence of methicillin-resistant *Staphylococcus aureus* nasal colonization among medical students in Jeddah, Saudi Arabia. Saudi Med J 2015; Vol. 36 (7): 807-812. DOI: 10.15537/smj.2015.7.11609.
6. A. I. Al-Omar, M. Arjunan and A. A. Al-Salamah (2013). Prevalence of Antibiotic Resistance and Minimum Inhibitory Concentrations of *Helicobacter pylori* Isolates from Clinical Gastritis Patients, Saudi Arabia. Global Journal of Gastroenterology and Hepatology, 2013; 1, 101-106. E-ISSN: 2308-6483/13.

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