
	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Teicoplanin Discs 30 µg (TEC30)

REF CT0647B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Teicoplanin Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs that contain specific amounts of the antimicrobial agent teicoplanin. The discs are labelled on both sides with details of the antimicrobial (TEC) and amount present (µg): TEC30.

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Teicoplanin AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Teicoplanin AST Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which teicoplanin has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- N/A

Gram-Negative

- N/A

Species with published breakpoints according to current CLSI M100 literature (all associated with 30 µg):

Gram-Positive

- *Enterococcus* spp.

Gram-Negative

- N/A

Species with published breakpoints according to current EUCAST literature (all associated with 30 µg):

Gram-Positive

- *Staphylococcus* spp. (*S. aureus*, Coagulase-negative staphylococci)
- *Enterococcus* spp.
- *Streptococcus* spp.
 - Group A, B, C and G
 - *Streptococcus pneumoniae*
 - Viridans group streptococci.

Gram-Negative

- N/A

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principles of Method

Teicoplanin AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology please refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Teicoplanin AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Teicoplanin AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Teicoplanin AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with CT0647B	
Component Description	Material Description
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table2. Description of Teicoplanin AST Disc Reagents	
Reagent	Description of Function
Teicoplanin	White to off-white crystalline powder. Antimicrobial agent. The drug mechanism of action involves the inhibition of peptidoglycan polymerisation therefore inhibition to cell wall synthesis resulting in cell death.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Teicoplanin AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Teicoplanin AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Raw data of zone sizes, in accordance with CLSI methodology, taken from the 3 most recent batches of Teicoplanin (30 µg) TEC30 Antimicrobial Susceptibility Test Discs (CT0647B).

CLSI Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	

3331956	<i>Staphylococcus aureus</i>	ATCC 25923	MHA	15	21	18	16	16	16	2 ± 0 (CV = 0)
Product Batch	Organism	ATCC number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3261723	<i>Staphylococcus aureus</i>	ATCC 25923	MHA	15	21	18	16	16	16	2 ± 0 (CV = 0)
Product Batch	Organism	ATCC number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
2475150	<i>Staphylococcus aureus</i>	ATCC 25923	MHA	15	21	18	16	16	16	2 ± 0 (CV = 0)

SD: Standard Deviation

Table 4. Raw data of zone sizes, in accordance with EUCAST methodology, taken from the 3 most recent batches of Teicoplanin (30 µg) TEC30 Antimicrobial Susceptibility Test Discs (CT0647B).

EUCAST Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3331956	<i>Enterococcus faecalis</i>	ATCC 29212	MHA	15	21	18	17	17	17	1 ± 0 (CV = 0)
	<i>Streptococcus pneumoniae</i>	ATCC 49619	MHA-F	18	24	21	19	19	19	2 ± 0 (CV = 0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3261723	<i>Enterococcus faecalis</i>	ATCC 29212	MHA	15	21	18	16	16	16	2 ± 0 (CV = 0)
	<i>Streptococcus pneumoniae</i>	ATCC 49619	MHA-F	18	24	21	20	20	20	1 ± 0 (CV = 0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	

2475150	<i>Enterococcus faecalis</i>	ATCC 29212	MHA	15	21	18	16	16	16	2 ± 0 (CV = 0)
	<i>Streptococcus pneumoniae</i>	ATCC 49619	MHA-F	18	24	21	20	20	20	1 ± 0 (CV = 0)

Clinical Performance Characteristics

A study collected biological samples from three hospitals in three different cities in Pakistan between February and November 2011.⁴ Only isolates containing *S. aureus* or *Enterococci* (identified by cultural, morphological, and biochemical properties) were tested for their susceptibility to teicoplanin by using Oxoid Teicoplanin (30 µg) TEC30AST Discs. All of the total 401 isolates recovered were sensitive to teicoplanin regardless of methicillin resistance or sensitivity in *S. aureus*. No out-of-specification results occurred whilst using Oxoid Teicoplanin AST Discs demonstrating the device's effectiveness for in vitro analysis of the *S. aureus* and *Enterococci* resistance profile.

Idelevich et al. (2016) aimed to test ADAGIO automated system for reading the inhibition zone diameters of fastidious bacteria for disc diffusion antimicrobial susceptibility testing.⁵ A range of clinical isolates were used in the testing making up a total of 144 clinical isolates. The study concluded that the ADAGIO system provided good categorical agreement with manual reading (0.8% very major errors and 3.0% major errors after visual adjustment). 100% of the bacterial species tested by Teicoplanin (30 µg) TEC30 AST Discs (*S. agalactiae*, *S. dysgalactiae*, *S. pyogenes*, *S. pneumoniae*, Viridans group streptococci) were shown to be susceptible to teicoplanin.

Kot et al., 2020, aimed to study the antimicrobial resistance profiles of 112 MRSA isolates collected from patients hospitalised during 2015 (18 isolates), 2016 (23 isolates), and 2017 (71 isolates). *S. aureus* isolates were confirmed by Gram staining, catalase activity, coagulase, and PCR analysis of the *nuc* gene. Teicoplanin (30 µg) TEC30 AST Discs among others from Oxoid were utilised for antimicrobial susceptibility testing. The study showed that the percentage of MRSA resistant to levofloxacin, ciprofloxacin, clindamycin and erythromycin significantly increased at the end of the study period ($p = 0.002$), however no (0%) isolates were resistant to teicoplanin. No out-of-specification results whilst using Teicoplanin (30 µg) TEC30AST Discs were reported, thus demonstrating the device's effectiveness in the in vitro analysis of *S. aureus* resistance profiles.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Ahmed, A., Hussain, S., Ijaz, T., & Hashemy, I. (2014). Susceptibility of methicillin-resistant staphylococcus aureus and enterococci to teicoplanin in Pakistan: The MRSET study. *Journal of the Pakistan Medical Association*, 64(3), 256–259.
- Idelevich, E. A., Becker, K., Schmitz, J., Knaack, D., Peters, G., & Köck, R. (2016b). Evaluation of an automated system for reading and interpreting disk diffusion antimicrobial susceptibility testing of fastidious bacteria. *PLoS ONE*, 11(7), 8–15. <https://doi.org/10.1371/journal.pone.0159183>
- Kot, B., Wierchowska, K., Piechota, M., & Gruzewska, A. (2019). Antimicrobial Resistance Patterns in Methicillin-Resistant *Staphylococcus aureus* from Patients Hospitalized during 2015-2017 in Hospitals in Poland. *Medical Principles and Practice*, 29(1), 61–68. <https://doi.org/10.1159/000501788>.

© 2021 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.


Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit

	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision Information

Version	Date of issue and modifications introduced
1.0	2022-11-01. New document

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
Literatura

- Clinical & Laboratory Standards Institute (CLSI). Proučite najnovije izdanje:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Proučite najnoviju verziju:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Ahmed, A., Hussain, S., Ijaz, T., & Hashemy, I. (2014). Susceptibility of methicillin-resistant staphylococcus aureus and enterococci to teicoplanin in Pakistan: The MRSET study. *Journal of the Pakistan Medical Association*, 64(3), 256–259.
- Idelevich, E. A., Becker, K., Schmitz, J., Knaack, D., Peters, G., & Köck, R. (2016b). Evaluation of an automated system for reading and interpreting disk diffusion antimicrobial susceptibility testing of fastidious bacteria. *PLoS ONE*, 11(7), 8–15. <https://doi.org/10.1371/journal.pone.0159183>
- Kot, B., Wierzbowska, K., Piechota, M., & Gruzewska, A. (2019). Antimicrobial Resistance Patterns in Methicillin-Resistant Staphylococcus aureus from Patients Hospitalized during 2015-2017 in Hospitals in Poland. *Medical Principles and Practice*, 29(1), 61–68. <https://doi.org/10.1159/000501788>.

© 2021 Thermo Fisher Scientific Inc. Sva su prava pridržana. Oznaka ATCC® zaštitni je znak Američke zbirke tipskih kultura (ATCC). Svi ostali zaštitni znakovi vlasništvo su društva Thermo Fisher Scientific Inc. i njegovih društava kćeri. Ove informacije nisu namijenjene poticanju upotrebe ovih proizvoda na bilo koji način koji bi mogao dovesti do povrede prava intelektualnog vlasništva drugih.

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Referencer

- Clinical & Laboratory Standards Institute (CLSI). Se den seneste udgave:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Se den seneste version:
 - EUCAST -agar-diffusionsmetode
 - Læsevejledning, EUCAST -agar-diffusionsmetode til antimikrobiel modtagelighedstestning
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Ahmed, A., Hussain, S., Ijaz, T., & Hashemy, I. (2014). Susceptibility of methicillin-resistant staphylococcus aureus and enterococci to teicoplanin in Pakistan: The MRSET study. *Journal of the Pakistan Medical Association*, 64(3), 256–259.
- Idelevich, E. A., Becker, K., Schmitz, J., Knaack, D., Peters, G., & Köck, R. (2016b). Evaluation of an automated system for reading and interpreting disk diffusion antimicrobial susceptibility testing of fastidious bacteria. *PLoS ONE*, 11(7), 8–15. <https://doi.org/10.1371/journal.pone.0159183>
- Kot, B., Wierzbowska, K., Piechota, M., & Gruzewska, A. (2019). Antimicrobial Resistance Patterns in Methicillin-Resistant Staphylococcus aureus from Patients Hospitalized during 2015-2017 in Hospitals in Poland. *Medical Principles and Practice*, 29(1), 61–68. <https://doi.org/10.1159/000501788>.

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Literatur

- Clinical & Laboratory Standards Institute (CLSI). Konsultieren Sie die neueste Ausgabe:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Konsultieren Sie die neueste Version:
 - EUCAST Disk Diffusion Method
 - Lesehilfe, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Teil 460 (2005).
- Ahmed, A., Hussain, S., Ijaz, T. & Hashemy, I. (2014). Susceptibility of methicillin-resistant staphylococcus aureus and enterococci to teicoplanin in Pakistan: The MRSET study. *Journal of the Pakistan Medical Association*, 64(3), 256–259.
- Idelevich, E. A., Becker, K., Schmitz, J., Knaack, D., Peters, G., & Köck, R. (2016b). Evaluation of an automated system for reading and interpreting disk diffusion antimicrobial susceptibility testing of fastidious bacteria. *PLoS ONE*, 11(7), 8–15. <https://doi.org/10.1371/journal.pone.0159183>
- Kot, B., Wierchowska, K., Piechota, M. & Gruzewska, A. (2019). Antimicrobial Resistance Patterns in Methicillin-Resistant Staphylococcus aureus from Patients Hospitalized during 2015-2017 in Hospitals in Poland. *Medical Principles and Practice*, 29(1), 61–68. <https://doi.org/10.1159/000501788>.

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1. Clinical & Laboratory Standards Institute (CLSI). Συμβουλευτείτε την τελευταία έκδοση:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Συμβουλευτείτε την τελευταία έκδοση:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Ahmed, A., Hussain, S., Ijaz, T. και Hashemy, I. (2014). Susceptibility of methicillin-resistant staphylococcus aureus and enterococci to teicoplanin in Pakistan: The MRSET study. *Journal of the Pakistan Medical Association*, 64(3), 256–259.
5. Idelevich, E. A., Becker, K., Schmitz, J., Knaack, D., Peters, G. και Köck, R. (2016b). Evaluation of an automated system for reading and interpreting disk diffusion antimicrobial susceptibility testing of fastidious bacteria. *PLoS ONE*, 11(7), 8–15. <https://doi.org/10.1371/journal.pone.0159183>
6. Kot, B., Wierzbowska, K., Piechota, M. και Gruzewska, A. (2019). Antimicrobial Resistance Patterns in Methicillin-Resistant Staphylococcus aureus from Patients Hospitalized during 2015-2017 in Hospitals in Poland. *Medical Principles and Practice*, 29(1), 61–68. <https://doi.org/10.1159/000501788>.

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Bibliografia

- Clinical & Laboratory Standards Institute (CLSI). Consultare l'ultima edizione:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- Comitato europeo sui test di suscettibilità antimicrobica (European Committee on Antimicrobial Susceptibility Testing, EUCAST). Consultare l'ultima versione:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Ahmed, A., Hussain, S., Ijaz, T. e Hashemy, I. (2014). Susceptibility of methicillin-resistant staphylococcus aureus and enterococchi to teicoplanin in Pakistan: The MRSET study. *Journal of the Pakistan Medical Association*, 64 (3), 256–259.
- Idevich, E. A., Becker, K., Schmitz, J., Knaack, D., Peters, G. e Köck, R. (2016b). Evaluation of an automated system for reading and interpreting disk diffusion antimicrobial susceptibility testing of fastidious bacteria. *PLoS ONE*, 11(7), 8–15. <https://doi.org/10.1371/journal.pone.0159183>
- Kot, B., Wierchowska, K., Piechota, M. e Gruzewska, A. (2019). Antimicrobial Resistance Patterns in Methicillin-Resistant Staphylococcus aureus from Patients Hospitalized during 2015-2017 in Hospitals in Poland. *Medical Principles and Practice*, 29 (1), 61–68. <https://doi.org/10.1159/000501788>.

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- Clinical & Laboratory Standards Institute (CLSI). Žr. naujausią leidimą:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Žr. naujausią versiją:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Ahmed, A., Hussain, S., Ijaz, T., & Hashemy, I. (2014). Susceptibility of methicillin-resistant staphylococcus aureus and enterococci to teicoplanin in Pakistan: The MRSET study. *Journal of the Pakistan Medical Association*, 64(3), 256–259.
- Idelevich, E. A., Becker, K., Schmitz, J., Knaack, D., Peters, G., & Köck, R. (2016b). Evaluation of an automated system for reading and interpreting disk diffusion antimicrobial susceptibility testing of fastidious bacteria. *PLoS ONE*, 11(7), 8–15. <https://doi.org/10.1371/journal.pone.0159183>
- Kot, B., Wierchowska, K., Piechota, M., & Gruzewska, A. (2019). Antimicrobial Resistance Patterns in Methicillin-Resistant Staphylococcus aureus from Patients Hospitalized during 2015-2017 in Hospitals in Poland. *Medical Principles and Practice*, 29(1), 61–68. <https://doi.org/10.1159/000501788>.

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5. Idelevich, E. A., Becker, K., Schmitz, J., Knaack, D., Peters, G. i Köck, R. (2016b). Evaluation of an automated system for reading and interpreting disk diffusion antimicrobial susceptibility testing of fastidious bacteria. *PLoS ONE*, 11(7), 8–15. <https://doi.org/10.1371/journal.pone.0159183>
6. Kot, B., Wierzchowska, K., Piechota, M. i Gruzewska, A. (2019). Antimicrobial Resistance Patterns in Methicillin-Resistant *Staphylococcus aureus* from Patients Hospitalized during 2015-2017 in Hospitals in Poland. *Medical Principles and Practice*, 29(1), 61–68. <https://doi.org/10.1159/000501788>.

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Referências

- Clinical & Laboratory Standards Institute (CLSI, Instituto de padrões clínicos e laboratoriais dos EUA). Consultar a edição mais recente:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST, Comité Europeu para o Teste à Suscetibilidade Antimicrobiana). Consultar a versão mais recente:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA, Agência dos alimentos e dos medicamentos dos EUA). CFR Title 21, Volume 5, Part 460 (2005).
- Ahmed, A., Hussain, S., Ijaz, T., & Hashemy, I. (2014). Susceptibility of methicillin-resistant staphylococcus aureus and enterococci to teicoplanin in Pakistan: The MRSET study. *Journal of the Pakistan Medical Association*, 64(3), 256–259.
- Idelevich, E. A., Becker, K., Schmitz, J., Knaack, D., Peters, G., & Köck, R. (2016b). Evaluation of an automated system for reading and interpreting disk diffusion antimicrobial susceptibility testing of fastidious bacteria. *PLoS ONE*, 11(7), 8–15. <https://doi.org/10.1371/journal.pone.0159183>
- Kot, B., Wierchowaska, K., Piechota, M., & Gruzewska, A. (2019). Antimicrobial Resistance Patterns in Methicillin-Resistant Staphylococcus aureus from Patients Hospitalized during 2015-2017 in Hospitals in Poland. *Medical Principles and Practice*, 29(1), 61–68. <https://doi.org/10.1159/000501788>.

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Odkazy

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Hlava 21, Zväzok 5, Časť 460 (2005).
- Ahmed, A., Hussain, S., Ijaz, T. a Hashemy, I. (2014). Susceptibility of methicillin-resistant staphylococcus aureus and enterococci to teicoplanin in Pakistan: The MRSET study. *Journal of the Pakistan Medical Association*, 64(3), 256 – 259.
- Idelevich, E. A., Becker, K., Schmitz, J., Knaack, D., Peters, G., & Köck, R. (2016b). Evaluation of an automated system for reading and interpreting disk diffusion antimicrobial susceptibility testing of fastidious bacteria. *PLoS ONE*, 11(7), 8 – 15. <https://doi.org/10.1371/journal.pone.0159183>
- Kot, B., Wierzychowska, K., Piechota, M., & Gruzewska, A. (2019). Antimicrobial Resistance Patterns in Methicillin-Resistant Staphylococcus aureus from Patients Hospitalized during 2015-2017 in Hospitals in Poland. *Medical Principles and Practice*, 29(1), 61 – 68. <https://doi.org/10.1159/000501788>.

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Referencias

- Clinical & Laboratory Standards Institute (CLSI). Consultar la última edición:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests.
 - M100, Performance Standards for Antimicrobial Susceptibility Testing.
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consultar la última versión:
 - EUCAST Disk Diffusion Method
 - «Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing».

3. Food and Drug Administration (FDA). Título 21 del CFR, volumen 5, parte 460 (2005).
4. Ahmed, A., Hussain, S., Ijaz, T. y Hashemy, I. (2014). «Susceptibility of methicillin-resistant staphylococcus aureus and enterococci to teicoplanin in Pakistan: The MRSET study». *Journal of the Pakistan Medical Association*, 64(3), 256–259.
5. Idelevich, E. A., Becker, K., Schmitz, J., Knaack, D., Peters, G. y Köck, R. (2016b). «Evaluation of an automated system for reading and interpreting disk diffusion antimicrobial susceptibility testing of fastidious bacteria». *PLoS ONE*, 11(7), 8–15. <https://doi.org/10.1371/journal.pone.0159183>
6. Kot, B., Wierchowska, K., Piechota, M. y Gruzewska, A. (2019). «Antimicrobial Resistance Patterns in Methicillin-Resistant Staphylococcus aureus from Patients Hospitalized during 2015-2017 in Hospitals in Poland». *Medical Principles and Practice*, 29(1), 61–68. <https://doi.org/10.1159/000501788>.

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