|                                                                                    |                                                                         |                 |
|------------------------------------------------------------------------------------|-------------------------------------------------------------------------|-----------------|
|   | <a href="https://www.thermofisher.com">https://www.thermofisher.com</a> |                 |
|  | EU                                                                      | +800 135 79 135 |
|                                                                                    | US                                                                      | 1 855 236 0910  |
|                                                                                    | CA                                                                      | 1 855 805 8539  |
|                                                                                    | ROW                                                                     | +31 20 794 7071 |

## Oxoid™ Tobramycin Discs 10 µg (TOB10)

**REF CT0056B**

### Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Tobramycin Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs that contain specific amounts of the antimicrobial agent tobramycin. The discs are labelled on both sides with details of the antimicrobial (TOB) and amount present (µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Tobramycin AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

#### Intended Use

Tobramycin Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which tobramycin has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling, and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,<sup>3</sup> CLSI M100<sup>1ac</sup> or EUCAST<sup>2</sup> breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

#### Gram-Positive

- *Staphylococcus aureus*

#### Gram-Negative

- Enterobacteriaceae
- *Pseudomonas aeruginosa*

Species with published breakpoints according to current CLSI M100 literature (all associated with 10 µg):

#### Gram-Positive

- N/A

#### Gram-Negative

- Enterobacterales
- *Pseudomonas aeruginosa*
- *Acinetobacter* spp.

Species with published breakpoints according to current EUCAST literature (all associated with 10 µg):

#### Gram-Positive

- *Staphylococcus* spp. (*Staphylococcus aureus* and CoNS)

#### Gram-Negative

- Enterobacterales (systemic infections and infections originating from the urinary tract)
- *Pseudomonas* spp. (systemic infections and infections originating from the urinary tract)
- *Acinetobacter* spp. (systemic infections and infections originating from the urinary tract)

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

#### Principle of Method

Tobramycin AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI<sup>1bc</sup>/EUCAST<sup>2ab</sup> methodology refer to the

relevant current standards. Tables showing CLSI<sup>1bc</sup>/EUCAST<sup>2</sup> compound/concentrations can be found in the documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

#### Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Tobramycin AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI<sup>1bc</sup> and/or EUCAST<sup>2</sup>.

#### Materials Provided

Tobramycin AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Tobramycin AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

| Table 1. Materials Provided with CT0056B     |                                                                                     | Table 2. Description of Tobramycin AST Disc Reagents |                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|----------------------------------------------|-------------------------------------------------------------------------------------|------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Component Description                        | Material Description                                                                | Reagent                                              | Description of Function                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Cartridge with spring, cap, and plunger (x5) | Assembly components and plastic cartridge containing 50x AST discs.                 | Tobramycin                                           | White to almost white powder. Tobramycin is a 4,6-disubstituted 2-deoxystreptamine (DOS) ring-containing aminoglycoside antibiotic with activity against various Gram-negative and some Gram-positive bacteria. In general, like other aminoglycosides, tobramycin is bactericidal and exhibits both immediate and delayed killing, which are attributed to different mechanisms, such as binding to bacterial cell membranes and inhibition of protein synthesis. |
| Desiccant tablet (x5)                        | Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge. |                                                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Foil                                         | Foil individually sealing each cartridge with its desiccant.                        |                                                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Susceptibility Test Discs (x250)             | Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack. |                                                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA<sup>3</sup>). The actual concentration is detailed on the Certificate of Analysis.

Oxoid Tobramycin AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

#### Shelf-life and Storage Conditions

Unopened cartridges of Tobramycin AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

## Analytical Performance Characteristics

**Table 1. Raw data of zone sizes, in accordance with CLSI methodology, taken from the 3 most recent batches of Tobramycin (10 µg) TOB10 Antimicrobial Susceptibility Test Discs (CT0056B).**

| CLSI Methodology |                               |              |            |             |             |                                      |              |    |    |                                                                                                                    |
|------------------|-------------------------------|--------------|------------|-------------|-------------|--------------------------------------|--------------|----|----|--------------------------------------------------------------------------------------------------------------------|
| Product Batch    | Organism                      | ATCC® number | Test Media | Lower Limit | Upper Limit | Calculated reference mid-point value | Reading (mm) |    |    | Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) |
|                  |                               |              |            |             |             |                                      | 1            | 2  | 3  |                                                                                                                    |
| 3304672          | <i>Escherichia coli</i>       | ATCC® 25922  | MHA        | 18          | 26          | 22                                   | 19           | 19 | 19 | -3 ± 0 (CV=0)                                                                                                      |
|                  | <i>Staphylococcus aureus</i>  | ATCC® 25923  | MHA        | 19          | 29          | 24                                   | 25           | 25 | 25 | 1 ± 0 (CV=0)                                                                                                       |
|                  | <i>Pseudomonas aeruginosa</i> | ATCC® 27853  | MHA        | 20          | 26          | 23                                   | 22           | 22 | 22 | -1 ± 0 (CV=0)                                                                                                      |
| Product Batch    | Organism                      | ATCC® number | Test Media | Lower Limit | Upper Limit | Calculated reference mid-point value | Reading (mm) |    |    |                                                                                                                    |
|                  |                               |              |            |             |             |                                      | 1            | 2  | 3  |                                                                                                                    |
| 3261712          | <i>Escherichia coli</i>       | ATCC® 25922  | MHA        | 18          | 26          | 22                                   | 23           | 23 | 23 | 1 ± 0 (CV=0)                                                                                                       |
|                  | <i>Staphylococcus aureus</i>  | ATCC® 25923  | MHA        | 19          | 29          | 24                                   | 24           | 24 | 24 | 0 ± 0 (CV=0)                                                                                                       |
|                  | <i>Pseudomonas aeruginosa</i> | ATCC® 27853  | MHA        | 20          | 26          | 23                                   | 22           | 22 | 22 | -1 ± 0 (CV=0)                                                                                                      |
| Product Batch    | Organism                      | ATCC® number | Test Media | Lower Limit | Upper Limit | Calculated reference mid-point value | Reading (mm) |    |    |                                                                                                                    |
|                  |                               |              |            |             |             |                                      | 1            | 2  | 3  |                                                                                                                    |
| 3235291          | <i>Escherichia coli</i>       | ATCC® 25922  | MHA        | 18          | 26          | 22                                   | 22           | 22 | 22 | 0 ± 0 (CV=0)                                                                                                       |
|                  | <i>Staphylococcus aureus</i>  | ATCC® 25923  | MHA        | 19          | 29          | 24                                   | 25           | 25 | 25 | 1 ± 0 (CV=0)                                                                                                       |
|                  | <i>Pseudomonas aeruginosa</i> | ATCC® 27853  | MHA        | 20          | 26          | 23                                   | 22           | 22 | 22 | -1 ± 0 (CV=0)                                                                                                      |

**Table 4. Raw data of zone sizes, in accordance with EUCAST methodology, taken from the 3 most recent batches of Tobramycin (10 µg) TOB10 Antimicrobial Susceptibility Test Discs (CT0056B).**

| EUCAST Methodology |                               |              |            |             |             |                                      |              |    |    |                                                                                                                    |
|--------------------|-------------------------------|--------------|------------|-------------|-------------|--------------------------------------|--------------|----|----|--------------------------------------------------------------------------------------------------------------------|
| Product Batch      | Organism                      | ATCC® number | Test Media | Lower Limit | Upper Limit | Calculated reference mid-point value | Reading (mm) |    |    | Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) |
|                    |                               |              |            |             |             |                                      | 1            | 2  | 3  |                                                                                                                    |
| 3304672            | <i>Escherichia coli</i>       | ATCC® 25922  | MHA        | 18          | 26          | 22                                   | 19           | 19 | 19 | -3 ± 0 (CV=0)                                                                                                      |
|                    | <i>Staphylococcus aureus</i>  | ATCC® 29213  | MHA        | 20          | 26          | 23                                   | 21           | 21 | 21 | -2 ± 0 (CV=0)                                                                                                      |
|                    | <i>Pseudomonas aeruginosa</i> | ATCC® 27853  | MHA        | 20          | 26          | 23                                   | 22           | 22 | 22 | -1 ± 0 (CV=0)                                                                                                      |
| Product Batch      | Organism                      | ATCC® number | Test Media | Lower Limit | Upper Limit | Calculated reference mid-point value | Reading (mm) |    |    | Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) |
|                    |                               |              |            |             |             |                                      | 1            | 2  | 3  |                                                                                                                    |
| 3261712            | <i>Escherichia coli</i>       | ATCC® 25922  | MHA        | 18          | 26          | 22                                   | 23           | 23 | 23 | 1 ± 0 (CV=0)                                                                                                       |
|                    | <i>Staphylococcus aureus</i>  | ATCC® 29213  | MHA        | 20          | 26          | 23                                   | 22           | 22 | 22 | -1 ± 0 (CV=0)                                                                                                      |
|                    | <i>Pseudomonas aeruginosa</i> | ATCC® 27853  | MHA        | 20          | 26          | 23                                   | 22           | 22 | 22 | -1 ± 0 (CV=0)                                                                                                      |
| Product Batch      | Organism                      | ATCC® number | Test Media | Lower Limit | Upper Limit | Calculated reference mid-point value | Reading (mm) |    |    | Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) |
|                    |                               |              |            |             |             |                                      | 1            | 2  | 3  |                                                                                                                    |
| 3235291            | <i>Escherichia coli</i>       | ATCC® 25922  | MHA        | 18          | 26          | 22                                   | 22           | 22 | 22 | 0 ± 0 (CV=0)                                                                                                       |
|                    | <i>Staphylococcus aureus</i>  | ATCC® 25923  | MHA        | 20          | 26          | 23                                   | 22           | 22 | 22 | -1 ± 0 (CV=0)                                                                                                      |
|                    | <i>Pseudomonas aeruginosa</i> | ATCC® 27853  | MHA        | 20          | 26          | 23                                   | 22           | 22 | 22 | -1 ± 0 (CV=0)                                                                                                      |

### Clinical Performance Characteristics

A study by Abduljabbar et al. (2020), aimed to detect the prevalence of multi-drug resistant and extended-spectrum β-lactamase-producing (ESBL) bacteria isolated from inpatients with burn infections and the antimicrobial sensitivity patterns of all bacterial isolates during three years.<sup>4</sup> Of the 295 swabs taken, 513 different bacterial strains were isolated and identified, 335 (65.3%) were Gram-negative and 178 (34.7%) were Gram-positive. Each isolate was tested for sensitivity against a range of antibiotics discs including Oxoid Tobramycin (10 µg) AST Discs. Results showed overall tobramycin resistance for all strains was 68.60%. Resistance to tobramycin was 78.8% in *Pseudomonas aeruginosa*, 70.7% in MRSA strains, 79.5% in *K. pneumoniae* strains,

55.5% in MSSA strains, 49% in *E. coli* strains, 59.3% of *A. baumannii* strains and 41.1% of *S. typhimurium* isolates. The study reported that the majority of isolates studied were multidrug resistant. Furthermore, it was reported that the elevated levels of resistance observed in the study could be due to biofilm formation by the pathogenic organisms, the overuse and unsuitability of different antimicrobials to treat burn infections has led to the emergence of resistance. In conclusion this article demonstrates that Oxoid Tobramycin (10 µg) AST Discs are effective in the clinical analysis of Gram-negative and Gram-positive clinical isolates given that the discs have displayed adequate susceptibility to some of the organisms tested when compared with other antimicrobials and they have detected emerging resistance in the organisms tested.

A 2019 study conducted by Ahman et al. aimed to identify all manufacturers of antimicrobial discs and evaluate the performance of a number of relevant antimicrobial discs.<sup>5</sup> To represent different antibiotic classes, 16 antimicrobial discs were selected including Oxoid Tobramycin TOB10 AST Discs. Discs were obtained on two occasions (in 2014 and 2017). Discs were tested in triplicate on Mueller Hinton Agar (MHA) against EUCAST quality control strains. Mean values were compared with targets and ranges in the EUCAST quality control tables. Oxoid Tobramycin discs demonstrated excellent and consistent quality in both years having a mean value within  $\pm 1$  mm of the target value, exceeding the performance of other discs.

### Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

### References

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
  - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests.
  - b. M100, Performance Standards for Antimicrobial Susceptibility Testing.
  - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
  - a. EUCAST Disk Diffusion Method.
  - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Abduljabbar, Ahmed, Jaloob Aljanaby, Israa Abduljabbar, 2020. "Prevalence of Aerobic Pathogenic Bacteria Isolated from Patients with Burn Infection and Their Antimicrobial Susceptibility Patterns in Al-Najaf City, Iraq- a Three-Year Cross-Sectional Study. [ Version 1; Peer Review: 2 Approved with Reservations]," 1–14.
5. Ahman, J., E. Matuschek, and G. Kahlmeter. 2019. "The Quality of Antimicrobial Discs from Nine Manufacturers—EUCAST Evaluations in 2014 and 2017." *Clinical Microbiology and Infection* 25 (3): 346–52.  
<https://doi.org/10.1016/j.cmi.2018.05.021>.

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

### Glossary of symbols

| Symbol/Label                                                                        | Meaning                                                                 |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------|
|  | Manufacturer                                                            |
|  | In Vitro Diagnostic Medical Device                                      |
|  | Temperature limit                                                       |
|  | Batch Code                                                              |
|  | Catalog Number                                                          |
|  | Do not re-use                                                           |
|  | Consult instructions for use or consult electronic instructions for use |

|                                                                                   |                                                                                             |
|-----------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|
|  | Contains sufficient for <n> tests                                                           |
|  | Use-by date                                                                                 |
|  | Do not use if package is damaged and<br>Consult instructions for use                        |
|  | Authorized representative in the<br>European Community/<br>European Union                   |
|  | Unique device identifier                                                                    |
|  | USA: Caution: Federal law restricts this<br>device to sale by or on order of a<br>Physician |
|  | European Conformity Mark                                                                    |
|  | UK Conformity Mark                                                                          |

#### Revision Information

| Version | Date of issue and modifications introduced |
|---------|--------------------------------------------|
| 1.0     | 2023-01-18                                 |

|                                                                                    |                                                                         |                 |
|------------------------------------------------------------------------------------|-------------------------------------------------------------------------|-----------------|
|   | <a href="https://www.thermofisher.com">https://www.thermofisher.com</a> |                 |
|  | EC                                                                      | +800 135 79 135 |
|                                                                                    | US                                                                      | 1 855 236 0910  |
|                                                                                    | CA                                                                      | 1 855 805 8539  |
|                                                                                    | ROW                                                                     | +31 20 794 7071 |

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### Références

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
  - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests.
  - b. M100, Performance Standards for Antimicrobial Susceptibility Testing.
  - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
  - a. EUCAST Disk Diffusion Method.
  - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Abduljabbar, Ahmed, Jaloob Aljanaby, Israa Abduljabbar, 2020. "Prevalence of Aerobic Pathogenic Bacteria Isolated from Patients with Burn Infection and Their Antimicrobial Susceptibility Patterns in Al-Najaf City, Iraq- a Three-Year Cross-Sectional Study. [ Version 1; Peer Review: 2 Approved with Reservations]," 1–14.
5. Åhman, J., E. Matuschek, and G. Kahlmeter. 2019. "The Quality of Antimicrobial Discs from Nine Manufacturers—EUCAST Evaluations in 2014 and 2017." *Clinical Microbiology and Infection* 25 (3): 346–52.  
<https://doi.org/10.1016/j.cmi.2018.05.021>.

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### Bibliografia

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
  - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests.
  - M100, Performance Standards for Antimicrobial Susceptibility Testing.
  - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
  - EUCAST Disk Diffusion Method.
  - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Abduljabbar, Ahmed, Jaloob Aljanaby, Israa Abduljabbar, 2020. "Prevalence of Aerobic Pathogenic Bacteria Isolated from Patients with Burn Infection and Their Antimicrobial Susceptibility Patterns in Al-Najaf City, Iraq- a Three-Year Cross-Sectional Study. [ Version 1; Peer Review: 2 Approved with Reservations]," 1–14.
- Ahman, J., E. Matuschek, and G. Kahlmeter. 2019. "The Quality of Antimicrobial Discs from Nine Manufacturers—EUCAST Evaluations in 2014 and 2017." *Clinical Microbiology and Infection* 25 (3): 346–52.  
<https://doi.org/10.1016/j.cmi.2018.05.021>.

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