Tryptone Soya Agar

EN

REF

CM0131B, CM0131K, CM0131R, CM0131T, CM0131Q, CM0131V, CM0131W, CM0131M

Intended Use

Tryptone Soya Agar (CM0131B, CM0131K, CM0131Q, CM0131R, CM0131T, CM0131V, CM0131W, CM0131M) is a highly nutritious general-purpose medium for the growth of fastidious and non-fastidious microorganisms from clinical samples (e.g. urine and wounds). The device can also be used for pharmaceutical, environmental, food and water samples. Tryptone Soya Agar (CM0131B, CM0131K, CM0131Q, CM0131R, CM0131T, CM0131V, CM0131W, CM0131M) is recommended as a reference medium when testing selective media to measure the degree of inhibition.

Tryptone Soya Agar (CM0131B, CM0131K, CM0131Q, CM0131R, CM0131T, CM0131V, CM0131W, CM0131M) are used in a diagnostic workflow to support the clinicians for the growth of a wide variety of fastidious and non-fastidious microorganisms from clinical samples of patients suspected of having microbial infections.

The devices are for professional use only, are not automated and nor are they companion diagnostics.

Summary and Explanation

Tryptone Soya Agar is a general-purpose, non-selective medium which will support the growth of a wide variety of organisms. It can be used for many purposes including the maintenance of stock cultures, as a reference medium when testing selective media to determine the degree of inhibition, isolation of pure cultures, or for general culture of microorganisms from a range of sample types.

In a clinical setting, Tryptone Soya Agar is used as part of a larger diagnostic workflow (often as the initial culture step) to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

CM0131B/K/R/T/Q/V/W/M is supplied as a powder which after preparation as listed in this IFU, will result in a solid growth medium.

Principle of Method

Casein and soy peptones provide nitrogen, amino acids, and peptides necessary for bacterial growth. Dextrose is a ready source of energy. Sodium chloride is a source of essential electrolytes and maintains osmotic equilibrium. Tryptone Soya Agar (TSA) may be supplemented with blood to facilitate the growth of more fastidious bacteria or antimicrobial agents to permit the selection of various microbial groups from mixed flora.

Typical Formula

	<u>grams per litre</u>
Pancreatic digest of casein	15.0
Enzymatic digest of soya bean	5.0
Sodium Chloride	5.0
Agar	15.0

Materials Provided

CM0131B: 500g of Tryptone Soya Agar base.
 CM0131K: 25kg of Tryptone Soya Agar base.
 CM0131M: 3 x 2.4kg of Tryptone Soya Agar base.
 CM0131Q: 6 x 4kg of Tryptone Soya Agar base.
 CM0131R: 2.5kg of Tryptone Soya Agar base.
 CM0131T: 5kg of Tryptone Soya Agar base.
 CM0131V: 13 x 2kg of Tryptone Soya Agar base.
 CM0131W: 3 x 4kg of Tryptone Soya Agar base.
 500g of Tryptone Soya Agar yields approximately 12.5L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dishes

Storage

- Store product in its original packaging between 10°C and 30°C.
 - Keep container tightly closed.
 - The product may be used until the expiry date stated on the label.
 - Protect from moisture.
 - Store away from light.
 - Allow reconstituted product to equilibrate to room temperature before use.
- Once reconstituted, store media between 2°C and 8°C.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.

- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established. Specimen Collection, Handling and Storage Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 01, ID 04, ID 07, ID 17, ID 12, B 1, B 2, B 5, B 11, B 57 and Q5.

Procedure

Suspend 40g in 1 litre of distilled water. Bring to the boil to dissolve completely. Sterilise by autoclaving at 121°C for 15 minutes. Allow to cool, mix well and pour into sterile containers. Once set, inoculate and streak the specimen onto the medium using a standard loop. Incubate plates for 18–24 hours at 30 - 35°C. Visually inspect plates to assess colony growth and colour under good lighting.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Positive Controls Incubation Conditions: 18-24 hours, 30 - 35°C Inoculum level: 10 – 100 cfu Colony count is $\geq 70\%$ of the control medium count	
<i>Staphylococcus aureus</i> ATCC® 9144™	0.5 – 1 mm, straw colonies
<i>Staphylococcus pyogenes</i> ATCC® 19615™	0.25 – 0.5 mm, pale straw colonies
<i>Streptococcus viridans</i> NCTC 1080	0.25-0.5 mm pale straw colonies
<i>Staphylococcus epidermidis</i> ATCC® 12228™	1 – 2 mm white/grey colonies
Enriched with 7% v/v horse blood Reactions after incubation at 37°C for 24 hours Medium is challenged with 10-100 cfu Colony count is $\geq 70\%$ of the control medium count	
<i>Staphylococcus pyogenes</i> ATCC® 19615™	0.25 – 0.5 mm, pale straw colonies, β haemolysis
<i>Streptococcus viridans</i> NCTC 1080	0.5-1 mm grey/green colonies, α haemolysis
<i>Streptococcus pneumoniae</i> ATCC® 6305™	0.5-1 mm grey/green colonies, α haemolysis
Reactions after incubation at 37°C for 48 hours under microaerophilic conditions Inoculum level: 10-100 cfu Colony count is $\geq 70\%$ of the control medium count	
<i>Haemophilus influenzae</i> ATCC® 19418™	Pinpoint-0.5mm colourless colonies

Zones of growth/no growth surrounding X, V and X+V factor discs (DD0003, DD0004 and DD0005) when plain plates are inoculated with the following organisms and incubated at 37°C for 18 hours:

	X	V	X+V
<i>Haemophilus parainfluenzae</i> ATCC® 33392™	0	≥20mm	≥20 mm
<i>Haemophilus influenzae</i> ATCC® 9334™	0	0	≥15 mm
<i>Haemophilus influenzae</i> ATCC® 19418™	0	0	≥15 mm
<i>Haemophilus influenzae</i> ATCC® 49247™	0	0	≥15 mm

Zones of inhibition with Bacitracin discs (DD0002) shall be 10-20mm when 7% v/v horse blood plates are inoculated with *Streptococcus pyogenes* ATCC® 19615 and incubated at 37°C for 18 hours.

Testing in accordance with ISO11133:2014

Inoculum level: 50 – 120 cfu

Colony count is ≥ 70% of the control medium count

Positive Controls	
Reactions after incubation at 30 ± 2°C for 24 ± 2 hours	
<i>Bacillus cereus</i> ATCC® 11778™	3 – 5 mm irregular, straw colonies
<i>Bacillus subtilis</i> ATCC® 6633™	2 – 4 mm irregular, straw colonies
<i>Escherichia coli</i> ATCC® 8739™	1-3 mm cream colonies
Reactions after incubation at 36 ± 2°C for 20 ± 2 hours	
<i>Escherichia coli</i> ATCC® 25922™	1-3 mm cream colonies
<i>Escherichia coli</i> ATCC® 11775™	1-3 mm cream colonies
<i>Escherichia coli</i> NCTC 13167	1-3 mm cream colonies
<i>Pseudomonas aeruginosa</i> ATCC® 10145™	1-4 mm straw colonies
<i>Enterococcus faecalis</i> ATCC® 29212™	0.5-2 mm straw colonies
Reactions after incubation at 37 ± 2°C for 24 ± 2 hours	
<i>Staphylococcus aureus</i> ATCC® 25923™	0.5-1 mm straw colonies
<i>Listeria monocytogenes</i> ATCC® 13932™	0.25 – 2 mm straw colonies
Reactions after incubation at 44 ± 2°C for 21 ± 3 hours	
<i>Escherichia coli</i> ATCC® 8739™	1-3 mm cream colonies
Reactions after anaerobic incubation at 44 ± 2°C for 21 ± 3 hours	
<i>Clostridium perfringens</i> ATCC® 13124™	1-2 mm straw colonies

Testing performed in accordance with current CLSI M22 A Enriched with 5% Sheep Blood

Positive Controls	
Reactions after incubation at 35 ± 2°C for 21 ± 3 hours Medium is challenged with 10 ² -10 ⁴ cfu	
<i>Streptococcus pyogenes</i> ATCC® 19615™	0.5-1 mm pale straw colonies, β haemolysis
<i>Streptococcus pneumoniae</i> ATCC® 6305™	0.5-2 mm grey/green colonies, α haemolysis
<i>Staphylococcus aureus</i> ATCC® 25923™	1-2 mm white/grey colonies
<i>Escherichia coli</i> ATCC® 25922™	1-2 mm straw colonies

Testing performed in accordance with current USP/EP/BP/JP

Positive Controls Reactions after incubation at 30 -35°C for 24 hours Medium is challenged with 10-100 cfu Colony count is ≥ 70% of the control medium count	
<i>Staphylococcus aureus</i> ATCC® 6538™	0.5-1 mm straw colonies
<i>Escherichia coli</i> ATCC® 8739™	1-3 mm cream colonies
<i>Bacillus subtilis</i> ATCC® 6633™	2-4 mm irregular, straw colonies
<i>Pseudomonas aeruginosa</i> ATCC® 9027™	1-4 mm straw colonies
<i>Salmonella</i> Typhimurium ATCC® 14028™	1-3 mm straw colonies
<i>Salmonella abony</i> NCTC 6017	1-3 mm straw colonies
Reactions after incubation at 30-35°C for 5 days Medium is challenged with 10-100 colony-forming units	
<i>Candida albicans</i> ATCC® 10231™	1-3 mm straw colonies
<i>Aspergillus brasiliensis</i> ATCC® 16404™	>10 mm colonies, white mycelia, with/without black spores

Limitations

The use of chocolate agar is preferable for species that require X and V factor for growth rather than blood agar or blood containing medium, because of risk of carryover of X factor. The X and V factor testing could also be done on a basic nutrient agar such as TSA, but the X and V factor discs need to be validated beforehand, in order to assess the presence of trace factors that could influence the results, usually during the identification of *Haemophilus influenzae* and *Haemophilus parainfluenzae*. Manufacturers' instructions should be followed when performing this test¹.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of fastidious and non-fastidious organisms is confirmed by the inclusion of well- characterised isolates in the QC processes performed as part of the manufacture of each batch of the devices. The precision of Tryptone Soya Agar (CM0131) was demonstrated by an overall pass rate of 100% obtained for the product over one month and a half of testing (12.MAY.2022 – 06.JUL.2022; 10 batches. This shows that the performance is reproducible.


Tryptone Soya Agar (CM0131) is tested in- house as part of the QC process since at least 1995. For target organisms, when using 10-100 cfu inoculum of *Streptococcus pyogenes* ATCC® 19615™, *Streptococcus viridans* NCTC 1080, *Staphylococcus aureus* ATCC® 9144™ and *Staphylococcus epidermidis* ATCC® 12228™ by incubating the medium at 30 - 35°C for 18 - 24 hours or at 37°C for 18 - 24 hours, the user can recover organisms equal or greater than 70% and colony size and morphology as listed in this document. Same percentage recovery and colony morphology conform to specifications is expected after enrichment of the medium with 7% v/v horse blood and incubation at 37°C for 24 hours of 10-100 cfu inoculum of *Streptococcus pyogenes* ATCC® 19615™, *Streptococcus viridans* NCTC 1080, *Streptococcus pneumoniae* ATCC® 6305™ and after incubation of 10 – 100 cfu of *Haemophilus influenzae* ATCC® 19418™ after enrichment of the medium with 7% v/v horse blood and under microaerophilic conditions at 37°C for 48 hours. The medium is also tested after enrichment with 5% sheep blood according to CLSI M22 A guidelines, and it is inoculated with 100 – 10000 cfu at 35 ± 2°C for 21 ± 3 hours for *Streptococcus pyogenes* ATCC® 19615™, *Streptococcus pneumoniae* ATCC® 6305™, *Staphylococcus aureus* ATCC® 25923™ and *Escherichia coli* ATCC® 25922™. Ultimately, the zones of growth/no growth for X, V, X+V Factor Diagnostic Discs are measured for *Haemophilus influenzae* ATCC® 9334™, *Haemophilus influenzae* ATCC® 19418™, *Haemophilus influenzae* ATCC® 49247™ and *Haemophilus parainfluenzae* ATCC® 33392™ after incubation at 37°C for 18 hours, while the zones of inhibition with Bacitracin Discs are measured after 7% v/v horse blood enrichment and inoculation with *Streptococcus pyogenes* ATCC® 19615™ incubated at 37°C for 18 hours.

Bibliography

1. Public Health England. (2021b). UK SMI ID 12: identification of Haemophilus species and the HACEK group of organisms - GOV.UK. <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>

Symbol Legend

Symbol	Definition
	Catalogue number

	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Consult instructions for use or consult electronic instructions for use
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2023 Thermo Fisher Scientific Inc. All rights reserved.
 ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
 All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke,
 Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of modifications introduced
3.0	2023-11-21

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©2023 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.