Tryptone Soya Broth

EN

[REF] CM0129B, CM0129K, CM0129R, CM0129T, CM0129V

Intended Use

Tryptone Soya Broth (CM0129B/K/R/T/V) devices are highly nutritious, general-purpose media for the growth of fastidious and non-fastidious microorganisms from a range of clinical samples (wounds, skin, respiratory secretions, genitourinal secretions etc.).

The devices are intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having microbial infections.

The devices are for professional use only, are not automated and nor are they companion diagnostics.

The device can also be used for testing non-clinical samples (food and pharma workflows).

Summary and Explanation

A fastidious organism is any organism that has complex or particular nutritional requirements. There are many different types of fastidious organisms which exhibit different requirements to ensure successful growth. For example: *Streptococcus pyogenes*, *Streptococcus pneumoniae*, and *Bacteroides fragilis*.

Non-fastidious microorganisms are organisms that grow and replicate without special nutritional supplements or conditions. They can synthesize all organic molecules required for their growth from simple sugars, amino acids, etc. Non-fastidious microorganisms can grow in nutrient sparse environments and can survive in a wider temperature range. This group of bacteria are involved in causing illnesses and infections in humans and includes *Escherichia coli*, *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Bacillus subtilis* and *Salmonella Typhimurium*^{1,2,3}.

Principle of Method

Growth of fastidious and non-fastidious microorganisms is achieved through the inclusion of a carefully formulated blend of peptones which provide necessary amino acids and nutrients to ensure growth of microorganisms. Glucose is the fermentable carbohydrate supplying energy for growth, whilst salt is required to maintain an osmotic equilibrium within the medium.

Typical Formula

	<u>grams per litre</u>
Pancreatic digest of casein	17.0
Enzymatic* digest of soya bean	3.0
Sodium chloride	5.0
Di-potassium hydrogen phosphate	2.5
Glucose	2.5

* contains papain

Materials Provided

CM0129B: 500g Tryptone Soya Broth

CM0129K: 25kg Tryptone Soya Broth

CM0129R: 2.5kg Tryptone Soya Broth

CM0129T: 5kg Tryptone Soya Broth

CM0129V: 8 x 3kg Tryptone Soya Broth

500g of Tryptone Soya Broth yields approximately 16.7L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers.
- Incubators.
- Quality control organisms.
- Sterile containers.

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store the prepared medium at room temperature away from light.

Warnings and Precautions

- If on skin wash with plenty of soap and water for at least 15 minutes. Get medical attention immediately if symptoms occur.
- If in eyes, rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Get medical attention.
- If inhaled, remove to fresh air. Get medical attention immediately if symptoms occur.
- If ingested, clean mouth with water and drink plenty of water afterwards. Get medical attention if symptoms occur.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to the nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimens should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 4, ID 7, ID 9, ID 16, ID 17, ID 22, ID 24, ID 25 and Q5.

Procedure

Add 30g to 1 litre of water (purified as required), mix well and distribute into final containers. Sterilise by autoclaving at 121°C for 15 minutes.

After incubation, take one loopful of broth from the edge of the surface of the fluid and inoculate selective media (e.g., Tryptone Soya Agar or Sabouraud Dextrose Agar)

Visually inspect plates to assess colony growth and morphology under good lighting.

Aerobic Cultivation

Tryptone Soya Broth may be used for the cultivation of aerobes and facultative anaerobes, including some fungi. Cultures should be examined at frequent intervals, as maximum growth is reached earlier than with less nutritious media and the phase of decline consequently begins sooner.

Anaerobic Cultivation

The addition of a small amount of agar renders the broth suitable for the cultivation of obligate anaerobes, such as *Clostridium* spp. For this purpose, the broth (with added agar) should be used soon after sterilisation, or heated and cooled just before inoculation.

Blood Culture

The superior growth-promoting properties of Tryptone Soya Broth make it especially useful for the isolation of organisms from blood or other normally sterile body fluids. Anticoagulants such as 'liquid' (sodium polyanethyl sulphonate) or sodium citrate may be added to the broth prior to sterilisation. Five to 10ml of blood may be added to 50ml of medium.

Interpretation

Turbidity indicates bacterial or fungal growth.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

Incubation Conditions: 18-24 hours at 30-35°C

Positive Controls	
Inoculum level: 10-100 cfu	
<i>Streptococcus pyogenes</i> ATCC ® 19615™	Turbid growth

Tested in accordance with current CLSI M22 A:

Incubation Conditions: 18-24 hours at 35 ± 2°C

Positive Controls	
Inoculum level: 10-100 cfu	
<i>Escherichia coli</i> ATCC ® 25922™	Turbid growth

<i>Staphylococcus aureus</i> ATCC ® 25923™	Turbid growth
---	---------------

Incubation Conditions: 5 days at 35 ± 2°C

Positive Controls Inoculum level: 10 ⁴ -10 ⁶ cfu	
<i>Streptococcus pneumoniae</i> ATCC ® 6305™	Turbid growth

Incubation Conditions: 5 days at 35 ± 2°C

Positive Controls Inoculum level: 10 ⁴ -10 ⁶ cfu	
<i>Streptococcus pneumoniae</i> ATCC ® 6305™	Turbid growth

Incubation Conditions: 5 days at 35 ± 2°C under anaerobic conditions

Positive Controls Inoculum level: 10 ⁴ -10 ⁶ cfu	
<i>Bacteroides fragilis</i> ATCC ® 25285™	Turbid growth

Tested in accordance with current USP/EP/BP/JP

Incubation Conditions: 24 hours at 30-35°C

Positive Controls Inoculum level: 10-100 cfu	
<i>Escherichia coli</i> ATCC ® 8739™	Turbid growth
<i>Staphylococcus aureus</i> ATCC ® 6538™	Turbid growth
<i>Pseudomonas aeruginosa</i> ATCC ® 9027™	Turbid growth
<i>Salmonella</i> abony NCTC 6017	Turbid growth
<i>Salmonella</i> typhimurium ATCC ® 14028™	Turbid growth

Incubation Conditions: 3 days at 30-35°C

Positive Controls Inoculum level: 10-100 cfu	
<i>Bacillus subtilis</i> ATCC ® 6633™	Flocculent/ surface growth
<i>Kocuria rhizophila</i> ATCC ® 9341™	Turbid growth

Incubation Conditions: 48 hours at 20-25°C

Positive Controls Inoculum level: 10-100 cfu	
<i>Bacillus subtilis</i> ATCC ® 6633™	Flocculent/ surface growth
<i>Candida albicans</i> ATCC ® 10231™	Flocculent/ surface growth

Incubation Conditions: 5 days at 20-25°C

Positive Controls Inoculum level: 10-100 cfu	
<i>Aspergillus brasiliensis</i> ATCC ® 16404™	White mycelia, black spores / no spores.

Limitations

Tryptone Soya Broth (CM0129B/K/R/T/V) devices are non-selective media intended to be used for the growth of some fastidious and non-fastidious microorganisms.

Due to variation in nutritional requirements some strains of the target organisms may be encountered that grow poorly or fail to grow in this medium.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Growth of fastidious and non-fastidious microorganisms is confirmed by the inclusion of well-characterised isolates in the QC processes performed as part of the manufacture of each batch of the devices. The precision of Tryptone Soya Broth EP/USP/BP/JP (CM0129B, CM0129K, CM0129R, CM0129T and CM0129V) was demonstrated by an overall pass rate of 100% obtained for the product over 2 months of testing (May 2022 - July 2022; 10 batches). This shows that the performance is reproducible.

Tryptone Soya Broth EP/USP/BP/JP (CM0129B, CM0129K, CM0129R, CM0129T and CM0129V) devices are tested in-house as part of the QC process since the products were launched in 1997 (CM0129B, CM0129R and CM0129T) and 1998 (CM0129K).

For target organisms, when using 10-100 cfu inoculum of *Streptococcus pyogenes* (ATCC® 19615™) and incubating the device at 30-35°C for 18-24 hours, the user can recover organisms with growth as listed in this document. For target organisms, when using 10⁴-10⁶ cfu inoculum of *Escherichia coli* (ATCC® 25922™) or *Staphylococcus aureus* (ATCC® 25923™) and incubating the device at 33-37°C for 18-24 hours, the user can recover organisms with growth as listed in this document. For target organisms, when using 10⁴-10⁶ cfu inoculum of *Streptococcus pneumoniae* (ATCC® 6305™) and incubating the device at 33-37°C for 5 days, the user can recover organisms with growth as listed in this document. For target organisms, when using 10⁴-10⁶ cfu inoculum of *Bacteroides fragilis* (ATCC® 25285™) and incubating the device at 33-37°C for 5 days under anaerobic conditions the user can recover organisms with growth as listed in this document

Tryptone Soya Broth EP/USP/BP/JP (CM0129B, CM0129K, CM0129R, CM0129T and CM0129V) devices are also tested in accordance with current USP/EP/BP/JP. For target organisms, when using 10-100 cfu inoculum of *Escherichia coli* (ATCC® 8739™), *Staphylococcus aureus* (ATCC® 6538™), *Pseudomonas aeruginosa* (ATCC® 9027™), *Salmonella abony* (NCTC 6017) or *Salmonella typhimurium* (ATCC® 14028™) and incubating the device at 30-35°C for 24 hours, the user can recover organisms with growth as listed in this document. For target organisms, when using 10-100 cfu inoculum of *Bacillus subtilis* (ATCC® 6633™) or *Kocuria rhizophila* (ATCC® 9341™) and incubating the device at 30-35°C for 3 days, the user can recover organisms with growth as listed in this document. For target organisms, when using 10-100 cfu inoculum of *Bacillus subtilis* (ATCC® 6633™) or *Candida albicans* (ATCC® 10231™) and incubating at 20-25°C for 48 hours, the user can recover organisms with growth as listed in this document. For target organisms, when using 10-100 cfu inoculum of *Aspergillus brasiliensis* (ATCC® 16404™) and incubating the device at 20-25°C for 5 days, the user can recover organisms with growth as listed in this document.

Bibliography

1. World Health Organization (WHO). 2017. 'Microbial Fact Sheets: Guidelines for drinking-water quality'. <https://www.who.int/publications/i/item/9789240045064>.
2. Public Health England. 2015a. 'Identification of *Enterobacteriaceae*'. UK Standards for Microbiology Investigations ID 16 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf.
3. Public Health England. 2015b. 'Identification of Vero cytotoxin-producing *Escherichia coli* including *Escherichia coli* O157'. UK Standards for Microbiology Investigations ID 22 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/436870/ID_22i4.pdf.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Consult instructions for use or consult electronic instructions for use
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/European Union
	European Conformity Assessment

	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of modifications introduced
3.0	2024-12-02 Update to the intended purpose to include sample types

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Bibliographie

1. World Health Organization (WHO). 2017. 'Microbial Fact Sheets: Guidelines for drinking-water quality'. <https://www.who.int/publications/i/item/9789240045064>.
2. Public Health England. 2015a. 'Identification of Enterobacteriaceae'. UK Standards for Microbiology Investigations ID 16 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf.
3. Public Health England. 2015b. 'Identification of Vero cytotoxin-producing Escherichia coli including Escherichia coli O157'. UK Standards for Microbiology Investigations ID 22 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/436870/ID_22i4.pdf.

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© 2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.