ESBL+AmpC screen disc kit

Disc tests for screening of ESBLs and/or AmpC producing organisms.

DESCRIPTION

Extended-spectrum β -lactamases (ESBLs) are enzymes hydrolyzing most penicillins and cephalosporins, including oxyimino- β -lactam compounds but not cephamycins and carbapenems. Most ESBLs belong to the Ambler class A of β -lactamases and are inhibited by β -lactamase inhibitors: clavulanic acid, sulbactam and tazobactam. ESBL production has been observed mostly in Enterobacteriaceae, particularly *Escherichia coli* and *Klebsiella pneumoniae*, but all other clinically-relevant Enterobacteriaceae species are also common ESBL-producers. In many areas, ESBL detection and characterization is recommended or mandatory for infection control purpose.

AmpC β -lactamases are enzymes encoded on the chromosome of many Enterobacteriaceae and a few other organisms. AmpC enzymes belong to class C, are active on penicillins but even more active on cephalosporins and can significantly hydrolyze cephamycins such as cefoxitin and cefotetan and oxyiminocephalosporins such as ceftazidime, cefotaxime, and ceftriaxone. Inhibitors of class A enzymes such as clavulanic acid, sulbactam, and tazobactam have much less effect on AmpC, which are strongly inhibited by boronic acid and cloxacillin.

CONTENTS OF THE PACKAGES

4 x 50 discs cartridges, each packaged in a "blister" with a dryer.

METHOD PRINCIPLE

Organisms suspected to be producers of either ESBL or AmpC enzymes may be tested with cefotaxime for evaluating the synergistic effects when combined with the following inhibitors:

- **Clavulanic acid** inhibits ESBL β -lactamases;
- **Cloxacillin** inhibits AmpC β -lactamases.

Combination Disc Tests (CDT) are carried out. Discs containing cefotaxime alone and in combination with clavulanic acid, cloxacillin and both of these inhibitors are applied. The inhibition zone around the cefotaxime disc combined with inhibitors is compared with the zone around the disc without the inhibitors.

GATHERING AND KEEPING SAMPLES

The colonies that are to be subjected to the susceptibility test are taken up by culture media that have been previously swabbed with the sample under examination.

TEST PROCEDURE

1. Using a fresh, pure culture prepare a suspension of the test organism equivalent to 0.5 McFarland.
2. Using a sterile cotton swab, spread the adjusted suspension over the entire area of a Mueller Hinton II agar plate.
3. Apply the discs on the inoculated plate, ensuring sufficient space between individual discs to allow for proper measurement of inhibition zones.
4. Incubate at $36\pm 1^\circ\text{C}$ for 18-24 hours.

EVALUATING THE RESULTS

At the end of the incubation period, measure the inhibition halos and interpret as indicated in the table below.

Interpretative Table. CDT method for screening of ESBL and/or AmpC.

Antibiotic Disc	Difference in zone diameters when compared with			β -lactamase
	Cefotaxime (CTX)	Cefotaxime+ Clavulanic acid (CTL)	Cefotaxime+ Cloxacillin (CTC)	
Cefotaxime+Clavulanic acid (CTL) or	≥ 5 mm	-	-	ESBL
Cefotaxime+Clavulanic acid+Cloxacillin (CTLC)	-	< 5 mm	≥ 5 mm	
Cefotaxime+Cloxacillin (CTC) or	≥ 5 mm	-	-	AmpC
Cefotaxime+Clavulanic acid+Cloxacillin (CTLC)	-	≥ 5 mm	< 5 mm	
Cefotaxime+Clavulanic acid (CTL) and	< 5 mm	-	-	ESBL+AmpC
Cefotaxime+Clavulanic acid+Cloxacillin (CTLC)	-	≥ 5 mm (AmpC)	≥ 5 mm (ESBL)	

If all zones are within 2 mm of each other, it means that the organism being tested is neither ESBL or AmpC producer.

QUALITY CONTROL

Appropriate strains for quality control of ESBL and AmpC detection tests:

Microorganism		ESBL phenotype	AmpC phenotype
<i>Klebsiella pneumoniae</i>	ATCC® 700603	+	-
<i>Enterobacater cloacae</i>	ATCC® BAA-1143	-	+
<i>Escherichia coli</i>	ATCC® 25922	-	-

LIMITS

Diffusion susceptibility tests use an *in vitro* technique and cannot therefore reproduce the extremely complex *in vivo* conditions. Nevertheless, it is a useful and important tool that helps the clinician choose the correct therapy. Many variable factors influence the final result of the diffusion susceptibility test. The main ones are: the culture medium used, impregnation of the discs, inoculation of the medium, temperature, time and incubation atmosphere of the plates, pre-incubation and pre-diffusion conditions, depth of the medium, etc.

PRECAUTIONS

The disc cannot be classified as being hazardous according to current legislation but fall within the specific field of application where a safety data sheet must be supplied because they can cause phenomena of sensitization in sensitive subjects if they come into contact with the skin.

The discs are disposable products. They are only for diagnostic *in vitro* use and are intended for professional use. They must be used in the laboratory by properly trained operators using approved aseptic and safety methods for pathogenic agents.

STORAGE

Store the unopened blister at -20°C to +8°C till the expiry date. Allow unopened cartridge to come to room temperature before removing it from the blister for minimising condensation on the discs. Leftover discs from an opened cartridge should be stored at 2-8°C for no more than 7 days. Return unused discs to the refrigerator as soon as the application of the discs has been completed. Dispose of expire discs.

ELIMINATING USED MATERIAL

After use, the discs and the material that comes into contact with the sample must be decontaminated and disposed of in accordance with current laboratory techniques for the decontamination and disposal of potentially infected material.

REFERENCES

- Jacoby G.A. AmpC β -Lactamases. Clin Microbiol Rev. 2009; 22(1): 161–182.
- EUCAST guidelines for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance. Version 1.0, 2013.

PRESENTATION

DESCRIPTION		PACKAGING	REF	TESTS
Cefotaxime 30 µg	CTX	1 x 50 Discs	99008	50
Cefotaxime 30 µg + Clavulanic acid	CTL	1 x 50 Discs		
Cefotaxime 30 µg + Cloxacillin	CTC	1 x 50 Discs		
Cefotaxime 30 µg + Clavulanic acid + Cloxacillin	CTLC	1 x 50 Discs		

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In Vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult accompanying documents	 Do not reuse

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