Kligler Iron Agar

REF CM0033B

EN

Intended Use

Kligler Iron Agar (CM0033B) is a differential medium for the identification of Enterobacteriaceae from clinical and non-clinical samples based on dextrose and lactose fermentation as well as hydrogen sulphide production.

Kligler Iron Agar (CM0033B) is used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The devices are for professional use only, are not automated, nor are they companion diagnostics.

Summary and Explanation

Members of the Enterobacteriaceae include *Escherichia coli*, *Proteus mirabilis*, *Salmonella typhimurium* and *Shigella flexneri*. Enterobacteriaceae are Gram-negative, straight rods that do not sporulate. Except for *Tatumella*, *Shigella*, and *Klebsiella* species, which are non-motile, some genera are motile via peritrichous flagella. Enterococci can grow and survive in harsh conditions, they can be found in the gastrointestinal tract of humans and animals, plants, soil, water, environment and fermented products¹.

Principle of Method

Phenol red is incorporated for detecting carbohydrate fermentation which can be seen as a colour change from orange/red to yellow in the medium. In the case of oxidative decarboxylation of the peptone, the pH of the media will rise and a colour change from orange/red to deep red will be seen. Production of hydrogen sulphide is indicated by a black colour in the butt of the tube.

Typical Formula

	<u>grams per litre</u>
'Lab-Lemco' powder	3.0
Yeast extract	3.0
Peptone	20.0
Sodium chloride	5.0
Lactose	10.0
Glucose	1.0
Iron (III) citrate	0.3
Sodium thiosulphate	0.3
Phenol red	0.05
Agar	12.0

Materials Provided

CM0033B: 500g of Kligler Iron Agar

500g of dehydrated Kligler Iron Agar yields approximately 9.1L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dish

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

- If inhaled: Remove to fresh air. Get medical attention immediately if symptoms occur.
- If ingested: Clean mouth with water and drink afterwards plenty of water. Get medical attention if symptoms occur.
- If in contact with skin: Wash off immediately with plenty of water for at least 15 minutes. Get medical attention immediately if symptoms occur.
- If in contact with eyes: Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Get medical attention.
- For in vitro diagnostic use only.
- For professional use only.

- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 16, B 41, S 7 and Q 5.

Procedure

Suspend 55g in 1 litre of distilled water. Bring to the boil to dissolve completely. Mix well and distribute into final containers. Sterilize by autoclaving at 121°C for 15 minutes. Allow to set as slopes with 2.5cm butts.

Interpretation

Once the medium is reconstituted:

The presence yellow butts and slopes indicates *Escherichia coli*. Black butts and red slopes indicates *Salmonella typhimurium* or *Proteus mirabilis*. Yellow butts and red slopes indicates *Shigella flexneri*.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18 hours @ 37°C

Inoculum level: 10⁴ – 10⁶ cfu

Control Strain	Butt	Slope	H ₂ S	Gas formation
<i>Escherichia coli</i> ATCC® 25922™	AG	A	-	+
<i>Klebsiella pneumoniae</i> ATCC® 13883™	AG	A	-	+
<i>Citrobacter freundii</i> ATCC® 8090™	AG	A	+	+
<i>Salmonella</i> Typhimurium ATCC® 14028™	AG	NC/Alk	+	+
<i>Salmonella</i> Enteritidis ATCC® 13076™	AG	NC/Alk	+	+
<i>Proteus mirabilis</i> ATCC® 29906™	A	NC/Alk	+	-
<i>Shigella sonnei</i> ATCC® 25931™	A	NC/Alk	-	-
<i>Shigella flexneri</i> ATCC® 12022™	A	NC/Alk	-	-
<i>Pseudomonas aeruginosa</i> ATCC® 27853™	Alk/NC	NC	-	-
<i>Alcaligenes faecalis</i> ATCC® 19018™	Alk	Alk	-	-

Key

AG = Acid (yellow) with gas formation

A = Acid (yellow)

NC = No change

Alk = Alkaline (red)

Hydrogen sulphide (H₂S)

Positive – Blackening

Negative = No blackening

Gas

Positive = Bubbles or splitting of agar

Negative = No bubbles or splitting of agar

Limitations

It is essential that Kligler Iron Agar is examined and reported at 18-24 hours. Early or late readings will give false results. Kligler Iron Agar will grow both oxidative and fermentative organisms. Confusion will result if care is not taken to distinguish between the two groups.

Always use a straight wire to inoculate the butt, to avoid splitting the agar with a loop.

Pure cultures are essential to avoid misinterpretation.

Do not use screw-capped tubes or bottles for Kligler Iron Agar medium. It is essential that air is freely available to growth on the slant.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of Enterobacteriaceae is confirmed by the inclusion of well-characterised isolates in the QC processes performed as part of the manufacture of each batch of the devices. The precision of Kligler Iron Agar (CM0033B) device was demonstrated by an overall pass rate of 100% obtained for the product over six years of testing (28.SEP.2016 – 25.FEB.2022; 10 batches). This shows that the performance is reproducible.


Kligler Iron Agar (CM0033B) device is tested in-house as part of the QC process since 1996. For target organisms, when using 10^4 – 10^6 cfu inoculum of *Escherichia coli* ATCC® 25922™, *Klebsiella pneumoniae* ATCC® 13883™, *Citrobacter freundii* ATCC® 8090™, *Salmonella* Typhimurium ATCC® 14028™, *Salmonella* Enteritidis ATCC® 13076™, *Proteus mirabilis* ATCC® 29906™, *Shigella sonnei* ATCC® 25931™, *Shigella flexneri* ATCC® 12022™, *Pseudomonas aeruginosa* ATCC® 27853™ or *Alcaligenes faecalis* ATCC® 19018™ and incubating the device at 37°C for 18 hours, the user can recover organisms with the typical reactions of *Enterobacteriaceae* species as listed in this document.

Bibliography

1. Public Health England. 2015. 'Identification of Enterobacteriaceae'. UK Standards for Microbiology Investigations. ID 16 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer

	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of modifications introduced
2.0	2023-12-13

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©2022 Thermo Fisher Scientific Inc. Toate drepturile rezervate.
Toate celelalte mărci comerciale sunt proprietatea Thermo Fisher Scientific Inc. și subsidiarelor acesteia.