C.L.E.D Medium

REF CM0301B, CM0301K, CM0301R, CM0301T & CM0301W

EN

Intended Use

C.L.E.D. Medium (CM0301B/K/R/W/T) devices are differential media for the isolation and enumeration of the common microorganisms causing urinary tract infections from urine samples.

The devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having urinary tract infections.

The devices are for professional use only, are not automated nor are they companion diagnostics.

Summary and Explanation

Urinary tract infections, often referred to as UTIs, are common infections that occur when bacteria invade the urinary tract by entering via the urethra. These infections can impact different parts of the urinary tract, but the most prevalent is cystitis (a bladder infection). The symptoms of cystitis can include frequent urination, or needing to urinate whilst having an empty bladder, a burning sensation or pain whilst urinating, blood in the urine and pressure or cramps in the lower abdominal area. The symptoms of kidney infections are very different. These can include chills, a fever, nausea and vomiting and pain in the lower back¹. There are many species of bacteria that cause UTIs, but most commonly are *Escherichia coli*, *Staphylococcus aureus* and *Proteus mirabilis*^{2,3}.

Due to the severity of UTIs in a vast number of patients, it is therefore highly important to be able to isolate and identify common UTI-causing microorganisms from urine samples and wound swabs. Early diagnosis is vital in the prevention and reduction of urinary tract infections.

Principle of Method

Peptones and Lab Lemco powder are present to supply the compounds and peptides required for the growth of bacteria and agar is the solidifying agent. Lactose provides a carbohydrate source. Bromothymol blue is a pH indicator which differentiates lactose fermenters (yellow) from non- fermenters. Cystine enhances the growth of cystine- dependent coliforms. Electrolytes are reduced in order to restrict the swarming of *Proteus* spp.

Typical Formula

	<u>grams per litre</u>
Peptone	4.0
'Lab-Lemco' Powder	3.0
Tryptone	4.0
Lactose	10.0
L-cysteine	0.128
Bromothymol blue	0.02
Agar	15.0

Materials Provided

CM0301B: 500g of dehydrated CLED Medium CM0301K: 25kg of dehydrated CLED Medium CM0301R: 2.5kg of dehydrated CLED Medium CM0301T: 5kg of dehydrated CLED Medium CM0301W: 2 x 7.24kg of dehydrated CLED Medium

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dishes

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 12°C.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.

- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local the recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B 59, B 60, ID 4 and Q 5.

Procedure

Suspend 36.2g in 1 litre of distilled water. Bring to the boil to dissolve completely. Sterilise by autoclaving at 121°C for 15 minutes. Cool to 50°C. Mix well and pour into sterile Petri dishes.

Interpretation

- Yellow colonies indicate *Staphylococcus aureus*
- Light yellow colonies indicate *Escherichia coli*
- Transparent blue colonies indicate *Proteus mirabilis*

Identification is presumptive and should be confirmed.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Positive Controls Inoculum level 10-100 cfu Colony counts should be ≥ 70% of the control medium Incubation Conditions: 18 hours @ 37°C	
<i>Enterococcus faecalis</i> ATCC® 19433™	0.25 - 1mm yellow colonies
<i>Proteus mirabilis</i> NCTC 10975™	0.25 – 1mm translucent blue colonies, no swarming
<i>Salmonella typhimurium</i> ATCC® 14028™	1 – 2mm blue colonies
Inoculum level 10 ⁴ – 10 ⁶ cfu Incubation conditions: 18 hours @ 37°C	
<i>Proteus mirabilis</i> NCTC 10975	Clear differentiation between lactose and non-lactose fermenters shall be seen.
<i>Escherichia coli</i> ATCC® 25922™	
Tested in accordance with current CLSI M22 A Inoculum level: 10-100 cfu Incubation conditions: 18-24 hours @ 35°C	
<i>Staphylococcus aureus</i> ATCC® 25923™	0.5-1mm yellow colonies
<i>Proteus vulgaris</i> ATCC® 8427™	0.25-1 mm translucent blue colonies, no swarming

<i>Escherichia coli</i> ATCC® 25922™	1-2 mm yellow colonies
Inoculum level 10 – 100 cfu Incubation Conditions: 48 hours @ 37°C	
<i>Lactobacillus fermentum</i> ATCC® 9338™	0.5 – 1 mm yellow colonies

Limitations

This medium is for laboratory use only and is not intended for the diagnosis of disease or other conditions. Identifications are presumptive and colonies should be identified using appropriate methods. A small number of atypical strains may give a weak reaction or fail to grow, especially when low numbers are present in the sample, for example due to a lack of electrolytes *Shigella* spp. do not usually grow.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of common UTI-causing microorganism strains is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the device. The precision of CLED medium was demonstrated by an overall pass rate of 100% obtained for the product over six months of testing (May – November 2021; 10 batches).

The devices are tested in-house as part of the QC process. For target organisms, when using 10-100 cfu inoculum of *Enterococcus faecalis*, *Proteus mirabilis* and *Salmonella typhimurium* and incubating the device at 37°C for 18 hours, the user can recover organisms with the appropriate morphology as listed in this document. For target organisms where the mixed medium is challenged with 10⁴-10⁶ cfu for


Escherichia coli and *Proteus mirabilis* and incubating the device at 37°C for 18 hours the user can recover organisms with colony size and morphology as listed in this document. For target organisms when using an inoculum level of 10- 100 cfu for *Escherichia coli*, *Staphylococcus aureus* and *Proteus vulgaris* and incubating the device at 35°C for 18-28 hours, in accordance with current CLSI M22 A, the user can recover organisms with the appropriate morphology listed in this document. For target organisms when using an inoculum level of 10-100 cfu for *Lactobacillus fermentum* ATCC® 9338™ and incubating the device at 37°C for 48 hours, the user can recover organisms with colony size and morphology that meets the defined acceptance criteria.

Bibliography

- 1) Centres for Disease Control and Prevention. 2021. 'Urinary Tract Infection'. <https://www.cdc.gov/antibiotic-use/uti.html>
- 2) Centres for Disease Control and Prevention. 2021. 'E.coli Prevention'. <https://www.cdc.gov/ecoli/ecoli-prevention.html>.
- 3) Centres for Disease Control and Prevention. 2021. 'Staphylococcus aureus in Healthcare Settings'. <https://www.cdc.gov/hai/organisms/staph.html>

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limitation (storage temp.)
	Use by (expiration date) YYYY- MM
	Keep away from sunlight
	Do not re-use
	Consult Instructions for Use
	Contains sufficient for <n> tests

	Do not use if packaging damaged
	Manufacturer
	Authorized representative in the European Community
	UK Conformity Assessment
	European Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures. CLSI is a trademark of the Clinical Laboratory and Standards Institute.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke,
Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of issue and modifications introduced
2.0	2023-11-16

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Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.