Schaedler Anaerobe Agar

EN

REF CM0437B & CM0437V

Intended Use

Schaedler Anaerobe Agar is a medium for the growth of aerobic and anaerobic microorganisms from faecal samples.

The devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The devices are for professional use only, are not automated, nor are they companion diagnostics.

Summary and Explanation

Anaerobic bacteria occur naturally and are the most common flora in the body. In their natural state they do not cause infection, however, when deep tissues become injured or exposed due to trauma or surgery, such as animal bites or root canal treatment, infections can happen. Anaerobic means "life without air" therefore oxygen is toxic to anaerobes, which can be explained by the absence of enzymes in the anaerobes of catalase, superoxide dismutase, and peroxidase enzymes. Anaerobic bacteria are important pathogens that can cause a variety of infections in humans. The site of anaerobic infection is commonly the site of normal colonization. Anaerobes are fastidious organisms and are difficult to grow if proper collection and culture methods are not used. The diagnosis requires clinical suspicion and proper microbiological identification. Significant disease-causing anaerobic microorganisms include, but are not limited to, *Bacteroides* species and *Clostridium* species¹.

Principle of Method

Schaedler Anaerobe Agar contains cysteine hydrochloride and glucose, as reducing substances, with the advantage that cysteine inhibits the growth of *Escherichia coli*. The inhibitory effect of cysteine on several enzymatic reactions of *Escherichia coli* *in vitro* has been reported².

Schaedler Anaerobe Agar is a suitable alternative to blood agar for the enumeration of clostridia. Strict anaerobic conditions for the successful recovery of obligate anaerobes when using this medium without the addition of blood has been reported.

Typical Formula

	<u>grams per litre</u>
Tryptone Soya Broth	10.0
Special peptone	5.0
Yeast extract	5.0
Glucose	5.0
Cysteine hydrochloride	0.4
Haemin	0.01
Tris buffer	0.75
Agar	13.5

Materials Provided

CM0437B: 500g of Schaedler Anaerobe Agar
CM0437V: 12 x 2.0kg of Schaedler Anaerobe Agar

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dish

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

- For *in vitro* diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use. Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.

- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any applicable federal, state or local regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material, following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between uses to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following recommended local guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 01 and S 7.

Procedure

Suspend 40g in 1 litre of distilled water. Bring to the boil to dissolve completely. Sterilise by autoclaving at 121°C for 15 minutes. Mix well and pour into sterile containers.

Interpretation

Once the medium is reconstituted:

1-2mm straw colonies will indicate *Clostridium tetani*, *Clostridium sporogenes* and *Bacteroides fragilis*. 1-3 mm straw colonies will indicate *Clostridium perfringens*. 0.5-2mm straw colonies will indicate *Clostridium sphenoides*, *Staphylococcus aureus* and *Enterococcus faecalis*.

Quality Control

It is the responsibility of the user to perform Quality Control testing, taking into account the intended use of the medium, and in accordance with any applicable local regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18 hours @ 37°C, anaerobic

Positive Controls Inoculum level: 10 – 100 cfu Colony count is ≥ 50% of the control medium count	
<i>Clostridium perfringens</i> NCTC 10613	1 – 3 mm straw colonies
<i>Clostridium sphenoides</i> ATCC® 19403™	0.5 - 2 mm straw colonies
<i>Clostridium tetani</i> ATCC® 9441™	1 – 2 mm straw colonies
<i>Clostridium sporogenes</i> ATCC® 19404™	1 – 2 mm straw colonies
<i>Bacteroides fragilis</i> ATCC® 23745™	1 – 2 mm straw colonies
<i>Staphylococcus aureus</i> ATCC® 25923™	0.5 – 2 mm straw colonies
<i>Enterococcus faecalis</i> ATCC® 19433™	0.5 – 2 mm straw colonies

Limitations

This medium is used to grow and isolate aerobic and anaerobic organisms from faecal samples (intestinal tract). Identification is presumptive and suspect colonies should be identified using appropriate methods. Not all strains of aerobic and anaerobic microorganisms will produce colonies of typical appearance.

Performance Characteristics


Accuracy has been demonstrated through review of the QC data. Correct detection of aerobic and anaerobic microorganism strains is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the device, which must meet the defined acceptance criteria. The precision of Schaedler Anaerobe Agar (applicable to CM0437B and V) was demonstrated by an overall pass rate of 98.77% obtained for the product over 26 years of testing (08.01.1996 – 15.07.2022, 731 batches). This shows that the performance is reproducible.

Schaedler Anaerobe Agar (applicable to CM0437B & V) has been tested in-house as part of the QC process since 1996. Target organisms, when using 10-100 cfu inoculum of *Clostridium perfringens* NCTC 10613, *Clostridium sphenoides* ATCC® 19403™, *Clostridium tetani* ATCC® 9441™, *Clostridium sporogenes* ATCC® 19404™, *Bacteroides fragilis* ATCC® 23745™, *Staphylococcus aureus* ATCC® 25923™ and *Enterococcus faecalis* ATCC® 19433™ and incubating anaerobically at 37°C for 18 hours, the user can recover organisms with colony morphology as listed in this document.

Bibliography

1. Noor, A. & S. Khetarpal. 2022. *Anaerobic Infections*. StatPearls Publishing, Florida, USA.
2. Kari C., Nagy Z., Kovacs P. and Hernadi F. (1971) J. Gen. Micro. 68. 349-356

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of modifications introduced
2.0	2023-12-19

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©2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.
Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Storbritannien