www.thermofisher.com

Sodium Thiosulphate	4.0
Ammonium iron (III) citrate	1.0
Brilliant Green	0.0125
Crystal Violet	0.009
Agar	15.0

MLCB Agar (Dehydrated)

REF CM0783B

EN

Intended Use

MLCB (mannitol lysine crystal violet brilliant green) Agar (CM0783B) is a selective medium for the isolation of *Salmonella* species (not *Salmonella* Typhi or *Salmonella* Paratyphi A) from faecal samples.

MLCB Agar (CM0783B) is intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having *Salmonella* infections.

The device is for professional use only, is not automated, nor is it a companion diagnostic.

Summary and Explanation

Salmonella species are rod shaped motile bacteria, with peritrichous flagella, 2 to 5µm in length and 0.7 to 1.5µm diameter, of the *Enterobacteriaceae* family¹. *Salmonella* species are broadly classified into two major groups that are *Salmonella enterica* and *Salmonella bongori* and the former subdivided into typhoidal and paratyphoid serovar. Almost 70% of the *Salmonella* infections in the UK are caused by the non-typhoidal serovar *Salmonella* Enteritidis¹. *Salmonella enterica* often causes gastroenteritis in humans and the mode of transmission can be food like undercooked meat, egg, dairy products, contaminated water, infected people and animals².

In 1967, a new *Salmonella* isolation medium called Brilliant Green Agar medium was described by Edel and Kampelmacher. The sensitivity of the medium in detected *Salmonella* species was evaluated using contaminated faecal samples from fattening calves. The medium offers greater inhibition of *Escherichia coli* and *Proteus* species than other formulations, due to the presence of Brilliant Green dye. Peptones within the medium provide nutritional sources for *Salmonella* species in addition to the presence of growth promoting B vitamins provided by yeast extract³. In 1968, the addition of Mannitol and Lysine by Inoue *et al.* allowed for greater distinction between strains by visual detection⁴.

Principle of Method

The medium contains a mixture of peptones, meat and yeast extracts to support good growth of salmonellae. The medium depends on the fermentation of mannitol and the decarboxylation of lysine together with an indicator system of thiosulphate and ferric ions to detect the production of hydrogen sulphide.

Salmonellae grow as large purple-black colonies due to hydrogen sulphide production. Most contaminating microorganisms that are able to grow develop as small colourless colonies.

Typical Formula

	grams per litre
Yeast Extract	5.0
Peptone	10.0
'Lab Lemco' Powder	2.0
Sodium Chloride	4.0
Mannitol	3.0
L-Lysine	5.0

Materials Provided

CM0783B: 500g of MLCB Agar.

500g of dehydrated MLCB Agar yields approximately 10.2L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers.
- Incubators.
- Quality control organisms.
- Petri dishes.

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

- Do not autoclave or overheat.
- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 24, S 7 and Q 5.

Procedure

Suspend 49.0g in 1 litre of distilled water. Mix and bring gently to the boil with frequent agitation to dissolve the medium completely. Cool to 50°C and pour approximately 20ml into sterile Petri dishes. DO NOT AUTOCLAVE OR OVERHEAT.

Interpretation

Examine for typical large purple-black colonies of H₂S positive *Salmonella*. Search carefully for H₂S negative strains that atypically grow as large mauve-grey colonies with a cratered centre. A proportion may show a black 'bull's-eye'.

Pick all colonies presumed to be *Salmonella* spp. and confirm by biochemical and serological testing.

Most contaminating organisms that are able to grow develop as small colourless colonies.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18 hours @ 37°C

Positive Controls	
Inoculum level: 50 – 200 cfu (colony forming units) Colony count is ≥ 70% of the control medium count	
<i>Salmonella</i> Virchow NCTC 5742	1-2 mm mauve colonies, black bull's-eye centre.
<i>Salmonella</i> Enteritidis ATCC® 13076™	1-2 mm mauve colonies, black bull's-eye centre.
<i>Salmonella</i> Typhimurium ATCC® 14028™	1-2 mm mauve colonies, black bull's-eye centre.
<i>Pseudomonas aeruginosa</i> ATCC® 25668™	1-2 mm mauve colonies, no blackening.
Inoculum level 10 ⁴ -10 ⁶ cfu	
<i>Salmonella gallinarum</i> NCTC 10532	0.5-1.5mm mauve colonies with slight or black centre
Negative Controls	
Inoculum level 10 ² -10 ⁴ cfu	
<i>Escherichia coli</i> ATCC® 25922™	No growth or pinpoint-0.5 mm colourless colonies
<i>Escherichia coli</i> ATCC® 11775™	No growth or pinpoint-0.5 mm colourless colonies

<i>Proteus mirabilis</i> ATCC® 29906™	No growth or pinpoint-0.5 mm colourless/brown colonies, no swarming
---------------------------------------	---

Limitations

Identifications are presumptive and colonies should be confirmed using appropriate methods.

The medium is non-selective but some strains with specific growth requirements may not grow.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of *Salmonella* species (not *Salmonella* Typhi or *Salmonella* Paratyphi A) is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the device. The precision of MLCB Agar (CM0783B) was demonstrated by an overall pass rate of 100% obtained for the product over one year of testing (05.02.2021 – 27.04.2022; 10 batches). This shows that the performance is reproducible.

MLCB Agar (CM0783B) is tested in-house as part of the QC process since the product was launched in 1996. For target organisms, when using 200 cfu inoculum of *Salmonella* Virchow NCTC 5742, *Salmonella* Enteritidis ATCC®13076™, *Salmonella* Typhimurium ATCC®14028™, *Pseudomonas aeruginosa* ATCC®25668™, *Salmonella* Gallinarum NCTC 10532, *Escherichia coli* ATCC®25922™, *Escherichia coli* ATCC®11775™, and *Proteus mirabilis* ATCC®29906™ and incubating the device at specified temperature and duration, the user can recover organisms with colony size and morphology as listed in this document.

Bibliography

- Public Health England. 2021. 'Standards for Microbiology Investigations Identification of *Salmonella* species.' Issue no: 4 | Issue date: 12.03.21. <https://www.gov.uk/government/publications/smi-id-24-identification-of-salmonella-species>
- Ford, L., Glass, K., Veitch, M., Wardell, R., Polkinghorne, B., Dobbins, T. & Kirk, M. D. (2016). Increasing incidence of *Salmonella* in Australia, 2000-2013. *PLoS one*, 11(10), e0163989.
- Edel, W., Kampelmacher, E. (1967). Studies on *Salmonella* infection in fattening calves. *Zentralbl Veterinarmed B* 14(2), 163-9.
- Takao Inoue et al. (1968). Proceedings of the Japanese Society of Veterinary Science. *Jap. J. Vet. Sci.* Number 169 (30).

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit

	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer- To indicate the entity importing the medical device into the locale. Applicable to the European Union.
Made in the United Kingdom	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.
 ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
 NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures.
 All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke,
 RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of modifications introduced
2.0	2023-06-21

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