Brilliance™ Candida Agar

EN

REF CM1002B

Intended Use

Brilliance Candida Agar Base (CM1002B) with Brilliance Candida Selective Supplement (SR0231E) added is a selective differential medium which enables isolation and identification of clinically important *Candida* spp. from clinical samples, and for identifying clinical isolates. Brilliance Candida Agar Base (CM1002B) with Brilliance Candida Selective Supplement (SR0231E) added can be used with a variety of sample types; the main sample types are sputum and urine, as well as bronchial secretion, genital, nasal, ear, throat, wound, skin, gastrostomy and central venous line site samples.

Brilliance Candida Agar Base (CM1002B) is used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having candidiasis. The device is for professional use only, is not automated, nor is it a companion diagnostic.

Summary and Explanation

Candida spp. are commensal yeasts that are found normally as part of the skin, oral, gut and vaginal flora of humans¹. *Candida* spp. are the causative agents of Candidiasis, a fungal infection which typically occurs in immunocompromised groups such as children or individuals with HIV⁺¹⁻¹². *Candida* spp. are also nosocomial pathogens, with invasive medical interventions such as catheterisation or immunosuppressant treatment contributes to the incidence of nosocomial candidiasis, such as urinary tract infections (UTIs) and bloodstream infections (candidaemia)⁸⁻¹³. Moreover, the emergence of antifungal resistance in *Candida* spp. complicates treatment and impacts patient outcomes, especially in immunocompromised or hospitalised individuals¹⁴⁻¹⁸.

Principle of Method

Differentiation of clinically important *Candida* spp. is achieved through the inclusion of two chromogens that are targeted by specific enzymes hexosaminidase and alkaline phosphatase. The action of these enzymes on the chromogens causes release of the coloured component inside the fungal cell, resulting in coloured colonies. The colour produced depends on which enzymes the organisms produce. The presence of hexosaminidase in *C. albicans*, and *Candida dubliniensis* results in green/blue colonies. The presence of hexosaminidase in *C. tropicalis*, and other metabolic reactions that cause localised drop in pH, results in dark blue colonies. The presence of alkaline phosphatase in *Candida krusei* results in brown or pink colonies. The presence of alkaline phosphatase in *Candida glabrata*, *Candida kefyr*, *Candida parapsilosis* and *Candida lusitanae* results in a variety of beige/brown/yellow colours due to the mixture of natural pigmentation and some alkaline phosphatase activity. Experienced users may be able to differentiate these species by colour and colony morphology. The opaque background allows identification of *Candida* spp., especially when mixed infections are present. The medium contains chloramphenicol, which inhibits bacterial growth.

Typical Formula

	grams per litre
Peptone	4.0
Chromogenic mix	13.6
Agar	13.6

Materials Provided

CM1002B: 500g of Brilliance Candida Agar Base

500g of Brilliance Candida Agar Base yields approximately 16.0L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dish
- Supplements (SR0231E)

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

Signal Word: Danger

Hazard Statements

H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled

MBD_BT_IFU-0572

Precautionary Statements

P261 - Avoid breathing dust/fume/gas/mist/vapors/spray
 P285 - In case of inadequate ventilation wear respiratory protection
 P342 + P311 - If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician
 P304 + P340 - IF INHALED: Remove person to fresh air and keep comfortable for breathing

For in vitro diagnostic use only.
 For professional use only.

Inspect the product packaging before first use.
 Do not use the product if there is any visible damage to the packaging (pot or cap).
 Do not use the product beyond the stated expiry date.
 Do not use the device if signs of contamination are present.
 It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
 Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 1, B 1, B 5, B 9, B 11, B 20, B 28, B 39, B 41, B 57 and S 6.

Procedure

Suspend 15.6 grams in 500 ml of distilled water and add the contents of one vial of *Brilliance™* Candida Selective Supplement (SR0231E), reconstituted as directed. Mix well and bring to the boil with frequent agitation. DO NOT AUTOCLAVE. Cool to 45°C, mix well and pour into sterile containers.

Interpretation

Once the medium is reconstituted, the presence of: Green colonies indicates *Candida albicans*.
 Dark blue colonies indicates *Candida tropicalis*.
 Dry irregular pink/brown colonies indicates *Candida krusei*.
 Beige/yellow colonies indicates *Candida kefyr* or *Candida glabrata*
 Yellow/ brown colonies indicates *Candida lusitaniae*
 Brown colonies indicates *Candida parapsilosis*

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 42-48 h @ 30°C

Positive Controls	
Inoculum level: 10 – 100 cfu Colony count is ≥ 70% of the control medium count	
<i>Candida albicans</i> ATCC® 10231™	1-2mm green colonies
<i>Candida albicans</i> ATCC® 18804™	1-2mm green colonies
<i>Candida albicans</i> ATCC® 2091™	0.5-1mm green colonies
<i>Candida tropicalis</i> ATCC® 750™	2-3mm dark blue colonies
<i>Candida krusei</i> ATCC® 6258™	5-10mm dry, irregular, pink/brown colonies
<i>Candida glabrata</i> NCPF® 3240	2-3mm beige/yellow colonies
<i>Candida lusitaniae</i> NCPF® 3516	1.5-2.5mm brown/yellow colonies

<i>Candida parapsilosis</i> ATCC® 22019	0.5-1mm brown colonies
Negative Controls Inoculum level: 10 ⁴ – 10 ⁶ cfu	
<i>Staphylococcus aureus</i> ATCC® 25923	No growth
<i>Escherichia coli</i> ATCC® 25922	No growth
<i>Pseudomonas aeruginosa</i> ATCC® 27853	No growth

Limitations

Identifications are presumptive and should be confirmed.

C. guilliermondii, *C. haemulonii* and *C. auris* may grow on the media, but identification of these organisms should be carried out using additional methods. *C. auris* colonies may appear as a variety of brown/yellow colours, *C. haemulonii* may appear as light brown to brown colonies, and *C. guilliermondii* may appear as dark purple colonies with brown-purple halos.

Candida dubliensis will produce colonies of similar colour to *Candida albicans* and additional testing may be required to differentiate between the two species. The green colour of *C. albicans* and *C. dubliensis* is caused by the same chromogenic reaction as the dark blue colour of *Candida tropicalis*. However, other reactions caused by the medium result in green colonies.

Candida glabrata, *Candida kefyr*, *C. parapsilosis* and *Candida lusitanae* appear as a variety of beige/brown/yellow colours due to the mixture of natural pigmentation and some alkaline phosphatase activity. Experienced users may be able to differentiate these species by colour and colony morphology.

Although most bacteria will be inhibited, those which are resistant to chloramphenicol may grow on the medium. Incubation other than as specified may affect the colour and growth of *Candida* spp.

Brilliance Candida Agar was found to correctly select against all non-target organisms, with the exception of *Saccharomyces cerevisiae* and *Aspergillus brasiliensis*, which grew with golden brown colonies and blue/green mycelium respectively. The mycelium of *A. brasiliensis* was found to be readily distinguishable from the *Candida* strains tested, however the non-target yeast *S. cerevisiae* produced growth highly similar to that of *C. glabrata*, *C. kefyr*, *C. parapsilosis* and *C. lusitanae*.

Analytical Performance

Two internal trials, as part of post-launch activities, were conducted at Thermo Fisher Scientific, Microbiology Division, Basingstoke, to evaluate the performance of *Brilliance Candida Agar*.

In the first study, 28 strains corresponding to ten *Candida* species, as well as 10 non-target organisms were tested. McFarland Solutions (1.0 for yeasts and moulds 0.5 for bacteria) were prepared from isolated colonies, with subsequent plates produced incubated in aerobic conditions at 30°C for 48 hours.

In the second study, 49 strains consisting of ten *Candida* species and two mixed *Candida* species cultures, each containing four strains from different target groups of *Candida*, were tested.¹⁹ Stock plates were produced on Sabouraud Dextrose Agar (PO0410B), from which 1.0 McFarland Solutions were prepared from isolated colonies; 10 µL were then streaked in duplicate onto half-plates of *Brilliance Candida Agar* and incubated in aerobic conditions at 30°C for 48 hours.

Brilliance Candida Agar isolated all *Candida* species tested and inhibited all bacterial growth. The results suggested that *Brilliance Candida Agar* was suitable for identifying several clinically important *Candida* spp.

Clinical Performance

Brilliance Candida Agar has been evaluated through a series of external trials conducted at UK hospitals, which compared and demonstrated the performance of the device in a clinical setting. The performance characteristics sensitivity and specificity have been assessed for this medium.

Brilliance Candida Agar has shown a high level of sensitivity and specificity and can be routinely used for the isolation and identification of clinically important *Candida* spp. from clinical samples including sputum and urine, as well as genital, nasal, ear, throat, wound, skin, gastrostomy and central venous line site samples, and for identifying clinical isolates, detecting clinically important *Candida* spp. following 48 hours incubation.

In one study conducted pre-launch at a UK hospital, six ATCC® isolates, 214 pure cultures of *Candida* spp. and other yeast species, and 14 mixed yeast cultures were cultured onto *Brilliance Candida Agar*. After incubation, growth and colony morphology on the plates was inspected.

Performance of *Brilliance Candida Agar*

Performance Characteristic	<i>Brilliance Candida Agar</i> (%)
----------------------------	------------------------------------

C. albicans (n=64) after 24 hours	
Sensitivity	100
Specificity	97.8
C. tropicalis (n=30) after 48 hours	
Sensitivity	100
Specificity	100
C. krusei (n=30) after 48 hours	
Sensitivity	96.7
Specificity	100

In a second study, the performance of *Brilliance* Candida Agar was evaluated using 319 clinical isolates from three hospitals. Samples were collected from various sources including sputum, genital, urine, ear/nose/throat, gastronomy, central venous sites, blood, and other miscellaneous sources (wound, skin, other).

138 strains were isolated directly from clinical specimens and the remaining 181 strains were sub-cultured from clinical isolates originally isolated on Sabouraud Dextrose Agar by suspending a few colonies in 0.85% saline and streaking 0.01 ml onto prepared plates.

Performance of *Brilliance* Candida Agar

Performance Characteristic	Incubation time	
	24 h	48 h
C. albicans primary isolation		
Sensitivity (%)	78.6 *	100
Specificity (%)	100*	100
C. albicans pure culture (n =181)		
Sensitivity (%)	100*	100
Specificity (%)	100*	100

*Not all results read after 24 h

Performance of *Brilliance* Candida Agar

Performance Characteristic	Primary isolation	Pure culture
C. tropicalis		
Sensitivity (%)	87.5	100
Specificity (%)	100	99.4
C. krusei		
Sensitivity (%)	100	100
Specificity (%)	100	100
Other <i>Candida</i> spp.		
Sensitivity (%)	100	98.4
Specificity (%)	100	100

The study showed that *Brilliance* Candida Agar has shown a high level of sensitivity and specificity and can be used for the isolation and identification of clinically important *Candida* spp. from clinical samples to aid clinicians in determining treatment options for patients suspected of having candidiasis.

Summary of results found in the studies evaluated in a literature review.

Study	Incubation time				
	24 hours		48 hours		72 hours
	Sensitivity	Specificity	Sensitivity	Specificity	Sensitivity
Scharmann <i>et al.</i> , 2020 ²⁰	32%*	69%	82%*	83%	ND
De Angelis <i>et al.</i> , ²¹	44.4%	ND	ND	ND	90.6%
Vecchione <i>et al.</i> , ²²	ND	ND	100%	100%	ND

ND – Not done

*Sensitivity calculated for *Candida albicans* identification

** Sensitivity calculated for all *Candida* spp


Bibliography

- Public Health England. 2015a. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.

2. Müller, F. M.C., A. H. Groll, and T. J. Walsh. 1999. 'Current Approaches to Diagnosis and Treatment of Fungal Infections in Children Infected with Human Immuno Deficiency Virus'. *European Journal of Pediatrics* 158 (3): 187–99. <https://doi.org/10.1007/s004310051047>.
3. Public Health England. 2013. 'Sexually Transmitted Infections'. UK Standards for Microbiology Investigations S 6 (1.2). <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>.
4. Public Health England. 2015b. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
5. Public Health England. 2016. 'Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
6. Public Health England. 2017a. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
7. Public Health England. 2017b. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
8. Public Health England. 2017c. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
9. Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
10. Public Health England. 2019a. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
11. Public Health England. 2019b. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
12. Sganga, Gabriele. 2009. 'Clinical Aspects of Invasive Candidiasis in the Surgical Patient'. *Drugs* 69: 29–32. <https://doi.org/10.2165/11315600-000000000-00000>.
13. Chen, Sharon C.A., D. Marriott, E. G. Playford, Q. Nguyen, D. Ellis, W. Meyer, T. C. Sorrell, et al. 2009. 'Candidaemia with Uncommon Candida Species: Predisposing Factors, Outcome, Antifungal Susceptibility, and Implications for Management'. *Clinical Microbiology and Infection* 15 (7): 662–69. <https://doi.org/10.1111/j.1469-0691.2009.02821.x>.
14. Farooqi, J. Q., K. Jabeen, N. Saeed, N. Iqbal, B. Malik, S. R. Lockhart, A. Zafar, M. E. Brandt, and R. Hasan. 2013. 'Invasive Candidiasis in Pakistan: Clinical Characteristics, Species Distribution and Antifungal Susceptibility'. *Journal of Medical Microbiology* 62 (2): 259–68. <https://doi.org/10.1099/jmm.0.048785-0>.
15. Pongrácz, Júlia, Emese Juhász, Miklós Iván, and Katalin Kristóf. 2015. 'Significance of Yeasts in Bloodstream Infection: Epidemiology and Predisposing Factors of Candidaemia in Adult Patients at a University Hospital (2010-2014)'. *Acta Microbiologica et Immunologica Hungarica* 62 (3): 317–29. <https://doi.org/10.1556/030.62.2015.3.9>.
16. Khadka, Sundar, Jeevan Bahadur Sherchand, Bharat Mani Pokhrel, Keshab Parajuli, Shyam Kumar Mishra, Sangita Sharma, Niranjan Shah, et al. 2017. 'Isolation, Speciation and Antifungal Susceptibility Testing of Candida Isolates from Various Clinical Specimens at a Tertiary Care Hospital, Nepal'. *BMC Research Notes* 10 (1): 1–5. <https://doi.org/10.1186/s13104-017-2547-3>.
17. Xiang, Guo Shi, Ling Bing Zeng, Jie Yu Zhang, Ying Ying, Xi Ren Deng, Xue Fei Hu, Yan Yan, Ke Hua Yu, Wei Wen Zou, and Xiao Tian Huang. 2018. 'In Vitro Antidrug Susceptibility Testing of Candida Species Isolated from Aseptic Body Fluids'. *Jundishapur Journal of Microbiology* 11 (8). <https://doi.org/10.5812/jjm.55547>.
18. Xiao, Meng, Xin Fan, Xin Hou, Sharon C.A. Chen, He Wang, Fanrong Kong, Zi Yong Sun, Yun Zhuo Chu, and Ying Chun Xu. 2018. 'Clinical Characteristics of the First Cases of Invasive Candidiasis in China Due to Pan-Echinocandin-Resistant Candida Tropicalis and Candida Glabrata Isolates with Delineation of Their Resistance Mechanisms'. *Infection and Drug Resistance* 11: 155–61. <https://doi.org/10.2147/IDR.S152785>.
19. Data held on file.
20. Scharmann, Ulrike, Lisa Kirchoff, Valerie le Saout Chapot, Jan Dziobaka, Hedda Luise Verhasselt, Raphael Stauf, Jan Buer, Joerg Steinmann, and Peter Michael Rath. 2020. 'Comparison of Four Commercially Available Chromogenic Media to Identify Candida Albicans and Other Medically Relevant Candida Species'. *Mycoses* 63: 823–31. <https://doi.org/10.1111/myc.13119>.
21. Angelis, Giulia De, Giulia Menchinelli, Riccardo Torelli, Elena De Carolis, Patrizia Posteraro, Maurizio Sanguinetti, and Brunella Posteraro. 2020. 'Different Detection Capabilities by Mycological Media for Candida Isolates from Mono- or Dual-Species Cultures'. *PLoS ONE* 15 (3): 1–9. <https://doi.org/10.1371/journal.pone.0226467>.
22. Vecchione, Alessandra, Walter Florio, Francesco Celandroni, Simona Barnini, Antonella Lupetti, and Emilia Ghelardi. 2017. 'Comparative Evaluation of Six Chromogenic Media for Presumptive Yeast Identification'. *Journal of Clinical Pathology* 70 (12): 1074–78. <https://doi.org/10.1136/jclinpath-2017-204396>.

Symbol Legend

Symbol	Definition
REF	Catalogue number
IVD	In Vitro Diagnostic Medical Device
LOT	Batch code

	Temperature limit
	Use-by date
	Keep away from sunlight
	Consult instructions for use or consult electronic instructions for use
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of modifications introduced
4.0	2024-10-28 Minor formatting changes Warnings and Precautions corrected to match SDS Updated limitations section Updated analytical/clinical performance section Updated bibliography section

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Literatura

1. Public Health England. 2015a. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
2. Müller, F. M.C., A. H. Groll, and T. J. Walsh. 1999. 'Current Approaches to Diagnosis and Treatment of Fungal Infections in Children Infected with Human Immuno Deficiency Virus'. *European Journal of Pediatrics* 158 (3): 187–99. <https://doi.org/10.1007/s004310051047>.
3. Public Health England. 2013. 'Sexually Transmitted Infections'. UK Standards for Microbiology Investigations S 6 (1.2). <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>.
4. Public Health England. 2015b. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
5. Public Health England. 2016. 'Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
6. Public Health England. 2017a. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
7. Public Health England. 2017b. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
8. Public Health England. 2017c. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
9. Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
10. Public Health England. 2019a. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
11. Public Health England. 2019b. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
12. Sganga, Gabriele. 2009. 'Clinical Aspects of Invasive Candidiasis in the Surgical Patient'. *Drugs* 69: 29–32. <https://doi.org/10.2165/11315600-000000000-00000>.
13. Chen, Sharon C.A., D. Marriott, E. G. Playford, Q. Nguyen, D. Ellis, W. Meyer, T. C. Sorrell et al. 2009. 'Candidaemia with Uncommon Candida Species: Predisposing Factors, Outcome, Antifungal Susceptibility, and Implications for Management'. *Clinical Microbiology and Infection* 15 (7): 662–69. <https://doi.org/10.1111/j.1469-0691.2009.02821.x>.
14. Farooqi, J. Q., K. Jabeen, N. Saeed, N. Iqbal, B. Malik, S. R. Lockhart, A. Zafar, M. E. Brandt, and R. Hasan. 2013. 'Invasive Candidiasis in Pakistan: Clinical Characteristics, Species Distribution and Antifungal Susceptibility'. *Journal of Medical Microbiology* 62 (2): 259–68. <https://doi.org/10.1099/jmm.0.048785-0>.
15. Pongrácz, Júlia, Emese Juhász, Miklós Iván, and Katalin Kristóf. 2015. 'Significance of Yeasts in Bloodstream Infection: Epidemiology and Predisposing Factors of Candidaemia in Adult Patients at a University Hospital (2010-2014)'. *Acta Microbiologica et Immunologica Hungarica* 62 (3): 317–29. <https://doi.org/10.1556/030.62.2015.3.9>.
16. Khadka, Sundar, Jeevan Bahadur Sherchand, Bharat Mani Pokhrel, Keshab Parajuli, Shyam Kumar Mishra, Sangita Sharma, Niranjan Shah, et al. 2017. 'Isolation, Speciation and Antifungal Susceptibility Testing of Candida Isolates from Various Clinical Specimens at a Tertiary Care Hospital, Nepal'. *BMC Research Notes* 10 (1): 1–5. <https://doi.org/10.1186/s13104-017-2547-3>.
17. Xiang, Guo Shi, Ling Bing Zeng, Jie Yu Zhang, Ying Ying, Xi Ren Deng, Xue Fei Hu, Yan Yan, Ke Hua Yu, Wei Wen Zou, and Xiao Tian Huang. 2018. 'In Vitro Antidrug Susceptibility Testing of Candida Species Isolated from Aseptic Body Fluids'. *Jundishapur Journal of Microbiology* 11 (8). <https://doi.org/10.5812/jjm.55547>.
18. Xiao, Meng, Xin Fan, Xin Hou, Sharon C.A. Chen, He Wang, Fanrong Kong, Zi Yong Sun, Yun Zhuo Chu, and Ying Chun Xu. 2018. 'Clinical Characteristics of the First Cases of Invasive Candidiasis in China Due to Pan-Echinocandin-Resistant Candida Tropicalis and Candida Glabrata Isolates with Delineation of Their Resistance Mechanisms'. *Infection and Drug Resistance* 11: 155–61. <https://doi.org/10.2147/IDR.S152785>.
19. Data held on file.
20. Scharmman, Ulrike, Lisa Kirchhoff, Valerie le Saout Chapot, Jan Dziobaka, Hedda Luise Verhasselt, Raphael Stauf, Jan Buer, Joerg Steinmann, and Peter Michael Rath. 2020. 'Comparison of Four Commercially Available Chromogenic Media to Identify Candida Albicans and Other Medically Relevant Candida Species'. *Mycoses* 63: 823–31. <https://doi.org/10.1111/myc.13119>.
21. Angelis, Giulia De, Giulia Menchinelli, Riccardo Torelli, Elena De Carolis, Patrizia Posteraro, Maurizio Sanguinetti, and Brunella Posteraro. 2020. 'Different Detection Capabilities by Mycological Media for Candida Isolates from Mono- or Dual-Species Cultures'. *PLoS ONE* 15 (3): 1–9. <https://doi.org/10.1371/journal.pone.0226467>.
22. Vecchione, Alessandra, Walter Florio, Francesco Celandroni, Simona Barnini, Antonella Lupetti, and Emilia Ghelardi. 2017. 'Comparative Evaluation of Six Chromogenic Media for Presumptive Yeast Identification'. *Journal of Clinical Pathology* 70 (12): 1074–78. <https://doi.org/10.1136/iclinpath-2017-204396>.

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** Sensitivitet beregnet for alle Candida spp

Litteratur

1. Public Health England. 2015a. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
2. Müller, F. M.C., A. H. Groll, and T. J. Walsh. 1999. 'Current Approaches to Diagnosis and Treatment of Fungal Infections in Children Infected with Human Immuno Deficiency Virus'. *European Journal of Pediatrics* 158 (3): 187–99. <https://doi.org/10.1007/s004310051047>.
3. Public Health England. 2013. 'Sexually Transmitted Infections'. UK Standards for Microbiology Investigations S 6 (1.2). <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>.
4. Public Health England. 2015b. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
5. Public Health England. 2016. 'Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
6. Public Health England. 2017a. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
7. Public Health England. 2017b. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
8. Public Health England. 2017c. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
9. Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
10. Public Health England. 2019a. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
11. Public Health England. 2019b. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
12. Sganga, Gabriele. 2009. 'Clinical Aspects of Invasive Candidiasis in the Surgical Patient'. *Drugs* 69: 29–32. <https://doi.org/10.2165/11315600-000000000-00000>.
13. Chen, Sharon C.A., D. Marriott, E. G. Playford, Q. Nguyen, D. Ellis, W. Meyer, T. C. Sorrell, et al. 2009. 'Candidaemia with Uncommon Candida Species: Predisposing Factors, Outcome, Antifungal Susceptibility, and Implications for Management'. *Clinical Microbiology and Infection* 15 (7): 662–69. <https://doi.org/10.1111/j.1469-0691.2009.02821.x>.
14. Farooqi, J. Q., K. Jabeen, N. Saeed, N. Iqbal, B. Malik, S. R. Lockhart, A. Zafar, M. E. Brandt, and R. Hasan. 2013. 'Invasive Candidiasis in Pakistan: Clinical Characteristics, Species Distribution and Antifungal Susceptibility'. *Journal of Medical Microbiology* 62 (2): 259–68. <https://doi.org/10.1099/jmm.0.048785-0>.
15. Pongrácz, Júlia, Emese Juhász, Miklós Iván, and Katalin Kristóf. 2015. 'Significance of Yeasts in Bloodstream Infection: Epidemiology and Predisposing Factors of Candidaemia in Adult Patients at a University Hospital (2010-2014)'. *Acta Microbiologica et Immunologica Hungarica* 62 (3): 317–29. <https://doi.org/10.1556/030.62.2015.3.9>.
16. Khadka, Sundar, Jeevan Bahadur Sherchand, Bharat Mani Pokhrel, Keshab Parajuli, Shyam Kumar Mishra, Sangita Sharma, Niranjana Shah, et al. 2017. 'Isolation, Speciation and Antifungal Susceptibility Testing of Candida Isolates from Various Clinical Specimens at a Tertiary Care Hospital, Nepal'. *BMC Research Notes* 10 (1): 1–5. <https://doi.org/10.1186/s13104-017-2547-3>.
17. Xiang, Guo Shi, Ling Bing Zeng, Jie Yu Zhang, Ying Ying, Xi Ren Deng, Xue Fei Hu, Yan Yan, Ke Hua Yu, Wei Wen Zou, and Xiao Tian Huang. 2018. 'In Vitro Antidrug Susceptibility Testing of Candida Species Isolated from Aseptic Body Fluids'. *Jundishapur Journal of Microbiology* 11 (8). <https://doi.org/10.5812/ijm.55547>.
18. Xiao, Meng, Xin Fan, Xin Hou, Sharon C.A. Chen, He Wang, Fanrong Kong, Zi Yong Sun, Yun Zhuo Chu, and Ying Chun Xu. 2018. 'Clinical Characteristics of the First Cases of Invasive Candidiasis in China Due to Pan-Echinocandin-Resistant Candida Tropicalis and Candida Glabrata Isolates with Delineation of Their Resistance Mechanisms'. *Infection and Drug Resistance* 11: 155–61. <https://doi.org/10.2147/IDR.S152785>.
19. Data i arkiv.
20. Scharmann, Ulrike, Lisa Kirchoff, Valerie le Saout Chapot, Jan Dziobaka, Hedda Luise Verhasselt, Raphael Stauf, Jan Buer, Joerg Steinmann, and Peter Michael Rath. 2020. 'Comparison of Four Commercially Available Chromogenic Media to Identify Candida Albicans and Other Medically Relevant Candida Species'. *Mycoses* 63: 823–31. <https://doi.org/10.1111/myc.13119>.
21. Angelis, Giulia De, Giulia Menchinelli, Riccardo Torelli, Elena De Carolis, Patrizia Posteraro, Maurizio Sanguinetti, and Brunella Posteraro. 2020. 'Different Detection Capabilities by Mycological Media for Candida Isolates from Mono- or Dual-Species Cultures'. *PLoS ONE* 15 (3): 1–9. <https://doi.org/10.1371/journal.pone.0226467>.
22. Vecchione, Alessandra, Walter Florio, Francesco Celandroni, Simona Barnini, Antonella Lupetti, and Emilia Ghelardi. 2017. 'Comparative Evaluation of Six Chromogenic Media for Presumptive Yeast Identification'. *Journal of Clinical Pathology* 70 (12): 1074–78. <https://doi.org/10.1136/jclinpath-2017-204396>.

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Literaturverzeichnis

- Public Health England. 2015a. „Investigation of Nasal Samples“. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
- Müller, F. M.C., A. H. Groll und T. J. Walsh. 1999. „Current Approaches to Diagnosis and Treatment of Fungal Infections in Children Infected with Human Immuno Deficiency Virus“. European Journal of Pediatrics 158 (3): 187–99. <https://doi.org/10.1007/s004310051047>.
- Public Health England. 2013. „Sexually Transmitted Infections“. UK Standards for Microbiology Investigations S 6 (1.2). <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>.
- Public Health England. 2015b. „Investigation of Throat Related Specimens.“ UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
- Public Health England. 2016. „Investigation of Dermatological Specimens for Superficial Mycoses“. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
- Public Health England. 2017a. „Genital Tract and Associated Specimens“. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
- Public Health England. 2017b. „Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture“. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
- Public Health England. 2017c. „Investigation of Intravascular Cannulae and Associated Specimens“. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
- Public Health England. 2018. „Swabs from Skin and Superficial Soft Tissue Infections“. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
- Public Health England. 2019a. „Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens“. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
- Public Health England. 2019b. „Investigation of Urine“. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
- Sganga, Gabriele. 2009. „Clinical Aspects of Invasive Candidiasis in the Surgical Patient“. Drugs 69: 29–32. <https://doi.org/10.2165/11315600-000000000-00000>.
- Chen, Sharon C.A., D. Marriott, E. G. Playford, Q. Nguyen, D. Ellis, W. Meyer, T. C. Sorrell et al. 2009. „Candidaemia with Uncommon Candida Species: Predisposing Factors, Outcome, Antifungal Susceptibility, and Implications for Management“. Clinical Microbiology and Infection 15 (7): 662–69. <https://doi.org/10.1111/j.1469-0691.2009.02821.x>.
- Farooqi, J. Q., K. Jabeen, N. Saeed, N. Iqbal, B. Malik, S. R. Lockhart, A. Zafar, M. E. Brandt und R. Hasan. 2013. „Invasive Candidiasis in Pakistan: Clinical Characteristics, Species Distribution and Antifungal Susceptibility“. Journal of Medical Microbiology 62 (2): 259–68. <https://doi.org/10.1099/jmm.0.048785-0>.
- Pongrácz, Júlia, Emese Juhász, Miklós Iván und Katalin Kristóf. 2015. ‘Significance of Yeasts in Bloodstream Infection: Epidemiology and Predisposing Factors of Candidaemia in Adult Patients at a University Hospital (2010-2014)’. Acta Microbiologica et Immunologica Hungarica 62 (3): 317–29. <https://doi.org/10.1556/030.62.2015.3.9>.
- Khadka, Sundar, Jeevan Bahadur Sherchand, Bharat Mani Pokhrel, Keshab Parajuli, Shyam Kumar Mishra, Sangita Sharma, Niranjana Shah et al. 2017. „Isolation, Speciation and Antifungal Susceptibility Testing of Candida Isolates from Various Clinical Specimens at a Tertiary Care Hospital, Nepal“. BMC Research Notes 10 (1): 1–5. <https://doi.org/10.1186/s13104-017-2547-3>.
- Xiang, Guo Shi, Ling Bing Zeng, Jie Yu Zhang, Ying Ying, Xi Ren Deng, Xue Fei Hu, Yan Yan, Ke Hua Yu, Wei Wen Zou, and Xiao Tian Huang. 2018. „In Vitro Antidrug Susceptibility Testing of Candida Species Isolated from Aseptic Body Fluids“. Jundishapur Journal of Microbiology 11 (8). <https://doi.org/10.5812/jjm.55547>.
- Xiao, Meng, Xin Fan, Xin Hou, Sharon C.A. Chen, He Wang, Fanrong Kong, Zi Yong Sun, Yun Zhuo Chu, and Ying Chun Xu. 2018. „Clinical Characteristics of the First Cases of Invasive Candidiasis in China Due to Pan-Echinocandin-Resistant

Candida Tropicalis and Candida Glabrata Isolates with Delineation of Their Resistance Mechanisms". Infection and Drug Resistance 11: 155–61. <https://doi.org/10.2147/IDR.S152785>.

19. In den Akten gespeicherte Daten.
20. Scharmann, Ulrike, Lisa Kirchhoff, Valerie le Saout Chapot, Jan Dziobaka, Hedda Luise Verhasselt, Raphael Stauf, Jan Buer, Joerg Steinmann und Peter Michael Rath. 2020. „Comparison of Four Commercially Available Chromogenic Media to Identify Candida Albicans and Other Medically Relevant Candida Species". Mycoses 63: 823–31. <https://doi.org/10.1111/myc.13119>.
21. Angelis, Giulia De, Giulia Menchinelli, Riccardo Torelli, Elena De Carolis, Patrizia Posteraro, Maurizio Sanguinetti und Brunella Posteraro. 2020. „Different Detection Capabilities by Mycological Media for Candida Isolates from Mono- or Dual-Species Cultures". PLoS ONE 15 (3): 1–9. <https://doi.org/10.1371/journal.pone.0226467>.
22. Vecchione, Alessandra, Walter Florio, Francesco Celandroni, Simona Barnini, Antonella Lupetti und Emilia Ghelardi. 2017. „Comparative Evaluation of Six Chromogenic Media for Presumptive Yeast Identification". Journal of Clinical Pathology 70 (12): 1074–78. <https://doi.org/10.1136/jclinpath-2017-204396>.

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- Public Health England. 2015a. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
- Müller, F. M.C., A. H. Groll, and T. J. Walsh. 1999. 'Current Approaches to Diagnosis and Treatment of Fungal Infections in Children Infected with Human Immuno Deficiency Virus'. European Journal of Pediatrics 158 (3): 187–99. <https://doi.org/10.1007/s004310051047>.
- Public Health England. 2013. 'Sexually Transmitted Infections'. UK Standards for Microbiology Investigations S 6 (1.2). <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>.
- Public Health England. 2015b. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
- Public Health England. 2016. 'Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
- Public Health England. 2017a. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
- Public Health England. 2017b. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
- Public Health England. 2017c. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
- Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
- Public Health England. 2019a. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
- Public Health England. 2019b. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
- Sganga, Gabriele. 2009. 'Clinical Aspects of Invasive Candidiasis in the Surgical Patient'. Drugs 69: 29–32. <https://doi.org/10.2165/11315600-000000000-00000>.
- Chen, Sharon C.A., D. Marriott, E. G. Playford, Q. Nguyen, D. Ellis, W. Meyer, T. C. Sorrell, et al. 2009. 'Candidaemia with Uncommon Candida Species: Predisposing Factors, Outcome, Antifungal Susceptibility, and Implications for Management'. Clinical Microbiology and Infection 15 (7): 662–69. <https://doi.org/10.1111/j.1469-0691.2009.02821.x>.
- Farooqi, J. Q., K. Jabeen, N. Saeed, N. Iqbal, B. Malik, S. R. Lockhart, A. Zafar, M. E. Brandt, and R. Hasan. 2013. 'Invasive Candidiasis in Pakistan: Clinical Characteristics, Species Distribution and Antifungal Susceptibility'. Journal of Medical Microbiology 62 (2): 259–68. <https://doi.org/10.1099/jmm.0.048785-0>.
- Pongrácz, Júlia, Emese Juhász, Miklós Iván, and Katalin Kristóf. 2015. 'Significance of Yeasts in Bloodstream Infection: Epidemiology and Predisposing Factors of Candidaemia in Adult Patients at a University Hospital (2010-2014)'. Acta Microbiologica et Immunologica Hungarica 62 (3): 317–29. <https://doi.org/10.1556/030.62.2015.3.9>.
- Khadka, Sundar, Jeevan Bahadur Sherchand, Bharat Mani Pokhrel, Keshab Parajuli, Shyam Kumar Mishra, Sangita Sharma, Niranjana Shah, et al. 2017. 'Isolation, Speciation and Antifungal Susceptibility Testing of Candida Isolates from Various Clinical Specimens at a Tertiary Care Hospital, Nepal'. BMC Research Notes 10 (1): 1–5. <https://doi.org/10.1186/s13104-017-2547-3>.
- Xiang, Guo Shi, Ling Bing Zeng, Jie Yu Zhang, Ying Ying, Xi Ren Deng, Xue Fei Hu, Yan Yan, Ke Hua Yu, Wei Wen Zou, and Xiao Tian Huang. 2018. 'In Vitro Antidrug Susceptibility Testing of Candida Species Isolated from Aseptic Body Fluids'. Jundishapur Journal of Microbiology 11 (8). <https://doi.org/10.5812/jjm.55547>.
- Xiao, Meng, Xin Fan, Xin Hou, Sharon C.A. Chen, He Wang, Fanrong Kong, Zi Yong Sun, Yun Zhuo Chu, and Ying Chun Xu. 2018. 'Clinical Characteristics of the First Cases of Invasive Candidiasis in China Due to Pan-Echinocandin-Resistant

Candida Tropicalis and Candida Glabrata Isolates with Delineation of Their Resistance Mechanisms'. Infection and Drug Resistance 11: 155–61. <https://doi.org/10.2147/IDR.S152785>.

19. Τα δεδομένα τηρούνται στο αρχείο.
20. Scharmann, Ulrike, Lisa Kirchhoff, Valerie le Saout Chapot, Jan Dziobaka, Hedda Luise Verhasselt, Raphael Stauf, Jan Buer, Joerg Steinmann, and Peter Michael Rath. 2020. 'Comparison of Four Commercially Available Chromogenic Media to Identify Candida Albicans and Other Medically Relevant Candida Species'. Mycoses 63: 823–31. <https://doi.org/10.1111/myc.13119>.
21. Angelis, Giulia De, Giulia Menchinelli, Riccardo Torelli, Elena De Carolis, Patrizia Posteraro, Maurizio Sanguinetti, and Brunella Posteraro. 2020. 'Different Detection Capabilities by Mycological Media for Candida Isolates from Mono- or Dual-Species Cultures'. PLoS ONE 15 (3): 1–9. <https://doi.org/10.1371/journal.pone.0226467>.
22. Vecchione, Alessandra, Walter Florio, Francesco Celandroni, Simona Barnini, Antonella Lupetti, and Emilia Ghelardi. 2017. 'Comparative Evaluation of Six Chromogenic Media for Presumptive Yeast Identification'. Journal of Clinical Pathology 70 (12): 1074–78. <https://doi.org/10.1136/jclinpath-2017-204396>.

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
- Public Health England. 2015a. «Investigation of Nasal Samples». UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
- Müller, F. M.C., A. H. Groll, and T. J. Walsh. 1999. «Current Approaches to Diagnosis and Treatment of Fungal Infections in Children Infected with Human Immuno Deficiency Virus». European Journal of Pediatrics 158 (3): 187–99. <https://doi.org/10.1007/s004310051047>.
- Public Health England. 2013. «Sexually Transmitted Infections». UK Standards for Microbiology Investigations S 6 (1.2). <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>.
- Public Health England. 2015b. «Investigation of Throat Related Specimens». UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
- Public Health England. 2016. «Investigation of Dermatological Specimens for Superficial Mycoses». UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
- Public Health England. 2017a. «Genital Tract and Associated Specimens». UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
- Public Health England. 2017b. «Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture». UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
- Public Health England. 2017c. «Investigation of Intravascular Cannulae and Associated Specimens». UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
- Public Health England. 2018. «Swabs from Skin and Superficial Soft Tissue Infections». UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
- Public Health England. 2019a. «Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens». UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
- Public Health England. 2019b. «Investigation of Urine». UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
- Sganga, Gabriele. 2009. «Clinical Aspects of Invasive Candidiasis in the Surgical Patient». Drugs 69: 29–32. <https://doi.org/10.2165/11315600-000000000-00000>.
- Chen, Sharon C.A., D. Marriott, E. G. Playford, Q. Nguyen, D. Ellis, W. Meyer, T. C. Sorrell, et al. 2009. «Candidaemia with Uncommon Candida Species: Predisposing Factors, Outcome, Antifungal Susceptibility, and Implications for Management». Clinical Microbiology and Infection 15 (7): 662–69. <https://doi.org/10.1111/j.1469-0691.2009.02821.x>.
- Farooqi, J. Q., K. Jabeen, N. Saeed, N. Iqbal, B. Malik, S. R. Lockhart, A. Zafar, M. E. Brandt, and R. Hasan. 2013. «Invasive Candidiasis in Pakistan: Clinical Characteristics, Species Distribution and Antifungal Susceptibility». Journal of Medical Microbiology 62 (2): 259–68. <https://doi.org/10.1099/jmm.0.048785-0>.
- Pongrácz, Júlia, Emese Juhász, Miklós Iván, and Katalin Kristóf. 2015. «Significance of Yeasts in Bloodstream Infection: Epidemiology and Predisposing Factors of Candidaemia in Adult Patients at a University Hospital (2010-2014)». Acta Microbiologica et Immunologica Hungarica 62 (3): 317–29. <https://doi.org/10.1556/030.62.2015.3.9>.
- Khadka, Sundar, Jeevan Bahadur Sherchand, Bharat Mani Pokhrel, Keshab Parajuli, Shyam Kumar Mishra, Sangita Sharma, Niranjana Shah, et al. 2017. «Isolation, Speciation and Antifungal Susceptibility Testing of Candida Isolates from Various Clinical Specimens at a Tertiary Care Hospital, Nepal». BMC Research Notes 10 (1): 1–5. <https://doi.org/10.1186/s13104-017-2547-3>.
- Xiang, Guo Shi, Ling Bing Zeng, Jie Yu Zhang, Ying Ying, Xi Ren Deng, Xue Fei Hu, Yan Yan, Ke Hua Yu, Wei Wen Zou, and Xiao Tian Huang. 2018. «In Vitro Antidrug Susceptibility Testing of Candida Species Isolated from Aseptic Body Fluids». Jundishapur Journal of Microbiology 11 (8). <https://doi.org/10.5812/jjm.55547>.
- Xiao, Meng, Xin Fan, Xin Hou, Sharon C.A. Chen, He Wang, Fanrong Kong, Zi Yong Sun, Yun Zhuo Chu, and Ying Chun Xu. 2018. «Clinical Characteristics of the First Cases of Invasive Candidiasis in China Due to Pan-Echinocandin-Resistant Candida Tropicalis and Candida Glabrata Isolates with Delineation of Their Resistance Mechanisms». Infection and Drug Resistance 11: 155–61. <https://doi.org/10.2147/IDR.S152785>.
- Data held on file.
- Scharmman, Ulrike, Lisa Kirchoff, Valerie le Saout Chapot, Jan Dziobaka, Hedda Luise Verhasselt, Raphael Stauf, Jan Buer, Joerg Steinmann, and Peter Michael Rath. 2020. «Comparison of Four Commercially Available Chromogenic Media to Identify Candida Albicans and Other Medically Relevant Candida Species». Mycoses 63: 823–31. <https://doi.org/10.1111/myc.13119>.
- Angelis, Giulia De, Giulia Menchinelli, Riccardo Torelli, Elena De Carolis, Patrizia Posteraro, Maurizio Sanguinetti, and Brunella Posteraro. 2020. «Different Detection Capabilities by Mycological Media for Candida Isolates from Mono- or Dual-Species Cultures». PLoS ONE 15 (3): 1–9. <https://doi.org/10.1371/journal.pone.0226467>.

22. Vecchione, Alessandra, Walter Florio, Francesco Celandroni, Simona Barnini, Antonella Lupetti, and Emilia Ghelardi. 2017. «Comparative Evaluation of Six Chromogenic Media for Presumptive Yeast Identification». *Journal of Clinical Pathology* 70 (12): 1074–78. <https://doi.org/10.1136/jclinpath-2017-204396>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Reino Unido

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Bibliograafia

- Public Health England, 2015a. Investigation of Nasal Samples. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
- Müller, F. M.C., A. H. Groll ja T. J. Walsh. 1999. Current Approaches to Diagnosis and Treatment of Fungal Infections in Children Infected with Human Immuno Deficiency Virus. *European Journal of Pediatrics* 158 (3): 187–199. <https://doi.org/10.1007/s004310051047>.
- Public Health England, 2013. Sexually Transmitted Infections. UK Standards for Microbiology Investigations S 6 (1.2). <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>.
- Public Health England, 2015b. Investigation of Throat Related Specimens. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
- Public Health England, 2016. Investigation of Dermatological Specimens for Superficial Mycoses. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
- Public Health England, 2017a. Genital Tract and Associated Specimens. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
- Public Health England, 2017b. Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
- Public Health England, 2017c. Investigation of Intravascular Cannulae and Associated Specimens. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
- Public Health England, 2018. Swabs from Skin and Superficial Soft Tissue Infections. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
- Public Health England, 2019a. Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
- Public Health England, 2019b. Investigation of Urine. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
- Sganga, Gabriele. 2009. Clinical Aspects of Invasive Candidiasis in the Surgical Patient. *Drugs* 69: 29–32. <https://doi.org/10.2165/11315600-000000000-00000>.
- Chen, Sharon C.A., D. Marriott, E. G. Playford, Q. Nguyen, D. Ellis, W. Meyer, T. C. Sorrell, et al., 2009. Candidaemia with Uncommon Candida Species: Predisposing Factors, Outcome, Antifungal Susceptibility, and Implications for Management. *Clinical Microbiology and Infection* 15 (7): 662–669. <https://doi.org/10.1111/j.1469-0691.2009.02821.x>.
- Farooqi, J. Q., K. Jabeen, N. Saeed, N. Iqbal, B. Malik, S. R. Lockhart, A. Zafar, M. E. Brandt ja R. Hasan. 2013. Invasive Candidiasis in Pakistan: Clinical Characteristics, Species Distribution and Antifungal Susceptibility. *Journal of Medical Microbiology* 62 (2): 259–268. <https://doi.org/10.1099/jmm.0.048785-0>.
- Pongrácz, Júlia, Emese Juhász, Miklós Iván ja Katalin Kristóf. 2015. Significance of Yeasts in Bloodstream Infection: Epidemiology and Predisposing Factors of Candidaemia in Adult Patients at a University Hospital (2010-2014). *Acta Microbiologica et Immunologica Hungarica* 62 (3): 317–329. <https://doi.org/10.1556/030.62.2015.3.9>.
- Khadka, Sundar, Jeevan Bahadur Sherchand, Bharat Mani Pokhrel, Keshab Parajuli, Shyam Kumar Mishra, Sangita Sharma, Niranjana Shah et al., 2017. Isolation, Speciation and Antifungal Susceptibility Testing of Candida Isolates from Various Clinical Specimens at a Tertiary Care Hospital, Nepal. *BMC Research Notes* 10 (1): 1–5. <https://doi.org/10.1186/s13104-017-2547-3>.
- Xiang, Guo Shi, Ling Bing Zeng, Jie Yu Zhang, Ying Ying, Xi Ren Deng, Xue Fei Hu, Yan Yan, Ke Hua Yu, Wei Wen Zou ja Xiao Tian Huang. 2018. In Vitro Antidrug Susceptibility Testing of Candida Species Isolated from Aseptic Body Fluids. *Jundishapur Journal of Microbiology* 11 (8). <https://doi.org/10.5812/jjm.55547>.
- Xiao, Meng, Xin Fan, Xin Hou, Sharon C.A. Chen, He Wang, Fanrong Kong, Zi Yong Sun, Yun Zhuo Chu ja Ying Chun Xu. 2018. Clinical Characteristics of the First Cases of Invasive Candidiasis in China Due to Pan-Echinocandin-Resistant Candida Tropicalis and Candida Glabrata Isolates with Delineation of Their Resistance Mechanisms. *Infection and Drug Resistance* 11: 155–161. <https://doi.org/10.2147/IDR.S152785>.
- Failis hoitavad andmed.
- Scharmman, Ulrike, Lisa Kirchoff, Valerie le Saout Chapot, Jan Dziobaka, Hedda Luise Verhasselt, Raphael Stauf, Jan Buer, Joerg Steinmann ja Peter Michael Rath. 2020. Comparison of Four Commercially Available Chromogenic Media to Identify Candida Albicans and Other Medically Relevant Candida Species. *Mycoses* 63: 823–831. <https://doi.org/10.1111/myc.13119>.
- Angelis, Giulia De, Giulia Menchinelli, Riccardo Torelli, Elena De Carolis, Patrizia Posteraro, Maurizio Sanguinetti ja Brunella Posteraro. 2020. Different Detection Capabilities by Mycological Media for Candida Isolates from Mono- or Dual-Species Cultures. *PLoS ONE* 15 (3): 1–9. <https://doi.org/10.1371/journal.pone.0226467>.
- Vecchione, Alessandra, Walter Florio, Francesco Celandroni, Simona Barnini, Antonella Lupetti ja Emilia Ghelardi. 2017. Comparative Evaluation of Six Chromogenic Media for Presumptive Yeast Identification. *Journal of Clinical Pathology* 70 (12): 1074–1078. <https://doi.org/10.1136/clinpath-2017-204396>.

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
Bibliographie

- Public Health England. 2015a. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
- Müller, F. M.C., A. H. Groll, and T. J. Walsh. 1999. 'Current Approaches to Diagnosis and Treatment of Fungal Infections in Children Infected with Human Immuno Deficiency Virus'. *European Journal of Pediatrics* 158 (3): 187–99. <https://doi.org/10.1007/s004310051047>.
- Public Health England. 2013. 'Sexually Transmitted Infections'. UK Standards for Microbiology Investigations S 6 (1.2). <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>.
- Public Health England. 2015b. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
- Public Health England. 2016. 'Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
- Public Health England. 2017a. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
- Public Health England. 2017b. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
- Public Health England. 2017c. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
- Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
- Public Health England. 2019a. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
- Public Health England. 2019b. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
- Sganga, Gabriele. 2009. 'Clinical Aspects of Invasive Candidiasis in the Surgical Patient'. *Drugs* 69: 29–32. <https://doi.org/10.2165/11315600-000000000-00000>.
- Chen, Sharon C.A., D. Marriott, E. G. Playford, Q. Nguyen, D. Ellis, W. Meyer, T. C. Sorrell, et al. 2009. 'Candidaemia with Uncommon Candida Species: Predisposing Factors, Outcome, Antifungal Susceptibility, and Implications for Management'. *Clinical Microbiology and Infection* 15 (7): 662–69. <https://doi.org/10.1111/j.1469-0691.2009.02821.x>.
- Farooqi, J. Q., K. Jabeen, N. Saeed, N. Iqbal, B. Malik, S. R. Lockhart, A. Zafar, M. E. Brandt, and R. Hasan. 2013. 'Invasive Candidiasis in Pakistan: Clinical Characteristics, Species Distribution and Antifungal Susceptibility'. *Journal of Medical Microbiology* 62 (2): 259–68. <https://doi.org/10.1099/jmm.0.048785-0>.
- Pongrácz, Júlia, Emese Juhász, Miklós Iván, and Katalin Kristóf. 2015. 'Significance of Yeasts in Bloodstream Infection: Epidemiology and Predisposing Factors of Candidaemia in Adult Patients at a University Hospital (2010-2014)'. *Acta Microbiologica et Immunologica Hungarica* 62 (3): 317–29. <https://doi.org/10.1556/030.62.2015.3.9>.
- Khadka, Sundar, Jeevan Bahadur Sherchand, Bharat Mani Pokhrel, Keshab Parajuli, Shyam Kumar Mishra, Sangita Sharma, Niranjana Shah, et al. 2017. 'Isolation, Speciation and Antifungal Susceptibility Testing of Candida Isolates from Various Clinical Specimens at a Tertiary Care Hospital, Nepal'. *BMC Research Notes* 10 (1): 1–5. <https://doi.org/10.1186/s13104-017-2547-3>.
- Xiang, Guo Shi, Ling Bing Zeng, Jie Yu Zhang, Ying Ying, Xi Ren Deng, Xue Fei Hu, Yan Yan, Ke Hua Yu, Wei Wen Zou, and Xiao Tian Huang. 2018. 'In Vitro Antidrug Susceptibility Testing of Candida Species Isolated from Aseptic Body Fluids'. *Jundishapur Journal of Microbiology* 11 (8). <https://doi.org/10.5812/ijm.55547>.
- Xiao, Meng, Xin Fan, Xin Hou, Sharon C.A. Chen, He Wang, Fanrong Kong, Zi Yong Sun, Yun Zhuo Chu, and Ying Chun Xu. 2018. 'Clinical Characteristics of the First Cases of Invasive Candidiasis in China Due to Pan-Echinocandin-Resistant *Candida Tropicalis* and *Candida Glabrata* Isolates with Delineation of Their Resistance Mechanisms'. *Infection and Drug Resistance* 11: 155–61. <https://doi.org/10.2147/IDR.S152785>.
- Données d'archive.
- Scharmann, Ulrike, Lisa Kirchoff, Valerie le Saout Chapot, Jan Dziobaka, Hedda Luise Verhasselt, Raphael Stauf, Jan Buer, Joerg Steinmann, and Peter Michael Rath. 2020. 'Comparison of Four Commercially Available Chromogenic Media to Identify *Candida Albicans* and Other Medically Relevant *Candida* Species'. *Mycoses* 63: 823–31. <https://doi.org/10.1111/myc.13119>.

21. Angelis, Giulia De, Giulia Menchinelli, Riccardo Torelli, Elena De Carolis, Patrizia Posteraro, Maurizio Sanguinetti, and Brunella Posteraro. 2020. 'Different Detection Capabilities by Mycological Media for Candida Isolates from Mono- or Dual-Species Cultures'. PLoS ONE 15 (3): 1–9. <https://doi.org/10.1371/journal.pone.0226467>.
22. Vecchione, Alessandra, Walter Florio, Francesco Celandroni, Simona Barnini, Antonella Lupetti, et Emilia Ghelardi. 2017. 'Comparative Evaluation of Six Chromogenic Media for Presumptive Yeast Identification'. Journal of Clinical Pathology 70 (12): 1074–78. <https://doi.org/10.1136/jclinpath-2017-204396>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Royaume-Uni

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Bibliografija

1. Public Health England. 2015a. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
2. Müller, F. M.C., A. H. Groll, and T. J. Walsh. 1999. 'Current Approaches to Diagnosis and Treatment of Fungal Infections in Children Infected with Human Immuno Deficiency Virus'. European Journal of Pediatrics 158 (3): 187–99. <https://doi.org/10.1007/s004310051047>.
3. Public Health England. 2013. 'Sexually Transmitted Infections'. UK Standards for Microbiology Investigations S 6 (1.2). <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>.
4. Public Health England. 2015b. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
5. Public Health England. 2016. 'Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
6. Public Health England. 2017a. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
7. Public Health England. 2017b. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
8. Public Health England. 2017c. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
9. Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
10. Public Health England. 2019a. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
11. Public Health England. 2019b. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
12. Sganga, Gabriele. 2009. 'Clinical Aspects of Invasive Candidiasis in the Surgical Patient'. Drugs 69: 29–32. <https://doi.org/10.2165/11315600-000000000-00000>.
13. Chen, Sharon C.A., D. Marriott, E. G. Playford, Q. Nguyen, D. Ellis, W. Meyer, T. C. Sorrell, et al. 2009. 'Candidaemia with Uncommon Candida Species: Predisposing Factors, Outcome, Antifungal Susceptibility, and Implications for Management'. Clinical Microbiology and Infection 15 (7): 662–69. <https://doi.org/10.1111/j.1469-0691.2009.02821.x>.
14. Farooqi, J. Q., K. Jabeen, N. Saeed, N. Iqbal, B. Malik, S. R. Lockhart, A. Zafar, M. E. Brandt, and R. Hasan. 2013. 'Invasive Candidiasis in Pakistan: Clinical Characteristics, Species Distribution and Antifungal Susceptibility'. Journal of Medical Microbiology 62 (2): 259–68. <https://doi.org/10.1099/jmm.0.048785-0>.
15. Pongrácz, Júlia, Emese Juhász, Miklós Iván, and Katalin Kristóf. 2015. 'Significance of Yeasts in Bloodstream Infection: Epidemiology and Predisposing Factors of Candidaemia in Adult Patients at a University Hospital (2010-2014)'. Acta Microbiologica et Immunologica Hungarica 62 (3): 317–29. <https://doi.org/10.1556/030.62.2015.3.9>.
16. Khadka, Sundar, Jeevan Bahadur Sherchand, Bharat Mani Pokhrel, Keshab Parajuli, Shyam Kumar Mishra, Sangita Sharma, Niranjana Shah, et al. 2017. 'Isolation, Speciation and Antifungal Susceptibility Testing of Candida Isolates from Various Clinical Specimens at a Tertiary Care Hospital, Nepal'. BMC Research Notes 10 (1): 1–5. <https://doi.org/10.1186/s13104-017-2547-3>.
17. Xiang, Guo Shi, Ling Bing Zeng, Jie Yu Zhang, Ying Ying, Xi Ren Deng, Xue Fei Hu, Yan Yan, Ke Hua Yu, Wei Wen Zou, and Xiao Tian Huang. 2018. 'In Vitro Antidrug Susceptibility Testing of Candida Species Isolated from Aseptic Body Fluids'. Jundishapur Journal of Microbiology 11 (8). <https://doi.org/10.5812/ijm.55547>.
18. Xiao, Meng, Xin Fan, Xin Hou, Sharon C.A. Chen, He Wang, Fanrong Kong, Zi Yong Sun, Yun Zhuo Chu, and Ying Chun Xu. 2018. 'Clinical Characteristics of the First Cases of Invasive Candidiasis in China Due to Pan-Echinocandin-Resistant Candida Tropicalis and Candida Glabrata Isolates with Delineation of Their Resistance Mechanisms'. Infection and Drug Resistance 11: 155–61. <https://doi.org/10.2147/IDR.S152785>.
19. Data held on file.
20. Scharmann, Ulrike, Lisa Kirchoff, Valerie le Saout Chapot, Jan Dziobaka, Hedda Luise Verhasselt, Raphael Stauf, Jan Buer, Joerg Steinmann, and Peter Michael Rath. 2020. 'Comparison of Four Commercially Available Chromogenic Media to Identify Candida Albicans and Other Medically Relevant Candida Species'. Mycoses 63: 823–31. <https://doi.org/10.1111/myc.13119>.
21. Angelis, Giulia De, Giulia Menchinelli, Riccardo Torelli, Elena De Carolis, Patrizia Posteraro, Maurizio Sanguinetti, and Brunella Posteraro. 2020. 'Different Detection Capabilities by Mycological Media for Candida Isolates from Mono- or Dual-Species Cultures'. PLoS ONE 15 (3): 1–9. <https://doi.org/10.1371/journal.pone.0226467>.
22. Vecchione, Alessandra, Walter Florio, Francesco Celandroni, Simona Barnini, Antonella Lupetti, and Emilia Ghelardi. 2017. 'Comparative Evaluation of Six Chromogenic Media for Presumptive Yeast Identification'. Journal of Clinical Pathology 70 (12): 1074–78. <https://doi.org/10.1136/jclinpath-2017-204396>.

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
Bibliografia

- Public Health England. 2015a. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
- Müller, F. M.C., A. H. Groll and T.J. Walsh. 1999. 'Current Approaches to Diagnosis and Treatment of Fungal Infections in Children Infected with Human Immuno Deficiency Virus'. *European Journal of Pediatrics* 158 (3): 187–99. <https://doi.org/10.1007/s004310051047>.
- Public Health England. 2013. 'Sexually Transmitted Infections'. UK Standards for Microbiology Investigations S 6 (1.2). <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>.
- Public Health England. 2015b. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
- Public Health England. 2016. 'Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
- Public Health England. 2017a. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
- Public Health England. 2017b. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
- Public Health England. 2017c. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
- Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
- Public Health England. 2019a. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
- Public Health England. 2019b. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
- Sganga, Gabriele. 2009. 'Clinical Aspects of Invasive Candidiasis in the Surgical Patient'. *Drugs* 69: 29–32. <https://doi.org/10.2165/11315600-000000000-00000>.
- Chen, Sharon C.A., D. Marriott, E. G. Playford, Q. Nguyen, D. Ellis, W. Meyer, T.C. Sorrell, et al. 2009. 'Candidaemia with Uncommon Candida Species: Predisposing Factors, Outcome, Antifungal Susceptibility, and Implications for Management'. *Clinical Microbiology and Infection* 15 (7): 662–69. <https://doi.org/10.1111/j.1469-0691.2009.02821.x>.
- Farooqi, J. Q., K. Jabeen, N. Saeed, N. Iqbal, B. Malik, S. R. Lockhart, A. Zafar, M. E. Brandt and R. Hasan. 2013. 'Invasive Candidiasis in Pakistan: Clinical Characteristics, Species Distribution and Antifungal Susceptibility'. *Journal of Medical Microbiology* 62 (2): 259–68. <https://doi.org/10.1099/jmm.0.048785-0>.
- Pongrácz, Júlia, Emese Juhász, Miklós Iván and Katalin Kristóf. 2015. 'Significance of Yeasts in Bloodstream Infection: Epidemiology and Predisposing Factors of Candidaemia in Adult Patients at a University Hospital (2010-2014)'. *Acta Microbiologica et Immunologica Hungarica* 62 (3): 317–29. <https://doi.org/10.1556/030.62.2015.3.9>.
- Khadka, Sundar, Jeevan Bahadur Sherchand, Bharat Mani Pokhrel, Keshab Parajuli, Shyam Kumar Mishra, Sangita Sharma, Niranjana Shah, et al. 2017. 'Isolation, Speciation and Antifungal Susceptibility Testing of Candida Isolates from Various Clinical Specimens at a Tertiary Care Hospital, Nepal'. *BMC Research Notes* 10 (1): 1–5. <https://doi.org/10.1186/s13104-017-2547-3>.
- Xiang, Guo Shi, Ling Bing Zeng, Jie Yu Zhang, Ying Ying, Xi Ren Deng, Xue Fei Hu, Yan Yan, Ke Hua Yu, Wei Wen Zou and Xiao Tian Huang. 2018. 'In Vitro Antidrug Susceptibility Testing of Candida Species Isolated from Aseptic Body Fluids'. *Jundishapur Journal of Microbiology* 11 (8). <https://doi.org/10.5812/jjm.55547>.
- Xiao, Meng, Xin Fan, Xin Hou, Sharon C. A. Chen, He Wang, Fanrong Kong, Zi Yong Sun, Yun Zhuo Chu and Ying Chun Xu. 2018. 'Clinical Characteristics of the First Cases of Invasive Candidiasis in China Due to Pan-Echinocandin-Resistant Candida Tropicalis and Candida Glabrata Isolates with Delineation of Their Resistance Mechanisms'. *Infection and Drug Resistance* 11: 155–61. <https://doi.org/10.2147/IDR.S152785>.
- Dati di archivio.
- Scharmann, Ulrike, Lisa Kirchoff, Valerie le Saout Chapot, Jan Dziobaka, Hedda Luise Verhasselt, Raphael Stauf, Jan Buer, Joerg Steinmann and Peter Michael Rath. 2020. 'Comparison of Four Commercially Available Chromogenic Media to Identify Candida Albicans and Other Medically Relevant Candida Species'. *Mycoses* 63: 823–31. <https://doi.org/10.1111/myc.13119>.
- Angelis, Giulia De, Giulia Menchinelli, Riccardo Torelli, Elena De Carolis, Patrizia Posteraro, Maurizio Sanguinetti and Brunella Posteraro. 2020. 'Different Detection Capabilities by Mycological Media for Candida Isolates from Mono- or Dual-Species Cultures'. *PLoS ONE* 15 (3): 1–9. <https://doi.org/10.1371/journal.pone.0226467>.

22. Vecchione, Alessandra, Walter Florio, Francesco Celandroni, Simona Barnini, Antonella Lupetti and Emilia Ghelardi. 2017. 'Comparative Evaluation of Six Chromogenic Media for Presumptive Yeast Identification'. Journal of Clinical Pathology 70 (12): 1074–78. <https://doi.org/10.1136/jclinpath-2017-204396>.

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
Bibliografia

- Public Health England. 2015a. „Investigation of Nasal Samples”. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
- Müller, F. M.C., A. H. Groll, and T. J. Walsh. 1999. „Current Approaches to Diagnosis and Treatment of Fungal Infections in Children Infected with Human Immuno Deficiency Virus”. *European Journal of Pediatrics* 158 (3): 187–99. <https://doi.org/10.1007/s004310051047>.
- Public Health England. 2013. „Sexually Transmitted Infections”. UK Standards for Microbiology Investigations S 6 (1.2). <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>.
- Public Health England. 2015b. „Investigation of Throat Related Specimens”. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
- Public Health England. 2016. „Investigation of Dermatological Specimens for Superficial Mycoses”. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
- Public Health England. 2017a. „Genital Tract and Associated Specimens”. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
- Public Health England. 2017b. „Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture”. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
- Public Health England. 2017c. „Investigation of Intravascular Cannulae and Associated Specimens”. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
- Public Health England. 2018. „Swabs from Skin and Superficial Soft Tissue Infections”. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
- Public Health England. 2019a. „Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens”. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
- Public Health England. 2019b. „Investigation of Urine”. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
- Sganga, Gabriele. 2009. „Clinical Aspects of Invasive Candidiasis in the Surgical Patient”. *Drugs* 69: 29–32. <https://doi.org/10.2165/11315600-000000000-00000>.
- Chen, Sharon C.A., D. Marriott, E. G. Playford, G. Nguyen, D. Ellis, W. Meyer, T. C. Sorrell, et al. 2009. „Candidaemia with Uncommon Candida Species: Predisposing Factors, Outcome, Antifungal Susceptibility, and Implications for Management”. *Clinical Microbiology and Infection* 15 (7): 662–69. <https://doi.org/10.1111/j.1469-0691.2009.02821.x>.
- Farooqi, J. Q., K. Jabeen, N. Saeed, N. Iqbal, B. Malik, S. R. Lockhart, A. Zafar, M. E. Brandt, and R. Hasan. 2013. „Invasive Candidiasis in Pakistan: Clinical Characteristics, Species Distribution and Antifungal Susceptibility”. *Journal of Medical Microbiology* 62 (2): 259–68. <https://doi.org/10.1099/jmm.0.048785-0>.
- Pongrácz, Júlia, Emese Juhász, Miklós Iván, and Katalin Kristóf. 2015. „Significance of Yeasts in Bloodstream Infection: Epidemiology and Predisposing Factors of Candidaemia in Adult Patients at a University Hospital (2010-2014)”. *Acta Microbiologica et Immunologica Hungarica* 62 (3): 317–29. <https://doi.org/10.1556/030.62.2015.3.9>.
- Khadka, Sundar, Jeevan Bahadur Sherchand, Bharat Mani Pokhrel, Keshab Parajuli, Shyam Kumar Mishra, Sangita Sharma, Niranjana Shah, et al. 2017. „Isolation, Speciation and Antifungal Susceptibility Testing of Candida Isolates from Various Clinical Specimens at a Tertiary Care Hospital, Nepal”. *BMC Research Notes* 10 (1): 1–5. <https://doi.org/10.1186/s13104-017-2547-3>.
- Xiang, Guo Shi, Ling Bing Zeng, Jie Yu Zhang, Ying Ying, Xi Ren Deng, Xue Fei Hu, Yan Yan, Ke Hua Yu, Wei Wen Zou, and Xiao Tian Huang. 2018. „In Vitro Antidrug Susceptibility Testing of Candida Species Isolated from Aseptic Body Fluids”. *Jundishapur Journal of Microbiology* 11 (8). <https://doi.org/10.5812/jjm.55547>.
- Xiao, Meng, Xin Fan, Xin Hou, Sharon C.A. Chen, He Wang, Fanrong Kong, Zi Yong Sun, Yun Zhuo Chu, and Ying Chun Xu. 2018. „Clinical Characteristics of the First Cases of Invasive Candidiasis in China Due to Pan-Echinocandin-Resistant *Candida Tropicalis* and *Candida Glabrata* Isolates with Delineation of Their Resistance Mechanisms”. *Infection and Drug Resistance* 11: 155–61. <https://doi.org/10.2147/IDR.S152785>.

19. Dane w posiadaniu firmy.
20. Scharmann, Ulrike, Lisa Kirchoff, Valerie le Saout Chapot, Jan Dziobaka, Hedda Luise Verhasselt, Raphael Stauf, Jan Buer, Joerg Steinmann, and Peter Michael Rath. 2020. „Comparison of Four Commercially Available Chromogenic Media to Identify Candida Albicans and Other Medically Relevant Candida Species”. *Mycoses* 63: 823–31. <https://doi.org/10.1111/myc.13119>.
21. Angelis, Giulia De, Giulia Menchinelli, Riccardo Torelli, Elena De Carolis, Patrizia Posteraro, Maurizio Sanguinetti, and Brunella Posteraro. 2020. „Different Detection Capabilities by Mycological Media for Candida Isolates from Mono- or Dual-Species Cultures”. *PLoS ONE* 15 (3): 1–9. <https://doi.org/10.1371/journal.pone.0226467>.
22. Vecchione, Alessandra, Walter Florio, Francesco Celandroni, Simona Barnini, Antonella Lupetti, and Emilia Ghelardi. 2017. „Comparative Evaluation of Six Chromogenic Media for Presumptive Yeast Identification”. *Journal of Clinical Pathology* 70 (12): 1074–78. <https://doi.org/10.1136/jclinpath-2017-204396>.

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Bibliografia

- Public Health England. 2015a. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
- Müller, F. M.C., A. H. Groll, and T. J. Walsh. 1999. 'Current Approaches to Diagnosis and Treatment of Fungal Infections in Children Infected with Human Immuno Deficiency Virus'. *European Journal of Pediatrics* 158 (3): 187–99. <https://doi.org/10.1007/s004310051047>.
- Public Health England. 2013. 'Sexually Transmitted Infections'. UK Standards for Microbiology Investigations S 6 (1.2). <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>.
- Public Health England. 2015b. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
- Public Health England. 2016. 'Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
- Public Health England. 2017a. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
- Public Health England. 2017b. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
- Public Health England. 2017c. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
- Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
- Public Health England. 2019a. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
- Public Health England. 2019b. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
- Sganga, Gabriele. 2009. 'Clinical Aspects of Invasive Candidiasis in the Surgical Patient'. *Drugs* 69: 29–32. <https://doi.org/10.2165/11315600-000000000-00000>.
- Chen, Sharon C.A., D. Marriott, E. G. Playford, Q. Nguyen, D. Ellis, W. Meyer, T. C. Sorrell, et al. 2009. 'Candidaemia with Uncommon Candida Species: Predisposing Factors, Outcome, Antifungal Susceptibility, and Implications for Management'. *Clinical Microbiology and Infection* 15 (7): 662–69. <https://doi.org/10.1111/j.1469-0691.2009.02821.x>.
- Farooqi, J. Q., K. Jabeen, N. Saeed, N. Iqbal, B. Malik, S. R. Lockhart, A. Zafar, M. E. Brandt, and R. Hasan. 2013. 'Invasive Candidiasis in Pakistan: Clinical Characteristics, Species Distribution and Antifungal Susceptibility'. *Journal of Medical Microbiology* 62 (2): 259–68. <https://doi.org/10.1099/jmm.0.048785-0>.
- Pongrácz, Júlia, Emese Juhász, Miklós Iván, and Katalin Kristóf. 2015. 'Significance of Yeasts in Bloodstream Infection: Epidemiology and Predisposing Factors of Candidaemia in Adult Patients at a University Hospital (2010-2014)'. *Acta Microbiologica et Immunologica Hungarica* 62 (3): 317–29. <https://doi.org/10.1556/030.62.2015.3.9>.
- Khadka, Sundar, Jeevan Bahadur Sherchand, Bharat Mani Pokhrel, Keshab Parajuli, Shyam Kumar Mishra, Sangita Sharma, Niranjana Shah, et al. 2017. 'Isolation, Speciation and Antifungal Susceptibility Testing of Candida Isolates from Various Clinical Specimens at a Tertiary Care Hospital, Nepal'. *BMC Research Notes* 10 (1): 1–5. <https://doi.org/10.1186/s13104-017-2547-3>.
- Xiang, Guo Shi, Ling Bing Zeng, Jie Yu Zhang, Ying Ying, Xi Ren Deng, Xue Fei Hu, Yan Yan, Ke Hua Yu, Wei Wen Zou, and Xiao Tian Huang. 2018. 'In Vitro Antidrug Susceptibility Testing of Candida Species Isolated from Aseptic Body Fluids'. *Jundishapur Journal of Microbiology* 11 (8). <https://doi.org/10.5812/jjm.55547>.
- Xiao, Meng, Xin Fan, Xin Hou, Sharon C.A. Chen, He Wang, Fanrong Kong, Zi Yong Sun, Yun Zhuo Chu, and Ying Chun Xu. 2018. 'Clinical Characteristics of the First Cases of Invasive Candidiasis in China Due to Pan-Echinocandin-Resistant Candida Tropicalis and Candida Glabrata Isolates with Delineation of Their Resistance Mechanisms'. *Infection and Drug Resistance* 11: 155–61. <https://doi.org/10.2147/IDR.S152785>.

19. Data held on file. (Dados conservados em arquivo.)
20. Scharmann, Ulrike, Lisa Kirchoff, Valerie le Saout Chapot, Jan Dziobaka, Hedda Luise Verhasselt, Raphael Stauf, Jan Buer, Joerg Steinmann, and Peter Michael Rath. 2020. 'Comparison of Four Commercially Available Chromogenic Media to Identify Candida Albicans and Other Medically Relevant Candida Species'. *Mycoses* 63: 823–31. <https://doi.org/10.1111/myc.13119>.
21. Angelis, Giulia De, Giulia Menchinelli, Riccardo Torelli, Elena De Carolis, Patrizia Posteraro, Maurizio Sanguinetti, and Brunella Posteraro. 2020. 'Different Detection Capabilities by Mycological Media for Candida Isolates from Mono- or Dual-Species Cultures'. *PLoS ONE* 15 (3): 1–9. <https://doi.org/10.1371/journal.pmed.0226467>
22. Vecchione, Alessandra, Walter Florio, Francesco Celandroni, Simona Barnini, Antonella Lupetti, and Emilia Ghelardi. 2017. 'Comparative Evaluation of Six Chromogenic Media for Presumptive Yeast Identification'. *Journal of Clinical Pathology* 70 (12): 1074–78. <https://doi.org/10.1136/jclinpath-2017-204396>.

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
Bibliografia

- Public Health England. 2015a. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
- Müller, F. M.C., A. H. Groll, a T. J. Walsh. 1999. 'Current Approaches to Diagnosis and Treatment of Fungal Infections in Children Infected with Human Immuno Deficiency Virus'. *European Journal of Pediatrics* 158 (3): 187–99. <https://doi.org/10.1007/s004310051047>.
- Public Health England. 2013. 'Sexually Transmitted Infections'. UK Standards for Microbiology Investigations S 6 (1.2). <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>.
- Public Health England. 2015b. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
- Public Health England. 2016. 'Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
- Public Health England. 2017a. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
- Public Health England. 2017b. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
- Public Health England. 2017c. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
- Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
- Public Health England. 2019a. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
- Public Health England. 2019b. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
- Sganga, Gabriele. 2009. 'Clinical Aspects of Invasive Candidiasis in the Surgical Patient'. *Drugs* 69: 29–32. <https://doi.org/10.2165/11315600-000000000-00000>.
- Chen, Sharon C.A., D. Marriott, E. G. Playford, Q. Nguyen, D. Ellis, W. Meyer, T. C. Sorrell, a kol. 2009. 'Candidaemia with Uncommon Candida Species: Predisposing Factors, Outcome, Antifungal Susceptibility, and Implications for Management'. *Clinical Microbiology and Infection* 15 (7): 662–69. <https://doi.org/10.1111/j.1469-0691.2009.02821.x>.
- Farooqi, J. Q., K. Jabeen, N. Saeed, N. Iqbal, B. Malik, S. R. Lockhart, A. Zafar, M. E. Brandt, a R. Hasan. 2013. 'Invasive Candidiasis in Pakistan: Clinical Characteristics, Species Distribution and Antifungal Susceptibility'. *Journal of Medical Microbiology* 62 (2): 259–68. <https://doi.org/10.1099/jmm.0.048785-0>.
- Pongrácz, Júlia, Emese Juhász, Miklós Iván, a Katalin Kristóf. 2015. 'Significance of Yeasts in Bloodstream Infection: Epidemiology and Predisposing Factors of Candidaemia in Adult Patients at a University Hospital (2010-2014)'. *Acta Microbiologica et Immunologica Hungarica* 62 (3): 317–29. <https://doi.org/10.1556/030.62.2015.3.9>.
- Khadka, Sundar, Jeevan Bahadur Sherchand, Bharat Mani Pokhrel, Keshab Parajuli, Shyam Kumar Mishra, Sangita Sharma, Niranjan Shah, a kol. 2017. 'Isolation, Speciation and Antifungal Susceptibility Testing of Candida Isolates from Various Clinical Specimens at a Tertiary Care Hospital, Nepal'. *BMC Research Notes* 10 (1): 1–5. <https://doi.org/10.1186/s13104-017-2547-3>.
- Xiang, Guo Shi, Ling Bing Zeng, Jie Yu Zhang, Ying Ying, Xi Ren Deng, Xue Fei Hu, Yan Yan, Ke Hua Yu, Wei Wen Zou, a Xiao Tian Huang. 2018. 'In Vitro Antidrug Susceptibility Testing of Candida Species Isolated from Aseptic Body Fluids'. *Jundishapur Journal of Microbiology* 11 (8). <https://doi.org/10.5812/jjm.55547>.
- Xiao, Meng, Xin Fan, Xin Hou, Sharon C.A. Chen, He Wang, Fanrong Kong, Zi Yong Sun, Yun Zhuo Chu, a Ying Chun Xu. 2018. 'Clinical Characteristics of the First Cases of Invasive Candidiasis in China Due to Pan-Echinocandin-Resistant Candida Tropicalis and Candida Glabrata Isolates with Delineation of Their Resistance Mechanisms'. *Infection and Drug Resistance* 11: 155–61. <https://doi.org/10.2147/IDR.S152785>.
- Údaje uchovávané v súbore.
- Scharmman, Ulrike, Lisa Kirchhoff, Valerie le Saout Chapot, Jan Dziobaka, Hedda Luise Verhasselt, Raphael Stauf, Jan Buer, Joerg Steinmann, a Peter Michael Rath. 2020. 'Comparison of Four Commercially Available Chromogenic Media to Identify Candida Albicans and Other Medically Relevant Candida Species'. *Mycoses* 63: 823–31. <https://doi.org/10.1111/myc.13119>.
- Angelis, Giulia De, Giulia Menchinelli, Riccardo Torelli, Elena De Carolis, Patrizia Posteraro, Maurizio Sanguinetti, a Brunella Posteraro. 2020. 'Different Detection Capabilities by Mycological Media for Candida Isolates from Mono- or Dual-Species Cultures'. *PLoS ONE* 15 (3): 1–9. <https://doi.org/10.1371/journal.pone.0226467>.
- Vecchione, Alessandra, Walter Florio, Francesco Celandroni, Simona Barnini, Antonella Lupetti, a Emilia Ghelardi. 2017. 'Comparative Evaluation of Six Chromogenic Media for Presumptive Yeast Identification'. *Journal of Clinical Pathology* 70 (12): 1074–78. <https://doi.org/10.1136/iclinpath-2017-204396>.

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Referenser

1. Public Health England. 2015a. 'Investigation of Nasal Samples'. UK Standards for Microbiology Investigations B 5 (7.1). <https://www.gov.uk/government/publications/smi-b-5-investigation-of-nose-swabs>.
2. Müller, F. M.C., A. H. Groll, and T. J. Walsh. 1999. 'Current Approaches to Diagnosis and Treatment of Fungal Infections in Children Infected with Human Immuno Deficiency Virus'. *European Journal of Pediatrics* 158 (3): 187–99. <https://doi.org/10.1007/s004310051047>.
3. Public Health England. 2013. 'Sexually Transmitted Infections'. UK Standards for Microbiology Investigations S 6 (1.2). <https://www.gov.uk/government/publications/smi-s-6-sexually-transmitted-infections>.
4. Public Health England. 2015b. 'Investigation of Throat Related Specimens'. UK Standards for Microbiology Investigations B 9 (9). <https://www.gov.uk/government/publications/smi-b-9-investigation-of-throat-swabs>.
5. Public Health England. 2016. 'Investigation of Dermatological Specimens for Superficial Mycoses'. UK Standards for Microbiology Investigation B 39 (3.1). <https://www.gov.uk/government/publications/smi-b-39-investigation-of-dermatological-specimens-for-superficial-mycoses>.
6. Public Health England. 2017a. 'Genital Tract and Associated Specimens'. UK Standards for Microbiology Investigation B 28 (4.6). <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
7. Public Health England. 2017b. 'Introduction to the Preliminary Identification of Medically Important Bacteria and Fungi from Culture'. UK Standards for Microbiology Investigations ID 1 (2). <https://www.gov.uk/government/publications/smi-id-1-preliminary-identification-of-medically-important-bacteria>.
8. Public Health England. 2017c. 'Investigation of Intravascular Cannulae and Associated Specimens'. UK Standards for Microbiology Investigations B 20 (6.1). <https://www.gov.uk/government/publications/smi-b-20-investigation-of-intravascular-cannulae-and-associated-specimens>.
9. Public Health England. 2018. 'Swabs from Skin and Superficial Soft Tissue Infections'. UK Standards for Microbiology Investigations B 11 (6.5). <https://www.gov.uk/government/publications/smi-b-11-investigation-of-skin-superficial-and-non-surgical-wound-swabs>.
10. Public Health England. 2019a. 'Investigation of Bronchoalveolar Lavage, Sputum and Associated Specimens'. UK Standards for Microbiology Investigations B 57 (3.5). <https://www.gov.uk/government/publications/smi-b-57-investigation-of-bronchoalveolar-lavage-sputum-and-associated-specimens>.
11. Public Health England. 2019b. 'Investigation of Urine'. UK Standards for Microbiology Investigation B 41 (8.7). <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>.
12. Sganga, Gabriele. (2009). 'Clinical Aspects of Invasive Candidiasis in the Surgical Patient'. *Drugs* 69: 29–32. <https://doi.org/10.2165/11315600-000000000-00000>.
13. Chen, Sharon C.A., D. Marriott, E. G. Playford, Q. Nguyen, D. Ellis, W. Meyer, T. C. Sorrell, et al. 2009. 'Candidaemia with Uncommon Candida Species: Predisposing Factors, Outcome, Antifungal Susceptibility, and Implications for Management'. *Clinical Microbiology and Infection* 15 (7): 662–69. <https://doi.org/10.1111/j.1469-0691.2009.02821.x>.
14. Farooqi, J. Q., K. Jabeen, N. Saeed, N. Iqbal, B. Malik, S. R. Lockhart, A. Zafar, M. E. Brandt, and R. Hasan. 2013. 'Invasive Candidiasis in Pakistan: Clinical Characteristics, Species Distribution and Antifungal Susceptibility'. *Journal of Medical Microbiology* 62 (2): 259–68. <https://doi.org/10.1099/jmm.0.048785-0>.
15. Pongrácz, Júlia, Emese Juhász, Miklós Iván, and Katalin Kristóf. 2015. 'Significance of Yeasts in Bloodstream Infection: Epidemiology and Predisposing Factors of Candidaemia in Adult Patients at a University Hospital (2010-2014)'. *Acta Microbiologica et Immunologica Hungarica* 62 (3): 317–29. <https://doi.org/10.1556/030.62.2015.3.9>.
16. Khadka, Sundar, Jeevan Bahadur Sherchand, Bharat Mani Pokhrel, Keshab Parajuli, Shyam Kumar Mishra, Sangita Sharma, Niranjana Shah, et al. 2017. 'Isolation, Speciation and Antifungal Susceptibility Testing of Candida Isolates from Various Clinical Specimens at a Tertiary Care Hospital, Nepal'. *BMC Research Notes* 10 (1): 1–5. <https://doi.org/10.1186/s13104-017-2547-3>.
17. Xiang, Guo Shi, Ling Bing Zeng, Jie Yu Zhang, Ying Ying, Xi Ren Deng, Xue Fei Hu, Yan Yan, Ke Hua Yu, Wei Wen Zou, and Xiao Tian Huang. 2018. 'In Vitro Antidrug Susceptibility Testing of Candida Species Isolated from Aseptic Body Fluids'. *Jundishapur Journal of Microbiology* 11 (8). <https://doi.org/10.5812/ijm.55547>.
18. Xiao, Meng, Xin Fan, Xin Hou, Sharon C.A. Chen, He Wang, Fanrong Kong, Zi Yong Sun, Yun Zhuo Chu, and Ying Chun Xu. 2018. 'Clinical Characteristics of the First Cases of Invasive Candidiasis in China Due to Pan-Echinocandin-Resistant Candida Tropicalis and Candida Glabrata Isolates with Delineation of Their Resistance Mechanisms'. *Infection and Drug Resistance* 11: 155–61. <https://doi.org/10.2147/IDR.S152785>.
19. Data held on file.
20. Scharmann, Ulrike, Lisa Kirchoff, Valerie le Saout Chapot, Jan Dziobaka, Hedda Luise Verhasselt, Raphael Stauf, Jan Buer, Joerg Steinmann, and Peter Michael Rath. 2020. 'Comparison of Four Commercially Available Chromogenic Media to Identify Candida Albicans and Other Medically Relevant Candida Species'. *Mycoses* 63: 823–31. <https://doi.org/10.1111/myc.13119>.
21. Angelis, Giulia De, Giulia Menchinelli, Riccardo Torelli, Elena De Carolis, Patrizia Posteraro, Maurizio Sanguinetti, and Brunella Posteraro. 2020. 'Different Detection Capabilities by Mycological Media for Candida Isolates from Mono- or Dual-Species Cultures'. *PLoS ONE* 15 (3): 1–9. <https://doi.org/10.1371/journal.pone.0226467>.
22. Vecchione, Alessandra, Walter Florio, Francesco Celandroni, Simona Barnini, Antonella Lupetti, and Emilia Ghelardi. 2017. 'Comparative Evaluation of Six Chromogenic Media for Presumptive Yeast Identification'. *Journal of Clinical Pathology* 70 (12): 1074–78. <https://doi.org/10.1136/jclinpath-2017-204396>.

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©2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.
Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Storbritannien