www.thermofisher.com

## Brain Heart Infusion Agar

REF

CM1136B

EN

### Intended Use

Brain Heart Infusion Agar (CM1136B) device is a highly nutritious, general purpose media for the growth of fastidious microorganisms from clinical samples (e.g. wounds, genital swabs).

Brain Heart Infusion Agar (CM1136B), when used with the addition of Chloramphenicol Selective Supplement (SR0078E and SR0078H), is intended to be used for the isolation of pathogenic fungi from clinical samples. Brain Heart Infusion Agar (CM1136B) devices and Brain Heart Infusion Agar (CM1136B) with the addition of Chloramphenicol Selective Supplement (SR0078E and SR0078H), are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having microbial infections.

The device is for professional use only, is not automated, and nor is it a companion diagnostic.

### Summary and Explanation

Fastidious aerobic and anaerobic microorganisms encompass a broad range of microorganisms, which can often be found as part of the normal skin and mucosal flora of humans. Such microorganisms may manifest and cause infection in the tissues of immunocompromised individuals or those undertaking a course of antibiotics. Infections with fastidious aerobic and anaerobic microorganisms can have a range of manifestations ranging from superficial to life-threatening. Prominent Gram-positive pathogens include *S. pyrogenes* and *S. pneumoniae*<sup>1</sup> and Gram-negative pathogens include *N. gonorrhoeae* and *N. meningitidis*<sup>2,3</sup>. These organisms are associated with a range of different infections, with some causing high morbidity and mortality in immunocompromised individuals.

Microorganisms exist in a wide range of ecological niches, including environmental bacteria, commensals, and human pathogens. Prominent Gram-positive human pathogens include *S. aureus* and *E. faecalis*<sup>1,4</sup>. Prominent Gram-negative human pathogens include *P. aeruginosa*<sup>5</sup>. These organisms are associated with a range of different infections, with some causing high morbidity and mortality in immunocompromised individuals.

### Principle of Method

Brain Heart Infusion Agar (CM1136B) device is a highly nutritious general-purpose media for the growth of fastidious microorganisms from clinical samples (e.g. wounds, genital swabs). These devices contain peptones and glucose, which provide the nutritional sources within the media. Sodium chloride maintains the osmotic equilibrium of the media, whilst disodium phosphate is added to maintain the pH. Agar is included to provide a solid growth medium.

### Typical Formula

	<u>grams per litre</u>
Brain infusion solids	12.5
Beef heart infusion solids	5.0
Proteose peptone	10.0
Sodium chloride	5.0
Glucose	2.0
Disodium phosphate	2.5
Agar	10.0

### Materials Provided

CM1136B: 500g of Brain Heart Infusion Agar

### Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dish
- Supplements (SR0050C, SR0051C)

### Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 10°C.

### Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.

- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product ([www.thermofisher.com](http://www.thermofisher.com)).

### Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

### Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B 29.

### Procedure

Suspend 47g in 1 litre of distilled water. Bring to the boil to dissolve completely. Sterilise by autoclaving at 121°C for 15 minutes. Mix well and pour into sterile containers. For blood Agar, cool to 50°C and enrich with 10% v/v sterile defibrinated blood.

### Interpretation

Once the medium is reconstituted as a plain plate, the presence of:

- Pale cream colonies indicate *Staphylococcus aureus* and *Streptococcus pneumoniae*.
- Irregular straw colonies indicate *Pseudomonas aeruginosa*.
- Pale straw colonies indicate *Enterococcus faecalis*, *Streptococcus pyogenes* and *Streptococcus pneumoniae*

Once the medium is reconstituted with 10% v/v horse blood, the presence of:

Pale straw colonies indicates *Streptococcus pyogenes*. Grey/ green colonies indicates *Streptococcus pneumoniae*.

Once the medium is reconstituted and incubated under microaerophilic conditions, the presence of grey/brown colonies indicates *Neisseria gonorrhoeae* or *Neisseria meningitidis*.

### Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 24 hours @ 37°C

<b>Positive Controls</b>	
Inoculum level: 10 – 100 cfu Colony count is ≥ 70% of the control medium count	
<i>Staphylococcus aureus</i> ATCC® 9144™	0.5 – 1 mm pale cream colonies
<i>Pseudomonas aeruginosa</i> ATCC® 27853™	1 – 3 mm irregular straw colonies
<i>Enterococcus faecalis</i> ATCC® 19433™	0.5 – 1 mm pale straw colonies
<i>Streptococcus pyogenes</i> ATCC® 19615™	0.5 – 1 mm pale straw colonies
<i>Streptococcus pneumoniae</i> ATCC® 6303™	1 – 2 mm pale straw colonies
<i>Streptococcus pneumoniae</i> ATCC® 6305™	0.5 – 1 mm pale cream colonies

Enriched with 10% v/v horse blood	
<i>Streptococcus pyogenes</i> ATCC® 19615™	0.25 – 1 mm pale straw colonies, β haemolysis
<i>Streptococcus pneumoniae</i> ATCC® 6303™	2 – 5 mm grey/ green colonies, α haemolysis
<i>Streptococcus pneumoniae</i> ATCC® 6305™	0.5 – 1 mm grey/ green colonies, α haemolysis
Incubation Conditions: 24 hours @ 37°C, microaerophilic conditions	
<i>Neisseria gonorrhoeae</i> ATCC® 19424™	1 – 3 mm grey/brown colonies
<i>Neisseria gonorrhoeae</i> NCTC 11148	1 – 3 mm grey/ brown colonies
<i>Neisseria meningitidis</i> ATCC® 13090™	1 – 3 mm grey/ brown colonies
<i>Neisseria meningitidis</i> ATCC® 13077™	1 – 3 mm grey/ brown colonies

### Limitations

Due to variation in nutritional requirements, some strains of the target organisms may be encountered that grow poorly or fail to grow on this medium. As with all media, atypical organisms may give anomalous reactions which may give a weak reaction and fail to grow, especially when low numbers are present in the sample. Haemolytic reactions will be strain dependent. All identifications are presumptive, therefore further biochemical, and serological tests must be carried out for further identification.

### Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of fastidious strains is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the device. The precision of Brain Heart Infusion Agar was demonstrated by an overall pass rate of 100% (CM1136B) obtained for the product over 13 years of testing (06.03.2009 – 18.03.2022; 96 batches). This shows that the performance is reproducible.



Brain Heart Infusion Agar devices are tested in-house as part of the QC process since the product was launched in 2010. For target organisms *Staphylococcus aureus* ATCC® 9144™, *Pseudomonas aeruginosa* ATCC® 27853™, *Enterococcus faecalis* ATCC® 19433™, *Streptococcus pyogenes* ATCC® 19615™, *Streptococcus pneumoniae* ATCC® 6303™, *Streptococcus pneumoniae* ATCC® 6305™ and incubating the device for 24 hours at 37°C the user can recover organisms with colony size and morphology as listed in this document. When the device is enriched with 10% v/v Horse blood, target organisms *Streptococcus pyogenes* ATCC® 19615™, *Streptococcus pneumoniae* ATCC® 6303™, *Streptococcus pneumoniae* ATCC® 6305™ and incubating the device for 24 hours at 37°C the user can recover organisms with colony size and morphology as listed in this document. For target organisms *Neisseria gonorrhoeae* ATCC® 19424™, *Neisseria gonorrhoeae* NCTC 11148, *Neisseria meningitidis* ATCC® 13090™ and *Neisseria meningitidis* ATCC® 13077™ and incubating the device for 24 hours at 37°C in microaerophilic conditions, the user can recover organisms with colony size and morphology as listed in this document.

### Bibliography

- Public Health England. 2021a. 'Identification of Streptococcus species, Enterococcus species and morphologically similar organisms.' ID 4 (4). <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
- Public Health England. 2014. 'Screening for *Neisseria meningitidis*.' B 51 (2). <https://www.gov.uk/government/publications/smi-b-51-screening-for-neisseria-meningitidis>.
- Public Health England. 2015a. 'Identification of *Neisseria* species.' ID 6 (3). <https://www.gov.uk/government/publications/smi-id-6-identification-of-neisseria-species>.
- Public Health England. 2020. 'Identification of *Staphylococcus* species, *Micrococcus* species and *Rothia* species.' ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
- Public Health England. 2015b. 'Identification of *Pseudomonas* species and other Non-Glucose Fermenters.' ID 17 (3). <https://www.gov.uk/government/publications/smi-id-17-identification-of-glucose-non-fermenting-gram-negative-rods>.

### Symbol Legend

Symbol	Definition
<b>REF</b>	Catalogue number
<b>IVD</b>	In Vitro Diagnostic Medical Device

	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

**Revision information**

Version	Date of modifications introduced
2.0	2023-12-14

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1. Public Health England. 2021a. 'Identification of *Streptococcus* species, *Enterococcus* species and morphologically similar organisms.' ID 4 (4). <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
2. Public Health England. 2014. 'Screening for *Neisseria meningitidis*.' B 51 (2). <https://www.gov.uk/government/publications/smi-b-51-screening-for-neisseria-meningitidis>.
3. Public Health England. 2015a. 'Identification of *Neisseria* species.' ID 6 (3). <https://www.gov.uk/government/publications/smi-id-6-identification-of-neisseria-species>.
4. Public Health England. 2020. 'Identification of *Staphylococcus* species, *Micrococcus* species and *Rothia* species.' ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
5. Public Health England. 2015b. 'Identification of *Pseudomonas* species and other Non-Glucose Fermenters.' ID 17 (3). <https://www.gov.uk/government/publications/smi-id-17-identification-of-glucose-non-fermenting-gram-negative-rods>.

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1. Public Health England. 2021a. 'Identification of Streptococcus species, Enterococcus species and morphologically similar organisms.' ID 4 (4). <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
2. Public Health England. 2014. 'Screening for *Neisseria meningitidis*.' B 51 (2). <https://www.gov.uk/government/publications/smi-b-51-screening-for-neisseria-meningitidis>.
3. Public Health England. 2015a. 'Identification of *Neisseria* species.' ID 6 (3). <https://www.gov.uk/government/publications/smi-id-6-identification-of-neisseria-species>.
4. Public Health England. 2020. 'Identification of *Staphylococcus* species, *Micrococcus* species and *Rothia* species.' ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
5. Public Health England. 2015b. 'Identification of *Pseudomonas* species and other Non-Glucose Fermenters.' ID 17 (3). <https://www.gov.uk/government/publications/smi-id-17-identification-of-glucose-non-fermenting-gram-negative-rods>.

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- Public Health England. 2021a. 'Identification of *Streptococcus* species, *Enterococcus* species and morphologically similar organisms.' ID 4 (4). <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
- Public Health England. 2014. 'Screening for *Neisseria meningitidis*.' B 51 (2). <https://www.gov.uk/government/publications/smi-b-51-screening-for-neisseria-meningitidis>.
- Public Health England. 2015a. 'Identification of *Neisseria* species.' ID 6 (3). <https://www.gov.uk/government/publications/smi-id-6-identification-of-neisseria-species>.
- Public Health England. 2020. 'Identification of *Staphylococcus* species, *Micrococcus* species and *Rothia* species.' ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
- Public Health England. 2015b. 'Identification of *Pseudomonas* species and other Non-Glucose Fermenters.' ID 17 (3). <https://www.gov.uk/government/publications/smi-id-17-identification-of-glucose-non-fermenting-gram-negative-rods>.

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1. Public Health England. 2021a. 'Identification of Streptococcus species, Enterococcus species and morphologically similar organisms.' ID 4 (4). <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
2. Public Health England. 2014. 'Screening for *Neisseria meningitidis*.' B 51 (2). <https://www.gov.uk/government/publications/smi-b-51-screening-for-neisseria-meningitidis>.
3. Public Health England 2015a. 'Identification of *Neisseria* species.' ID 6 (3). <https://www.gov.uk/government/publications/smi-id-6-identification-of-neisseria-species>.
4. Public Health England. 2020. 'Identification of *Staphylococcus* species, *Micrococcus* species and *Rothia* species.' ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
5. Public Health England. 2015b. 'Identification of *Pseudomonas* species and other Non-Glucose Fermenters.' ID 17 (3). <https://www.gov.uk/government/publications/smi-id-17-identification-of-glucose-non-fermenting-gram-negative-rods>.

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### Bibliografia

1. Public Health England. 2021a. "Identification of *Streptococcus* species, *Enterococcus* species and morphologically similar organisms." ID 4 (4). <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
2. Public Health England. 2014. "Screening for *Neisseria meningitidis*." B 51 (2). <https://www.gov.uk/government/publications/smi-b-51-screening-for-neisseria-meningitidis>.
3. Public Health England 2015a. "Identification of *Neisseria* species." ID 6 (3). <https://www.gov.uk/government/publications/smi-id-6-identification-of-neisseria-species>.
4. Public Health England. 2020. "Identification of *Staphylococcus* species, *Micrococcus* species and *Rothia* species." ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
5. Public Health England. 2015b. "Identification of *Pseudomonas* species and other Non-Glucose Fermenters." ID 17 (3). <https://www.gov.uk/government/publications/smi-id-17-identification-of-glucose-non-fermenting-gram-negative-rods>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Regno Unito

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1. Public Health England. 2021a. 'Identification of *Streptococcus* species, *Enterococcus* species and morphologically similar organisms.' ID 4 (4). <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
2. Public Health England. 2014. 'Screening for *Neisseria meningitidis*.' B 51 (2). <https://www.gov.uk/government/publications/smi-b-51-screening-for-neisseria-meningitidis>.
3. Public Health England 2015a. 'Identification of *Neisseria* species.' ID 6 (3). <https://www.gov.uk/government/publications/smi-id-6-identification-of-neisseria-species>.
4. Public Health England. 2020. 'Identification of *Staphylococcus* species, *Micrococcus* species and *Rothia* species.' ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
5. Public Health England. 2015b. 'Identification of *Pseudomonas* species and other Non-Glucose Fermenters.' ID 17 (3). <https://www.gov.uk/government/publications/smi-id-17-identification-of-glucose-non-fermenting-gram-negative-rods>.

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
### Bibliografie

1. Public Health England. 2021a. „Identification of *Streptococcus* species, *Enterococcus* species and morphologically similar organisms.” ID 4 (4). <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
2. Public Health England. 2014. „Screening for *Neisseria meningitidis*.” B 51 (2). <https://www.gov.uk/government/publications/smi-b-51-screening-for-neisseria-meningitidis>.
3. Public Health England. 2015a. „Identification of *Neisseria* species.” ID 6 (3). <https://www.gov.uk/government/publications/smi-id-6-identification-of-neisseria-species>.
4. Public Health England. 2020. „Identification of *Staphylococcus* species, *Micrococcus* species and *Rothia* species.” ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
5. Public Health England. 2015b. „Identification of *Pseudomonas* species and other Non-Glucose Fermenters.” ID 17(3). <https://www.gov.uk/government/publications/smi-id-17-identification-of-pseudomonas-species-and-other-non-glucose-fermenters>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Regatul Unit

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Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.
Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Storbritannien