	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Cefepime Discs 30µg (FEP30)

REF CT0771B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Cefepime Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent cefepime. The discs are labelled on both sides with details of the antimicrobial (FEP) and amount present (µg): FEP30 (30µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Cefepime AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which Cefepime has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100,^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA tables:

Gram-Negative

- Enterobacterales
- *Pseudomonas aeruginosa*

Gram-Positive

- *Streptococcus* spp. β-Hemolytic Group
- *Streptococcus* spp. Viridans Group

Species with published breakpoints according to current CLSI tables:

Gram-Negative

- Enterobacterales
- *Pseudomonas aeruginosa*
- *Acinetobacter* spp.
- *Haemophilus influenzae* and *parainfluenzae*
- *Neisseria gonorrhoeae*

Gram-Positive

- *Streptococcus* spp. β-Hemolytic Group
- *Streptococcus* spp. Viridans Group

Species with published breakpoints according to current EUCAST tables:

Gram-Negative

- Enterobacterales
- *Pseudomonas* spp.
- *Haemophilus influenzae*
- *Moraxella catarrhalis*
- *Aeromonas* spp.

Gram-Positive

- Viridans group *streptococci*

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Cefepime AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology, refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Cefepime AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST.²

Materials Provided

Cefepime AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. They are supplied in cartridges of 50 discs; each pack contains 5 cartridges. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Component Description	Material Description
Cartridge with spring, cap and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Reagent	Description of Function
Cefepime	White to almost white crystalline powder. Cefepime is a fourth-generation cephalosporin antibiotic used in the treatment of various bacterial infections caused by susceptible bacteria, such as pneumonia, urinary tract infections, and skin infections. Cephalosporins are bactericidal and have the same mode of action as other beta-lactam antibiotics (such as penicillin). Cephalosporins disrupt the synthesis of the peptidoglycan layer of bacterial cell walls.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Cefepime AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Cefepime AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Raw data of zone sizes, in accordance with CLSI methodology, for from the 3 most recent batches of Cefepime (30 µg) FEP30 Antimicrobial Susceptibility Test Discs (CT0771B).

CLSI Methodology										
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated Reference Mid-Point Value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3278313	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	23	29	26	26	26	26	0 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC® 25922	MHA	31	37	34	32	32	32	-2 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	25	31	28	28	28	28	0 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49247	HTM	25	31	28	28	28	28	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA+5% SB	28	35	31.5	34	34	34	2.5 ± 0 (CV=0)
	<i>Neisseria gonorrhoeae</i>	ATCC® 49226	GCV	37	46	41.5	41	42	41	-0.17 ± 0.58 (CV=-3.460) (0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated Reference Mid-point Value	Reading (mm)			Average (mm) + SD
							1	2	3	
3266014	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	23	29	26	26	26	26	0 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC® 25922	MHA	31	37	34	32	32	32	-2 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	25	31	28	28	28	28	0 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49247	HTM	25	31	28	28	28	28	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA+5% SB	28	35	31.5	31	31	31	-0.5 ± 0 (CV=0)
	<i>Neisseria gonorrhoeae</i>	ATCC® 49226	GCV	37	46	41.5	40	40	40	-1.5 ± 0 (CV=0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated Reference Mid-point Value	Reading (mm)			Average (mm) + SD
							1	2	3	
3247846	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	23	29	26	27	27	27	1 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC® 25922	MHA	31	37	34	35	35	35	1 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	25	31	28	29	29	29	1 ± 0 (CV=0)

	<i>Haemophilus influenzae</i>	ATCC® 49247	HTM	25	31	28	28	28	28	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA+5% SB	28	35	31.5	33	33	33	1.5 ± 0 (CV=0)
	<i>Neisseria gonorrhoeae</i>	ATCC® 49226	GCV	37	46	41.5	40	40	40	-1.5 ± 0 (CV=0)

Table 4. Raw data of zone sizes, in accordance with EUCAST methodology, taken from the 3 most recent batches Cefepime (30 µg) FEP30 Antimicrobial Susceptibility Test Discs (CT0771B).

EUCAST Methodology										
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated Reference Mid-Point Value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3278313	<i>Escherichia coli</i>	ATCC® 25922	MHA	31	37	34	32	32	32	-2 ± 0 (CV=0.0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	25	31	28	28	28	28	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA-F	31	37	34	34	35	35	0.67 ± 0.58 (CV=0.87)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated Reference Mid-Point Value	Reading (mm)			Average (mm) + SD
3266014	<i>Escherichia coli</i>	ATCC® 25922	MHA	31	37	34	32	32	32	-2 ± 0 (CV=0.0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	25	31	28	28	28	28	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA-F	31	37	34	33	33	33	-1 ± 0 (CV=0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated Reference Mid-Point Value	Reading (mm)			Average (mm) + SD
3247846	<i>Escherichia coli</i>	ATCC® 25922	MHA	31	37	34	35	35	35	1 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	25	31	28	29	29	29	1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA-F	31	37	34	35	35	35	1 ± 0 (CV=0)

Clinical Performance Characteristics

A study was conducted to evaluate the resistance activity of *Pseudomonas aeruginosa* obtained from clinical specimens against different generations of cephalosporin.⁴ One hundred clinical isolates of *Pseudomonas aeruginosa* were collected from different clinical sites including burns, wounds, and ears. They were tested for sensitivity to antibiotics by means of disc diffusion susceptibility tests using the Kirby-Bauer method using prepared antibiotic discs containing different µg of various antibiotics including Cefepime (30µg).

The majority of the isolates (n = 47) were recovered from infected wounds. *Pseudomonas aeruginosa* showed absolute resistance to all antibiotics used except Ceftazidime and Cefepime, which are third and fourth generation cephalosporins, respectively. Cefepime had the highest susceptibility rate of 90%, minimal resistance of 6%, and intermediate of 4%. From the isolates on wound swabs, Cefepime showed a significant difference with a mean value of 1.17. No significant difference was recorded in the mean effect of the second-generation cephalosporin. The third generation, Cefpodoxime, Cefixime and Ceftazidime, showed significant differences with mean values of 2.96 and 15.67 and 9.94, respectively. This study provided evidence that fourth generation is still the most effective of all generations.

A cross-sectional hospital-based study was carried out to determine the prevalence of urinary tract infections (UTIs) and antimicrobial sensitivity patterns of uropathogens.⁵ Patients (20 to 90 years old) had symptomatic UTIs. A total of 500 midstream urine samples were collected and cultured on Blood agar, MacConkey agar, and CLED medium (Oxoid). The significant growth bacteria were processed for identification and antibiotics susceptibility testing was performed by Kirby Bauer's disc diffusion technique. The zone of inhibition was calculated and interpreted using the standard chart and the organisms were reported as susceptible, intermediate, or resistant. Different antibiotic discs (Oxoid) were used in this study including Cefepime (FEP). *E. coli* ATCC 252395 was used as the control strain.

A total of 11.6% urine samples showed significant growth of uropathogenic bacteria. Out of total positive cultures, 20.7% were Gram-positive bacteria and 79.3 were Gram-negative bacteria. The most prevalent uropathogenic bacteria were *E. coli* (41.4%), which showed the highest sensitivity (91.7%) to cefepime. *P. aeruginosa* showed 80% sensitivity to cefepime. Cefepime showed high susceptibility in the clinical analysis of *E. coli* and *P. aeruginosa* with no "out of specification" or inconsistent results reported.

This study aimed to determine the bacterial profile of bloodstream infections and their antibiotic susceptibility pattern. 822 blood culture samples were screened.⁶ Positive cultures were examined, preliminary signs of bacterial growth were detected, and organisms identified. Antimicrobial susceptibility testing was carried out (Kirby-Bauer disc diffusion method /Mueller-Hinton agar/CLSI guidelines). Based on the zone of inhibition obtained, isolates were classified into sensitive, intermediate, and resistant patterns. For each group of organisms, a separate set of antimicrobials were used including Cefepime (30 mcg). 105 (12.77 %) samples cultured positive; 717 (87.22 %) were negative. *E. coli* were 79% sensitive to Cefepime, and *Staphylococcus aureus* were 100% sensitive to Cefepime. This correlates with the sensitivity study done in India (Vasudeva et al. 2016).

There was variation in the antibiotic sensitivity rate of organisms isolated in this study when compared to past studies. This may be because sensitivity of organisms to antibiotics is variable and depends upon prevalence of strains, antibiotics use, and resistance patterns in a particular area. Cefepime showed high susceptibility in the clinical analysis of *E. coli* and 100% susceptibility of *Staphylococcus* spp. with no "out of specification" or inconsistent results reported.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Oladipo, E.K, J.K Oloke, I.O Omomowo, A.O Oyeniran, E.H Awoyelu, and S.O Ogundele. (2015). "Incidence of Cephalosporin Resistance among Clinical Isolates of Pseudomonas Aeruginosa in Ibadan, South-Western Nigeria." *International Journal of Medicine and Biomedical Research* 4 (3): 135–41. <https://doi.org/10.14194/ijmbr.4.3.4>.
- Ullah, Anwar. (2018). "Prevalence of Symptomatic Urinary Tract Infections and Antimicrobial Susceptibility Patterns of Isolated Uropathogens in Kohat Region of Pakistan." *MOJ Biology and Medicine* 3 (2). <https://doi.org/10.15406/mojbm.2018.03.00082>.
- Iqbal, Hasina, Nabeela Mahboob, Mushtaque Ahmed, Kazi Zulfiquer Mamun, Abdur Rahim, and Khan Abul Kalam Azad. (2020). "Bacteriological Profile and Antimicrobial Susceptibility Patterns of Blood Culture Isolates among Febrile Patients in a Tertiary Care Hospital, Bangladesh." *Journal of Dhaka Medical College* 27 (2): 114–22. <https://doi.org/10.3329/jdmc.v27i2.45820>.

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code

	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision Information

Version	Date of issue and modifications introduced
1.0	2023-01-13

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Oladipo, E.K, J.K Oloke, I.O Omomowo, A.O Oyeniran, E.H Awoyelu, and S.O Ogundele. (2015). "Incidence of Cephalosporin Resistance among Clinical Isolates of Pseudomonas Aeruginosa in Ibadan, South-Western Nigeria." *International Journal of Medicine and Biomedical Research* 4 (3): 135–41. <https://doi.org/10.14194/ijmbr.4.3.4>.
- Ullah, Anwar. (2018). "Prevalence of Symptomatic Urinary Tract Infections and Antimicrobial Susceptibility Patterns of Isolated Uropathogens in Kohat Region of Pakistan." *MOJ Biology and Medicine* 3 (2). <https://doi.org/10.15406/mojbm.2018.03.00082>.
- Iqbal, Hasina, Nabeela Mahboob, Mushtaque Ahmed, Kazi Zulfiqar Mamun, Abdur Rahim, and Khan Abul Kalam Azad. (2020). "Bacteriological Profile and Antimicrobial Susceptibility Patterns of Blood Culture Isolates among Febrile Patients in a Tertiary Care Hospital, Bangladesh." *Journal of Dhaka Medical College* 27 (2): 114–22. <https://doi.org/10.3329/jdmc.v27i2.45820>.

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Referencije

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Oladipo, E.K, J.K Oloke, I.O Omomowo, A.O Oyeniran, E.H Awoyelu, and S.O Ogundele. (2015). “Incidence of Cephalosporin Resistance among Clinical Isolates of *Pseudomonas Aeruginosa* in Ibadan, South-Western Nigeria.” *International Journal of Medicine and Biomedical Research* 4 (3): 135–41. <https://doi.org/10.14194/ijmbr.4.3.4>.
- Ullah, Anwar. (2018). “Prevalence of Symptomatic Urinary Tract Infections and Antimicrobial Susceptibility Patterns of Isolated Uropathogens in Kohat Region of Pakistan.” *MOJ Biology and Medicine* 3 (2). <https://doi.org/10.15406/mojbm.2018.03.00082>.
- Iqbal, Hasina, Nabeela Mahboob, Mushtaque Ahmed, Kazi Zulfiquer Mamun, Abdur Rahim, and Khan Abul Kalam Azad. (2020). “Bacteriological Profile and Antimicrobial Susceptibility Patterns of Blood Culture Isolates among Febrile Patients in a Tertiary Care Hospital, Bangladesh.” *Journal of Dhaka Medical College* 27 (2): 114–22. <https://doi.org/10.3329/jdmc.v27i2.45820>.

© 2023. Thermo Fisher Scientific Inc. Sva prava pridržana. ATCC® je zaštitni znak društva ATCC. Svi ostali zaštitni znakovi vlasništvo su društva Thermo Fisher Scientific Inc. i njegovih podružnica. Ove informacije nisu namijenjene za poticanje na uporabu tih proizvoda na način kojim se mogu povrijediti prava intelektualnog vlasništva drugih društava.

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Reference

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Oladipo, E.K, J.K Oloke, I.O Omomowo, A.O Oyeneran, E.H Awoyelu, and S.O Ogundele. (2015). “Incidence of Cephalosporin Resistance among Clinical Isolates of Pseudomonas Aeruginosa in Ibadan, South-Western Nigeria.” *International Journal of Medicine and Biomedical Research* 4 (3): 135–41. <https://doi.org/10.14194/ijmbr.4.3.4>.
- Ullah, Anwar. (2018). “Prevalence of Symptomatic Urinary Tract Infections and Antimicrobial Susceptibility Patterns of Isolated Uropathogens in Kohat Region of Pakistan.” *MOJ Biology and Medicine* 3 (2). <https://doi.org/10.15406/mojbm.2018.03.00082>.
- Iqbal, Hasina, Nabeela Mahboob, Mushtaque Ahmed, Kazi Zulfiqur Mamun, Abdur Rahim, and Khan Abul Kalam Azad. (2020). “Bacteriological Profile and Antimicrobial Susceptibility Patterns of Blood Culture Isolates among Febrile Patients in a Tertiary Care Hospital, Bangladesh.” *Journal of Dhaka Medical College* 27 (2): 114–22. <https://doi.org/10.3329/jdmc.v27i2.45820>.

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Litteratur

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Oladipo, E.K, J.K Oloke, I.O Omomowo, A.O Oyeniran, E.H Awoyelu, and S.O Ogundele. (2015). "Incidence of Cephalosporin Resistance among Clinical Isolates of Pseudomonas Aeruginosa in Ibadan, South-Western Nigeria." *International Journal of Medicine and Biomedical Research* 4 (3): 135–41. <https://doi.org/10.14194/ijmbr.4.3.4>.
5. Ullah, Anwar. (2018). "Prevalence of Symptomatic Urinary Tract Infections and Antimicrobial Susceptibility Patterns of Isolated Uropathogens in Kohat Region of Pakistan." *MOJ Biology and Medicine* 3 (2). <https://doi.org/10.15406/mojbm.2018.03.00082>.
6. Iqbal, Hasina, Nabeela Mahboob, Mushtaque Ahmed, Kazi Zulfiqar Mamun, Abdur Rahim, and Khan Abul Kalam Azad. (2020). "Bacteriological Profile and Antimicrobial Susceptibility Patterns of Blood Culture Isolates among Febrile Patients in a Tertiary Care Hospital, Bangladesh." *Journal of Dhaka Medical College* 27 (2): 114–22. <https://doi.org/10.3329/jdmc.v27i2.45820>.

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Viitteet

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Oladipo, E.K, J.K Oloke, I.O Omomowo, A.O Oyeniran, E.H Awoyelu, and S.O Ogundele. (2015). "Incidence of Cephalosporin Resistance among Clinical Isolates of Pseudomonas Aeruginosa in Ibadan, South-Western Nigeria." *International Journal of Medicine and Biomedical Research* 4 (3): 135–41. <https://doi.org/10.14194/ijmbr.4.3.4>.
5. Ullah, Anwar. (2018). "Prevalence of Symptomatic Urinary Tract Infections and Antimicrobial Susceptibility Patterns of Isolated Uropathogens in Kohat Region of Pakistan." *MOJ Biology and Medicine* 3 (2). <https://doi.org/10.15406/mojbm.2018.03.00082>.
6. Iqbal, Hasina, Nabeela Mahboob, Mushtaque Ahmed, Kazi Zulfiqer Mamun, Abdur Rahim, and Khan Abul Kalam Azad. (2020). "Bacteriological Profile and Antimicrobial Susceptibility Patterns of Blood Culture Isolates among Febrile Patients in a Tertiary Care Hospital, Bangladesh." *Journal of Dhaka Medical College* 27 (2): 114–22. <https://doi.org/10.3329/jdmc.v27i2.45820>.

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Références

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Oladipo, E.K, J.K Oloke, I.O Omomowo, A.O Oyeniran, E.H Awoyelu, and S.O Ogundele. (2015). "Incidence of Cephalosporin Resistance among Clinical Isolates of *Pseudomonas Aeruginosa* in Ibadan, South-Western Nigeria." *International Journal of Medicine and Biomedical Research* 4 (3): 135–41. <https://doi.org/10.14194/ijmbr.4.3.4>.
- Ullah, Anwar. (2018). "Prevalence of Symptomatic Urinary Tract Infections and Antimicrobial Susceptibility Patterns of Isolated Uropathogens in Kohat Region of Pakistan." *MOJ Biology and Medicine* 3 (2). <https://doi.org/10.15406/mojbm.2018.03.00082>.
- Iqbal, Hasina, Nabeela Mahboob, Mushtaque Ahmed, Kazi Zulfiquer Mamun, Abdur Rahim, and Khan Abul Kalam Azad. (2020). "Bacteriological Profile and Antimicrobial Susceptibility Patterns of Blood Culture Isolates among Febrile Patients in a Tertiary Care Hospital, Bangladesh." *Journal of Dhaka Medical College* 27 (2): 114–22. <https://doi.org/10.3329/jdmc.v27i2.45820>.

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- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Oladipo, E.K, J.K Oloke, I.O Omomowo, A.O Oyeniran, E.H Awoyelu, and S.O Ogundele. (2015). “Incidence of Cephalosporin Resistance among Clinical Isolates of Pseudomonas Aeruginosa in Ibadan, South-Western Nigeria.” International Journal of Medicine and Biomedical Research 4 (3): 135–41. <https://doi.org/10.14194/ijmbr.4.3.4>.
- Ullah, Anwar. (2018). “Prevalence of Symptomatic Urinary Tract Infections and Antimicrobial Susceptibility Patterns of Isolated Uropathogens in Kohat Region of Pakistan.” MOJ Biology and Medicine 3 (2). <https://doi.org/10.15406/mojbm.2018.03.00082>.
- Iqbal, Hasina, Nabeela Mahboob, Mushtaq Ahmed, Kazi Zulfiquer Mamun, Abdur Rahim, and Khan Abul Kalam Azad. (2020). “Bacteriological Profile and Antimicrobial Susceptibility Patterns of Blood Culture Isolates among Febrile Patients in a Tertiary Care Hospital, Bangladesh.” Journal of Dhaka Medical College 27 (2): 114–22. <https://doi.org/10.3329/jdmc.v27i2.45820>.

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Bibliografia

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Oladipo, E.K, J.K Oloke, I.O Omomowo, A.O Oyeniran, E.H Awoyelu, and S.O Ogundele. (2015). "Incidence of Cephalosporin Resistance among Clinical Isolates of Pseudomonas Aeruginosa in Ibadan, South-Western Nigeria." *International Journal of Medicine and Biomedical Research* 4 (3): 135–41. <https://doi.org/10.14194/ijmbr.4.3.4>.
5. Ullah, Anwar. (2018). "Prevalence of Symptomatic Urinary Tract Infections and Antimicrobial Susceptibility Patterns of Isolated Uropathogens in Kohat Region of Pakistan." *MOJ Biology and Medicine* 3 (2). <https://doi.org/10.15406/mojbm.2018.03.00082>.
6. Iqbal, Hasina, Nabeela Mahboob, Mushtaque Ahmed, Kazi Zulfiquer Mamun, Abdur Rahim, and Khan Abul Kalam Azad. (2020). "Bacteriological Profile and Antimicrobial Susceptibility Patterns of Blood Culture Isolates among Febrile Patients in a Tertiary Care Hospital, Bangladesh." *Journal of Dhaka Medical College* 27 (2): 114–22. <https://doi.org/10.3329/jdmc.v27i2.45820>.

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Atsauces

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Oladipo, E.K, J.K Oloke, I.O Omomowo, A.O Oyeniran, E.H Awoyelu, and S.O Ogundele. (2015). "Incidence of Cephalosporin Resistance among Clinical Isolates of Pseudomonas Aeruginosa in Ibadan, South-Western Nigeria." *International Journal of Medicine and Biomedical Research* 4 (3): 135–41. <https://doi.org/10.14194/ijmbr.4.3.4>.
5. Ullah, Anwar. (2018). "Prevalence of Symptomatic Urinary Tract Infections and Antimicrobial Susceptibility Patterns of Isolated Uropathogens in Kohat Region of Pakistan." *MOJ Biology and Medicine* 3 (2). <https://doi.org/10.15406/mojbm.2018.03.00082>.
6. Iqbal, Hasina, Nabeela Mahboob, Mushtaque Ahmed, Kazi Zulfiqur Mamun, Abdur Rahim, and Khan Abul Kalam Azad. (2020). "Bacteriological Profile and Antimicrobial Susceptibility Patterns of Blood Culture Isolates among Febrile Patients in a Tertiary Care Hospital, Bangladesh." *Journal of Dhaka Medical College* 27 (2): 114–22. <https://doi.org/10.3329/jdmc.v27i2.45820>.

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1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Oladipo, E.K, J.K Oloke, I.O Omomowo, A.O Oyeniran, E.H Awoyelu, and S.O Ogundele. (2015). "Incidence of Cephalosporin Resistance among Clinical Isolates of *Pseudomonas Aeruginosa* in Ibadan, South-Western Nigeria." *International Journal of Medicine and Biomedical Research* 4 (3): 135–41. <https://doi.org/10.14194/ijmbr.4.3.4>.
5. Ullah, Anwar. (2018). "Prevalence of Symptomatic Urinary Tract Infections and Antimicrobial Susceptibility Patterns of Isolated Uropathogens in Kohat Region of Pakistan." *MOJ Biology and Medicine* 3 (2). <https://doi.org/10.15406/mojbm.2018.03.00082>.
6. Iqbal, Hasina, Nabeela Mahboob, Mushtaque Ahmed, Kazi Zulfiqer Mamun, Abdur Rahim, and Khan Abul Kalam Azad. (2020). "Bacteriological Profile and Antimicrobial Susceptibility Patterns of Blood Culture Isolates among Febrile Patients in a Tertiary Care Hospital, Bangladesh." *Journal of Dhaka Medical College* 27 (2): 114–22. <https://doi.org/10.3329/jdmc.v27i2.45820>.

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Referências

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Oladipo, E.K, J.K Oloke, I.O Omomowo, A.O Oyeniran, E.H Awoyelu, and S.O Ogundele. (2015). “Incidence of Cephalosporin Resistance among Clinical Isolates of Pseudomonas Aeruginosa in Ibadan, South-Western Nigeria.” *International Journal of Medicine and Biomedical Research* 4 (3): 135–41. <https://doi.org/10.14194/ijmbr.4.3.4>.
5. Ullah, Anwar. (2018). “Prevalence of Symptomatic Urinary Tract Infections and Antimicrobial Susceptibility Patterns of Isolated Uropathogens in Kohat Region of Pakistan.” *MOJ Biology and Medicine* 3 (2). <https://doi.org/10.15406/mojbm.2018.03.00082>.
6. Iqbal, Hasina, Nabeela Mahboob, Mushtaque Ahmed, Kazi Zulfiquer Mamun, Abdur Rahim, and Khan Abul Kalam Azad. (2020). “Bacteriological Profile and Antimicrobial Susceptibility Patterns of Blood Culture Isolates among Febrile Patients in a Tertiary Care Hospital, Bangladesh.” *Journal of Dhaka Medical College* 27 (2): 114–22. <https://doi.org/10.3329/jdmc.v27i2.45820>.

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Referințe

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Oladipo, E.K, J.K Oloke, I.O Omomowo, A.O Oyeniran, E.H Awoyele, and S.O Ogundele. (2015). “Incidence of Cephalosporin Resistance among Clinical Isolates of Pseudomonas Aeruginosa in Ibadan, South-Western Nigeria.” *International Journal of Medicine and Biomedical Research* 4 (3): 135–41. <https://doi.org/10.14194/ijmbr.4.3.4>.
- Ullah, Anwar. (2018). “Prevalence of Symptomatic Urinary Tract Infections and Antimicrobial Susceptibility Patterns of Isolated Uropathogens in Kohat Region of Pakistan.” *MOJ Biology and Medicine* 3 (2). <https://doi.org/10.15406/mojbm.2018.03.00082>.
- Iqbal, Hasina, Nabeela Mahboob, Mushtaque Ahmed, Kazi Zulfiquer Mamun, Abdur Rahim, and Khan Abul Kalam Azad. (2020). “Bacteriological Profile and Antimicrobial Susceptibility Patterns of Blood Culture Isolates among Febrile Patients in a Tertiary Care Hospital, Bangladesh.” *Journal of Dhaka Medical College* 27 (2): 114–22. <https://doi.org/10.3329/jdmc.v27i2.45820>.

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Referencie

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Oladipo, E.K, J.K Oloke, I.O Omomowo, A.O Oyeneran, E.H Awoyelu, and S.O Ogundele. (2015). "Incidence of Cephalosporin Resistance among Clinical Isolates of Pseudomonas Aeruginosa in Ibadan, South-Western Nigeria." International Journal of Medicine and Biomedical Research 4 (3): 135–41. <https://doi.org/10.14194/ijmbr.4.3.4>.
5. Ullah, Anwar. (2018). "Prevalence of Symptomatic Urinary Tract Infections and Antimicrobial Susceptibility Patterns of Isolated Uropathogens in Kohat Region of Pakistan." MOJ Biology and Medicine 3 (2). <https://doi.org/10.15406/mojbm.2018.03.00082>.
6. Iqbal, Hasina, Nabeela Mahboob, Mushtaque Ahmed, Kazi Zulfiqar Mamun, Abdur Rahim, and Khan Abul Kalam Azad. (2020). "Bacteriological Profile and Antimicrobial Susceptibility Patterns of Blood Culture Isolates among Febrile Patients in a Tertiary Care Hospital, Bangladesh." Journal of Dhaka Medical College 27 (2): 114–22. <https://doi.org/10.3329/jdmc.v27i2.45820>.

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Referenser

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Oladipo, E.K, J.K Oloke, I.O Omomowo, A.O Oyeniran, E.H Awoyelu, and S.O Ogundele. (2015). "Incidence of Cephalosporin Resistance among Clinical Isolates of Pseudomonas Aeruginosa in Ibadan, South-Western Nigeria." *International Journal of Medicine and Biomedical Research* 4 (3): 135–41. <https://doi.org/10.14194/ijmbr.4.3.4>.
- Ullah, Anwar. (2018). "Prevalence of Symptomatic Urinary Tract Infections and Antimicrobial Susceptibility Patterns of Isolated Uropathogens in Kohat Region of Pakistan." *MOJ Biology and Medicine* 3 (2). <https://doi.org/10.15406/mojbm.2018.03.00082>.
- Iqbal, Hasina, Nabeela Mahboob, Mushtaque Ahmed, Kazi Zulfiqer Mamun, Abdur Rahim, and Khan Abul Kalam Azad. (2020). "Bacteriological Profile and Antimicrobial Susceptibility Patterns of Blood Culture Isolates among Febrile Patients in a Tertiary Care Hospital, Bangladesh." *Journal of Dhaka Medical College* 27 (2): 114–22. <https://doi.org/10.3329/jdmc.v27i2.45820>.

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