Key Code TSMX5916D
www.oxoid.com/ifu

Europe + 800 135 79 135
CA 1 855 805 8539

US 1 855 236 0910
ROW +31 20 794 7071

Oxoid Biochemical Identification System (O.B.I.S.) PYR

REF ID0580M  60

INTENDED USE

O.B.I.S. PYR is a rapid, qualitative, colorimetric test for the detection of PYRase activity in streptococci and *Citrobacter* spp. grown on plated media. The device is used in a diagnostic workflow to aid clinicians in the treatment options for patients suspected of having bacterial infections.

The device is not automated, is for professional use only and is not a companion diagnostic.

PRINCIPLE OF THE TEST

enterococci from other streptococcal groups including Group D streptococci (previously called faecal streptococci)^{1,2}, thus offering a rapid diagnostic alternative to time-consuming culture methods. In addition, the same enzymatic activity may be used to aid differentiation of *Salmonella* spp. from *Citrobacter* spp. and other Enterobacteriaceae^{3,4}.

Some PYRase tests incorporate other chromogenic substrates, such as β -naphthylamides. These deliver similar diagnostic performance, but should be handled and disposed of with care because of their potential toxicity⁵.

The O.B.I.S. PYR test uses Test Cards impregnated with L-pyrogutamic acid 7-amido-4-methylcoumarin and a developing solution containing dimethylaminocinnamaldehyde. The enzymatic hydrolysis of this substrate by enterococci, Group A streptococci and *Citrobacter* spp. produces a purple colour following the addition of the PYR Developing Solution⁶.

COMPONENTS OF THE KIT (ID0580M)

Each O.B.I.S. PYR kit contains the following reagents, sufficient for 60 tests:

- | | |
|----------------|--|
| ID0581M | PYR Test Cards: 1 Pouch containing 10 cards and a moisture absorbent sachet. There are 6 test reaction areas on each card. 60 tests in total.

Each Test Card is impregnated with L-Pyrogutamic acid 7-amido-4-methylcoumarin. |
| ID0200M | PYR Developing Solution – 1 dropper bottle containing 7 ml of 1M hydrochloric acid and 0.5% w/v dimethylaminocinnamaldehyde |
| ID0100M | Buffer – 1 dropper bottle containing 7 ml of PBS (Phosphate Buffered Saline) |

Instruction Leaflet

Materials required but not provided

Plastic microbiological inoculating loops.

PRECAUTIONS

This product is for *in vitro* diagnostic use only.

Do not use the O.B.I.S. PYR reagents beyond the expiry date.

Specimen material may contain pathogenic organisms, handle with appropriate precautions.

The PYR Developing Solution (ID0200M) contains a weak acid.

Avoid direct contact by wearing suitable protective equipment. If the material comes into contact with the skin, mucous membranes or eyes immediately wash the area by rinsing with plenty of water.

Used Test Cards and inoculating loops must be disposed of as biohazardous waste and incinerated or autoclaved for 15 minutes at 121°C.

Please refer to the Safety Data Sheet (SDS) on company website and product labelling for information on potentially hazardous components.

Directions should be read and followed carefully.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction do not use device.

STORAGE AND OPENING

This kit must be stored at 2–8°C. Allow the pouches to reach room temperature before opening to prevent the formation of condensation on the Test Cards.

Open the pouches by cutting with scissors between the seal at the end of the pouch and the zip lock seal.

Once opened, remove the number of Test Cards required for testing within next 10 minutes and reseal the pouch immediately.

If fewer tests are required, cut the Test Card and return the unused portions to the pouch. Do not return used portions to the pouch because they will be contaminated under these conditions.

When stored as indicated, reagents will retain their activity until the expiry date shown on the box.

QUALITY CONTROL PROCEDURE

Each day the kit is used the following control procedures should be performed:

- Positive Control** – Use a known PYRase positive strain such as *Enterococcus faecalis* ATCC™ 29212 (Thermo Scientific™ Culti-Loops™ R4607030), *Streptococcus pyogenes* ATCC™ 19615 (Thermo Scientific™ Culti-Loops™ R4607000) or *Citrobacter freundii* ATCC™ 8090 (Thermo Scientific™ Culti-Loops™ R4601800). Follow the method given in the test procedure. Ensure that a purple colour forms within 20 seconds.
- Negative Control** – Use a known PYRase negative strain such as *Streptococcus agalactiae* ATCC™ 13813 or *Salmonella enteritidis* ATCC™ 13706 (Thermo Scientific™ Culti-Loops™ R4608200). Follow the method given in the test procedure. Ensure that no purple colour forms within 20 seconds.

Do not use the reagents if reactions with control organisms are incorrect.

SPECIMEN COLLECTION AND PREPARATION

For details of specimen collection and treatment a standard textbook should be consulted⁷.

When identifying enterococci and Group A streptococci: Fresh primary or secondary cultures grown overnight on nonselective media such as blood agar give best results. Colonies tested must be Gram-positive cocci and catalase-negative. In case of insufficient growth, a subculture should be performed.

When identifying *Citrobacter* spp. from *Salmonella* spp. and

Enterobacteriaceae:

Colonies from non-selective media or selective media (such as XLD Medium, MLCB Agar, MacConkey Agar, Desoxycholate Citrate Agar, Salmonella Shigella Agar, Brilliant Green Agar or Hektoen Enteric Agar (modified test procedure required – see below)) may be used. Plates should be tested a maximum of one hour after being removed from the incubator.

Colonies should be Gram-negative, oxidase-negative and urease-negative⁷.

TEST PROCEDURE

- Apply one suspect colony (0.5 mm or larger) onto the test area using a plastic loop (enough to make a visible smear).
- Moisten test area with 1 drop of Buffer (ID0100M).
- Incubate the inoculated Test Card at room temperature (15–30°C) for 5 minutes.
- Dispense 1 drop of Developing Solution (ID0200M) onto the test area. Development of a vivid purple colour on and around the colonies within 20 seconds confirms PYRase activity.

READING AND INTERPRETATION OF RESULTS

Positive Result

A positive result is indicated by the development of a vivid purple colour in the inoculated portion of the test area within a 20 second period following addition of Developing Solution.

Negative Result

A negative result is indicated by lack of colour development in the inoculated portion of the test area within the 20 second period following addition of the Developing Solution.

MODIFIED TEST PROCEDURE – HEKTOEN ENTERIC AGAR

- Apply one suspect colony (0.5 mm or larger), to the test area using a plastic loop (enough to make a visible smear).
- Moisten the test area with one drop of Buffer Solution (ID0100M).
- Place the inoculated card into a 75 mm x 100 mm plastic sleeve (available from Oxoid on request).
- Incubate at 37°C $\pm$ 2°C for 30 minutes.
- Dispense one drop of Developing Solution (ID0200M) onto the test area. Development of a vivid purple colour on and around the smeared colony within 20 seconds confirms PYRase activity. (The purple colour may be slightly obscured by the black pigmentation of ferrous sulfide if a black colony is selected.)
- A negative reaction is denoted by no change in colour of the inoculum smear.
- Development of a greenish colour within one minute is considered a negative PYR reaction and is due to the presence of indole. These organisms are not *Salmonella* spp.

Interpretation Chart

EXPECTED RESULTS		
	Reactivity with	
Organism	Lancefield Group	O.B.I.S. Group
<i>Streptococcus pyogenes</i>	A	+
<i>Streptococcus agalactiae</i>	B	–
<i>Streptococcus</i> Group C	C	–
<i>Streptococcus bovis</i>	D	–
<i>Enterococcus faecalis</i>	D	+
<i>Enterococcus faecium</i>	D	+
<i>Streptococcus</i> Group F	F	–
<i>Streptococcus</i> Group G	G	–
<i>Citrobacter</i>	–	+
Salmonellae	–	–
<i>Escherichia coli</i>	–	–

LIMITATIONS OF THE TEST

O.B.I.S. PYR is intended for the detection of PYRase activity in Gram-positive, catalase-negative cocci and *Citrobacter* spp.

Less commonly encountered isolates of lactococci and aerococci may be PYRase positive. The confirmation of enterococci or Group A streptococci can be achieved by serological grouping with a suitable test, e.g. Oxoid Streptococcal Grouping Kit (DR0595A).

Some strains of Enterobacter cloacae are PYRase negative.

The reactions with O.B.I.S. PYR are a marker for enzyme activity and atypical strains may occasionally occur.

A slight colour change may develop with negative reactions. This is restricted to the immediate site of inoculation. Refer to the positive and negative controls to aid interpretation.

Incubation beyond 20 seconds following addition of the PYR Developing Solution may produce non-specific colour reactions. It is important therefore that the test is read as indicated..

REFERENCES

- Ellner, P. D., Williams, D. A., Hosmer, M. E. and Cohenford, M. (1985). Preliminary evaluation of a rapid colorimetric method for the presumptive identification of Group A streptococci and enterococci. J. Clin. Microbiol. 22. 880–881.
- Facklam, R. R., Padula, J. F., Thacker, L. G., Wortham, E. C. and Sconyers, B. J. (1974). Presumptive identification of Group A, B and D streptococci. Appl. Microbiol. 27. 107–109.
- Chagla, A. H., Borczyk, A. A., Aldom, J. E., Dalla Rosa, S. and Cole, D. D. (1993). Evaluation of the L-pyrrolidonyl-beta-naphthylamide hydrolysis test for differentiation of members of the families Enterobacteriaceae and Vibrionaceae. J. Clin. Microbiol. 31. 1946–1948.
- Bennett, A. R., MacPhee, S., Betts, R. and Post, D. (1999). Use of the enzyme pyrrolidonyl peptidase (Oxoid PYR test) to distinguish *Citrobacter* from *Salmonella* isolated from foods. Letters in Applied Microbiology. 28. 175–178.
- Bascomb, S. and Manafi, M. (1988). Use of enzyme tests in the characterisation and identification of aerobic and facultatively anaerobic Gram-Positive cocci. Clinical Microbiology Reviews 11. 318–340.
- Data on file at Oxoid.

7. Carroll K. C. et al. Manual of Clinical Microbiology, 12th Edition, 2019. ASM Press, Washington, D.C. USA.

SYMBOL LEGEND	
	Catalogue Number
	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Contains sufficient for <N> tests
	Do not use if package is damaged
	Do not re-use
	Batch Code (Lot Number)
	Use By (Expiration Date)
	Importer
	Unique Device Identifier
	Authorized representative in the European Community
	UK Conformity Assessed
	European Conformity Assessment
	Manufacturer

IFU X5916D Revised May 2025 Printed in the UK

 Oxoid Ltd
Wade Road, Basingstoke
Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.
©2025 Thermo Fisher Scientific Inc. All rights reserved.
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid™ Biochemical Identification System (O.B.I.S.) PYR

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Bennett, A. R., MacPhee, S., Betts, R. and Post, D. (1999). Use of the enzyme pyrrolidonyl peptidase (Oxid PYR test) to distinguish *Citrobacter* from *Salmonella* isolated from foods. Letters in Applied Microbiology. 28. 175–178.
5.

Bascomb, S. and Manafi, M. (1988). Use of enzyme tests in the characterisation and identification of aerobic and facultatively anaerobic Gram-Positive cocci. Clinical Microbiology Reviews 11. 318–340.
6.

Data on file at Oxoid.
7.

Carroll K. C. et al. Manual of Clinical Microbiology, 12th Edition, 2019. ASM Press, Washington, D.C. USA.

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IFU X5916D květen 2025

 Oxoid Ltd, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

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REFERENZEN

- Ellner, P. D., Williams, D. A., Hosmer, M. E. und Cohenford, M. (1985). Preliminary evaluation of a rapid colorimetric method for the presumptive identification of Group A streptococci and enterococci. J. Clin. Microbiol. 22. 880–881.
- Facklam, R. R., Padula, J. F., Thacker, L. G., Wortham, E. C. und Sconyers, B. J. (1974). Presumptive identification of Group A, B und D streptococci. Appl. Microbiol. 27. 107–109.
- Chagla, A. H., Borczyk, A. A., Aldom, J. E., Dalla Rosa, S. und Cole, D. D. (1993). Evaluation of the L-pyrrolidonyl-beta-naphthylamide hydrolysis test for differentiation of members of the families Enterobacteriaceae und Vibrionaceae. J. Clin. Microbiol. 31. 1946–1948.
- Bennett, A. R., MacPhee, S., Betts, R. und Post, D. (1999). Use of the enzyme pyrrolidonyl peptidase (Oxoid PYR test) to distinguish *Citrobacter* from *Salmonella* isolated from foods. Letters in Applied Microbiology. 28. 175–178.
- Bascomb, S. und Manafi, M. (1988). Use of enzyme tests in the characterisation and identification of aerobic und facultatively anaerobic Gram-Positive cocci. Clinical Microbiology Reviews 11. 318–340.
- Archivdaten bei Oxoid.
- Carroll K. C. et al. Manual of Clinical Microbiology, 12th Edition, 2019. ASM Press, Washington, D.C. USA.

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REFERENCIAS

- Ellner, P. D., Williams, D. A., Hosmer, M. E. and Cohenford, M. (1985). Preliminary evaluation of a rapid colorimetric method for the presumptive identification of Group A streptococci and enterococci. J. Clin. Microbiol. 22. 880–881.
- Facklam, R. R., Padula, J. F., Thacker, L. G., Wortham, E. C. and Sconyers, B. J. (1974). Presumptive identification of Group A, B and D streptococci. Appl. Microbiol. 27. 107–109.
- Chagla, A. H., Borczyk, A. A., Aldom, J. E., Dalla Rosa, S. and Cole, D. D. (1993). Evaluation of the L-pyrrolidonyl-beta-naphthylamide hydrolysis test for differentiation of members of the families Enterobacteriaceae and Vibrionaceae. J. Clin. Microbiol. 31. 1946–1948.
- Bennett, A. R., MacPhee, S., Betts, R. and Post, D. (1999). Use of the enzyme pyrrolidonyl peptidase (Oxoid PYR test) to distinguish *Citrobacter* from *Salmonella* isolated from foods. Letters in Applied Microbiology 28. 175–178.
- Bascomb, S. and Manafi, M. (1988). Use of enzyme tests in the characterisation and identification of aerobic and facultatively anaerobic Gram-Positive cocci. Clinical Microbiology Reviews 11. 318–340.
- Data on file at Oxoid.
- Carroll K. C. et al. Manual of Clinical Microbiology, 12th Edition, 2019. ASM Press, Washington, D.C. EE. UU.

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IFU X5916D Mai 2025

 Oxoid Ltd, Wade Road, Basingstoke, Hampshire, RG24 8PW, Reino Unido

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RÉFÉRENCES

- Ellner, P. D., Williams, D. A., Hosmer, M. E. and Cohenford, M. (1985). Preliminary evaluation of a rapid colorimetric method for the presumptive identification of Group A streptococci and enterococci. J. Clin. Microbiol. 22. 880–881.
- Facklam, R. R., Padula, J. F., Thacker, L. G., Wortham, E. C. and Sconyers, B. J. (1974). Presumptive identification of Group A, B and D streptococci. Appl. Microbiol. 27. 107–109.
- Chagla, A. H., Borczyk, A. A., Aldom, J. E., Dalla Rosa, S. and Cole, D. D. (1993). Evaluation of the L-pyrrolidonyl-beta-naphthylamide hydrolysis test for differentiation of members of the families Enterobacteriaceae and Vibrionaceae. J. Clin. Microbiol. 31. 1946–1948.
- Bennett, A. R., MacPhee, S., Betts, R. and Post, D. (1999). Use of the enzyme pyrrolidonyl peptidase (Oxoid PYR test) to distinguish *Citrobacter* from *Salmonella* isolated from foods. Letters in Applied Microbiology 28. 175–178.
- Bascomb, S. and Manafi, M. (1988). Use of enzyme tests in the characterisation and identification of aerobic and facultatively anaerobic Gram-Positive cocci. Clinical Microbiology Reviews 11. 318–340.
- Data on file at Oxoid.
- Carroll K. C. et al. Manual of Clinical Microbiology, 12th Edition, 2019. ASM Press, Washington, D.C. ÉTATS-UNIS.

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IFU X5916D Mai 2025

 Oxoid Ltd, Wade Road, Basingstoke, Hampshire, RG24 8PW, Royaume-Uni

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- Ellner, P. D., Williams, D. A., Hosmer, M. E. and Cohenford, M. (1985). Preliminary evaluation of a rapid colorimetric method for the presumptive identification of Group A streptococci and enterococci. J. Clin. Microbiol. 22. 880–881.
- Facklam, R. R., Padula, J. F., Thacker, L. G., Wortham, E. C. and Sconyers, B. J. (1974). Presumptive identification of Group A, B and D streptococci. Appl. Microbiol. 27. 107–109.
- Chagla, A. H., Borczyk, A. A., Aldom, J. E., Dalla Rosa, S. and Cole, D. D. (1993). Evaluation of the L-pyrrolidonyl-beta-naphthylamide hydrolysis test for differentiation of members of the families Enterobacteriaceae and Vibrionaceae. J. Clin. Microbiol. 31. 1946–1948.
- Bennett, A. R., MacPhee, S., Betts, R. and Post, D. (1999). Use of the enzyme pyrrolidonyl peptidase (Oxoid PYR test) to distinguish *Citrobacter* from *Salmonella* isolated from foods. Letters in Applied Microbiology. 28. 175–178.
- Bascomb, S. and Manafi, M. (1988). Use of enzyme tests in the characterisation and identification of aerobic and facultatively anaerobic Gram-Positive cocci. Clinical Microbiology Reviews 11. 318–340.
- Data on file at Oxoid.
- Carroll K. C. et al. Manual of Clinical Microbiology, 12th Edition, 2019. ASM Press, Washington, D.C. USA.

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BIBLIOGRAFIA

- Ellner, P. D., Williams, D. A., Hosmer, M. E. and Cohenford, M. (1985). Preliminary evaluation of a rapid colorimetric method for the presumptive identification of Group A streptococci and enterococci. J. Clin. Microbiol. 22. 880–881.
- Facklam, R. R., Padula, J. F., Thacker, L. G., Wortham, E. C. and Sconyers, B. J. (1974). Presumptive identification of Group A, B and D streptococci. Appl. Microbiol. 27. 107–109.
- Chagla, A. H., Borczyk, A. A., Aldom, J. E., Dalla Rosa, S. and Cole, D. D. (1993). Evaluation of the L-pyrrolidonyl-beta-naphthylamide hydrolysis test for differentiation of members of the families Enterobacteriaceae and Vibrionaceae. J. Clin. Microbiol. 31. 1946–1948.
- Bennett, A. R., MacPhee, S., Betts, R. and Post, D. (1999). Use of the enzyme pyrrolidonyl peptidase (Oxoid PYR test) to distinguish *Citrobacter* from *Salmonella* isolated from foods. Letters in Applied Microbiology. 28. 175–178.
- Bascomb, S. and Manafi, M. (1988). Use of enzyme tests in the characterisation and identification of aerobic and facultatively anaerobic Gram-Positive cocci. Clinical Microbiology Reviews 11. 318–340.
- Data in archivio presso Oxoid.
- Carroll K. C. et al. Manual of Clinical Microbiology, 12th Edition, 2019. ASM Press, Washington, D.C. USA.

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REFERÊNCIAS

- Ellner, P. D., Williams, D. A., Hosmer, M. E. and Cohenford, M. (1985). Preliminary evaluation of a rapid colorimetric method for the presumptive identification of Group A streptococci and enterococci. J. Clin. Microbiol. 22. 880–881.
- Facklam, R. R., Padula, J. F., Thacker, L. G., Wortham, E. C. and Sconyers, B. J. (1974). Presumptive identification of Group A, B and D streptococci. Appl. Microbiol. 27. 107–109.
- Chagla, A. H., Borczyk, A. A., Aldom, J. E., Dalla Rosa, S. and Cole, D. D. (1993). Evaluation of the L-pyrrolidonyl-beta-naphthylamide hydrolysis test for differentiation of members of the families Enterobacteriaceae and Vibrionaceae. J. Clin. Microbiol. 31. 1946–1948.
- Bennett, A. R., MacPhee, S., Betts, R. and Post, D. (1999). Use of the enzyme pyrrolidonyl peptidase (Oxoid PYR test) to distinguish *Citrobacter* from *Salmonella* isolated from foods. Letters in Applied Microbiology. 28. 175–178.
- Bascomb, S. and Manafi, M. (1988). Use of enzyme tests in the characterisation and identification of aerobic and facultatively anaerobic Gram-Positive cocci. Clinical Microbiology Reviews 11. 318–340.
- Dados em arquivo na Oxoid.
- Carroll K. C. et al. Manual of Clinical Microbiology, 12th Edition, 2019. ASM Press, Washington, D.C. USA.

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REFERINȚE

- Ellner, P. D., Williams, D. A., Hosmer, M. E. and Cohenford, M. (1985). Preliminary evaluation of a rapid colorimetric method for the presumptive identification of Group A streptococci and enterococci. J. Clin. Microbiol. 22. 880–881.
- Facklam, R. R., Padula, J. F., Thacker, L. G., Wortham, E. C. and Sconyers, B. J. (1974). Presumptive identification of Group A, B and D streptococci. Appl. Microbiol. 27. 107–109.
- Chagla, A. H., Borczyk, A. A., Aldom, J. E., Dalla Rosa, S. and Cole, D. D. (1993). Evaluation of the L-pyrrolidonyl-beta-naphthylamide hydrolysis test for differentiation of members of the families Enterobacteriaceae and Vibrionaceae. J. Clin. Microbiol. 31. 1946–1948.
- Bennett, A. R., MacPhee, S., Betts, R. and Post, D. (1999). Use of the enzyme pyrrolidonyl peptidase (Oxoid PYR test) to distinguish *Citrobacter* from *Salmonella* isolated from foods. Letters in Applied Microbiology. 28. 175–178.
- Bascomb, S. and Manafi, M. (1988). Use of enzyme tests in the characterisation and identification of aerobic and facultatively anaerobic Gram-Positive cocci. Clinical Microbiology Reviews 11. 318–340.
- Date din dosar de la Oxoid.
- Carroll K. C. et al. Manual of Clinical Microbiology, 12th Edition, 2019. ASM Press, Washington, D.C. USA.

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REFERENCIE

- Ellner, P. D., Williams, D. A., Hosmer, M. E. and Cohenford, M. (1985). Preliminary evaluation of a rapid colorimetric method for the presumptive identification of Group A streptococci and enterococci. J. Clin. Microbiol. 22. 880–881.
- Facklam, R. R., Padula, J. F., Thacker, L. G., Wortham, E. C. and Sconyers, B. J. (1974). Presumptive identification of Group A, B and D streptococci. Appl. Microbiol. 27. 107–109.
- Chagla, A. H., Borczyk, A. A., Aldom, J. E., Dalla Rosa, S. and Cole, D. D. (1993). Evaluation of the L-pyrrolidonyl-beta-naphthylamide hydrolysis test for differentiation of members of the families Enterobacteriaceae and Vibrionaceae. J. Clin. Microbiol. 31. 1946–1948.
- Bennett, A. R., MacPhee, S., Betts, R. and Post, D. (1999). Use of the enzyme pyrrolidonyl peptidase (Oxoid PYR test) to distinguish Citrobacter from Salmonella isolated from foods. Letters in Applied Microbiology. 28. 175–178.
- Bascomb, S. and Manafi, M. (1988). Use of enzyme tests in the characterisation and identification of aerobic and facultatively anaerobic Gram-Positive cocci. Clinical Microbiology Reviews 11. 318–340.
- Údaje v súboroch v spoločnosti Oxoid.
- Carroll K. C. et al. Manual of Clinical Microbiology, 12th Edition, 2019. ASM Press, Washington, D.C. USA.

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REFERENSER

- Ellner, P. D., Williams, D. A., Hosmer, M. E. and Cohenford, M. (1985). Preliminary evaluation of a rapid colorimetric method for the presumptive identification of Group A streptococci and enterococci. J. Clin. Microbiol. 22. 880–881.
- Facklam, R. R., Padula, J. F., Thacker, L. G., Wortham, E. C. and Sconyers, B. J. (1974). Presumptive identification of Group A, B and D streptococci. Appl. Microbiol. 27. 107–109.
- Chagla, A. H., Borczyk, A. A., Aldom, J. E., Dalla Rosa, S. and Cole, D. D. (1993). Evaluation of the L-pyrrolidonyl-beta-naphthylamide hydrolysis test for differentiation of members of the families Enterobacteriaceae and Vibrionaceae. J. Clin. Microbiol. 31. 1946–1948.
- Bennett, A. R., MacPhee, S., Betts, R. and Post, D. (1999). Use of the enzyme pyrrolidonyl peptidase (Oxoid PYR test) to distinguish *Citrobacter* from *Salmonella* isolated from foods. Letters in Applied Microbiology. 28. 175–178.
- Bascomb, S. and Manafi, M. (1988). Use of enzyme tests in the characterisation and identification of aerobic and facultatively anaerobic Gram-Positive cocci. Clinical Microbiology Reviews 11. 318–340.
- Arkiverade data hos Oxoid.
- Carroll K. C. et al. Manual of Clinical Microbiology, 12th Edition, 2019. ASM Press, Washington, D.C. USA.

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©2025 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.