Urea Broth

Liquid medium for detection of urea-utilizing organisms from clinical and nonclinical samples.

DESCRIPTION

Urea Broth is a liquid medium used for the differentiation of microorganisms, especially within the Enterobacteriaceae, on the basis of urease production.

This medium is particularly recommended for differentiating members of the genus *Proteus* from those of *Salmonella* and *Shigella* in the diagnosis of enteric infections.

TYPICAL FORMULA

	(g/l)
Yeast Extract	0.1
Urea	20.0
Disodium Phosphate	9.5
Monopotassium Phosphate	9.1
Phenol Red	0.01
Final pH 6.8 ± 0.2 at 25°C	

METHOD PRINCIPLE

Yeast extract provides amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Urea is a source of nitrogen for those organisms producing urease. Sodium phosphate and potassium phosphate are the buffering agents. Phenol red is the pH indicator.

Organisms which possess the urease enzyme hydrolyze urea and produce ammonia which is detected by a change of the medium to pink-red. Anyway, this highly buffered medium usually reacts only to the high outputs of ammonia by *Proteus*, *Morganella* and *Providencia* in the first 24 hours of incubation.

PREPARATION

Dehydrated medium Suspend 38.7 g of the powder in 1 liter of distilled or deionized water. Mix thoroughly until completely dissolved. Filter sterilize. DO NOT BOIL OR AUTOCLAVE THE MEDIUM.

TEST PROCEDURE

Inoculate the medium heavily with a pure culture of the microorganism under investigation. Incubate at 35 ± 2°C for 8-48 hours.

INTERPRETING RESULTS

Turbidity indicates microbial growth. Observe for a pink to red color development indicating a positive reaction. A negative reaction shows no color change.

APPEARANCE

Dehydrated medium: free-flowing, homogeneous, light orange to light pink.
Prepared medium: clear, yellow-orange.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store tubes at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.
Medium in tubes: 1 year.

QUALITY CONTROL

Tubes are inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: ≤ 100 CFU

Incubation conditions: aerobically at $35 \pm 2^\circ\text{C}$ for 8-48 hours.

QC Table.

Microorganism		Growth	Urease Reaction
<i>Proteus mirabilis</i>	ATCC® 25933	Good	Positive
<i>Salmonella</i> Typhimurium	ATCC® 14028	Good	Negative
<i>Escherichia coli</i>	ATCC® 25922	Good	Negative

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

1. MacFaddin, J.F. (2000) Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Lippincott Williams & Wilkins, Philadelphia, PA.
2. W.B. Christensen (1946) Urea decomposition as a means of differentiating *Proteus* and paracolon cultures from each other and from *Salmonella* and *Shigella* types. J. Bacteriol., 52, 461.
3. C.A. Stuart, E. Van Stratum, R. Rustigian (1945) Further studies on urease production by *Proteus* and related organisms. J. Bacteriol., 49, 437.
4. R. Rustigian, C.A. Stuart (1941) Decomposition of urea by *Proteus*. Proc. Soc. Exptl. Biol. Med., 47, 108-112.

PRESENTATION		Contents	Ref.
Urea Broth	Tubes	20 x 10 ml tubes	24416
Urea Broth	Tubes	20 x 3 ml tubes	24516
Urea Broth	Dehydrated medium	500 g of powder	610311
Urea Broth	Dehydrated medium	100 g of powder	620311

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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