1. INTENDED USE

BactiCard™ Neisseria is a rapid, qualitative, enzymatic test for the presumptive identification of *Neisseria* species grown on plated media by the detection of β-galactosidase, indoxyl butyrate esterase, γ-glutamyl aminopeptidase and prolyl aminopeptidase. The device is used in a diagnostic workflow to aid clinicians in the treatment options for patients suspected of having bacterial infections.

The device is not automated, is for professional use only and is not a companion diagnostic.

2. SUMMARY AND EXPLANATION

Chromogenic substrates were utilized by D’Amato et al. in 1978 to assay for enzymes produced by *Neisseria* spp.<sup>1</sup> Distinct enzymatic patterns were observed that have since been employed in various methods for differentiating these organisms.<sup>2-4</sup> In 1986, Sperry et al. reported an increase in the detection of prolylaminopeptidase-positive strains of *Neisseria meningitidis* which they attributed to the greater sensitivity of the filter paper methodology.<sup>5</sup>

3. PRINCIPLE

BactiCard Neisseria utilizes four chromogenic substrates impregnated on individual test circles to detect preformed enzymes associated with pathogenic *Neisseria*. Specific color reactions presumptively identify *Neisseria gonorrhoeae*, *Neisseria meningitidis*, *Neisseria lactamica*, and *Moraxella catarrhalis*. The enzymes β-galactosidase (BGAL) and indoxyl butyrate esterase hydrolyze the substrates 5-bromo-4-chloro-3-indolyl-β-D-galactoside and 5-bromo-4-chloro-3-indolyl butyrate (IB) (respectively) to form a blue-green complex in the presence of oxygen. Gamma-glutamyl aminopeptidase (GLUT) hydrolyzes gamma-glutamyl-β-naphthylamide and prolylaminopeptidase (PRO) hydrolyzes L proline-β-naphthylamide and reacts with the BactiCard Neisseria Color Developer to form a red color.

4. REAGENTS AND MATERIALS SUPPLIED

BactiCard Neisseria Cards

**BactiCard Neisseria Test Cards** (25 – 5 sets of 5)  
Reactive Ingredients:  
5-Bromo-4-chloro-3-indolyl butyrate (IB)  
5-Bromo-4-chloro-3-indolyl-β-D-galactoside (BGAL)  
Gamma-glutamyl-β-naphthylamide (GLUT)  
L-Proline-β-naphthylamide (PRO)

BactiCard Neisseria Rehydrating Fluid

**Rehydrating Fluid** (1 bottle, 6.0 ml) with preservative

BactiCard Neisseria Color Developer

**Color Developer** (1 bottle, 6.0 ml)  
Reactive Ingredient:  
0.02% p-Dimethylaminocinnamaldehyde

Applicator Sticks

**Applicator Sticks** (1 vial)  
**Instructions for Use (IFU)**

5. PRECAUTIONS

This product is for *in vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after use. Directions should be read and followed carefully.

Caution!

1. Color Developer is toxic and harmful by inhalation, contact with skin, or if swallowed. May impair fertility or cause harm to unborn child.

2. Please refer to the Safety Data Sheet (SDS) on company website and product labelling for information on potentially hazardous components.

Directions should be read and followed carefully.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction do not use device.

Contents

Acetic acid 64-19-7  
Dimethyl sulfoxide 67-68-5  
2-Methoxyethanol 109-86-4  
Formamide 75-12-7  
Ethyl alcohol 64-17-5

DANGER

US

US & EU

|           |                                                                                                 |
|-----------|-------------------------------------------------------------------------------------------------|
| H360FD    | May damage fertility. May damage the unborn child                                               |
| H335      | May cause respiratory irritation                                                                |
| H336      | May cause drowsiness or dizziness                                                               |
| H350      | May cause cancer                                                                                |
| H373      | May cause damage to organs through prolonged or repeated exposure                               |
| P201      | Obtain special instructions before use.                                                         |
| P281      | Use personal protective equipment as required.                                                  |
| P202      | Do not handle until all safety precautions have been read and understood                        |
| P260      | Do not breathe dust/fume/gas/mist/vapors/spray                                                  |
| P271      | Use only outdoors or in a well-ventilated area                                                  |
| P304+P340 | IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing |
| 308+P313  | IF exposed or concerned: Get medical advice/attention.                                          |
| P405      | Store locked up                                                                                 |
| P403+P233 | Store in a well-ventilated place. Keep container tightly closed                                 |
| P501      | Dispose of contents/container to an approved waste disposal plant                               |

Emergency Telephone Number

INFOTRAC - 24 Hour Number: 1-800-535-5053

Outside of the United States, call 24 Hour Number: 001-352-323-3500 (Call Collect)

6. STORAGE

This product is ready for use and no further preparation is necessary. Store product in its original container at 2-8°C until used. Do not freeze or overheat. A precipitate formed at storage temperatures in the Color Developer is common and does not affect test performance. Allow product to equilibrate to room temperature before use. Do not incubate prior to use. Once opened, maintain unused test cards at 2-8°C in the foil pouch provided. Seal to protect from moisture.

7. PRODUCT DETERIORATION

This product should not be used if (1) the color has changed, (2) the expiration date has passed, or (3) there are other signs of deterioration.

8. SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected and handled following recommended guidelines.<sup>6</sup>

9. MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Oxidase reagent, (7) Gram stain reagents, (8) Microscope slides, (9) Microscope, and (10) Catalase reagent.

10. PROCEDURE

Test isolates acceptable for testing with BactiCard Neisseria include oxidase-positive, gram-negative diplococci with Gram stain morphology consistent with *N. gonorrhoeae*, *N. lactamica*, and *N. meningitidis*. Remove test isolates from pure, 18-48 hour cultures growing on a selective medium, such as Thayer Martin, Martin Lewis, or New York City Medium. Consult appropriate references for additional information<sup>6</sup>

1. Label BactiCard Neisseria with the specimen identification.
2. Inoculate only one isolate per test card.
3. Add one (1) drop of BactiCard Neisseria Rehydrating Fluid to each test circle.
4. Inoculate each test circle with a visible inoculum of the test isolate using an applicator stick (provided).
5. Incubate at room temperature (RT) for two (2) minutes.
6. Observe the IB test circle for a blue-green color indicating a positive test. If positive, do not proceed beyond this point. Consult EXPECTED VALUES for identification.
7. If the IB test is negative, continue incubation at RT for 13 additional minutes (15 minutes total test time).
8. Observe the BGAL test circle for a blue-green color indicating a positive test. If positive, do not proceed beyond this point. Consult EXPECTED VALUES for identification.
9. If the BGAL test is negative (no blue-green color), proceed by adding one (1) drop of BactiCard Neisseria Color Developer to the GLUT and PRO test circles.
10. Observe GLUT and PRO test circles for development of a pink to red color at 30 seconds. Consult EXPECTED VALUES for identification.

11. INTERPRETATION

IB Test:

Positive Test - Blue-green color development within 2 minutes

Negative Test - No color change within 2 minutes

BGAL Test:

Positive Test - Blue-green color development within 15 minutes

Negative Test - No color change

GLUT Test:

Positive Test - Pink to red color development within 30 seconds after the addition of Color Developer

Negative Test - No color change

PRO Test:

Positive Test - Pink to red color development within 30 seconds after the addition of Color Developer

Negative Test - No color change

12. EXPECTED VALUES

| Organism                                   | IB | BGAL | GLUT | PRO |
|--------------------------------------------|----|------|------|-----|
| <i>Moraxella catarrhalis</i>               | +  | *    | *    | *   |
| <i>Neisseria lactamica</i>                 | –  | +    | *    | *   |
| <i>Neisseria gonorrhoeae</i> <sup>a</sup>  | –  | –    | –    | +   |
| <i>Neisseria meningitidis</i> <sup>b</sup> | –  | –    | +    | V   |

V = Variable, \* = not applicable

<sup>a</sup> PRO-negative strains of *N. gonorrhoeae* have been reported.<sup>7</sup>

<sup>b</sup> GLUT-negative strains of *N. meningitidis* have been reported.<sup>11</sup>

13. QUALITY CONTROL

All lot numbers of BactiCard Neisseria have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory procedures. If aberrant quality control results are noted, patient results should not be reported. Quality control testing should be run with each shipment and new lot number received.

| CONTROL <sup>a</sup>                         | RESULTS |      |      |     |
|----------------------------------------------|---------|------|------|-----|
|                                              | IB      | BGAL | GLUT | PRO |
| <i>Moraxella catarrhalis</i><br>ATCC® 25238  | +       | *    | *    | *   |
| <i>Neisseria lactamica</i><br>ATCC® 23970    | –       | +    | *    | *   |
| <i>Neisseria gonorrhoeae</i><br>ATCC® 43069  | –       | –    | –    | +   |
| <i>Neisseria meningitidis</i><br>ATCC® 13090 | –       | –    | +    | V   |

<sup>a</sup>Test incubated under conditions specified in the PROCEDURE section.

V = Variable, \* = not applicable

14. PERFORMANCE CHARACTERISTICS<sup>13</sup>

BactiCard *Neisseria* was evaluated with 449 isolates which included *N. gonorrhoeae* (217), *N. meningitidis* (98), *N. lactamica* (42), and *M. catarrhalis* (92). Greater than 99% agreement was obtained when compared to conventional biochemical methods. To determine the age of cultures to be used with the test kit, 53 stock isolates including *N. gonorrhoeae* (24), *N. lactamica* (4), *N. meningitidis* (9), and *M. catarrhalis* (16) were evaluated at 24, 48, and 72 hours. Due to weak PRO reactions with *N. gonorrhoeae* and false-negative IB reactions with *M. catarrhalis*, the use of 72 hour cultures is not recommended. Parallel tests using 24 and 48 hour cultures showed comparable results with regards to reaction rate and intensity. Greater than 96% agreement was obtained when compared to conventional biochemical systems.

15. LIMITATIONS

- All isolates from medical legal cases should be confirmed with at least two methodologies.
- Kingella* spp. are gram-negative, coccobacillary organisms which may be isolated on Thayer Martin medium. The superoxol test (30% hydrogen peroxide) can be used to differentiate *N. gonorrhoeae* (positive) and *Kingella* spp. (negative).<sup>8,12</sup>
- Neisseria cinerea*, if recovered on selective media, may be confused with *N. gonorrhoeae* because it can be PRO-positive. Although isolation of *N. cinerea* from a genital site is rare and most strains are inhibited on selective media, such isolates should be differentiated from *N. gonorrhoeae* by subculture to nutrient agar or tryptic soy agar; *N. cinerea* will grow, *N. gonorrhoeae* will not.<sup>6,9</sup>
- Some PRO-positive strains of *Neisseria subflava* biovar *perflava* may grow on selective media. Careful attention to colony morphology (i.e., pigmentation, colony size, and consistency) is necessary to differentiate pathogenic *Neisseria* spp. from nonpathogenic species.<sup>10</sup>
- Isolates that have atypical characteristics or are recovered from extragenital sites require further testing with other methods.<sup>6</sup>
- Saprophytic *Neisseria* spp. may appear on selective media, especially from nonurogenital specimens, after 36-48 hours incubation.<sup>6</sup>
- Organisms recovered from primary isolation on nonselective media, such as chocolate agar, should only be tested after subculture to appropriate selective media. Most *Neisseria* spp. (*N. sicca*, *N. mucosa*, *N. cinerea*) possess the enzyme prolylaminopeptidase and may be confused with *N. gonorrhoeae* if not tested properly. Further biochemical and serological testing may be required.
- Certain auxotypes of gonococci grow slowly, produce smaller than normal colonies on agar plates, and are difficult to identify biochemically. Those strains that require arginine, hypoxanthine, and uracil (AHU-dependent auxotypes) are more likely to cause asymptomatic urethral infections and disseminated gonococcal infections.<sup>6</sup>
- Interpretation of the IB test past the 2-minute incubation period may result in false-positive reactions.
- PRO-negative strains of *N. gonorrhoeae* have been reported. Therefore, gram-negative, oxidase-positive diplococci isolated on selective media that are PRO-negative cannot be confirmed negative for *N. gonorrhoeae* without additional testing.<sup>7</sup>
- GLUT-negative strains of *N. meningitidis* have been reported.<sup>11</sup> PRO-positive strains of *N. meningitidis* that are GLUT-negative may be misidentified as *N. gonorrhoeae* with BactiCard *Neisseria*. Additional testing, such as carbohydrate acidification (i.e., maltose and glucose), is required to definitively identify PRO-positive, GLUT-negative isolates that are otherwise characteristic of *N. meningitidis* or *N. gonorrhoeae*.

16. BIBLIOGRAPHY

1. D'Amato, R.F., L.A. Eriquez, K.M. Tomfohrde, and E. Singerman. 1978. J. Clin. Microbiol. 7:77-81.

2. Yajko, D.M., A. Chu, and W.K. Hadley. 1984. J. Clin. Microbiol. 19:380-382.

3. Philip, A. and G.C. Carton. 1985. J. Clin. Microbiol. 22:101-104.

4. Hosmer, M.E., M.A. Cohenford, and P.D. Ellner. 1986. J. Clin. Microbiol. 24:141-142.

5. Sperry, J.F., M.A. Cohenford, P. Campognone, W. Lawton, and D.O. Chee. 1986. J. Clin. Microbiol. 24:145.

6. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry, and M.A. Pfaller. 2007. Manual of Clinical Microbiology. 9th ed. ASM Press, Washington, D.C.

7. Blackmore, T., G. Hererra, S. Shi, P. Bridgewater, L. Wheeler, and J. Byrne. 2005. J. Clin. Microbiol. 43:4189-4190.

8. Brown, J.D. and K. Thomas. 1985. J. Clin. Microbiol. 21:857-858.

9. Knapp, J.S., P.A. Totten, M.H. Mulks, and B.H. Minshew. 1984. J. Clin. Microbiol. 19:63-67.

10. Janda, W.M. and V. Sobieski. 1988. Eur. J. Clin. Microbiol. Infect. Dis. 7:25-29.

11. Takahashi, H., H. Tanaka, H. Inouy, T. Kuroki, Y. Watanabe, S. Yamai, and H. Watanabe. 2002. J. Clin. Microbiol. 40:3035-3037.

12. Saginur, R., B. Clecner, J. Portnoy, and J. Mendelson. 1982. J. Clin. Microbiol. 15:475-477.

13. Data on file.

17. PACKAGING

 R21110 .....  25

18. SYMBOL LEGEND

|                                                                                     |                                                     |
|-------------------------------------------------------------------------------------|-----------------------------------------------------|
|    | Catalogue Number                                    |
|   | In Vitro Diagnostic Medical Device                  |
|  | Batch Code (Lot Number)                             |
|  | Consult Instructions for Use (IFU)                  |
|  | Temperature Limitations (Storage temp.)             |
|  | Contains sufficient for <N> tests                   |
|  | Do not use if package is damaged                    |
|  | Do not re-use                                       |
|  | Use By (Expiration Date)                            |
|  | Importer                                            |
|  | Authorized representative in the European Community |
|  | UK Conformity Assessed                              |
|  | European Conformity Mark                            |
|  | Manufacturer                                        |

©2022 Thermo Fisher Scientific Inc. All rights reserved.

ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

 12076 Santa Fe Trail Drive, Lenexa, KS 66215, USA

[www.thermofisher.com/microbiology](http://www.thermofisher.com/microbiology)

Tel: (800) 255-6730 • International: (913) 888-0939

[www.oxid.com/IFU](http://www.oxid.com/IFU)

Europe +800 135 79 135 • US 1 855 2360 190

CA 1 855 805 8539 • ROW +31 20 794 7071

| Version  | Date of modifications introduced              |
|----------|-----------------------------------------------|
| IFU21110 | August 2023<br>Correction of hazard pictogram |

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© 2022 Thermo Fisher Scientific Inc. Všechna práva vyhrazena.