Oxoid™ Meropenem Discs 10 µg (MEM10)

REF CT0774B

EN

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the Oxoid™ AST generic instructions for use supplied with the product and available online.

Thermo Scientific™ Oxoid™ Meropenem (10 µg) MEM10 Antimicrobial Susceptibility Test (AST) Discs are 6mm paper discs that contain specific amounts of the antimicrobial agent meropenem. The discs are labelled on both sides with details of the antimicrobial (MEM) and the drug content (10µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack.

Oxoid Meropenem AST discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Oxoid Meropenem AST discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which Meropenem has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate, or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA, CLSI M100,^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- N/A

Gram-Negative

- *Enterobacteriaceae*
- *Pseudomonas aeruginosa*
- *Acinetobacter ssp*
- *Haemophilus influenzae*
- *Haemophilus parainfluenzae*
- *Neisseria meningitidis*

Species with published breakpoints according to current CLSI literature:

Gram-Positive

- N/A

Gram-Negative

- *Enterobacterales*
- *Pseudomonas aeruginosa*
- *Acinetobacter spp*
- *Burkholderia cepacia complex*
- *Haemophilus influenza*
- *Haemophilus parainfluenzae*
- *Neisseria meningitidis*

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- *Listeria monocytogenes*
- *Aerococcus sanguinicola* and *urinae*

Gram-Negative

- *Enterobacterales*
- *Pseudomonas spp*
- *Acinetobacter spp*
- *Haemophilus influenzae*
- *Moraxella catarrhalis*
- *Kingella Kingae*
- *Burkholderia pseudomallei*

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Oxoid Meropenem AST discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in the documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Oxoid Meropenem AST discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Oxoid Meropenem AST discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Meropenem AST discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with Meropenem AST Discs	
Component Description	Material Description
Cartridge with spring, cap and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs
Desiccant tablet (x5)	Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge, 5 cartridges per pack.

Table 2. Description of Meropenem AST Disc Reagent	
Reagent	Description of Function
Meropenem	Antimicrobial agent. White powder. Meropenem is a broad-spectrum carbapenem that inhibits cell wall synthesis in Gram-positive and Gram-negative bacteria. It penetrates cell walls and binds to penicillin-binding protein targets.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Oxoid Meropenem AST discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Oxoid Meropenem AST discs have a shelf-life of 12 months if stored under the recommended conditions. Unopened cartridges must be stored at -20° to 8°C, until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque air-tight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of

condensation. Once opened from their desiccant-containing packaging, discs should be used within 3 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Raw Data of Zone Sizes in Accordance with CLSI Methodology Taken from the Three Most Recent Batches of Oxoid Meropenem AST Discs.

CLSI Methodology										
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated Reference Mid-Point Value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3179516	<i>Escherichia coli</i>	ATCC® 25922	MHA	28	35	31.5	30	30	30	1.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	29	37	33	35	35	35	2 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	27	33	30	31	31	31	1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49247	HTM	20	28	24	22	22	22	2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA+ 5%SB	28	35	31.5	32	32	32	0.5 ± 0%
3179515	<i>Escherichia coli</i>	ATCC® 25922	MHA	28	35	31.5	31	31	31	0.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	29	37	33	33	33	33	0 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	27	33	30	32	32	32	2 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49247	HTM	20	28	24	21	21	21	3 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA+ 5%SB	28	35	31.5	31	31	31	0.5 ± 0 (CV=0)

Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated Reference Mid-Point Value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3174317	<i>Escherichia coli</i>	ATCC® 25922	MHA	28	35	31.5	30	30	30	1.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	29	37	33	36	36	36	3 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	27	33	30	32	32	32	2 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49247	HTM	20	28	24	20	20	20	4 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA+ 5%SB	28	35	31.5	32	32	32	0.5 ± 0 (CV=0)

Table 4. Raw Data of Zone Sizes in Accordance with EUCAST Methodology Taken from the Three Most Recent Batches of Oxoid Meropenem AST Discs.

EUCAST Methodology										
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3179516	<i>Escherichia coli</i>	ATCC® 29922	MHA	28	35	31.5	30	30	30	1.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	27	33	30	31	31	31	1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49766	MHA-F	27	35	31	31	31	31	1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA-F	30	38	34	35	35	36	3 ± 0 (CV=0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated Reference Mid-Point Value	Reading (mm)			Average (mm) + SD
							1	2	3	
3179515	<i>Escherichia coli</i>	ATCC® 29922	MHA	28	35	31.5	31	31	31	1 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	27	33	30	32	32	32	2 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49766	MHA-F	27	35	31	33	33	33	1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA-F	30	38	34	32	32	33	1.67 ± 0 (CV=0.58)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated Reference Mid-Point Value	Reading (mm)			Average (mm) + SD
							1	2	3	
3174317	<i>Escherichia coli</i>	ATCC® 29922	MHA	28	35	31.5	30	30	30	1.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	27	33	30	32	32	32	2 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49766	MHA-F	27	35	31	31	31	31	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA-F	30	38	34	36	36	36	2 ± 0 (CV=0)

Clinical Performance Characteristics

51 clinical isolates were included in AST susceptibility testing of Oxoid discs via disk diffusion according to CLSI and EUCAST breakpoints guidelines.⁴ Oxoid Meropenem AST Discs were shown to be effective in the clinical analysis of *K. pneumoniae*, given that it detected resistance and was able to identify specific resistance genotypes (KPC/VIM production).

Another study was performed to determine bacterial susceptibility of non-fermentative Gram-negative organisms to various antibacterial agents.⁵ Isolates identified from this study included *Pseudomonas aeruginosa* (n=1790), *Acinetobacter* spp. (n=1365), and *Burkholderia* spp. (n=59). This study provided evidence that Oxoid Meropenem AST Discs perform effectively in clinical settings.

A third study investigated antimicrobial resistance and molecular epidemiology of extended-spectrum β-lactamase (ESBL)-producing *Escherichia coli* isolates.⁶ Antimicrobial susceptibility testing was carried out and interpreted using the disc diffusion method and guidelines in accordance with CLSI (2014 Version). The *E. coli* isolates were tested for susceptibility to 17 different antibiotics from Oxoid, UK. Oxoid Meropenem AST Discs were shown to perform effectively in clinical settings.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.


References

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Vading, M. S. (2011). Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing *Klebsiella Pneumoniae* with the EUCAST and CLSI Breakpoint Systems, *Clinical Microbiology and Infection*, P668-674.
5. Fu, W. Z. (2003). The Susceptibility of Non-Fermentative Gram-Negative Bacilli to Cefperazone and Sulbactam Compared with Other Antibacterial Agents, *International Journal of Antimicrobial Agents*, 22 (4): 44444–48. [https://doi.org/10.1016/S0924-8579\(03\)00109-2](https://doi.org/10.1016/S0924-8579(03)00109-2).
6. Miao, Z. S. (2017). Antimicrobial Resistance and Molecular Epidemiology of ESBL-Producing *Escherichia Coli* Isolated from Outpatients in Town Hospitals of Shandong Province, China, *Frontiers in Microbiology*.

Glossary of Symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests

	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark
	Period-after-opening – the number of days the product can be used after opening
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

Revision information

Version	Date of issue and modifications introduced
2.0	2023-07-21

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- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Vading, M. S. (2011). Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing *Klebsiella Pneumoniae* with the EUCAST and CLSI Breakpoint Systems, *Clinical Microbiology and Infection*, P668-674.
- Fu, W. Z. (2003). The Susceptibility of Non-Fermentative Gram-Negative Bacilli to Cefperazone and Sulbactam Compared with Other Antibacterial Agents, *International Journal of Antimicrobial Agents*, 22 (4): 4444–48. [https://doi.org/10.1016/S0924-8579\(03\)00109-2](https://doi.org/10.1016/S0924-8579(03)00109-2).
- Miao, Z. S. (2017). Antimicrobial Resistance and Molecular Epidemiology of ESBL-Producing *Escherichia Coli* Isolated from Outpatients in Town Hospitals of Shandong Province, China, *Frontiers in Microbiology*.

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Referencije

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Vading, M. S. (2011). Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing *Klebsiella Pneumoniae* with the EUCAST and CLSI Breakpoint Systems, *Clinical Microbiology and Infection*, P668-674.
- Fu, W. Z. (2003). The Susceptibility of Non-Fermentative Gram-Negative Bacilli to Ceftazidime and Sulbactam Compared with Other Antibacterial Agents, *International Journal of Antimicrobial Agents*, 22 (4): 4444–48. [https://doi.org/10.1016/S0924-8579\(03\)00109-2](https://doi.org/10.1016/S0924-8579(03)00109-2).
- Miao, Z. S. (2017). Antimicrobial Resistance and Molecular Epidemiology of ESBL-Producing *Escherichia Coli* Isolated from Outpatients in Town Hospitals of Shandong Province, China, *Frontiers in Microbiology*.

© 2023 Thermo Fisher Scientific Inc. Sva prava pridržana. ATCC® je zaštitni znak društva ATCC. Svi ostali zaštitni znakovi vlasništvo su društva Thermo Fisher Scientific Inc. i njegovih podružnica. Ove informacije nisu namijenjene za poticanje na uporabu tih proizvoda na način kojim se mogu povrijediti prava intelektualnog vlasništva drugih društava.

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Reference

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Vading, M. S. (2011). Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing *Klebsiella Pneumoniae* with the EUCAST and CLSI Breakpoint Systems, *Clinical Microbiology and Infection*, P668-674.
- Fu, W. Z. (2003). The Susceptibility of Non-Fermentative Gram-Negative Bacilli to Ceftazidime and Sulbactam Compared with Other Antibacterial Agents, *International Journal of Antimicrobial Agents*, 22 (4): 44444–48. [https://doi.org/10.1016/S0924-8579\(03\)00109-2](https://doi.org/10.1016/S0924-8579(03)00109-2).
- Miao, Z. S. (2017). Antimicrobial Resistance and Molecular Epidemiology of ESBL-Producing *Escherichia Coli* Isolated from Outpatients in Town Hospitals of Shandong Province, China, *Frontiers in Microbiology*.

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Viited

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Vading, M. S. (2011). Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing *Klebsiella pneumoniae* with the EUCAST and CLSI Breakpoint Systems, *Clinical Microbiology and Infection*, P668-674.
- Fu, W. Z. (2003). The Susceptibility of Non-Fermentative Gram-Negative Bacilli to Cefperazone and Sulbactam Compared with Other Antibacterial Agents, *International Journal of Antimicrobial Agents*, 22 (4): 44444–48. [https://doi.org/10.1016/S0924-8579\(03\)00109-2](https://doi.org/10.1016/S0924-8579(03)00109-2).
- Miao, Z. S. (2017). Antimicrobial Resistance and Molecular Epidemiology of ESBL-Producing *Escherichia coli* Isolated from Outpatients in Town Hospitals of Shandong Province, China, *Frontiers in Microbiology*.

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Références

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (-)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Vading, M. S. (2011). Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing *Klebsiella pneumoniae* with the EUCAST and CLSI Breakpoint Systems, *Clinical Microbiology and Infection*, P668-674.
- Fu, W. Z. (2003). The Susceptibility of Non-Fermentative Gram-Negative Bacilli to Ceftazidime and Sulbactam Compared with Other Antibacterial Agents, *International Journal of Antimicrobial Agents*, 22 (4): 4444-48. [https://doi.org/10.1016/S0924-8579\(03\)00109-2](https://doi.org/10.1016/S0924-8579(03)00109-2).
- Miao, Z. S. (2017). Antimicrobial Resistance and Molecular Epidemiology of ESBL-Producing *Escherichia coli* Isolated from Outpatients in Town Hospitals of Shandong Province, China, *Frontiers in Microbiology*.

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Referenzen

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Vading, M. S. (2011). Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing *Klebsiella pneumoniae* with the EUCAST and CLSI Breakpoint Systems, *Clinical Microbiology and Infection*, P668-674.
- Fu, W. Z. (2003). The Susceptibility of Non-Fermentative Gram-Negative Bacilli to Ceftazidime and Sulbactam Compared with Other Antibacterial Agents, *International Journal of Antimicrobial Agents*, 22 (4): 4444–48. [https://doi.org/10.1016/S0924-8579\(03\)00109-2](https://doi.org/10.1016/S0924-8579(03)00109-2).
- Miao, Z. S. (2017). Antimicrobial Resistance and Molecular Epidemiology of ESBL-Producing *Escherichia coli* Isolated from Outpatients in Town Hospitals of Shandong Province, China, *Frontiers in Microbiology*.

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Bibliografia

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (-)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Vading, M. S. (2011). Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing *Klebsiella Pneumoniae* with the EUCAST and CLSI Breakpoint Systems, *Clinical Microbiology and Infection*, P668-674.
5. Fu, W. Z. (2003). The Susceptibility of Non-Fermentative Gram-Negative Bacilli to Ceftazidime and Sulbactam Compared with Other Antibacterial Agents, *International Journal of Antimicrobial Agents*, 22 (4): 4444-48. [https://doi.org/10.1016/S0924-8579\(03\)00109-2](https://doi.org/10.1016/S0924-8579(03)00109-2).
6. Miao, Z. S. (2017). Antimicrobial Resistance and Molecular Epidemiology of ESBL-Producing *Escherichia Coli* Isolated from Outpatients in Town Hospitals of Shandong Province, China, *Frontiers in Microbiology*.

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Referências

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (-)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Vading, M. S. (2011). Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing *Klebsiella Pneumoniae* with the EUCAST and CLSI Breakpoint Systems, *Clinical Microbiology and Infection*, P668-674.
5. Fu, W. Z. (2003). The Susceptibility of Non-Fermentative Gram-Negative Bacilli to Cefperazone and Sulbactam Compared with Other Antibacterial Agents, *International Journal of Antimicrobial Agents*, 22 (4): 44444–48. [https://doi.org/10.1016/S0924-8579\(03\)00109-2](https://doi.org/10.1016/S0924-8579(03)00109-2).
6. Miao, Z. S. (2017). Antimicrobial Resistance and Molecular Epidemiology of ESBL-Producing *Escherichia Coli* Isolated from Outpatients in Town Hospitals of Shandong Province, China, *Frontiers in Microbiology*.

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Referințe

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Vading, M. S. (2011). Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing *Klebsiella Pneumoniae* with the EUCAST and CLSI Breakpoint Systems, *Clinical Microbiology and Infection*, P668-674.
5. Fu, W. Z. (2003). The Susceptibility of Non-Fermentative Gram-Negative Bacilli to Cefperazone and Sulbactam Compared with Other Antibacterial Agents, *International Journal of Antimicrobial Agents*, 22 (4): 4444–48. [https://doi.org/10.1016/S0924-8579\(03\)00109-2](https://doi.org/10.1016/S0924-8579(03)00109-2).
6. Miao, Z. S. (2017). Antimicrobial Resistance and Molecular Epidemiology of ESBL-Producing *Escherichia Coli* Isolated from Outpatients in Town Hospitals of Shandong Province, China, *Frontiers in Microbiology*.

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Referencias

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Vading, M. S. (2011). Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing *Klebsiella Pneumoniae* with the EUCAST and CLSI Breakpoint Systems, *Clinical Microbiology and Infection*, P668-674.
5. Fu, W. Z. (2003). The Susceptibility of Non-Fermentative Gram-Negative Bacilli to Ceftazidime and Sulbactam Compared with Other Antibacterial Agents, *International Journal of Antimicrobial Agents*, 22 (4): 4444–48. [https://doi.org/10.1016/S0924-8579\(03\)00109-2](https://doi.org/10.1016/S0924-8579(03)00109-2).
6. Miao, Z. S. (2017). Antimicrobial Resistance and Molecular Epidemiology of ESBL-Producing *Escherichia Coli* Isolated from Outpatients in Town Hospitals of Shandong Province, China, *Frontiers in Microbiology*.

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Referenser

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Vading, M. S. (2011). Comparison of Disk Diffusion, Etest and VITEK2 for Detection of Carbapenemase-Producing *Klebsiella Pneumoniae* with the EUCAST and CLSI Breakpoint Systems, *Clinical Microbiology and Infection*, P668-674.
- Fu, W. Z. (2003). The Susceptibility of Non-Fermentative Gram-Negative Bacilli to Cefperazone and Sulbactam Compared with Other Antibacterial Agents, *International Journal of Antimicrobial Agents*, 22 (4): 44444–48. [https://doi.org/10.1016/S0924-8579\(03\)00109-2](https://doi.org/10.1016/S0924-8579(03)00109-2).
- Miao, Z. S. (2017). Antimicrobial Resistance and Molecular Epidemiology of ESBL-Producing *Escherichia Coli* Isolated from Outpatients in Town Hospitals of Shandong Province, China, *Frontiers in Microbiology*.

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