www.thermofisher.com

Thermo
SCIENTIFIC

Typical Formula

	grams per litre
'Lab-Lemco' powder	1.0
Yeast extract	2.0
Peptone	5.0
Sodium chloride	5.0

Nutrient Broth

EN

REF CM0001B & CM0001R

Intended Use

Nutrient Broth (CM0001B/R) is a general-purpose medium for the cultivation of fastidious and non-fastidious microorganisms which can be enriched for special purposes.

Nutrient Broth (CM0001B/R) device is used in a diagnostic workflow to support the clinicians for the growth of a wide variety of fastidious and non-fastidious microorganisms from clinical samples of patients suspected of having microbial infections.

The devices are for professional use only, are not automated, nor are they companion diagnostics.

Summary and Explanation

Bacteria can be classified as fastidious or non-fastidious depending on their nutritional and growth requirements: while non-fastidious bacteria do not normally require special nutritional supplements or conditions because they can synthesize all the necessary organic molecules. Fastidious bacteria will not grow without specific conditions or factors present because they have very complex and restricted nutritional and/or environmental requirements. Both groups of bacteria include medically important species, and therefore it is important to have a highly nutritious and versatile type of medium that can be enriched with blood or incubated at different temperatures and oxygen conditions in order to target the clinical species of interest, either fastidious or non-fastidious.

Non-fastidious microorganisms are organisms that grow and replicate without special nutritional supplements or conditions. They can synthesize all organic molecules required for their growth. Non-fastidious microorganisms can grow in nutrient sparse environments and can survive in a wider temperature range. This groups of bacteria are involved in causing illnesses and infections in humans and includes *E. coli* and *S. aureus*^{1,2,3,4,5}.

The first step in processing a clinical specimen in order to isolate organisms in the laboratory is the inoculation of culture medium. Cell culture is the most commonly used method for cultivating and isolating intracellular bacteria from clinical samples. There are two possible approaches. First, culture conditions can be standardised to allow for the efficient growth of a diverse range of bacteria. The second strategy consists isolating a previously known bacterium and exposing it to specially adapted culture conditions in order to optimise bacterial growth. Non-selective culture media contain no inhibitors and should allow the growth of the majority of microorganisms found in clinical samples.

Principle of Method

Lab-Lemco beef extract is combined with peptone and sodium chloride to form the basic bouillon described by Loeffler and other early bacteriologists. Yeast extract is added to provide vitamins and minerals to help speed the growth of most organisms.

Materials Provided

CM0001B: 500g of Nutrient Broth
CM0001R: 2.5kg of Nutrient Broth

500g of dehydrated Nutrient Broth yields approximately 38.5L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dish
- Defibrinated Horse Blood (SR0050B/C/D/E/F/R)
- Defibrinated Sheep Blood (SR0051B/C/D/E)

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 25°C.

Warnings and Precautions

- If inhaled: Remove to fresh air. Get medical attention immediately if symptoms occur.
- If ingested: clean mouth with water and drink afterwards plenty of water. Get medical attention if symptoms occur.
- If contact with skin: Wash off immediately with plenty of water for at least 15 minutes. Get medical attention immediately if symptoms occur.
- If contact with eyes: Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Get medical attention.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) Q 5.

Procedure

Dissolve 13g in 1 litre of distilled water. Mix well and distribute into final containers. Sterilize by autoclaving at 121°C for 15 minutes.

Interpretation

Turbidity indicates bacterial growth.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18 hours @ 37°C

Positive Controls	
Inoculum level: 10 – 100 cfu	
<i>Enterococcus faecalis</i> ATCC® 19433™	Turbid growth
<i>Streptococcus pyogenes</i> ATCC® 19615™	Turbid growth
<i>Staphylococcus aureus</i> ATCC® 25923™	Turbid growth
<i>Escherichia coli</i> ATCC® 25922™	Turbid growth
<i>Bacillus subtilis</i> ATCC® 6633™	Turbid growth
<i>Pseudomonas aeruginosa</i> ATCC® 27853™	Turbid growth

Limitations

Nutrient Broth (CM0001B/R) are non-selective media intended to be used for the maintenance of microorganisms isolated from other culture media products previously inoculated using clinical samples, and there are no expectations of specimens negative for the assay but maintain non-fastidious microorganisms in a condition that can mimic the test condition.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of non-fastidious organisms is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the devices. The precision of Nutrient Broth (CM0001B/R) devices was demonstrated by an overall pass rate of 100% obtained for the product over one year of testing (10.NOV.2020 – 17.AUG.2022; 10 batches). This shows that the performance is reproducible.

Nutrient Broth (CM0001B/R) device is tested in-house as part of the QC process since the products were launched in 1996. For target organisms, when using 10-100 cfu inoculum of *Enterococcus faecalis* ATCC® 19433™, *Streptococcus pyogenes* ATCC® 19615™, *Staphylococcus aureus* ATCC® 25923™, *Escherichia coli* ATCC® 25922™, *Bacillus subtilis* ATCC® 6633™ or *Pseudomonas aeruginosa* ATCC® 27853™ and incubating the device at 37°C for 18 hours, the user can expect growth as listed in this document.

Bibliography

- World Health Organization (WHO). 2017. 'Microbial Fact Sheets: Guidelines for drinking-water quality'. https://www.who.int/water_sanitation_health/publications/gdw-q4-with-add1-chap11.pdf.
- Public Health England. 2015. 'Identification of Enterobacteriaceae'. UK Standards for Microbiology Investigations ID 16 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf
- Public Health England. 2015. 'Identification of Vero cytotoxin-producing Escherichia coli including Escherichia coli O157'. UK Standards for Microbiology Investigations ID 22 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/436870/ID_22i4.pdf
- Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species'. UK Standards for Microbiology Investigations ID 07 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887570/UK_SMI_ID_07_i4.pdf
- Barton, Michelle, Michael Hawkes, Dorothy Moore, John Conly, Lindsay Nicolle, Upton Allen, Nora Boyd, et al. 2006. 'Guidelines for the Prevention and Management of Community-Associated Methicillin-Resistant Staphylococcus aureus: A Perspective for Canadian Health Care Practitioners.' Canadian Journal of Infectious Diseases and Medical Microbiology 17 (c): 4C-24C. <https://doi.org/10.1155/2006/971352>

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use

	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

ATCC Licensed Derivative[®]

©2022 Thermo Fisher Scientific Inc. All rights reserved.
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke,
RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of modifications introduced
2.0	2023-06-13

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Typisk formel

	Gram pr. liter
"Lab-Lemco"-pulver	1,0
Gærekstrakt	2,0
Pepton	5,0
Natriumchlorid	5,0

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Bibliografi

- Verdenssundhedsorganisationen (WHO). 2017. "Microbial Fact Sheets: Guidelines for drinking-water quality". https://www.who.int/water_sanitation_health/publications/gdw-q4-with-add1-chap11.pdf.
- Public Health England. 2015. "Identification of Enterobacteriaceae". UK Standards for Microbiology Investigations ID 16 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf
- Public Health England. 2015. "Identification of Vero cytotoxin-producing *Escherichia coli* including *Escherichia coli* O157". UK Standards for Microbiology Investigations ID 22 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/436870/ID_22i4.pdf
- Public Health England. 2020. "Identification of *Staphylococcus* species, *Micrococcus* species and *Rothia* species". UK Standards for Microbiology Investigations ID 07 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887570/UK_SMI_ID_07_i4.pdf
- Barton, Michelle, Michael Hawkes, Dorothy Moore, John Conly, Lindsay Nicolle, Upton Allen, Nora Boyd et al. 2006. "Guidelines for the Prevention and Management of Community-Associated Methicillin-Resistant *Staphylococcus aureus*: A Perspective for Canadian Health Care Practitioners." *Canadian Journal of Infectious Diseases and Medical Microbiology* 17 (c): 4C-24C. <https://doi.org/10.1155/2006/971352>

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Bibliographie

1. World Health Organization (WHO). 2017. 'Microbial Fact Sheets: Guidelines for drinking-water quality'. https://www.who.int/water_sanitation_health/publications/gdw-q4-with-add1-chap11.pdf.
2. Public Health England. 2015. 'Identification of Enterobacteriaceae'. UK Standards for Microbiology Investigations ID 16 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf
3. Public Health England. 2015. 'Identification of Vero cytotoxin-producing Escherichia coli including Escherichia coli O157'. UK Standards for Microbiology Investigations ID 22 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/436870/ID_22i4.pdf
4. Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species'. UK Standards for Microbiology Investigations ID 07 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887570/UK_SMI_ID_07_i4.pdf
5. Barton, Michelle, Michael Hawkes, Dorothy Moore, John Conly, Lindsay Nicolle, Upton Allen, Nora Boyd, et al. 2006. 'Guidelines for the Prevention and Management of Community-Associated Methicillin-Resistant Staphylococcus aureus: A Perspective for Canadian Health Care Practitioners.' Canadian Journal of Infectious Diseases and Medical Microbiology 17 (c): 4C-24C. <https://doi.org/10.1155/2006/971352>

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Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, UK

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
Bibliographie

1. Organisation mondiale de la santé (OMS). 2017. « Microbial Fact Sheets: Guidelines for drinking-water quality ». https://www.who.int/water_sanitation_health/publications/gdw-q4-with-add1-chap11.pdf.
2. Public Health England. 2015. « Identification of Enterobacteriaceae ». UK Standards for Microbiology Investigations ID 16 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf
3. Public Health England. 2015. « Identification of Vero cytotoxin-producing Escherichia coli including Escherichia coli O157 ». UK Standards for Microbiology Investigations ID 22 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/436870/ID_22i4.pdf
4. Public Health England. 2020. « Identification of Staphylococcus species, Micrococcus species and Rothia species ». UK Standards for Microbiology Investigations ID 07 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887570/UK_SMI_ID_07i4.pdf
5. Barton, Michelle, Michael Hawkes, Dorothy Moore, John Conly, Lindsay Nicolle, Upton Allen, Nora Boyd, et al. 2006. « Guidelines for the Prevention and Management of Community-Associated Methicillin-Resistant Staphylococcus aureus: A Perspective for Canadian Health Care Practitioners ». Canadian Journal of Infectious Diseases and Medical Microbiology 17 (c): 4C-24C. <https://doi.org/10.1155/2006/971352>

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Oxoid Limited, Wade Road, Basingstoke,
RG24 8PW, Royaume-Uni

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Bibliografia

- World Health Organization (WHO). 2017. "Microbial Fact Sheets: Guidelines for drinking-water quality". https://www.who.int/water_sanitation_health/publications/gdw-q4-with-add1-chap11.pdf.
- Public Health England. 2015. "Identification of Enterobacteriaceae". UK Standards for Microbiology Investigations ID 16 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf
- Public Health England. 2015. "Identification of Vero cytotoxin-producing Escherichia coli including Escherichia coli O157". UK Standards for Microbiology Investigations ID 22 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/436870/ID_22i4.pdf
- Public Health England. 2020. "Identification of Staphylococcus species, Micrococcus species and Rothia species". UK Standards for Microbiology Investigations ID 07 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887570/UK_SMI_ID_07_i4.pdf
- Barton, Michelle, Michael Hawkes, Dorothy Moore, John Conly, Lindsay Nicolle, Upton Allen, Nora Boyd, et al. 2006. "Guidelines for the Prevention and Management of Community-Associated Methicillin-Resistant Staphylococcus aureus: A Perspective for Canadian Health Care Practitioners." Canadian Journal of Infectious Diseases and Medical Microbiology 17 (c): 4C-24C. <https://doi.org/10.1155/2006/971352>

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Literatūra

1. Pasaulio sveikatos organizacija (PSO). 2017. 'Microbial Fact Sheets: Guidelines for drinking-water quality'. https://www.who.int/water_sanitation_health/publications/gdw-q4-with-add1-chap11.pdf.
2. Public Health England. 2015. 'Identification of Enterobacteriaceae'. UK Standards for Microbiology Investigations ID 16 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf
3. Public Health England. 2015. 'Identification of Vero cytotoxin-producing Escherichia coli including Escherichia coli O157'. UK Standards for Microbiology Investigations ID 22 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/436870/ID_22i4.pdf
4. Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species'. UK Standards for Microbiology Investigations ID 07 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887570/UK_SMI_ID_07_i4.pdf
5. Barton, Michelle, Michael Hawkes, Dorothy Moore, John Conly, Lindsay Nicolle, Upton Allen, Nora Boyd, et al. 2006. 'Guidelines for the Prevention and Management of Community-Associated Methicillin-Resistant Staphylococcus aureus: A Perspective for Canadian Health Care Practitioners.' Canadian Journal of Infectious Diseases and Medical Microbiology 17 (c): 4C-24C. <https://doi.org/10.1155/2006/971352>

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Bibliografia

- World Health Organization (WHO). 2017. „Microbial Fact Sheets: Guidelines for drinking-water quality”. https://www.who.int/water_sanitation_health/publications/gdw-q4-with-add1-chap11.pdf.
- Public Health England. 2015. „Identification of Enterobacteriaceae”. UK Standards for Microbiology Investigations ID 16 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf
- Public Health England. 2015. „Identification of Vero cytotoxin-producing Escherichia coli including Escherichia coli O157”. UK Standards for Microbiology Investigations ID 22 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/436870/ID_22i4.pdf
- Public Health England. 2020. „Identification of Staphylococcus species, Micrococcus species and Rothia species”. UK Standards for Microbiology Investigations ID 07 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887570/UK_SMI_ID_07_i4.pdf
- Barton, Michelle, Michael Hawkes, Dorothy Moore, John Conly, Lindsay Nicolle, Upton Allen, Nora Boyd, et al. 2006. „Guidelines for the Prevention and Management of Community-Associated Methicillin-Resistant Staphylococcus aureus: A Perspective for Canadian Health Care Practitioners.” Canadian Journal of Infectious Diseases and Medical Microbiology 17 (c): 4C-24C. <https://doi.org/10.1155/2006/971352>

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Bibliografie

1. World Health Organization (WHO). 2017. „Microbial Fact Sheets: Guidelines for drinking-water quality”. https://www.who.int/water_sanitation_health/publications/gdw-q4-with-add1-chap11.pdf.
2. Public Health England. 2015. „Identification of Enterobacteriaceae.” UK Standards for Microbiology Investigations ID 16 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf
3. Public Health England. 2015. „Identification of Vero cytotoxin-producing Escherichia coli including Escherichia coli O157”. UK Standards for Microbiology Investigations ID 22 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/436870/ID_22i4.pdf.
4. Public Health England. 2020. „Identification of Staphylococcus species, Micrococcus species and Rothia species.” UK Standards for Microbiology Investigations ID 07 (4). https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887570/UK_SMI_ID_07_i4.pdf
5. Barton, Michelle, Michael Hawkes, Dorothy Moore, John Conly, Lindsay Nicolle, Upton Allen, Nora Boyd, et al. 2006. „Guidelines for the Prevention and Management of Community-Associated Methicillin-Resistant Staphylococcus aureus: A Perspective for Canadian Health Care Practitioners.” Canadian Journal of Infectious Diseases and Medical Microbiology 17 (c): 4C-24C. <https://doi.org/10.1155/2006/971352>

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