www.thermofisher.com

MacConkey Agar

EN

REF CM0007A CM0007B, CM0007R, CM0007T, CM0007K

Intended Use

MacConkey Agar (CM0007A/B/K/R/T) devices are differential media for the isolation and differentiation of lactose and non-lactose fermenting coliforms and intestinal bacteria pathogens from a wide range of clinical samples direct from the patient, such as urine and faeces, or indirectly from swabs taken from wounds or infections. MacConkey Agar devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The devices are for professional use only, are not automated, nor are they companion diagnostics.

Summary and Explanation

Due to its ability to support the growth of pathogenic Gram- positive cocci (e.g., staphylococci and enterococci) as well as Enterobacteriaceae, MacConkey Agar is particularly recommended for the cultivation of pathogens which may be present in a variety of specimens such as urine, faeces or wound swabs. It provides several diagnostic indications in addition to bile tolerance, such as colony morphology and chromogenesis.

Enterococci can grow and survive in harsh conditions, they can be found in the gastrointestinal tract of humans and animals, plants, soil, water, environment and fermented products. *Enterococcus faecalis* and *Enterococcus faecium* are the most common enterococci isolated from human infections. Enterococci demonstrates both intrinsic and acquired resistance to a number of antibiotics¹.

Staphylococci are gram positive and the most frequently associated species with human infection. *Staphylococcus aureus* can be associated with severe infection therefore it is important to distinguish it from the opportunistic coagulase negative staphylococci².

Members of the Enterobacteriaceae include *Escherichia*, *Proteus*, *Klebsiella*, *Salmonella* and *Shigella*. Within the *Escherichia* spp. there are five known species to cause human disease. There are four species of *Proteus*, of which three cause disease. The genus *Klebsiella* contains six species, of which four are related to humans and include *K. pneumoniae* that can cause bacteraemia and hepatic infections and have been isolated from a number of unusual infections, including endocarditis, primary gas-containing mediastinal abscess, peritonitis, acute cholecystitis, crepitant myonecrosis, pyomyositis, necrotising fasciitis, psoas muscle abscess, fascial space infections of the head and neck, and septic arthritis³.

Principle of Method

The formulation of MacConkey Agar (CM0007A/B/K/R/T) does not contain crystal violet, unlike other MacConkey Agar formulations such as MacConkey Agar No. 3 (CM0115B/T/K/Q/V/W/R) and Sorbitol MacConkey Agar (CM0813B/K/R), which differ in their intended purpose. The formulation of MacConkey Agar (CM0007A/B/K/R/T) is less selective and permits the growth of *Staphylococcus* spp. and *Enterococcus* spp. MacConkey Agar (CM0007A/B/K/R/T) devices contain peptone to supply nutrients for growth, agar as the solidifying agent and sodium chloride to maintain the osmotic equilibrium. Bile salts are inhibitory to non-intestinal bacteria and help prevent *Proteus* spp. swarming. Lactose is added as a carbon source. Differentiation of lactose fermenters is achieved by the combination of lactose and the indicator dye, neutral red. Lactose fermenters produce acid, reducing the pH and causing the dye to change from yellow to red. This results in lactose fermenting bacteria appearing as red-pink colonies which may be surrounded by zones of precipitated bile salts caused by the action of the acid produced by lactose fermentation. Non-lactose fermenters such as *Salmonella* spp., appear as straw or colourless colonies.

Typical Formula

	grams per litre
Peptone	20.0
Lactose	10.0
Bile salts	5.0
Sodium chloride	5.0
Neutral red indicator	0.075
Agar	12.0

Materials Provided

(500g makes 9.6 litres)

CM0007A: 100g of dehydrated MacConkey Agar

CM0007B: 500g of dehydrated MacConkey Agar

CM0007K: 25kg of dehydrated MacConkey Agar

CM0007R: 2.5kg of dehydrated MacConkey Agar

CM0007T: 5kg of dehydrated MacConkey Agar

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri Dishes

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse thoroughly with water, also under the eyelids.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B 1, B 4, B 11, B 14, B 17, B 25, B 26, B 37 B 57 and B 60.

Procedure

- Suspend 52g per 1 litre of distilled water.
- Bring to the boil to dissolve completely.
- Sterilize by autoclaving at 121°C for 15 minutes.
- Cool to 50°C. Mix well and pour into sterile petri dishes.
- Store using suggested guidelines until use.

Interpretation

Lactose fermenting bacteria appear as red-pink colonies, which may be surrounded by zones of precipitated bile salts.

Non-lactose fermenters such as *Salmonella* spp., appear as straw or colourless colonies.

Identification is presumptive and should be confirmed.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18 - 24 h @ 37°C aerobic

Positive Controls Colony count is ≥ 70% of the control medium count Inoculum level: 10 - 100 cfu	
<i>Escherichia coli</i> ATCC 25922	1-2mm pink/red colonies
<i>Klebsiella pneumoniae</i> ATCC® 29665™	1-3mm red/pink, mucoid colonies

<i>Enterococcus faecalis</i> ATCC® 29212™	0.5-1mm red colonies
<i>Staphylococcus aureus</i> ATCC® 25923™	0.5-1mm pink colonies
<i>Salmonella typhimurium</i> ATCC® 14028™	1-3mm straw colonies
<i>Salmonella nottingham</i> NCTC 7832	1-3mm straw colonies
<i>Shigella flexneri</i> ATCC® 12022™	1-4mm straw colonies
<i>Pseudomonas aeruginosa</i> ATCC® 27853™	1-3mm straw colonies
<i>Proteus mirabilis</i> ATCC® 29906™	1-2mm straw colonies, no swarming

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of lactose and non-lactose fermenting coliform organisms and intestinal bacteria pathogens is confirmed by the inclusion of a well- characterised isolate in the QC processes performed as part of the manufacture of each batch of the device. The precision of MacConkey Agar (CM0007A/B/R/T/K) was demonstrated by an overall pass rate of 100% over 2 years of testing (October 2019 - December 2021). This shows that the performance is reproducible.

The devices are tested in-house as part of the QC process. For target organisms, when using 10-100 cfu inoculum of, *Escherichia coli*, *Klebsiella pneumoniae*, *Enterococcus faecalis*, *Staphylococcus aureus*, *Salmonella typhimurium*, *Salmonella nottingham*, *Shigella flexneri*, *Pseudomonas aeruginosa*, and *Proteus mirabilis* and incubating the device at 37°C for 18-24 hours, the user can recover organisms with colony size and morphology as detailed in this document.

Limitations


Streptococci and other organisms requiring specific growth factors may grow poorly or not at all on this medium, therefore the specimen should also be cultivated on a non- selective medium such as blood agar. Some strains of *Proteus* spp. may be able to swarm. Prolonged incubation may lead to confusion of results.

Bibliography

1. Standing Committee of Analysts. The Microbiology of Drinking Water (2009) Part 4.
2. Hoogendijk J. L., *Antonie van Leeuwenhoek* 1962, 28 (1) 315-320.
3. MacConkey A. T., 1905 *J Hyg*, Jul; 5(3):333-379.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Consult instructions for use
	Do not use if packaging damaged and consult instructions for use

	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of issue and modifications introduced
2.0	2023-12-12

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Všetky ostatné ochranné známky sú vlastníctvom Thermo Fisher Scientific Inc. a jej dcérskych spoločností.
Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, SPOJENÉ KRÁĽOVSTVO