www.thermofisher.com

Peptone Water

REF CM0009B

Intended Use

Peptone water (CM0009B) is a medium for the dilution of microorganisms isolated from clinical and non-clinical samples. Can also be supplemented with carbohydrates and indicator for identification purposes.

The device is intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The device is for professional use only, is not automated, nor is it a companion diagnostic.

Summary and Explanation

Pathogenic aerobic and anaerobic microorganisms encompass a broad range of manifestations, which can often be found as part of the normal skin and mucosal flora of humans¹. *Escherichia coli* (*E. coli*) are non-spore forming bacteria that are able to grow in aerobic and anaerobic conditions². Whilst *Staphylococcus aureus* is also a nosocomial pathogen, with invasive medical interventions such as catheterization or immunosuppressant treatment, contributing to the incidence of peripheral intravenous catheter (PIVC) failure and catheter-related bloodstream infections (CRBSI)¹.

Infections from fastidious aerobic and anaerobic microorganisms can have a range of manifestations ranging from superficial to life-threatening^{1,2}. Intestinal pathogens like *Escherichia coli* increase their ability to persist and cause disease by utilizing various mechanisms to outcompete the host intestinal microbiota³. Additionally, infection with *Staphylococcus aureus*, a common commensal of the human upper respiratory tract and skin microbiomes, can cause superficial skin infections such as cellulitis and abscesses to life-threatening illnesses such as pneumonia, toxic shock syndrome and sepsis¹.

Principle of Method

Peptone is added to provide nutrients essential for bacterial growth and sodium chloride helps to maintain osmotic equilibrium.

Peptone Water, adjusted to pH 8.4, is suitable for the cultivation and enrichment of *Vibrio cholerae* from infected material.

Peptone Water may be modified for use in carbohydrate fermentation tests by the addition of Andrade's indicator. The indicator which is pink at pH 5.0 and yellow at pH 8.0 is prepared by adding sodium hydroxide to acid fuchsin until it becomes yellow. When added to Peptone Water it is colourless to slightly pink. Filter-sterilised 'sugar' solutions are added to the base medium after sterilisation. These solutions are usually at 10% w/v concentrations, and it is important to allow for dilution of the Peptone Water when making up the initial volume of medium.

Typical Formula

	grams per litre
Peptone	10.0
Sodium chloride	5.0

Materials Provided

CM0009B: 500g of Peptone Water

500g of Peptone Water yields approximately 33.3L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dish

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 25°C.

Warnings and Precautions

- If inhaled: remove to fresh air. Get medical attention immediately if symptoms occur.
- If ingested: clean mouth with water and drink afterwards plenty of water. Get medical attention if symptoms occur.
- In contact with skin: wash off immediately with plenty of water for at least 15 minutes. Get medical attention.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 1, ID 7, [ID 16](#) and Q 5.

Procedure

Dissolve 15g in 1 litre of distilled water. Mix well and distribute into final containers. Sterilize by autoclaving at 121°C 15 minutes. When sterile solutions are to be added after autoclaving, reduce the volume of water for reconstitution by an equal amount.

Interpretation

Once the medium is reconstituted:
Turbidity indicates bacterial growth

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Inoculate the media with 40µl of a 4-hour broth culture. After incubation add a few drops of Kovacs' reagent to each tube to assess indole production. Indole positive reaction is indicated by a red layer after the addition of Kovacs' reagent.

Incubation Conditions: 18 hours @ 37°C and 44°C

Positive Controls		
Inoculum level: 10 ⁴ – 10 ⁶ cfu		
	37°C	44°C
<i>Escherichia coli</i> ATCC® 25922™	Turbid growth; indole positive	Turbid growth; indole positive
<i>Escherichia coli</i> ATCC® 8739™	Turbid growth; indole positive	Turbid growth; indole positive
<i>Klebsiella pneumoniae</i> ATCC® 13883™	Turbid growth; indole negative	Turbid growth; indole negative
Add 1ml of the following 10% w/v sterile carbohydrate solutions to 9ml volumes of medium in bottles with Durham tubes:- Glucose, Mannitol, Lactose and Sucrose. Inoculate the enriched media and a control with 40µl of a 4 hour broth culture of <i>Salmonella typhimurium</i> ATCC® 14028. After incubation add 1ml of 0.03% w/v phenol red solution to each tube to demonstrate acid production. Acid production is indicated by a colour change to yellow		
Reactions after incubation at 37°C for 18 hours		
Glucose	Turbid growth, acid and gas	
Mannitol	Turbid growth, acid and gas	
Lactose	Turbid growth, no reaction and no gas	
Sucrose	Turbid growth, no reaction and no gas	
Control	Turbid growth, no reaction and no gas	
Inoculate the enriched media and a control with 40µl of a 4 hour broth culture of <i>Staphylococcus aureus</i> ATCC® 25923. After incubation add 1ml of 0.03% w/v phenol red solution to each tube to demonstrate acid production. Acid production is indicated by a colour change to yellow		

Reactions after incubation at 37°C for 18 hours	
Glucose	Turbid growth, acid and gas
Mannitol	Turbid growth, acid and gas
Lactose	Turbid growth, no reaction and no gas
Sucrose	Turbid growth, no reaction and no gas
Control	Turbid growth, no reaction and no gas

Limitations

Some sugar solutions may affect the pH of the Peptone Water; check and correct if so.

Sub-cultures may be necessary to ensure purity of the inoculant. Mixed or contaminated cultures will give false reactions

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of isolation of common microorganisms is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the devices. The precision of Peptone Water (CM0009B) devices was demonstrated by an overall pass rate of 100% obtained for the product over one year of testing (25.JUN.2020 – 24.JUN.2022; 10 batches). This shows that the performance is reproducible.

Peptone Water (CM0009B) device is tested in-house as part of the QC process since the products were launched in 1996. For target organisms, when using 10⁴ – 10⁶ cfu inoculum of *Escherichia coli* ATCC® 25922™, *Escherichia coli* ATCC® 8739™, *Klebsiella pneumoniae* ATCC® 13883™ and incubating the device at 37°C and 44°C for 18 hours, the user can recover organisms with growth as listed in this document.

Bibliography

- Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species'. ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>
- Public Health England. 2015. 'SMI ID 16: identification of Enterobacteriaceae'. ID 17 (3). <https://www.gov.uk/government/publications/smi-id-16-identification-of-enterobacteriaceae>
- Ormsby, M. J., Logan, M., Johnson, S. A., McIntosh, A., Fallata, G., Papadopolou, R., ... & Wall, D. M. (2019). Inflammation associated ethanolamine facilitates infection by Crohn's disease-linked adherent-invasive *Escherichia coli*. *EBioMedicine*, 43, 325-332.

Symbol Legend

Symbol	Definition
REF	Catalogue number
IVD	In Vitro Diagnostic Medical Device
LOT	Batch code

	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
Made in the United Kingdom	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke,
RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of modifications introduced
1.0	2022-11-22 Original Document

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© 2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.