www.thermofisher.com

Agar

Urea Agar Base

REF CM0053B

Intended Use

Urea Agar Base (CM0053B) is a non-selective medium intended for use in culture media. Urea Agar Base (CM0053B), when supplemented with Urea 40% Solution (SR0020K) for the preparation of Christensen's medium, is intended to be used to detect rapid urease activity of the *Proteae* and non-rapid urease activity of some *Enterobacteriaceae* in faecal and urine samples.

The device is intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having *Enterobacteriaceae* infections.

The device is for professional use only, is not automated and nor is it a companion diagnostic.

Summary and Explanation

Urinary tract infections, often referred to as UTIs, are common infections that occur when bacteria invade the urinary tract by entering via the urethra. These infections can impact different parts of the urinary tract, but the most prevalent is cystitis (a bladder infection). UTIs occur more frequently in women, as the urethra is shorter, and therefore closer to the rectum. This increases the likelihood of infection as it is easier for bacteria to enter the urinary tract. Kidney infections are another form of UTI. These are more serious than cystitis, but far less common. The symptoms of cystitis can include frequent urination, or needing to urinate whilst having an empty bladder, a burning sensation or pain whilst urinating, blood in the urine and pressure or cramps in the lower abdominal area. The symptoms of kidney infections are very different. These can include chills, a fever, nausea and vomiting and pain in the lower back¹. There are many species of bacteria that cause UTIs, but the three covered in this file are *Escherichia coli*, *Staphylococcus aureus* and *Proteus mirabilis*.

Due to the severity of UTIs in a vast number of patients, it is therefore highly important to be able to isolate and identify common UTI-causing microorganisms from urine samples. Early diagnosis is vital in the prevention and reduction of urinary tract infections.

Principle of Method

Urea Agar is a lightly buffered medium containing urea and phenol red, a pH indicator. When organisms utilize urea, ammonia is formed which turns the medium alkaline. The indicator, phenol red, changes the medium colour from pale-yellow to pink red in an alkaline environment. Gelatin peptone promotes rapid growth of enteric gram-negative bacilli, permitting a decrease in incubation time. Dextrose stimulates urease activity in those organisms which hydrolyze urea slowly.

Typical Formula

	grams per litre
Peptone	1.0
Glucose	1.0
Sodium chloride	5.0
Di-sodium phosphate	1.2
Potassium dihydrogen phosphate	0.8
Phenol red	0.012

Materials Provided

CM0053B: 500g of Urea Agar Base

500g of dehydrated Urea Agar Base yields approximately 20.8L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dish
- Supplements (SR0020K)

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

Warning:

- Harmful if swallowed.
- Causes severe skin burns and eye damage.
- Wear protective gloves/ protective clothing/ eye protection/ face protection.
- IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
- IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- Immediately call a POISON CENTRE or doctor/physician.
- IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.
- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any

federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI).

Procedure

Suspend 2.4g in 95ml of distilled water. Bring to the boil to dissolve completely. Sterilize by autoclaving at 115°C for 20 minutes. Cool to 50°C and aseptically add the contents of 1 vial of Urea 40% Solution (SR0020K). Mix well, aseptically dispense 10ml amounts into sterile containers and allow to set in the slope position.

Interpretation

Once the medium is reconstituted and challenged with 10^4 – 10^6 cfu and incubated at 30°C for 6 hours:

Weak urease positive, pale pink slope indicates

Morganella morganii

Urease positive, pink slope indicates *Proteus mirabilis*

Once the medium is reconstituted and challenged with 10^4 – 10^6 cfu and incubated at 30°C for 24 hours:

Urease positive, pink slope indicates *Morganella morganii*

Weak urease positive, pale pink slope indicates *Klebsiella pneumoniae*.

Urease negative, no colour change indicates *Shigella sonnei*, *Enterobacter aerogenes*, *Escherichia coli* (ATCC® 25922™ & ATCC® 11775™) and *Salmonella Nottingham*

Once the medium is reconstituted and challenged with 10^4 – 10^6 cfu and incubated at 30°C for up to 72 hours:

Urease negative, no colour change indicates *Candida albicans*

Weak urease positive, pale pink slope indicates

Cryptococcus albidus.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Positive Controls	
Inoculum level: 10^4 – 10^6 cfu Incubation Conditions: 30°C for 6 hours:	
<i>Morganella morganii</i> ATCC® 25830™	Weak urease +ve, pale pink slope, or negative

<i>Proteus mirabilis</i> ATCC® 29906™	Urease +ve, pink slope
Inoculum level: 10^4 – 10^6 cfu Incubation Conditions: 30°C for 24 hours:	
<i>Morganella morganii</i> ATCC® 25830™	Urease +ve, pink slope
<i>Klebsiella pneumoniae</i> ATCC® 13883™	Weak urease +ve, pale pink slope
<i>Shigella sonnei</i> ATCC® 25931™	Urease -ve, no colour change
<i>Enterobacter aerogenes</i> ATCC® 13048™	Urease -ve, no colour change
<i>Escherichia coli</i> ATCC® 25922™	Urease -ve, no colour change
<i>Escherichia coli</i> ATCC® 11775™	Urease -ve, no colour change
<i>Salmonella nottingham</i> NCTC 7832	Urease -ve, no colour change
Inoculum level: 10^4 – 10^6 cfu Incubation Conditions: 30°C for up to 72 hours:	
<i>Candida albicans</i> ATCC® 10231™	Urease -ve, no colour change
<i>Cryptococcus albidus</i> ATCC® 34140™	Weak urease +ve, pale pink slope

Limitations

Identifications are presumptive and colonies should be confirmed using appropriate methods.

Urea test media rely on demonstration of alkalinity and are not specific for detection of urease activity. Peptones in the media may be hydrolyzed releasing amino acid residues, raising pH, and resulting in false-positive reactions. A control test using the same test medium without urea can be used to facilitate interpretation of questionable reactions.

For Urea Agar results to be available for rapid urease-positive organisms, the slant should be read within the first 2-6 hours after inoculation/incubation. Other members of the *Enterobacteriaceae* such as *Klebsiella* and *Enterobacter*, exhibit a delayed urease reaction which may require 24-48 hours of incubation.

Performance Characteristics


Accuracy has been demonstrated through review of the QC data. Correct detection of *Enterobacteriaceae* is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the devices. The precision of Urea Agar Base (CM0053B) devices was demonstrated by an overall pass rate of 100% obtained for the product over two years of testing (04.JUN.2020 – 06.JUN.2022; 10 batches). This shows that the performance is reproducible.

Urea Agar Base (CM0053B) device is tested in-house as part of the QC process, as the products have been on the market since at least 1997. For target organisms, when using 10^4 – 10^6 cfu inoculum of *Morganella morganii* or *Proteus mirabilis* and incubating the device at 30°C for 6 hours, the user can expect urease reactions as listed in this document.

Bibliography

- Centers for Disease Control and Prevention. 2021. 'Urinary Tract Infection'. <https://www.cdc.gov/antibiotic-use/uti.html>.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.
 ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
 NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures.
 All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke,
 RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of modifications introduced
1.0	2022-09-29 Original Document

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