Mannitol Salt Agar

EN

REF CB0085G, CM0085B, CM0085K, CM0085R, CM0085T & CM0085W

Intended Use

Mannitol Salt Agar (CB0085G, CM0085B/K/R/T/W) devices are selective media for the isolation of pathogenic staphylococci from clinical samples, including wound, faecal and respiratory tract samples.

Mannitol Salt Agar (CB0085G, CM0085B/K/R/T/W) devices are intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having staphylococcal infections.

The devices are for professional use only, are not automated and nor are they companion diagnostics.

Summary and Explanation

Staphylococci are Gram-positive cocci, widespread in nature and generally non-pathogenic^{1,2}. However, some staphylococci can cause disease in humans^{1,2}. The staphylococcus most frequently associated with disease in humans is *Staphylococcus aureus*¹, an opportunistic pathogen commonly found as a commensal of the human upper respiratory tract and skin microbiomes². *S. aureus* enters the body through cuts or abrasions where it can cause a spectrum of disease, from minor infections such as cellulitis and abscesses to life-threatening illnesses such as pneumonia, toxic shock syndrome and sepsis².

Principle of Method

Mannitol Salt Agar devices contain Lab Lemco and peptone as a source of nitrogen and other nutrients. Mannitol fermentation aids the differentiation of *Staphylococcus* spp. and phenol red is the pH indicator; it is yellow at acid pH and red/pink at alkaline pH. High salt content will inhibit most other organisms.

Typical Formula

	grams per litre
'Lab-Lemco' powder	1.0
Peptone	10.0
Mannitol	10.0
Sodium chloride	75.0
Phenol red	0.025
Agar	15.0

Materials Provided

CB0085G: 10kg of dehydrated Mannitol Salt Agar

CM0085B: 500g of dehydrated Mannitol Salt Agar CM0085K: 25kg of dehydrated Mannitol Salt Agar CM0085R: 2.5kg of dehydrated

Mannitol Salt Agar CM0085T: 5kg of dehydrated Mannitol Salt Agar CM0085W: 11.1kg of dehydrated Mannitol Salt Agar

500g yields approximately 4.5L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dish

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.

- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 7, B17, B29, B57 and S7.

Procedure

Suspend 111g in 1 litre of distilled water. Bring to the boil to dissolve completely. Sterilise by autoclaving at 121°C for 15 minutes. Cool to 50°C. Mix well and pour into sterile containers.

Interpretation

Once the medium is reconstituted, yellow colonies with yellow halos indicate *Staphylococcus aureus*. Pink colonies and pink medium indicate *Staphylococcus epidermidis*.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 36 hours @ 37°C

Positive Controls Inoculum level: 10-100 cfu Colony count is ≥ 70% of the control medium count	
<i>Staphylococcus aureus</i> ATCC® 25923™	0.5-2 mm yellow colonies with yellow halo
<i>Staphylococcus aureus</i> ATCC® 6538™	0.5-2 mm yellow colonies with yellow halo
<i>Staphylococcus aureus</i> MRSA NCTC 10442	0.75-4 mm yellow colonies with yellow halo
<i>Staphylococcus epidermidis</i> ATCC® 12228™	0.5-2 mm pink colonies and pink medium
Negative Controls Inoculum level: 10 ⁴ -10 ⁶ cfu	
<i>Escherichia coli</i> ATCC® 25922™	No growth
Inoculum level: 10 ³ -10 ⁴ cfu	
<i>Escherichia coli</i> ATCC® 8739™	No growth or pinpoint-0.5 mm yellow colonies
Inoculum level: 10 ² -10 ³ cfu Colony count is ≥ 40% of the control medium count	
<i>Proteus mirabilis</i> ATCC® 29906™	Pinpoint to 2 mm pink colonies, no swarming
Testing performed in accordance with current CLSI M22 A Incubation Conditions: 36 hours @ 35°C Inoculum level 10-100 cfu Colony count is ≥ 40% of the control medium count	

<i>Proteus vulgaris</i> ATCC® 8427™	Pinpoint - 2 mm pink colonies, no swarming
-------------------------------------	--

Limitations

Identifications are presumptive and colonies must be confirmed using appropriate biochemical/serological methods. It is recommended that strains are sub-cultured to a non-selective medium containing less salt before carrying out identification tests, for example the tube coagulase test. Some strains of *S. aureus* may exhibit a delayed mannitol fermentation reaction; negative plates should be re-incubated for a further 24 hours. Extended incubation may result in accompanying flora growing through. *Proteus* species grow as pink colonies but will not swarm. Enterococci may grow as pin-point colonies. Strains of *S. aureus* with particular growth requirements may not grow on the medium and those with atypical enzyme patterns may give anomalous results.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of pathogenic staphylococci is confirmed by the inclusion of well-characterised isolates in the QC processes performed as part of the manufacture of each batch of the devices. The precision of Mannitol Salt Agar (CB0085G, CM0085B/K/R/T/W) was demonstrated by an overall pass rate of 100% obtained for the product over one year of testing (2021 – 2022; 10 batches). This shows that the performance is reproducible.

Mannitol Salt Agar (CB0085G, CM0085B/K/R/T/W) devices are tested in-house as part of the QC process since the products were launched in 1997 (CM0085B), 1998 (CM0085T), 2000 (CM0085R), 2004 (CB0085G) and 2009 (CM0085K). For target organisms, when using 10-100 cfu inoculum of *Staphylococcus aureus* (ATCC® 25923™; ATCC® 6538™), *Staphylococcus aureus* MRSA (NCTC 10442™) and *Staphylococcus epidermidis* (ATCC® 12228™), and incubating the device at 37°C for 36 hours, the user can recover organisms with colony size and morphology as listed in this document. For non-target organisms, when using 10²-10³ cfu inoculum of *Proteus mirabilis* (ATCC® 29906™), 10⁴-10⁶ cfu inoculum of *Escherichia coli* (ATCC® 25922™) or 10³-10⁴ cfu inoculum of *Escherichia coli* (ATCC® 8739™) and incubating the device at 37°C for 36 hours, the user can expect no growth or colony size and morphology as listed in this document. When using 10-100 cfu inoculum of *Proteus vulgaris* (ATCC® 8427™) and incubating the device at 35°C for 36 hours, the user can recover organisms with colony size and morphology as listed in this document.

Bibliography

1. Public Health England. 2020a. "Identification of Staphylococcus Species, Micrococcus Species and Rothia Species." UK Standards for Microbiology Investigations. ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
2. Foster, Timothy J. 1996. "Staphylococcus." In Medical Microbiology. 4th Edition., 4th ed. Galveston.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests

	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures.

CLSI is a trademark of the Clinical Laboratory and Standards Institute.

All other marks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

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© 2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.