Charcoal Agar Base (Dehydrated)

REF CM0119B

EN

Intended Use

Charcoal Agar Base (CM0119B) device is non-selective general-purpose dehydrated culture medium for the cultivation and isolation of fastidious organisms especially *Bordetella pertussis*, *Bordetella parapertussis* and *Haemophilus influenzae* from nasopharyngeal/pernasal samples and for subcultures.

Charcoal Agar device is used in a diagnostic workflow to aid clinicians in determining potential treatment option for patients suspected of having bacterial infections, the device is for professional use only, is not automated nor is a companion diagnostic

Summary and Explanation

Pertussis infection is severe and highly contagious respiratory illness also known as whooping cough. Gram-negative bacteria *Bordetella pertussis* and *Bordetella parapertussis* infections especially in infants and young children can lead to different complications such as pneumonia, weight loss, dehydration, and seizures¹. In older children and adults, this infection is often mild or subclinical and thus it is frequently unrecognized². Mooi et al. reported that resurgence in pertussis infection is associated with increased toxin production by *Bordetella pertussis* strains. *B. pertussis* infection is the cause of 13–20 % of prolonged cough in adults and adolescents. This infection often remains unrecognized, providing a reservoir of infection for unvaccinated or incompletely vaccinated infants³. *B. pertussis* is responsible for producing frequent outbreaks around the world every 3 to 5 years. For example, in 2008, 195,000 deaths occurred worldwide².

Haemophilus influenzae is a Gram-negative coccobacillus that can cause respiratory tract infection for example pneumonia, epiglottitis, sinusitis in both adults and children. There are two main types of *H. influenzae* on basis of the capsule. One type is encapsulated *H. influenzae* and the other is nontypeable *H. influenzae* (NTHi). The NTHi cannot be typed by the serum agglutination test. Detection of *H. influenzae* by bacterial culture is cheap and facilitates the detection of antibiotic resistance⁴. The nasopharynx is the niche for *Haemophilus influenzae* including other respiratory bacteria. *Haemophilus influenzae* is the main bacterial pathogen associated with acute lower respiratory infections (ALRIs) in children below 5 years of age⁵.

Principle of Method

'Lab-Lemco' powder and peptone provide the nutritional sources within the medium. Starch and charcoal bacteriological are neutralizers to toxins that effect *Bordetella pertussis*. Sodium chloride supplies electrolytes and maintains the osmotic equilibrium. Nicotinic acid is an essential growth factor required by *Bordetella pertussis*. Agar is solidifying agent.

Typical Formula

	grams per litre
'Lab-Lemco' powder	10.0
Peptone	10.0
Soluble starch	10.0
Activated carbon	4.0
Sodium chloride	5.0
Nicotinic acid	0.001
Agar	12.0

Materials Provided

500g of dehydrated Charcoal Agar Base that yields approximately 9.8L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Selective supplements (SR0082E)
- Petri dish

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 8°C.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.

- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B 6 and ID 5.

Procedure

Suspend 51g in 1 litre of distilled water. Bring to the boil to dissolve completely. Sterilise by autoclaving at 121°C for 15 minutes. Cool to 50°C, add 10% of defibrinated blood and mix gently. The medium is made selective for the isolation of *Bordetella pertussis* and *Bordetella parapertussis* by the addition of Bordetella Selective Supplement SR0082E.

To one vial add 2ml of sterile distilled water and dissolve the contents completely. Add this solution to 500ml of sterile, molten Charcoal Agar, cooled to 50°C, together with 10% v/v defibrinated horse blood SR0050. Mix well before pouring into sterile Petri dishes. For *Haemophilus influenzae*, omit the selective agents and convert to 'chocolate' agar

Interpretation

Once the media is reconstituted, *Bordetella pertussis* (NCTC 8605 and ATCC® 12742™), *Bordetella parapertussis* ATCC® 12742™ and *Haemophilus influenzae* ATCC® 35056™ will show as 0.25 - 1mm grey colonies. *Bordetella pertussis* ATCC® 8467™ will show as pinpoint – 1 mm grey colonies.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: ≥5 days @ 37°C

Positive Controls	
Inoculum level 10-100 cfu Colony count is ≥ 50% of the control medium count	
<i>Bordetella pertussis</i> NCTC 8605	0.25-1mm grey colonies
<i>Bordetella pertussis</i> ATCC® 8467™	Pinpoint – 1mm grey colonies
<i>Bordetella pertussis</i> ATCC® 12742™	0.25-1mm grey colonies
<i>Bordetella parapertussis</i> NCTC 10521	0.25-1mm grey colonies
<i>Haemophilus influenzae</i> ATCC® 35056™	0.25-1.5mm grey colonies

Limitations

Stuart's transport medium or similar formulation media should not be used for Bordetella-containing specimens. Two pernasal swabs should be taken from each patient, one through each nostril. Make sure charcoal remains in suspension when dispensing the medium by gently swirling the flask, Lysed horse blood is used in transport medium but whole blood is used in isolation medium. Most naso- pharyngeal flora is inhibited by cephalixin but *Pseudomonas aeruginosa* and some fungi may grow through. Amphotericin B can be added (12µg/ml) as an antifungal agent to prevent the growth of filamentous fungi, However, this level of amphotericin B can be inhibitory to *Bordetella pertussis* and should not be used routinely.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of *Bordetella pertussis*, *Bordetella parapertussis* and *Haemophilus influenzae* is confirmed by the inclusion of well-characterised isolates in the QC processes performed as part of the manufacture of each batch of the device. The precision of Charcoal Agar Base (CM0119B) was demonstrated by an overall pass rate of 100% obtained for the product over 2 years of testing (February 2019 – October 2021; 7 batches). This shows that the performance is reproducible.

The devices are tested in-house as part of the QC process. For target organisms, when using 10-100 cfu inoculum of *Bordetella pertussis*, *Bordetella parapertussis* and *Haemophilus influenzae* and incubating the device at 37°C for up to 5 days, the user can recover organisms with the appropriate colony size and morphology as listed in this document.

Bibliography

1. Relich, R. F., Leber, A., Young, S., Schutzbank, T., Dunn, R., Farhang, J., & Uphoff, T. S. 2019. "Multicenter Clinical Evaluation of the Automated Aries Bordetella Assay. ." *Journal of clinical microbiology*.
2. Dominguez, D. C. 2020. "A Profile of the Simplex™ Bordetella Direct Assay for the Detection and Differentiation of Bordetella Pertussis and Bordetella Parapertussis in Nasopharyngeal Swabs." *Expert Review of Molecular Diagnostics* 889-894.
3. Grogan JA, Logan C, O'Leary J, Rush R, O'Sullivan N. 2011. " Real-time PCR-based detection of Bordetella pertussis and Bordetella parapertussis in an Irish paediatric population. ." *J Med Microbiol.* 722-729. Group, Pneumonia Etiology Research for Child Health (PERCH) Study. 2019. "Causes of severe pneumonia requiring hospital admission in children without HIV infection from Africa and Asia: the PERCH multi-country case-control study. ." *Lancet (London, England)*.
4. Fan X, Liu X, Ji L, Cai D, Jiang J, Zhu J, Sun A, Yan J. 2018. "Epidemiological analysis and rapid detection by one-step multiplex PCR assay of Haemophilus influenzae in children with respiratory tract infections in Zhejiang Province, China. ." *BMC Infect Dis.* 414.
5. Muhandule Birindwa, A., Gonzales-Siles, L., Nordén, R., Geravandi, S., Tumusifu Manegabe, J., Morisho, L., Saili Mushobekwa, S., Andersson, R., & Skovbjerg, S. 2020. "High bacterial and viral load in the upper respiratory tract of children in the Democratic Republic of the Congo." *PLoS One*.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Consult instructions for use or consult electronic instructions for use
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment

	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of issue and modifications introduced
2.0	2023-12-13

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Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.