www.thermofisher.com

Thermo
SCIENTIFIC

- Quality control organisms and reference media (CM0007 and CM0071).
- Suitable tubes for creating agar slopes.

Triple Sugar Iron Agar (ISO)

EN

REF CM0277B

Intended Use

Triple Sugar Iron Agar (CM0277B) is a medium used for the differentiation of Enterobacteriaceae from clinical samples such as faecal, pus and exudates, by carbohydrate fermentations and hydrogen sulphide production.

The device is used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The device is for professional use only, is not automated, nor is it a companion diagnostic.

Summary and Explanation

A composite medium for the differentiation of Enterobacteriaceae according to their ability to ferment lactose, sucrose and glucose, and to produce hydrogen sulphide.

In addition to performing most of the functions of Kligler Iron Agar, the sucrose content of Triple Sugar Iron Agar permits the recognition and exclusion of sucrose-fermenting species. These organisms may ferment lactose slowly or not at all during the incubation period, but they readily metabolise sucrose. Some *Proteus* and other species may give similar reactions to salmonellae and shigellae and it is necessary to distinguish them by their ability to hydrolyse urea. For this reason, Triple Sugar Iron Agar should be used in parallel with Urea Broth or Urea Agar.

Principle of Method

Phenol red is incorporated for detecting carbohydrate fermentation which can be seen as a colour change from orange/red to yellow in the medium. In the case of oxidative decarboxylation of the peptone, the pH of the media will rise and a colour change from orange/red to deep red will be seen. Production of hydrogen sulfide is indicated by a black colour in the butt of the tube.

Typical Formula

	grams per litre
'Lab Lemco' Powder	3.0
Yeast Extract	3.0
Peptone	20.0
Sodium Chloride	5.0
Lactose	10.0
Sucrose	10.0
Glucose	1.0
Ferric Citrate	0.3
Sodium Thiosulphate	0.3
Phenol Red	0.024
Agar	12.0

Materials Provided

CM0277B: 500g of Triple Sugar Iron Agar.

500g of dehydrated Triple Sugar Iron Agar yields approximately 7.7L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers.
- Incubators.

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store medium between 2°C and 25°C.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For *in vitro* diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between uses to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimens should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B 14, S 7 and Q 5.

Procedure

Suspend 65g in 1 litre of distilled water. Bring to the boil to dissolve completely. Mix well and distribute into tubes. Sterilise by autoclaving at 121°C for 15 minutes. Allow the medium to set in sloped form with a butt about 2.5 cm in depth.

- Pick a single colony from the surface of a selective plating medium and smear a MacConkey Agar (CM0007) plate. Incubate for 18 hours at 37°C and inoculate two separate tubes of media from one single isolated colony:
 - Triple Sugar Iron Agar - smear the slope and stab the butt.
 - Urea Broth Base (CM0071) with added Urea Solution (SR0020).
- Incubate at 35°C.
- Examine the Urea Broth tube after 5 hours and again after 18 hours incubation. Discard tubes showing a red or pink coloration, which is due to urea hydrolysis by *Proteus* or other organisms.
- Where there is no urea hydrolysis, examine the Triple Sugar Iron Agar tubes after 18 hours.

Interpretation

Typical reactions observed on CM0277B are:

Organism	Butt	Slope	H ₂ S
<i>Enterobacter aerogenes</i>	AG	A	-
<i>Enterobacter cloacae</i>	AG	A	-
<i>Escherichia coli</i>	AG	A	-
<i>Proteus vulgaris</i>	AG	A	+
<i>Morganella morganii</i>	A or AG	NC or ALK	-
<i>Shigella dysenteriae</i>	A	NC or ALK	-
<i>Shigella sonnei</i>	A	NC or ALK	-
<i>Salmonella</i> Typhi	A	NC or ALK	+
<i>Salmonella paratyphi</i>	AG	NC or ALK	-
<i>Salmonella</i> Enteritidis	AG	NC or ALK	+
<i>Salmonella</i> Typhimurium	AG	NC or ALK	+

AG = acid (yellow) and gas formation

A = acid (yellow)

NC = no change

ALK = alkaline (red)

+ = hydrogen sulphide (black)

- = no hydrogen sulphide (no black)

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Organism	Butt	Slope	H ₂ S
<i>Escherichia coli</i> ATCC® 8739™	AG	A	-
<i>Escherichia coli</i> ATCC® 11775™	AG	A	-
<i>Shigella sonnei</i> ATCC® 25931™	A	NC	-
<i>Salmonella enteritidis</i> ATCC® 13076™	AG	NC	+
<i>Salmonella typhimurium</i> ATCC® 14028™	AG	NC	+
<i>Salmonella abony</i> NCTC 6017	AG	NC	+
<i>Salmonella nottingham</i> NCTC 7832	AG	NC	+
<i>Proteus hauseri</i> ATCC® 13315™	A	A	+
<i>Enterobacter aerogenes</i> ATCC® 13048™	AG	A	-
<i>Pseudomonas aeruginosa</i> ATCC® 9027™	ALK	ALK	-
<i>Staphylococcus aureus</i> ATCC® 6538™	A	A	-

Limitations

Identifications are presumptive and colonies should be confirmed using appropriate methods.

It is now thought that Triple Sugar Iron Agar is not suitable for the detection of hydrogen sulphide production by sucrose-fermenting organisms, such as some *Citrobacter* and *Proteus* species, in which the sucrose fermentation masks the hydrogen sulphide indicator in the medium¹.

Performance Characteristics


Accuracy has been demonstrated through review of internal QC data. Testing of well-characterised isolates in the QC processes is performed as part of the manufacture of each batch of the devices.

Bibliography

- Bulmash J. M. and Fulton M. (1966) *J. Bact.* 88. 1813.

Symbol Legend

Symbol	Definition
REF	Catalogue number
IVD	In Vitro Diagnostic Medical Device
LOT	Batch code

	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke,
RG24 8PW, United kingdom

For technical assistance please contact your local distributor.

Revision information

Version	Date of modifications introduced
2.0	2023-04-03

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© 2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc.