Mueller Hinton Agar Base

REF CM0337B, CM0337K, CM0337M, CM0337R, CM0337T & CM0337W

EN

Intended Use

Mueller Hinton Agar (CM0337B/K/M/R/T/W) devices are antimicrobial susceptibility testing agars recommended for disc diffusion and also recommended for agar dilution MIC testing.

Mueller Hinton Agar (CM0337B/K/M/R/T/W) devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The devices are for professional use only, are not automated and nor are they companion diagnostics.

Summary and Explanation

The rapid global emergence of antimicrobial resistance threatens the effective treatment of patients with infections. A number of bacteria have been designated as urgent, serious, and worrying threats¹.

Assays for detection of resistance can be performed using genotypic or phenotypic techniques. The minimum inhibitory concentration (MIC) serves as the 'gold standard' phenotypic resistance determination technique, with the organism classified as susceptible, intermediate, or resistant to the antimicrobial agent based on clinical MIC breakpoints. MIC assays, however, can be labour-intensive, time-consuming and require experience to interpret. Disc diffusion is widely used in routine clinical microbiology laboratories as it is low cost, simple to perform, versatile and almost all antimicrobial agents and pathogens can be tested. Disc diffusion is accurate and reproducible when carried out in strict accordance with guidelines/recommendations (such as those published by the European Committee on Antimicrobial Susceptibility Testing (EUCAST) or Clinical and Laboratory Standards Institute (CLSI)). Mueller Hinton media have been recommended in CLSI standards for more than 25 years and the agar is now also recommended by the Europe-wide EUCAST standard, as well as ISO. Mueller Hinton Agar (CM0337B/K/M/R/T/W) devices are therefore key tools for laboratories to support correct MIC and disc diffusion testing of microorganisms.

Principle of Method

Mueller Hinton Agar (CM0337B/K/M/R/T/W) devices contain dehydrated infusion from beef and casein hydrolysate to provide nitrogen, vitamins, carbon, and amino acids. Soluble starch is added to absorb toxic metabolites and its hydrolysis yields dextrose, which serves as an energy source. Agar is added as the solidifying agent.

Mueller Hinton Agar (CM0337B/K/M/R/T/W) devices are stable, reproducible, and standardized solid media recommended for testing the susceptibility of bacteria to antimicrobial agents using disc diffusion techniques based on the Kirby-Bauer method, such as those recommended by CLSI or EUCAST, or by MIC agar dilution.

Un-supplemented Mueller Hinton Agar (CM0337B/K/M/R/T/W) devices are recommended by CLSI and EUCAST for antimicrobial susceptibility testing of non-fastidious strains. Supplements (such as defibrinated blood and nicotinamide adenine dinucleotide (NAD)) can be added to the base medium (CM0337B/K/M/R/T/W) to make it suitable for fastidious strains.

Typical Formula

	<u>grams per litre</u>
Beef, dehydrated infusion from	300.0
Casein hydrolysate	17.5
Soluble starch	1.5
Agar	17.0

Materials Provided

CM0337B: 500g of Mueller Hinton Agar base.
 CM0337K: 25kg of Mueller Hinton Agar base.
 CM0337M: 3 x 2.28kg of Mueller Hinton Agar base.
 CM0337R: 2.5kg of Mueller Hinton Agar base.
 CM0337T: 5kg of Mueller Hinton Agar base.
 CM0337W: 3 x 3.8kg of Mueller Hinton Agar base.

500g of dehydrated Mueller Hinton Agar yields approximately 13.1L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Diluents
- Sterile containers
- Antimicrobial discs / agents

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.
- Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

- Ensure adequate ventilation. Wear personal protective equipment / face protection. Avoid contact with skin, eyes or clothing.
- Avoid ingestion and inhalation. Avoid dust formation.
- If on skin wash off immediately with plenty of water for at least 15 minutes. Get medical attention immediately if symptoms occur.
- If in eyes, rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes
If inhaled, remove to fresh air. Get medical attention immediately if symptoms occur.
If ingested, clean mouth with water and drink afterwards plenty of water. Get medical attention if symptoms occur.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.
- Antimicrobial susceptibility testing must be carried out in strict accordance with current guidelines/recommendations in order to maintain acceptable levels of accuracy.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Microorganisms isolated from clinical samples are tested on Mueller Hinton Agar (CM0337) devices. Clinical samples are not inoculated directly onto the medium.

Procedure

Suspend 38g in 1 litre of distilled water. Bring to the boil to dissolve completely. Sterilise by autoclaving at 121°C for 15 minutes. Allow to cool, mix well and pour into sterile containers.

Interpretation

Antibiotic susceptibility tests should be performed and interpreted in accordance with current guidelines (such as those published by EUCAST or CLSI).

Quality Control

It is the responsibility of the user to perform Quality Control testing, taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Organism	Antimicrobials Tested (Disc content in µg)	Disc diffusion zone diameter (mm)
<i>Staphylococcus aureus</i> ATCC® 29213™	Penicillin (P1)	12 - 18
	Cefoxitin (FOX30)	24 - 30
	Ciprofloxacin (CIP5)	21 - 27
	Erythromycin (E15)	23 - 29
	Gentamicin (CN10)	19 - 25
	Linezolid (LZD10)	21 - 27

	Tetracycline (TE30)	23 - 31
<i>Staphylococcus aureus</i> ATCC®25923™	Erythromycin (E15)	22 - 30
	Tetracycline (TE30)	24 - 30
	Ciprofloxacin (CIP5)	22 - 30
	Amoxicillin/clavulanic acid (AMC30)	28 - 36
	Ampicillin-sulbactam (SAM20)	29 - 37
	Linezolid (LZD30)	25 - 32
	Cefoxitin (FOX30)	23 - 29
	Gentamicin (CN10)	19 - 27
	Penicillin (P10)	26 - 37
<i>Escherichia coli</i> ATCC® 25922™	Ampicillin (AMP10)	15 - 22
	Chloramphenicol (C30)	21 - 27
	Gentamicin (CN10)	19 - 26
	Trimethoprim/sulfamethoxazole (SXT25)	23 - 29
	Amoxicillin/clavulanic acid (AMC30)	18 - 24
	Ciprofloxacin (CIP5)	30 - 40
	Sulfisoxazole (SF300)	15 - 23
	Tetracycline (TE30)	18 - 25
	Cefotaxime (CTX5)	25 - 31
Tigecycline (TGC15)	20 - 27	
<i>Escherichia coli</i> ATCC®35218™	Amoxicillin/clavulanic acid (AMC30)	17 - 22
<i>Pseudomonas aeruginosa</i> ATCC® 27853™	Gentamicin (CN10)	16 - 21
	Imipenem (IPM10)	20 - 28
	Aztreonam (ATM30)	23 - 29
	Ciprofloxacin (CIP5)	25 - 33
	Tobramycin (TOB10)	19 - 25
	Ceftazidime (CAZ10)	21 - 27
	Piperacillin/ Tazobactam (TZP36)	23 - 29
<i>Staphylococcus aureus</i> ATCC®43300™	Cefoxitin (FOX30)	≤ 21
<i>Staphylococcus aureus</i> NCTC12493™	Cefoxitin (FOX30)	14 - 20
<i>Enterococcus faecalis</i> ATCC®33186	Sulfamethoxazole/trimethoprim (SXT25)	≤ 20
<i>Enterococcus faecalis</i> ATCC®29212™	Ampicillin (AMP2)	15 - 21
	Imipenem (IPM10)	24 - 30

	Linezolid (LZD10)	19 - 25
	Nitrofurantoin (F100)	18 - 24
	Trimethoprim (W5)	24 - 32
	Sulfamethoxazole/trimethoprim (SXT25)	26 - 34
	Vancomycin (VA5)	10 - 16
<i>Streptococcus pneumoniae</i>* ATCC®49619™	Vancomycin (VA5)	17 - 23
<i>Haemophilus influenzae</i>* ATCC®49247™	Ampicillin (AMP2)	6 - 12
<i>Haemophilus influenzae</i>* ATCC®49766™	Cefixime (CFM5)	29 - 35

Limitations

Some strains with specific growth requirements may not grow.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct MIC or disc diffusion results are confirmed by the inclusion of well-characterised isolates and antimicrobials in the QC processes performed as part of the manufacture of each batch of the devices. The precision of Mueller Hinton Agar (CM337) devices was demonstrated by an overall pass rate of 100 % obtained for the product over one month of testing (17.06.2022 – 26.07.2022; 10 batches). This shows that the performance is reproducible.

These media devices have been tested in-house as part of the QC process since they were first manufactured; Mueller Hinton Agar (CM0337B/K/M/R/T/W) devices have been on the market since 1997. When using optimum inoculum levels and incubation conditions, the user can achieve zone sizes within the range stated in the specification.

Bibliography

- Centers for Disease Control and Prevention (CDC). 2013. 'Office of Infectious Disease. Antibiotic resistance threats in the United States'. <http://www.cdc.gov/drugresistance/threat-report-2013>.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use

	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

©2023 Thermo Fisher Scientific Inc. All rights reserved.

ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of modifications introduced
2.0	2022-11-15

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© 2023 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.