Hektoen Enteric Agar

[REF] CM0419K, CM0419W, CM0419B, CM0419R, CM0419T, CM0419V

Intended Use

Hektoen Enteric Agar is a differential medium for the isolation of *Shigella* and *Salmonella* species from faecal samples. Can also be used for testing food samples.

The device is for professional use only, is not automated, and nor is it a companion diagnostic.

Summary and Explanation

Hektoen Enteric Agar was developed by King and Metzger¹ at the Hektoen Institute. Taylor and Schelhaut² found the medium to be of value for the differentiation of salmonellae and shigellae and for better growth of the latter. Hektoen Enteric Agar is one of the culture media recommended in the Bacteriological Analytical Manual³ for the isolation of salmonella from foods. It is also recommended in ISO 21567⁴ for the isolation of shigella from foods.

Principle of Method

This medium contains greater quantities of peptone to offset the inhibitory effects of bile salts on *Shigella* species. The additional carbohydrates (sucrose and salicin) give better differentiation than lactose alone and the lower toxicity of the double indicator, acid fuchsin and bromothymol blue, improves recovery. The increased lactose content helps early recognition of slow lactose fermenting organisms. Thiosulphate and ferric citrate are present to detect hydrogen sulphide producing organisms. Bile salts are present to inhibit Gram-positive bacteria.

Typical Formula

	grams per litre
Proteose peptone	12.0
Yeast extract	3.0
Lactose	12.0
Sucrose	12.0
Salicin	2.0
Bile salts No.3	9.0
Sodium chloride	5.0
Sodium thiosulphate	5.0
Ammonium ferric citrate	1.5
Acid fuchsin	0.1
Bromothymol blue	0.065
Agar	14.0

Physical Characteristics of Reconstituted Media

Colour	Olive yellow to ochre yellow
Clarity	Transparent
Fill weight	17g ± 5%
pH	7.5 ± 0.2

Materials Provided

- **CM0419B:** 500g dehydrated culture media
- **CM0419R:** 2.5Kg dehydrated culture media
- **CM0419T:** 5Kg dehydrated culture media
- **CM0419W:** 7.6Kg dehydrated culture media
- **CM0419V:** 7 x 3.8Kg dehydrated culture media
- **CM0419K:** 25Kg dehydrated culture media

Materials Required but Not Supplied

- (1) Inoculating loops, swabs, collection containers
- (2) Incubators
- (3) Quality control organisms
- (4) Petri dishes

Storage

This product is ready to use and no further preparation is necessary. Store product in its original packaging at 2–12°C until used. Store away from light. Allow product to equilibrate to room temperature before use. Do not incubate prior to use.

Warnings and Precautions

This product is for *in vitro* diagnostic use and should only be used by trained individuals. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products. It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully.

Inspect the product packaging before first use. Do not use the product if there is any visible damage to the packaging or plates. Do not use the product beyond the stated expiry date. Do not use the device if signs of contamination are present. Do not use the device if the colour has changed or there are other signs of deterioration.

For professional use only. Refer to the Safety Data Sheet (SDS) at www.thermofisher.com

Specimen Collection, Handling and Storage

Specimens should be collected and handled following the recommended guidelines.

Procedure

Refer to local protocols and guidelines.

Suspend 76g of the medium in 1 litre of distilled water and soak for 10 minutes. Heat gently and allow to boil for a few seconds to dissolve the agar. DO NOT AUTOCLAVE. Cool to 50°C and pour plates.

The procedure followed may depend upon the specimen type.

Clinical samples

- (1) Inoculate the medium with fresh stool samples suspended in e.g. Ringers solution or inoculate directly with a swab.
- (2) Streak the inoculum to obtain well separated colonies.
- (3) Incubate aerobically for 18 to 24 hours at 36°C ± 1°C.
- (4) Examine plates for typical colonies.

Food samples

Isolation of salmonella and shigella may be performed after enrichment in selective broths. Refer to the appropriate standard for the methodology.

Interpretation of Results

Shigella - green, moist, raised colonies.

Salmonella – blue-green, transparent colonies with or without black centres.

Coliforms (rapid lactose/sucrose/salicin fermenters) – salmon pink to orange colonies surrounded by a zone of bile precipitation.

Quality Control

This medium can be tested with the following strains:

Incubation Conditions: 18 – 24 h @ 36 ± 1°C, aerobic

Positive Controls	
Inoculum of 50 – 120 colony forming units (cfu).	
<i>Salmonella</i> Typhimurium ATCC®14028	2 – 3 mm, transparent colonies, black centres.
Colony count shall be greater than or equal to 50% of the control medium (Tryptone Soya Agar).	
Inoculum of 10 ³ – 10 ⁴ cfu.	

<i>Shigella flexneri</i> ATCC®12022	Good growth, green shiny colonies.
Negative Controls	
Inoculum of $\geq 10^4$ cfu.	
<i>Enterococcus faecalis</i> ATCC®29212	Complete inhibition (≤ 10 cfu).
Inoculum of 50 – 120 cfu	
<i>Escherichia coli</i> ATCC®25922	No growth or 0.5 – 2 mm orange colonies.

Note:

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

Performance

Performance was evaluated using 38 bacterial strains including the following: *Salmonella* spp., *Shigella* spp., *E. coli*, *Proteus* spp., *Pseudomonas aeruginosa*, *Citrobacter* spp., *Enterococcus faecalis* and *Providencia* spp. All organisms gave expected growth characteristics according to the current product specification.

Limitations

All identifications are presumptive and should be confirmed using appropriate methods. Due to variation in nutritional requirements some strains may grow poorly or fail to grow on this medium. A second enteric medium may be inoculated in parallel to improve the chances of recovery of the target organism from a mixed population. *Proteus* spp. may resemble salmonella or shigella species.

Bibliography

- (1) King S. and Metzger W. I. (1968) *Appl. Microbiol.* 16. 577-561.
- (2) Taylor W. I. and Schelhaut D. (1971) *Appl. Microbiol.* 21. 32-37.
- (3) US Food and Drug Administration. *Bacteriological Analytical Manual (BAM)* Chapter 5 *Salmonella* 07/20/2020. On line. AOAC International, Gaithersburg, Md
- (4) ISO 21567:2004 *Microbiology of food and animal feeding stuffs - Horizontal method for the detection of Shigella spp.*

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use

	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Manufacturer
	Authorized representative in the European Community/European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Made in the United Kingdom

The ATCC Licensed Derivative® Emblem, the ATCC Licensed Derivative® word mark, and the ATCC catalog marks are trademarks of ATCC. Thermo Fisher Scientific Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

© 2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. ATCC® is a trademark of American Type Culture Collection. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Oxoid Deutschland GmbH, Am Lippeglaeis 4-8, 46483 Wesel, Germany

For technical assistance please contact your local distributor.

Version	Date of issue and modifications introduced
1.0	2022-11-08. New document

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Oxoid Deutschland GmbH, Am Lippegelacis 4-8, 46483