www.thermofisher.com

X.L.D. Medium

EN

REF CM0469B/R/T/W*

*This instructions for use (IFU) document is intended to be read in conjunction with the IFU for X.L.D. Medium (product code: PO0164A)

Intended Use

The XLD Agar is a selective medium for the isolation of *Salmonella* and *Shigella* species from faecal samples. The devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections and are for professional use only and not intended for self-testing. The devices are not automated nor are a companion diagnostic.

Summary and Explanation

Xylose Lysine Desoxycholate (XLD) Agar was originally formulated by Taylor¹ for the isolation and identification of *shigellae* from stool specimens. It has since been found to be a satisfactory medium for the isolation and presumptive identification of both *salmonellae* and *shigellae*².

The recovery of *Shigella* spp. has previously been neglected despite the high incidence of shigellosis due to inadequate isolation media³.

The sensitivity and selectivity of XLD Agar exceeds that of the traditional plating media e.g. Eosin Methylene Blue, *Salmonella-Shigella* and Bismuth Sulphite agars, which tend to suppress the growth of *shigellae*.

Many favourable comparisons between XLD Agar and these other media have been recorded in the literature^{2,4,5,6,7}.

The recovery of *salmonellae* and *shigellae* is not obscured by profuse growth of other species³, therefore XLD Agar is ideal for the screening of samples containing mixed flora and suspected of harboring enteric pathogens such as medical specimens or food products.

Principle of Method

The XLD Agar relies on xylose fermentation, lysine decarboxylation (XLD) and the production of hydrogen sulphide for the primary differentiation of *shigellae* and *salmonellae* from non-pathogenic bacteria.

Rapid xylose fermentation is almost universal amongst enteric bacteria, except for members of the *Shigella*, *Providencia* and *Edwardsiella* genera. Xylose is present so that *Shigella* spp. may be identified by a negative reaction. *Salmonellae* are differentiated from non-pathogenic xylose fermenters by the incorporation of lysine in the medium. *Salmonellae* exhaust the xylose and decarboxylate the lysine, thus altering the pH to alkaline and mimicking the *shigella* reaction. Phenol red is included as a pH indicator. The presence of *salmonellae* and *Edwardsiella* spp. is differentiated from that of *shigellae* by a hydrogen sulphide indicator.

The high acid level produced by fermentation of lactose and sucrose prevents lysine positive coliforms from reverting the pH to an alkaline value and non-pathogenic hydrogen sulphide producers do not decarboxylate lysine. The acid level also prevents blackening by these microorganisms until after the 24-hour examination for pathogens.

Sodium desoxycholate is incorporated as an inhibitor in the medium. The concentration used inhibits Gram-positive bacteria and allows for the inhibition of coliforms without

Thermo

SCIENTIFIC

decreasing the ability to support *shigellae* and *salmonellae*. Yeast extract is present as a nutrient source.

Typical Formula

	grams per litre
Yeast extract	3
L-Lysine HCl	5
Xylose	3.75
Lactose	7.5
Sucrose	7.5
Sodium deoxycholate	1
Sodium chloride	5
Sodium thiosulphate	6.8
Ferric ammonium citrate	0.8
Phenol red	0.08
Agar	12.5

Materials Provided

CM0469B – 500g of dehydrated X.L.D. powder that yields approximately 9.4L after reconstitution

CM0469R – 2.5kg of dehydrated X.L.D. powder that yields approximately 47L after reconstitution

CM0469T – 5KG of dehydrated X.L.D. powder that yields approximately 94L after reconstitution

CM0469W – 2x 5.3kg of dehydrated X.L.D. powder that yields approximately 200L after reconstitution

Materials Required but Not Supplied

- (1) Inoculating loops, swabs, collection containers
- (2) Incubators
- (3) Quality control organisms
- (4) Petri dish

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.

- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B 29 (Public Health England, 2020).

Procedure

Suspend 53g in 1 litre of distilled water. With frequent agitation, bring to the boil to dissolve completely. Cool to 50°C. Mix well and pour into sterile Petri dishes.

Interpretation

The presence of colony growth indicates the presence of Gram-negative enteric pathogens.

Degradation of xylose, lactose and sucrose to acid causes phenol red indicator to change its colour to yellow.

Bacteria that decarboxylate lysine to cadaverine can be recognized by the appearance of a red colouration around the colonies due to an increase in pH.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 21-27 h @ 37 ± 2°C aerobic

Inoculum level: 10 ³ -10 ⁶ cfu	
<i>Salmonella abony</i> NCTC 6017	1-3 mm red colonies, black centre
<i>Salmonella</i> Enteritidis ATCC® 13076™	1-2 mm red colonies, black centre
<i>Salmonella</i> Typhimurium ATCC® 14028™	1-2 mm red colonies, black centre
<i>Salmonella</i> Virchow NCTC 5742	1-2 mm red colonies, black centre

<i>Salmonella arizonae</i> ATCC® 13314™	1-3 mm red colonies, black centre
<i>Salmonella Nottingham</i> NCTC 7832	1-3 mm red colonies, black centre
<i>Shigella sonnei</i> ATCC® 9290™	0.5-7 mm irregular/smooth red colonies
<i>Shigella flexneri</i> ATCC® 12022™	0.5-2 mm irregular, red colonies
Inoculum level: 10 ⁴ -10 ⁶ cfu	
<i>Shigella sonnei</i> ATCC® 25931™	0.5-7 mm irregular/smooth red colonies
Inoculation with pure cultures Inoculum level: 10-100 cfu <i>Pseudomonas aeruginosa</i> ATCC® 9027™ colony count is ≤ 90% of the control medium count.	
<i>Pseudomonas aeruginosa</i> ATCC® 9027™	No growth or 0.5-2 mm red colonies
<i>Proteus mirabilis</i> ATCC® 12453™	0.5-2 mm orange/red colonies, with or without black centre, no swarming
<i>Proteus mirabilis</i> ATCC® 29906™	0.5-2 mm orange/red colonies, with or without black centre, no swarming
<i>Serratia marcescens</i> ATCC® 8100™	1-2 mm orange/yellow colonies
<i>Citrobacter freundii</i> ATCC® 8090™	0.5-2 mm yellow colonies
<i>Klebsiella pneumoniae</i> ATCC® 29665™	2-4 mm yellow, mucoid colonies
Negative Controls Inoculum level: 10 ⁴ -10 ⁶ cfu	
<i>Staphylococcus aureus</i> ATCC® 6538™	No growth
<i>Escherichia coli</i> ATCC® 11775™	No growth or 0.5-4mm yellow colonies

Testing performed in accordance with ISO11133:2014

Incubation Conditions: 21-27 h @ 37 ± 2°C aerobic. Inoculum level: 50-120 cfu Colony count is ≤ 90% of the control medium count.	
<i>Salmonella</i> Enteritidis ATCC® 13076™	1-3 mm red colonies, black centre
<i>Salmonella</i> Typhimurium ATCC® 14028™	1-3 mm red colonies, black centre

Inoculum level: 10 ⁴ -10 ⁵ cfu	
<i>Escherichia coli</i> ATCC® 8739™	No growth or 0.5-4 mm yellow colonies
<i>Escherichia coli</i> ATCC® 25922™	No growth or 0.5-4 mm yellow colonies
Inoculum level: 10 ⁴ -10 ⁶ cfu	
<i>Enterococcus faecalis</i> ATCC® 29212™	No growth
<i>Enterococcus faecalis</i> ATCC® 19433™	No growth

Testing performed in accordance with current CLSI M22 A

Incubation Conditions: 18-24 h @ 35°C Medium is challenged with 10-100 cfu Colony count is ≥70% of the control medium count.	
<i>Shigella flexneri</i> ATCC® 12022™	0.5-2 mm irregular, red colonies
<i>Salmonella</i> Typhimurium ATCC® 14028™	1-2mm red colonies, black centre
Inoculum level: 10 ⁴ -10 ⁶ cfu	
<i>Enterococcus faecalis</i> ATCC® 29212™	No growth
Inoculum level: 10 ⁴ -10 ⁶ cfu Inhibited strains shall produce no growth or at least a 1 log (10) reduction	
<i>Escherichia coli</i> ATCC® 22922™	No growth or 0.5-4mm yellow colonies

Limitations

Identifications are presumptive and should be confirmed using appropriate biochemical/serological methods. *Pseudomonads* may grow on X.L.D. medium and may appear as red colonies, i.e. false positives, although the oxidase reaction can be used to differentiate them. Some strains of *Proteus* spp. may also appear as false positives and some may develop black centres.

Other non-target organisms may also be able to grow. Some serotypes such as *Salmonella Gallinarum*., *Salmonella Paratyphi* and *Salmonella Pullorum* may not develop black centres.

Lactose positive *salmonella* may appear as yellow colonies with or without black centres.

Due to variation in nutritional requirements or sensitivity to selective agents some strains of the target organisms may be encountered that grow poorly or fail to grow on this medium.

Performance Characteristics

20 *Salmonella* species were spiked into minced beef, minced turkey, cottage cheese, whole liquid egg and ready prepared lasagne at a final dilution of 1-100 cfu/ml.

Performance	Overall	Minced Beef	Minced Turkey	Cottage Cheese	Whole Liquid Egg	Pre-prepared Lasagne
Sensitivity	99%	95%	100%	100%	100%	100%

Performance was also evaluated using 38 bacterial strains including the following; *Salmonella* spp., *Shigella* spp., *E. coli*, *Proteus* spp., *Pseudomonas aeruginosa*, *Citrobacter* spp., *Enterococcus faecalis* and *Providencia* spp. All organisms gave expected growth characteristics according to the current product specification.

Bibliography


1. Taylor W.I. (1965) Am. J. Clin. Path. 44(4): 471-475.
2. McCarthy M.D. (1966) N.Z. J. Med. Lab. Technol. 20. 127-131
3. Isenberg H.D., Kominos S. and Siegel M. (1969) Appl. Microbiol. 18. 656-659.
4. Taylor W. I. and Harris B. (1965) Am. J. Clin. Path. 44. 476-479.
5. Taylor W.I. and Schelhart D. (1968) Appl. Microbiol. 16. 1387-1392.
6. Taylor W. I. and Schelhart D. (1969) Appl. Microbiol. 18. 393-395.
7. Dunn C. and Martin W.J. (1971) Appl. Microbiol. 22. 17-22

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment

	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

For technical assistance please contact your local distributor.

Revision information

Version	Date of issue and modifications introduced
2.0	2023-06-09

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©2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.