
	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Doxycycline Discs 30 µg (DO30)

REF CT0018B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Oxoid Doxycycline Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs that contain specific amounts of the antimicrobial agent doxycycline. The discs are labelled on both sides with details of the antimicrobial (DO) and amount present (µg): DO30.

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Oxoid Doxycycline AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Doxycycline Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which doxycycline has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling, and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of FDA,³ CLSI M100,^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- *Streptococcus pneumoniae*

Gram-Negative

- Enterobacteriaceae
- *Acinetobacter* spp.

Species with published breakpoints according to current CLSI M100 literature (all associated with 30 µg):

Gram-Positive

- *Streptococcus pneumoniae*
- *Staphylococcus* spp.
- *Enterococcus* spp.

Gram-Negative

- Enterobacterales
- *Acinetobacter* spp.

Species with published breakpoints according to current EUCAST literature (all associated with 30 µg):

Gram-Positive

- N/A

Gram-Negative

- N/A

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Doxycycline AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Doxycycline AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Doxycycline AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Doxycycline AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with CT0018B	
Component Description	Material Description
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Pale beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Doxycycline AST Disc Reagents	
Reagent	Description of Function
Doxycycline	Yellow crystalline powder. Antimicrobial agent. The drug mechanism of action involves inhibition of translation by binding to the 16S rRNA portion of the ribosome 9, thus, preventing binding of tRNA to the RNA-30S bacterial ribosomal subunit, which is necessary for the delivery of amino acids for protein synthesis. As a result of the above actions, the initiation of protein synthesis by polyribosome formation is blocked. This stops the replication of bacteria and produces a bacteriostatic effect.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Doxycycline AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Doxycycline AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Raw data of zone sizes, in accordance with CLSI methodology, taken from the 3 most recent batches of Doxycycline (30 µg) DO30 Antimicrobial Susceptibility Test Discs (CT0018B).

CLSI Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3289651	<i>Escherichia coli</i>	ATCC 25922	MHA	18	24	21	20	20	20	-1 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC 25923	MHA	23	29	26	28	28	28	2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC 49619	MHA +5% SB	25	34	29.5	31	31	31	1.5 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3195822	<i>Escherichia coli</i>	ATCC 25922	MHA	18	24	21	20	20	20	-1 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC 25923	MHA	23	29	26	26	26	26	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC 49619	MHA +5% SB	25	34	29.5	32	32	32	2.5 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3138678	<i>Escherichia coli</i>	ATCC 25922	MHA	18	24	21	22	22	22	1 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC 25923	MHA	23	29	26	27	27	27	1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC 49619	MHA +5% SB	25	34	29.5	30	30	30	0.5 ± 0 (CV=0)

Clinical Performance Characteristics

In a 2011 study by Tripodi, Fortunato and Cuccurullo, the objective was to assess the efficacy of multiple antibiotics (alone or in combination) and their potential interactions, and compare the in vitro results with the clinical outcome of a cohort of heart transplant (HTx) recipients with pulmonary nocardiosis (PN).⁴ Over a 20-year period, Nocardia isolates were identified from 12 HTx recipients with PN. The in vitro activity of various antibacterial agents including doxycycline, alone or in combination, was assessed using disk-diffusion, minimal inhibitory concentration (MIC), and time-kill methodology. The in vitro results were compared with the clinical outcome of the patients. All isolates were susceptible to amikacin, netilmicin, and linezolid. Gentamicin, tobramycin, and streptomycin were active against 8, 7, and 2 isolates, respectively. Doxycycline and imipenem were active

against 8 of the 9 isolates. Notably, meropenem was active against only 2 isolates. Cotrimoxazole was active against 7 isolates. In conclusion, this study demonstrated that Oxoid Doxycycline AST Discs were effective in the clinical analysis of *Nocardia* spp., given the susceptibilities observed and that no out of specification or inconsistent results were reported.

A study by Reta, Wubie and Mekuria from 2017, aimed to determine the nasal colonization, associated factors, and antimicrobial susceptibility patterns of *Staphylococcus aureus* isolates among pre-school children.⁵ A total of 400 nasal swab samples were collected; 52 *S. aureus* isolates were detected from the samples; none of those was classified as MRSA. The isolates were then tested for susceptibility against several antibiotics including Oxoid Doxycycline (30 µg) DO30 AST Discs, which demonstrated 69.2% sensitivity. The findings of this investigation demonstrate that the Oxoid Doxycycline AST Discs are an effective tool to determine the antimicrobial susceptibility pattern of *S. aureus* isolates responsible for nasal colonization within a population, and they constitute a valid instrument to detect the development of any resistance, as well as prevent treatments failures.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

1. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
2. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version.
3. EUCAST Disk Diffusion Method Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing.
4. Tripodi, M, R Fortunato, and S Cuccurullo. 2011. "In Vitro Activity of Multiple Antibiotic Combinations against *Nocardia* : Relationship with a Short-Term Treatment Strategy in Heart Transplant Recipients with Pulmonary Nocardiosis," no. 2: 335–43. <https://doi.org/10.1111/j.1399-3062.2010.00588.x>.
5. Reta, Alemayehu, Moges Wubie, and Getnet Mekuria. 2017. "Nasal Colonization and Antimicrobial Susceptibility Pattern of *Staphylococcus Aureus* among Pre - School Children in Ethiopia." *BMC Research Notes*, 1–7. <https://doi.org/10.1186/s13104-017-3079-6>

© 2021 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of Symbols


Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date

	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision Information

Version	Date of issue and modifications introduced
1.0	2022-11-02. New document

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Literatura

- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Clinical & Laboratory Standards Institute (CLSI). Proučite najnovije izdanje
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Proučite najnovije izdanje.
- EUCAST Disk Diffusion Method Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing.
- Tripodi, M, R Fortunato, and S Cuccurullo. 2011. "In Vitro Activity of Multiple Antibiotic Combinations against Nocardia : Relationship with a Short-Term Treatment Strategy in Heart Transplant Recipients with Pulmonary Nocardiosis," no. 2: 335–43. <https://doi.org/10.1111/j.1399-3062.2010.00588.x>.
- Reta, Alemayehu, Moges Wubie, and Getnet Mekuria. 2017. "Nasal Colonization and Antimicrobial Susceptibility Pattern of Staphylococcus Aureus among Pre - School Children in Ethiopia." *BMC Research Notes*, 1–7. <https://doi.org/10.1186/s13104-017-3079-6>

© 2021 Thermo Fisher Scientific Inc. Sva su prava pridržana. Oznaka ATCC® zaštitni je znak Američke zbirke tipskih kultura (ATCC). Svi ostali zaštitni znakovi vlasništvo su društva Thermo Fisher Scientific Inc. i njegovih društava kćeri. Ove informacije nisu namijenjene poticanju upotrebe ovih proizvoda na bilo koji način koji bi mogao dovesti do povrede prava intelektualnog vlasništva drugih.

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Referencer

1. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
2. Clinical & Laboratory Standards Institute (CLSI). Se den seneste udgave
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Se den seneste version.
3. EUCAST Disk Diffusion Method Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing.
4. Tripodi, M, R Fortunato, and S Cuccurullo. 2011. "In Vitro Activity of Multiple Antibiotic Combinations against *Nocardia* : Relationship with a Short-Term Treatment Strategy in Heart Transplant Recipients with Pulmonary Nocardiosis," no. 2: 335–43. <https://doi.org/10.1111/j.1399-3062.2010.00588.x>.
5. Reta, Alemayehu, Moges Wubie, and Getnet Mekuria. 2017. "Nasal Colonization and Antimicrobial Susceptibility Pattern of *Staphylococcus Aureus* among Pre - School Children in Ethiopia." *BMC Research Notes*, 1–7. <https://doi.org/10.1186/s13104-017-3079-6>

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Literatur

1. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Teil 460 (2005).
2. Clinical & Laboratory Standards Institute (CLSI). Konsultieren Sie die neueste Ausgabe
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Konsultieren Sie die neueste Version.
3. EUCAST Disk Diffusion Method Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing.
4. Tripodi, M, R Fortunato und S Cuccurullo. 2011. „In Vitro Activity of Multiple Antibiotic Combinations against *Nocardia* : Relationship with a Short-Term Treatment Strategy in Heart Transplant Recipients with Pulmonary Nocardiosis“, Nr. 2: 335–43. <https://doi.org/10.1111/j.1399-3062.2010.00588.x>.
5. Reta, Alemayehu, Moges Wubie und Getnet Mekuria. 2017. „Nasal Colonization and Antimicrobial Susceptibility Pattern of *Staphylococcus Aureus* among Pre - School Children in Ethiopia.“ *BMC Research Notes*, 1–7. <https://doi.org/10.1186/s13104-017-3079-6>

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- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Clinical & Laboratory Standards Institute (CLSI). Συμβουλευτείτε την τελευταία έκδοση
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Συμβουλευτείτε την τελευταία έκδοση.
- EUCAST Disk Diffusion Method Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing.
- Tripodi, M, R Fortunato και S Cuccurullo. 2011. "In Vitro Activity of Multiple Antibiotic Combinations against Nocardia : Relationship with a Short-Term Treatment Strategy in Heart Transplant Recipients with Pulmonary Nocardiosis," no. 2: 335–43. <https://doi.org/10.1111/j.1399-3062.2010.00588.x>.
- Reta, Alemayehu, Moges Wubie και Getnet Mekuria. 2017. "Nasal Colonization and Antimicrobial Susceptibility Pattern of Staphylococcus Aureus among Pre - School Children in Ethiopia." *BMC Research Notes*, 1–7. <https://doi.org/10.1186/s13104-017-3079-6>

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Irodalomjegyzék

- Food and Drug Administration (Élelmiszer- és Gyógyszerügyi Hatóság, FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Clinical Laboratory Standards Institute (Klinikai és Laboratóriumi Szabványügyi Intézet, CLSI). Tekintse meg a legújabb kiadást:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria. European Committee on Antimicrobial Susceptibility Testing (Az antimikrobiális érzékenység tesztelésével foglalkozó európai bizottság, EUCAST). Tekintse meg a legújabb kiadást.
- EUCAST Disk Diffusion Method Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing.
- Tripodi, M, R Fortunato, and S Cuccurullo. 2011. "In Vitro Activity of Multiple Antibiotic Combinations against *Nocardia* : Relationship with a Short-Term Treatment Strategy in Heart Transplant Recipients with Pulmonary Nocardiosis," no. 2: 335–43. <https://doi.org/10.1111/j.1399-3062.2010.00588.x>.
- Reta, Alemayehu, Moges Wubie, and Getnet Mekuria. 2017. "Nasal Colonization and Antimicrobial Susceptibility Pattern of *Staphylococcus Aureus* among Pre - School Children in Ethiopia." *BMC Research Notes*, 1–7. <https://doi.org/10.1186/s13104-017-3079-6>

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Bibliografia

1. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
2. Clinical & Laboratory Standards Institute (CLSI). Consultare l'ultima edizione
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria. Comitato europeo sui test di suscettibilità antimicrobica (European Committee on Antimicrobial Susceptibility Testing, EUCAST). Consultare l'ultima versione.
3. EUCAST Disk Diffusion Method Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing.
4. Tripodi, M, R Fortunato e S Cuccurullo. 2011. "In Vitro Activity of Multiple Antibiotic Combinations against *Nocardia* : Relationship with a Short-Term Treatment Strategy in Heart Transplant Recipients with Pulmonary Nocardiosis," no. 2: 335-43. <https://doi.org/10.1111/j.1399-3062.2010.00588.x>.
5. Reta, Alemayehu, Moges Wubie e Getnet Mekuria. 2017. "Nasal Colonization and Antimicrobial Susceptibility Pattern of *Staphylococcus Aureus* among Pre - School Children in Ethiopia." *BMC Research Notes*, 1-7. <https://doi.org/10.1186/s13104-017-3079-6>

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1. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
2. Clinical & Laboratory Standards Institute (CLSI). Žr. naujausią leidimą
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Žr. naujausią versiją.
3. EUCAST Disk Diffusion Method Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing.
4. Tripodi, M, R Fortunato, and S Cuccurullo. 2011. "In Vitro Activity of Multiple Antibiotic Combinations against *Nocardia*: Relationship with a Short-Term Treatment Strategy in Heart Transplant Recipients with Pulmonary Nocardiosis," no. 2: 335–43. <https://doi.org/10.1111/j.1399-3062.2010.00588.x>.
5. Reta, Alemayehu, Moges Wubie, and Getnet Mekuria. 2017. "Nasal Colonization and Antimicrobial Susceptibility Pattern of *Staphylococcus Aureus* among Pre - School Children in Ethiopia." *BMC Research Notes*, 1–7. <https://doi.org/10.1186/s13104-017-3079-6>

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1. Agencja ds. Żywności i Leków (FDA, Food and Drug Administration). CFR Title 21, Volume 5, Part 460 (2005).
2. Instytut Norm Klinicznych i Laboratoryjnych (CLSI, Clinical Laboratory Standards Institute). Sprawdź najnowsze wydanie
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria. Europejski Komitet ds. Oznaczania Lekowrażliwości (EUCAST, European Committee on Antimicrobial Susceptibility Testing). Sprawdź najnowsze wydanie.
3. EUCAST Disk Diffusion Method Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing.
4. Tripodi, M, R Fortunato i S Cuccurullo. 2011. „In Vitro Activity of Multiple Antibiotic Combinations against *Nocardia*: Relationship with a Short-Term Treatment Strategy in Heart Transplant Recipients with Pulmonary Nocardiosis”, no. 2: 335–43. <https://doi.org/10.1111/j.1399-3062.2010.00588.x>.
5. Reta, Alemayehu, Moges Wubie i Getnet Mekuria. 2017. „Nasal Colonization and Antimicrobial Susceptibility Pattern of *Staphylococcus Aureus* among Pre - School Children in Ethiopia”. *BMC Research Notes*, 1–7. <https://doi.org/10.1186/s13104-017-3079-6>

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Referências

1. Food and Drug Administration (FDA, Agência dos alimentos e dos medicamentos dos EUA). CFR Title 21, Volume 5, Part 460 (2005).
2. Clinical & Laboratory Standards Institute (CLSI, Instituto de padrões clínicos e laboratoriais dos EUA). Consultar a edição mais recente
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria. European Committee on Antimicrobial Susceptibility Testing (EUCAST, Comité Europeu para o Teste à Suscetibilidade Antimicrobiana). Consultar a versão mais recente.
3. EUCAST Disk Diffusion Method Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing.
4. Tripodi, M, R Fortunato, and S Cuccurullo. 2011. "In Vitro Activity of Multiple Antibiotic Combinations against *Nocardia*: Relationship with a Short-Term Treatment Strategy in Heart Transplant Recipients with Pulmonary Nocardiosis," no. 2: 335–43. <https://doi.org/10.1111/j.1399-3062.2010.00588.x>.
5. Reta, Alemayehu, Moges Wubie, and Getnet Mekuria. 2017. "Nasal Colonization and Antimicrobial Susceptibility Pattern of *Staphylococcus Aureus* among Pre - School Children in Ethiopia." *BMC Research Notes*, 1–7. <https://doi.org/10.1186/s13104-017-3079-6>

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Odkazy

- Food and Drug Administration (FDA). CFR Hlava 21, Zväzok 5, Časť 460 (2005).
- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version.
- EUCAST Disk Diffusion Method Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing.

4. Tripodi, M, R Fortunato, a S Cuccurullo. 2011. „In Vitro Activity of Multiple Antibiotic Combinations against Nocardia : Relationship with a Short-Term Treatment Strategy in Heart Transplant Recipients with Pulmonary Nocardiosis,“ č. 2: 335–43. <https://doi.org/10.1111/j.1399-3062.2010.00588.x>.
5. Reta, Alemayehu, Moges Wubie, a Getnet Mekuria. 2017. „Nasal Colonization and Antimicrobial Susceptibility Pattern of Staphylococcus Aureus among Pre - School Children in Ethiopia.“ *BMC Research Notes*, 1–7. <https://doi.org/10.1186/s13104-017-3079-6>

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Referencias

1. Food and Drug Administration (FDA). Título 21 del CFR, volumen 5, parte 460 (2005).
2. Clinical & Laboratory Standards Institute (CLSI). Consultar la última edición.
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests.
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing.
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consultar la última versión.
3. Método de difusión en disco EUCAST, «Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing».
4. Tripodi, M, R Fortunato y S Cuccurullo. 2011. «In Vitro Activity of Multiple Antibiotic Combinations against *Nocardia* : Relationship with a Short-Term Treatment Strategy in Heart Transplant Recipients with Pulmonary Nocardiosis», n.º 2: 335-43. <https://doi.org/10.1111/j.1399-3062.2010.00588.x>.
5. Reta, Alemayehu, Moges Wubie y Getnet Mekuria. 2017. «Nasal Colonization and Antimicrobial Susceptibility Pattern of *Staphylococcus Aureus* among Pre - School Children in Ethiopia». *BMC Research Notes*, 1-7. <https://doi.org/10.1186/s13104-017-3079-6>.

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