	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Nalidixic Acid Discs 30 µg (NA30)

REF CT0031B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Nalidixic Acid Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent nalidixic. The discs are labelled on both sides with details of the antimicrobial (NA) and amount present (µg): NA30 (30µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Nalidixic Acid AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which nalidixic Acid has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- N/A

Gram-Negative

- N/A

Species with published breakpoints according to current CLSI M100:

Gram-Negative

- Enterobacterales

Gram-Positive

- N/A

Species with published breakpoints according to current EUCAST literature:

Gram-Negative

- *Haemophilus influenzae* (screen only)
- *Moraxella catarrhalis* (screen only)
- *Pasteurella multocida* (screen only)

Gram-Positive

- N/A

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Nalidixic Acid AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface.

The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Nalidixic Acid AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Nalidixic Acid AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Nalidixic Acid AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with CT0031B	
Component Description	Material Description
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Pale beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Nalidixic Acid AST Disc Reagents	
Reagent	Description of Function
Nalidixic Acid	Nalidixic acid is a quinolone antibiotic used to treat urinary tract infections. It is a synthetic 1,8-naphthyridine antimicrobial agent with a limited bactericidal spectrum. It works as an inhibitor of the A subunit of bacterial DNA gyrase.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Nalidixic Acid AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Nalidixic Acid AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at -20°C to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 1. Raw data of zone sizes, in accordance with EUCAST methodology, for 3 batches of Nalidixic Acid (30 µg) NA30 Antimicrobial Susceptibility Test Discs (CT0031B).

Antimicrobial Susceptibility Test Discs (OT0001B)										
EUCAST Methodology										
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3258254	<i>Escherichia coli</i>	ATCC® 25922	MHA	22	28	25	26	26	26	1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49766	MHA-F	26	32	29	29	29	29	0 ± 0 (CV=0)
Product	Organism						Reading (mm)			Average (mm) + SD

		ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	1	2	3	
3221272	<i>Escherichia coli</i>	ATCC® 25922	MHA	22	28	25	26	26	26	1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49766	MHA-F	26	32	29	30	30	30	1 ± 0 (CV=0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3184452	<i>Escherichia coli</i>	ATCC® 25922	MHA	22	28	25	26	26	27	1.33 ± 0.58 (CV=0.43)
	<i>Haemophilus influenzae</i>	ATCC® 49766	MHA-F	26	32	29	29	29	28	-0.33 ± 0.58 (CV=-1.73)

Table 2. Raw data of zone sizes, in accordance with CLSI methodology, taken from the 3 most recent batches of Nalidixic Acid (30 µg) NA30 Antimicrobial Susceptibility Test Discs (CT0031B).

CLSI Methodology										
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid- point values (mm) + SD (Coefficient of variance)
							1	2	3	
3258254	<i>Escherichia coli</i>	ATCC® 25922	MHA	22	28	25	26	26	26	1 ± 0 (CV=0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3221272	<i>Escherichia coli</i>	ATCC® 25922	MHA	22	28	25	26	26	26	1 ± 0 (CV=0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3184452	<i>Escherichia coli</i>	ATCC® 25922	MHA	22	28	25	26	26	27	1.33 ± 0.58 (CV=0.43)

SD: Standard Deviation

Clinical Performance Characteristics

A study was conducted to compare Oxoid Nalidixic Acid AST discs to those from other manufacturers who produce discs with a similar content of antibiotic.⁴ Discs were tested with an organism of known sensitivity and the microorganisms were chosen to cover all the degrees of sensitivity. For each antibiotic disc two to four organisms were used. The results showed that Oxoid Nalidixic Acid AST Discs have the antimicrobial content that is stated for their discs and demonstrated that they are an effective method for in vitro analysis of *E. coli* and *K. aerogenes*.

A 2004 study collected 21 *Salmonella typhi* and *Salmonella paratyphi* blood cultures from patients.⁵ Susceptibility testing was performed via disc diffusion method using Oxoid Nalidixic Acid AST Discs and other antibiotics. MICs for nalidixic acid and ciprofloxacin were determined using the agar dilution method. The study concluded that reduced fluoroquinolone susceptibility among *Salmonella* spp. isolates impacts clinical efficacy among enteric fever patients and demonstrated that the device is an effective method for in vitro analysis of *Salmonella* spp. and the detection of fluoroquinolone resistance.

Twenty-eight *Y. enterocolitica* 4/O:3 strains were collected from clinical stool samples and were studied for antimicrobial resistance through disc diffusion and broth microdilution methods performed in accordance with Clinical and Laboratory Standards Institute (CLSI).⁶ The disc diffusion test was performed on Mueller-Hinton agar (MHA) plates for 10 antibiotics, including Oxoid Nalidixic Acid AST Discs. Broth microdilution was then performed to determine MICs for the sample antibiotics.

All strains of *Y. enterocolitica* were susceptible to nalidixic acid, making it one of the most active antimicrobials against this microorganism. Oxoid Nalidixic Acid AST Discs demonstrated that they were effective for in vitro analysis of *Y. enterocolitica*.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005)
4. D. F. J. Brown and D. Kothari (1975). Comparison of antibiotic discs from different sources (Brown1 and Kothari 1975). *Journal of Clinical Pathology*, 1975 Oct; 28(10): 779–783. DOI: <https://doi.org/10.1136/jcp.28.10.779>
5. R. Slinger, Marc Desjardins, A. E. McCarthy, K. Ramotar, et al (2004). Suboptimal clinical response to ciprofloxacin in patients with enteric fever due to *Salmonella* spp. with reduced fluoroquinolone susceptibility: a case series. *BMC Infect Dis* 4, 36 (2004). DOI: <https://doi.org/10.1186/1471-2334-4-36>
6. C. Meyer, A. Stolle, M. Fredriksson – Ahomaa (2011). Comparison of Broth Microdilution and Disk Diffusion Test for Antimicrobial Resistance Testing in *Yersinia enterocolitica* 4/O:3 Strains. *Microbial Drug Resistance*, Volume 17, Number 3, 2011

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number

	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision Information

Version	Date of issue and modifications introduced
1.0	2023-01-17

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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Referencije

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005)
- D. F. J. Brown and D. Kothari (1975). Comparison of antibiotic discs from different sources (Brown1 and Kothari 1975). *Journal of Clinical Pathology*, 1975 Oct; 28(10): 779–783. DOI: <https://doi.org/10.1136/jcp.28.10.779>
- R. Slinger, Marc Desjardins, A. E. McCarthy, K. Ramotar, et al (2004). Suboptimal clinical response to ciprofloxacin in patients with enteric fever due to *Salmonella spp.* with reduced fluoroquinolone susceptibility: a case series. *BMC Infect Dis* 4, 36 (2004). DOI: <https://doi.org/10.1186/1471-2334-4-36>
- C. Meyer, A. Stolle, M. Fredriksson – Ahomaa (2011). Comparison of Broth Microdilution and Disk Diffusion Test for Antimicrobial Resistance Testing in *Yersinia enterocolitica* 4/O:3 Strains. *Microbial Drug Resistance*, Volume 17, Number 3, 2011

© 2023. Thermo Fisher Scientific Inc. Sva prava pridržana. ATCC® je zaštitni znak društva ATCC. Svi ostali zaštitni znakovi vlasništvo su društva Thermo Fisher Scientific Inc. i njegovih podružnica. Ove informacije nisu namijenjene za poticanje na uporabu tih proizvoda na način kojim se mogu povrijediti prava intelektualnog vlasništva drugih društava.

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Litteratur

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005)
4. D. F. J. Brown and D. Kothari (1975). Comparison of antibiotic discs from different sources (Brown and Kothari 1975). *Journal of Clinical Pathology*, 1975 Oct; 28(10): 779–783. DOI: <https://doi.org/10.1136/jcp.28.10.779>
5. R. Slinger, Marc Desjardins, A. E. McCarthy, K. Ramotar, et al (2004). Suboptimal clinical response to ciprofloxacin in patients with enteric fever due to *Salmonella* spp. with reduced fluoroquinolone susceptibility: a case series. *BMC Infect Dis* 4, 36 (2004). DOI: <https://doi.org/10.1186/1471-2334-4-36>
6. C. Meyer, A. Stolle, M. Fredriksson – Ahomaa (2011). Comparison of Broth Microdilution and Disk Diffusion Test for Antimicrobial Resistance Testing in *Yersinia enterocolitica* 4/O:3 Strains. *Microbial Drug Resistance*, Volume 17, Number 3, 2011

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Références

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005)
- D. F. J. Brown and D. Kothari (1975). Comparison of antibiotic discs from different sources (Brown1 and Kothari 1975). *Journal of Clinical Pathology*, 1975 Oct; 28(10): 779–783. DOI: <https://doi.org/10.1136/jcp.28.10.779>
- R. Slinger, Marc Desjardins, A. E. McCarthy, K. Ramotar, et al (2004). Suboptimal clinical response to ciprofloxacin in patients with enteric fever due to *Salmonella* spp. with reduced fluoroquinolone susceptibility: a case series. *BMC Infect Dis* 4, 36 (2004). DOI: <https://doi.org/10.1186/1471-2334-4-36>
- C. Meyer, A. Stolle, M. Fredriksson – Ahomaa (2011). Comparison of Broth Microdilution and Disk Diffusion Test for Antimicrobial Resistance Testing in *Yersinia enterocolitica* 4/O:3 Strains. *Microbial Drug Resistance*, Volume 17, Number 3, 2011

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- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005)
- D. F. J. Brown and D. Kothari (1975). Comparison of antibiotic discs from different sources (Brown1 and Kothari 1975). *Journal of Clinical Pathology*, 1975 Oct; 28(10): 779–783. DOI: <https://doi.org/10.1136/jcp.28.10.779>

5. R. Slinger, Marc Desjardins, A. E. McCarthy, K. Ramotar, et al (2004). Suboptimal clinical response to ciprofloxacin in patients with enteric fever due to *Salmonella* spp. with reduced fluoroquinolone susceptibility: a case series. *BMC Infect Dis* 4, 36 (2004). DOI: <https://doi.org/10.1186/1471-2334-4-36>
6. C. Meyer, A. Stolle, M. Fredriksson – Ahomaa (2011). Comparison of Broth Microdilution and Disk Diffusion Test for Antimicrobial Resistance Testing in *Yersinia enterocolitica* 4/O:3 Strains. *Microbial Drug Resistance*, Volume 17, Number 3, 2011

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Bibliografia

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005)
- D. F. J. Brown and D. Kothari (1975). Comparison of antibiotic discs from different sources (Brown and Kothari 1975). *Journal of Clinical Pathology*, 1975 Oct; 28(10): 779–783. DOI: <https://doi.org/10.1136/jcp.28.10.779>
- R. Slinger, Marc Desjardins, A. E. McCarthy, K. Ramotar, et al (2004). Suboptimal clinical response to ciprofloxacin in patients with enteric fever due to *Salmonella* spp. with reduced fluoroquinolone susceptibility: a case series. *BMC Infect Dis* 4, 36 (2004). DOI: <https://doi.org/10.1186/1471-2334-4-36>
- C. Meyer, A. Stolle, M. Fredriksson – Ahomaa (2011). Comparison of Broth Microdilution and Disk Diffusion Test for Antimicrobial Resistance Testing in *Yersinia enterocolitica* 4/O:3 Strains. *Microbial Drug Resistance*, Volume 17, Number 3, 2011

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