	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Sulphamethoxazole/Trimethoprim 19:1 Discs

25 µg (SXT25)

REF CT0052B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Sulphamethoxazole/Trimethoprim 19:1 (25 µg) SXT25 Antimicrobial Susceptibility Test Discs (AST) are 6mm discs which contain specific amounts of the antimicrobial agents sulphamethoxazole and trimethoprim. The discs are labelled on both sides with details of the antimicrobial (SXT) and amount present (25µg).

Discs are supplied in cartridges containing 50 discs each, and there are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Sulphamethoxazole/Trimethoprim Antimicrobial Susceptibility Test (AST) Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Sulphamethoxazole/Trimethoprim Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which sulphamethoxazole and trimethoprim have been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- *Streptococcus pneumoniae*

Gram-Negative

- *Enterobacteriaceae*
- *Haemophilus influenzae* and *parainfluenzae*

Species with published breakpoints according to current CLSI literature:

Gram-Positive

- *Staphylococcus* spp.
- *Streptococcus pneumoniae*

Gram-Negative

- *Enterobacterales*
- *Acinetobacter* spp.
- *Burkholderia cepacia*
- *Stenotrophomonas maltophilia*
- *Haemophilus influenzae* and *parainfluenzae*
- *Neisseria meningitidis*

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- *Staphylococcus* spp.
- *Streptococcus* groups A, B, C and G
- *Streptococcus pneumoniae*
- *Listeria monocytogenes*

Gram-Negative

- *Enterobacterales*
- *Stenotrophomonas maltophilia*
- *Acinetobacter* spp.
- *Haemophilus influenzae*
- *Moraxella catarrhalis*
- *Pasteurella multocida*
- *Kingella kingae*
- *Aeromonas* spp.
- *Achromobacter xylosoxidans*
- *Burkholderia pseudomallei*

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principles of Method

Sulphamethoxazole/Trimethoprim AST Discs are used in the semi-quantitative agar diffusion test method for *in vitro* susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Sulphamethoxazole/Trimethoprim AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Sulphamethoxazole/Trimethoprim AST Discs consist of 6mm diameter paper discs impregnated with a specific antimicrobial concentration. The discs are marked on both sides with an alpha-numeric code identifying the antimicrobial agent and concentration. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to table 2.

Component Description	Material Description
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Pale beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Reagent	Description of Function
Sulfamethoxazole	White crystalline powder. Antimicrobial agent. A sulfonamide bacteriostatic antibiotic that interferes with folic acid synthesis, in susceptible bacteria.
Trimethoprim	White or off-white powder. Antimicrobial agent that inhibits dihydrofolic acid reduction to tetrahydrofolic acid resulting in sequential inhibition of enzymes of the folic acid pathway.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Sulphamethoxazole/Trimethoprim AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Sulphamethoxazole/Trimethoprim AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation.

Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Raw data of zone sizes, in accordance with CLSI methodology, taken from the 3 most recent batches of Sulphamethoxazole/Trimethoprim 19:1 (25 µg) SXT25 Antimicrobial Susceptibility Test Discs (CT0052B).

CLSI Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3179494	<i>Escherichia coli</i>	25922	MHA	23	29	26	27	27	27	1 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	24	32	28	32	32	32	4 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	24	32	28	27	27	27	-1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA+5%SB	20	28	24	23	23	23	-1 ± 0 (CV=0)
3179495	<i>Escherichia coli</i>	25922	MHA	23	29	26	26	26	26	0 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	24	32	28	30	30	30	2 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	24	32	28	28	28	28	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA+5%SB	20	28	24	24	24	24	0 ± 0 (CV=0)
3134308	<i>Escherichia coli</i>	25922	MHA	23	29	26	26	26	26	0 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	24	32	28	27	27	27	-1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	24	32	28	26	26	26	-2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA+5%SB	20	28	24	24	24	24	0 ± 0 (CV=0)

Table 4. Raw data of zone sizes, in accordance with EUCAST methodology, taken from the 3 most recent batches of Sulphamethoxazole/Trimethoprim 19:1 (25 µg) SXT25 Antimicrobial Susceptibility Test Discs (CT0052B).

EUCAST Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values. (mm) + SD (Coefficient of variance)
							1	2	3	
3179494	<i>Escherichia coli</i>	25922	MHA	23	29	26	27	27	27	1 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	29213	MHA	26	32	29	31	31	31	2 ± 0 (CV=0)
	<i>Enterococcus faecalis</i>	29212	MHA	26	34	30	32	32	32	2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MH + 5% HB+ 20 mg/L NAD in CO ₂	18	26	22	24	24	24	2 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49766	MH + 5% HB+ 20 mg/L NAD in CO ₂	27	35	31	31	31	31	0 ± 0 (CV=0)
3179495	<i>Escherichia coli</i>	25922	MHA	23	29	26	26	26	26	0 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	29213	MHA	26	32	29	31	31	31	2 ± 0 (CV=0)
	<i>Enterococcus faecalis</i>	29212	MHA	26	34	30	31	31	31	1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MH + 5% HB+ 20 mg/L NAD in CO ₂	18	26	22	23	23	23	1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49766	MH + 5% HB+ 20 mg/L NAD in CO ₂	27	35	31	32	32	32	1 ± 0 (CV=0)
3134308	<i>Escherichia coli</i>	25922	MHA	23	29	26	26	26	26	0 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	29213	MHA	26	32	29	28	28	28	-1 ± 0 (CV=0)
	<i>Enterococcus faecalis</i>	29212	MHA	26	34	30	30	30	30	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MH + 5% HB+ 20 mg/L NAD in CO ₂	18	26	22	25	25	25	3 ± 0 (CV=0)

	<i>Haemophilus influenzae</i>	49766	MH + 5% HB+ 20 mg/L NAD in CO ₂	27	35	31	30	30	30	-1 ± 0 (CV=0)
--	-------------------------------	-------	--	----	----	----	----	----	----	---------------

Clinical Performance Characteristics

A study by Hegde et al. (2019) described the epidemiology and antimicrobial susceptibility patterns of culture-confirmed *Shigella* infections in facility-based surveillance sites in Guatemala.⁴ From 2007 to 2012, stool samples from 5399 patients yielded 261 cases of *Shigella* spp., with the most frequently isolated serotypes being *Shigella flexneri* (59%) and *Shigella sonnei* (36%). AST was performed, at the hospital laboratories, on laboratory confirmed *Shigella* isolates (254) using the disc diffusion method with a range of antibiotics, including Oxoid Sulphamethoxazole/Trimethoprim SXT25 AST Discs, and using CLSI breakpoints. Multidrug resistance (MDR) was defined as resistance to at least ampicillin and trimethoprim sulfamethoxazole, an important MDR phenotype for *Shigella*. Overall, the study found that antimicrobial resistance to commonly used antimicrobials was high, as was the proportion that was resistant to at least three antimicrobial classes. The resistance to trimethoprim sulfamethoxazole was detected in 83% of the isolates. 120 cases were MDR (51%), which is markedly higher than in the United States (16%). Overall, the study concluded that although resistance to both ampicillin and trimethoprim sulfamethoxazole (MDR resistance) increased over the surveillance period, continued surveillance of these dynamics is critical for understanding both prevention and treatment. Therefore, as the study did not report any out of scope results and/or inconsistencies with using AST discs, they are still effective tools in diagnostic workflow to aid clinicians in treatment options.

A study by Ahmad et al., in 2015, determined the frequency and resistance pattern of *E. coli* in patients of CA-UTIs in a Northern part of Pakistan.⁵ Urine specimens were collected from patients who were clinically diagnosed with CA-UTIs. In vitro AST of *E. coli* were performed with standardized commercial susceptibility discs, which included Oxoid Sulphamethoxazole/Trimethoprim SXT25 AST Discs. The sensitivity of antibiotics was studied on Mueller-Hinton agar (MHA). Out of the 50 specimens, 20 showed pure growth of *E. coli* on culture. Overall, *E. coli* showed a high rate of resistance towards ampicillin (90%), tetracycline (70%), erythromycin (70%), trimethoprim-sulfamethoxazole (55%), ciprofloxacin (50%) and sparfloxacin (45%). Overall, the study concluded that the rate of resistance of *E. coli* against commonly used antibacterials was quite high and a majority of the strains showed multidrug resistance. However, the use of AST discs, including Oxoid Sulphamethoxazole/Trimethoprim SXT25 AST Discs, showed consistent results, and no out-of-specification results were documented.

A study was conducted to determine the antibiotic susceptibility pattern of uropathogenic *E. coli*, prevalence of the extended-spectrum β -lactamase (ESBL) gene group and class 1 integrons.⁶ A total of 500 clinical isolates of uropathogenic *E. coli* were isolated from different hospitals in Karachi. AST was performed by Kirby-Bauer methods on MHA according to the guidelines of CLSI. *E. coli* ATCC® 25922 was used as a control in the study. A range of antibiotics were tested, including cotrimoxazole (trimethoprim/sulfamethoxazole). Out of 500 isolates, 105 were identified as MDR *E. coli*. MDR *E. coli* exhibited an increasing trend of resistance; more than 50% resistance was observed in antibiotics, including cotrimoxazole. As uropathogenic ESBL-producing *E. coli* is an emerging health issue spreading worldwide, therefore the study emphasised the importance of exploring the epidemiology of resistance, as well as the application of recommendations for proper control of antimicrobial usage and infection control measures. Overall, the study concluded that active surveillance and antimicrobial stewardship are highly recommended, and therefore, as the study did not report any out of scope results and/or inconsistencies with using AST discs, they are still effective tools in diagnostic workflow to aid clinicians in treatment options.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.


References

1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Hegde, Sonia, Stephen R. Benoit, Wences Arvelo, Kim Lindblade, Beatriz López, John P. McCracken, Chris Bernart, Aleida Roldan, and Joe P. Bryan. 2019. 'Burden of Laboratory-Confirmed Shigellosis Infections in Guatemala 2007-2012: Results from a Population-Based Surveillance System'. *BMC Public Health* 19 (Suppl 3): 474. <https://doi.org/10.1186/s12889-019-6780-7>
5. Ahmad W, Jamshed F, Ahmad W. Frequency of Escherichia Coli in Patients with Community Acquired Urinary Tract Infection and Their Resistance Pattern Against Some Commonly Used Anti Bacterials. *Journal of Ayub Medical College Abbottabad*. 2015 Apr-Jun;27(2):333-7. PMID: 26411110
6. Zeeshan Khan, Fouzia, Tehseen Nawaz, Zulfiqar Ali Mirani, Saeed Khan, Yasir Raza, and Shahana Urooj Kazmi. 2018. 'Study of Class 1 Integrons in Multidrug-Resistant Uropathogenic Escherichia Coli Isolated from Different Hospitals in Karachi'. *Iranian Journal of Basic Medical Sciences* 21 (10): 1079–82. <https://doi.org/10.22038/IJBMS.2018.28807.6966>.

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision Information

Version	Date of issue and modifications introduced
1.0	2023-01-18

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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Reference

1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Hegde, Sonia, Stephen R. Benoit, Wences Arvelo, Kim Lindblade, Beatriz López, John P. McCracken, Chris Bernart, Aleida Roldan, and Joe P. Bryan. 2019. 'Burden of Laboratory-Confirmed Shigellosis Infections in Guatemala 2007-2012: Results from a Population-Based Surveillance System'. *BMC Public Health* 19 (Suppl 3): 474. <https://doi.org/10.1186/s12889-019-6780-7>
5. Ahmad W, Jamshed F, Ahmad W. Frequency of Escherichia Coli in Patients with Community Acquired Urinary Tract Infection and Their Resistance Pattern Against Some Commonly Used Anti Bacterials. *Journal of Ayub Medical College Abbottabad*. 2015 Apr-Jun;27(2):333-7. PMID: 26411110
6. Zeeshan Khan, Fouzia, Tehseen Nawaz, Zulfiqar Ali Mirani, Saeed Khan, Yasir Raza, and Shahana Urooj Kazmi. 2018. 'Study of Class 1 Integrons in Multidrug-Resistant Uropathogenic Escherichia Coli Isolated from Different Hospitals in Karachi'. *Iranian Journal of Basic Medical Sciences* 21 (10): 1079–82. <https://doi.org/10.22038/IJBMS.2018.28807.6966>

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Referencie

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Hegde, Sonia, Stephen R. Benoit, Wences Arvelo, Kim Lindblade, Beatriz López, John P. McCracken, Chris Bernart, Aleida Roldan, and Joe P. Bryan. 2019. 'Burden of Laboratory-Confirmed Shigellosis Infections in Guatemala 2007-2012: Results from a Population-Based Surveillance System'. BMC Public Health 19 (Suppl 3): 474. <https://doi.org/10.1186/s12889-019-6780-7>
- Ahmad W, Jamshed F, Ahmad W. Frequency of Escherichia Coli in Patients with Community Acquired Urinary Tract Infection and Their Resistance Pattern Against Some Commonly Used Anti Bacterials. Journal of Ayub Medical College Abbottabad. 2015 Apr-Jun;27(2):333-7. PMID: 26411110
- Zeeshan Khan, Fouzia, Tehseen Nawaz, Zulfiqar Ali Mirani, Saeed Khan, Yasir Raza, and Shahana Urooj Kazmi. 2018. 'Study of Class 1 Integrons in Multidrug-Resistant Uropathogenic Escherichia Coli Isolated from Different Hospitals in Karachi'. Iranian Journal of Basic Medical Sciences 21 (10): 1079–82. <https://doi.org/10.22038/IJBMS.2018.28807.6966>.

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Referencias

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Hegde, Sonia, Stephen R. Benoit, Wences Arvelo, Kim Lindblade, Beatriz López, John P. McCracken, Chris Bernart, Aleida Roldan, and Joe P. Bryan. 2019. 'Burden of Laboratory-Confirmed Shigellosis Infections in Guatemala 2007-2012: Results from a Population-Based Surveillance System'. *BMC Public Health* 19 (Suppl 3): 474. <https://doi.org/10.1186/s12889-019-6780-7>
- Ahmad W, Jamshed F, Ahmad W. Frequency of Escherichia Coli in Patients with Community Acquired Urinary Tract Infection and Their Resistance Pattern Against Some Commonly Used Anti Bacterials. *Journal of Ayub Medical College Abbottabad*. 2015 Apr-Jun;27(2):333-7. PMID: 26411110
- Zeeshan Khan, Fouzia, Tehseen Nawaz, Zulfiqar Ali Mirani, Saeed Khan, Yasir Raza, and Shahana Urooj Kazmi. 2018. 'Study of Class 1 Integrons in Multidrug-Resistant Uropathogenic Escherichia Coli Isolated from Different Hospitals in Karachi'. *Iranian Journal of Basic Medical Sciences* 21 (10): 1079–82. <https://doi.org/10.22038/IJBMS.2018.28807.6966>.

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Oxoid Limited, Wade Road, Basingstoke,
Hampshire, RG24 8PW, Reino Unido

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