	https://www.thermofisher.com	
		
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Amikacin Discs 30 µg (AK30)

REF CT0107B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Amikacin Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent Amikacin. The discs are labelled on both sides with details of antimicrobial (AK) and amount present (µg): AK30 (30ug).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Amikacin AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disk should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for *in vitro* susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which Amikacin has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA³, CLSI M100^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Negative

- Enterobacteriaceae
- *Pseudomonas aeruginosa*
- *Acinetobacter* spp.

Gram-Positive

- *Staphylococcus aureus*

Species with published breakpoints according to current CLSI literature:

Gram-Negative

- Enterobacterales
- *Pseudomonas aeruginosa*
- *Acinetobacter* spp.

Species with published breakpoints according to current EUCAST literature:

Gram-Negative

- Enterobacterales (systemic infections and infections originating from the urinary tract)
- *Pseudomonas* spp. (systemic infections and infections originating from the urinary tract)
- *Acinetobacter* spp. (systemic infections and infections originating from the urinary tract)

Gram-Positive

- *Staphylococcus* spp. (*S. aureus* and coagulase-negative staphylococci)

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Amikacin AST Discs are used in the semi-quantitative agar diffusion test method for *in vitro* susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation

referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Amikacin AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Amikacin AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Amikacin AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

desiccant.

Table 1. Materials Provided with CT0107B	
Component Description	Material Description
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Pale beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Amikacin AST Disc Reagents	
Reagent	Description of Function
Amikacin	Antimicrobial agent. White crystalline powder. Mechanism of action: binds to bacterial 30S ribosomal subunits, interferes with mRNA binding and tRNA acceptor sites, interfering with bacterial growth

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Amikacin AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Amikacin AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this document.

Analytical Performance Characteristics

Table 3. Raw data of zone sizes, in accordance with CLSI methodology, taken from the 3 most recent batches of Amikacin AK30 (30µg) Antimicrobial Susceptibility Test Discs (CT0107B).

CLSI Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3116894	<i>Escherichia coli</i>	ATCC 25922	MHA	19	26	22.5	19	19	19	-3.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC 25923	MHA	20	26	23	23	23	23	0 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC 27853	MHA	18	26	22	21	21	20	-1.33 ± 0.58 (CV=-0.43)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3104599	<i>Escherichia coli</i>	ATCC 25923	MHA	19	26	22.5	24	24	24	-1.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC 25922	MHA	20	26	23	21	21	21	-2 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC 49766	HTM	18	26	22	23	23	23	1 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3001687	<i>Escherichia coli</i>	ATCC 25923	MHA	19	26	22.5	22	22	22	-0.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC 25922	MHA	20	26	23	21	21	21	-2 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC 49766	HTM	18	26	22	23	23	22	0.67 ± 0.58 (CV=0.87)

Table 1. Raw data of zone sizes, in accordance with EUCAST methodology, taken from the 3 most recent batches of Amikacin AK30 (30µg) Antimicrobial Susceptibility Test Discs (CT0107B).

EUCAST Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3116894	<i>Staphylococcus aureus</i>	ATCC 25923	MHA	18	24	21	22	22	22	1 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC 25922	MHA	19	26	22.5	19	19	19	-3.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC 27853	MHA	18	26	22	21	21	20	-1.33 ± 0.58 (CV=-0.43)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3104599	<i>Escherichia coli</i>	ATCC 25923	MHA	19	26	21	21	21	21	0 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC 25922	MHA	20	26	22.5	24	24	24	1.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC 49766	HTM	18	26	22	23	23	23	1 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3001687	<i>Escherichia coli</i>	ATCC 25923	MHA	19	26	21	23	23	23	2 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC 25922	MHA	20	26	22.5	22	22	22	-0.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC 49766	HTM	18	26	22	23	23	22	0.67 ± 0.58 (CV=0.87)

Clinical Performance Characteristics

A study by Acar, Levent, and Atalan 2015, aimed to analyse interrelated *Salmonella enteritidis* strains isolated in Turkey using antimicrobial susceptibility testing, plasmid profile analysis and pulsed-field gel electrophoresis.⁴ During the study period, from 2000 to 2010, a total of 55 *S. enteritidis* strains, isolated from faeces and environmental sources, were selected from the collection of the Public Health Institution of Turkey and were studied. Antimicrobial susceptibility testing was carried out using the standard disk diffusion method and was interpreted in accordance with CLSI guidelines. Each isolate was tested for sensitivity against a selection of 20 antibiotics from Oxoid including, Oxoid Amikacin AST Discs. *E. coli* ATCC 25922 was used as the quality control strain. Results showed all strains were 100% susceptible to amikacin (Acar, Levent, and Atalan 2015). This provides evidence that Amikacin discs perform effectively in a clinical setting and are effective against clinical Gram-negative isolates.

A study by Alabi et al. 2017, aimed to screen for presence of transferable resistance genes and mobile genetic elements (MGEs) such as integron among multi-drug resistant (MDR) *P. mirabilis*.⁵ During the study period of 12 months (January to December 2011), a total of 108 non-duplicated clinical isolates of *Proteus mirabilis* were randomly collected from Microbiology units of five tertiary hospitals in Southwest Nigeria. Antimicrobial susceptibility testing was performed using the disc diffusion method in accordance with CLSI guidelines against a panel of twelve antibiotics including, Thermo Scientific™ Oxoid™ Amikacin (30 µg)

AK30 Antimicrobial Susceptibility Test Discs (CT0107B). Results showed amikacin resistance was 4.6%, the second lowest of the range (Table 2.) (Alabi et al. 2017). This provides evidence that Amikacin discs perform effectively and consistently in a clinical setting and are effective against clinical Gram-negative isolates. Furthermore, no out of specification or inconsistent results were reported in the study.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Acar, Sümeyra, Belkıs Levent, and Ekrem Atalan. 2015. "A Molecular Epidemiological Investigation of Multistate Outbreaks of Salmonella Enteritidis from Clinical and Environmental Samples in Turkey, 2000–2010." *Turkish Journal of Medical Sciences* 45 (1): 76–83. <https://doi.org/10.3906/sag-1401-33>.
5. Alabi, Olumuyiwa Samuel, Nuno Mendonça, Olufemi Ezekiel Adeleke, and Gabriela Jorge da Silva. 2017. "Molecular Screening of Antibiotic-Resistant Determinants among Multidrug-Resistant Clinical Isolates of *Proteus Mirabilis* from SouthWest Nigeria." *African Health Sciences* 17 (2): 356–65. <https://doi.org/10.4314/ahs.v17i2.9>.

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of Symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use

	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision Information

Version	Date of issue and modifications introduced
1.0	2023-01-12.

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Acar, Sümeyra, Belkis Levent, and Ekrem Atalan. 2015. "A Molecular Epidemiological Investigation of Multistate Outbreaks of Salmonella Enteritidis from Clinical and Environmental Samples in Turkey, 2000–2010." *Turkish Journal of Medical Sciences* 45 (1): 76–83. <https://doi.org/10.3906/sag-1401-33>.
5. Alabi, Olumuyiwa Samuel, Nuno Mendonça, Olufemi Ezekiel Adeleke, and Gabriela Jorge da Silva. 2017. "Molecular Screening of Antibiotic-Resistant Determinants among Multidrug-Resistant Clinical Isolates of Proteus Mirabilis from SouthWest Nigeria." *African Health Sciences* 17 (2): 356–65. <https://doi.org/10.4314/ahs.v17i2.9>.

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Referencije

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Acar, Sümeyra, Belkis Levent, and Ekrem Atalan. 2015. "A Molecular Epidemiological Investigation of Multistate Outbreaks of Salmonella Enteritidis from Clinical and Environmental Samples in Turkey, 2000–2010." *Turkish Journal of Medical Sciences* 45 (1): 76–83. <https://doi.org/10.3906/sag-1401-33>.
5. Alabi, Olumuyiwa Samuel, Nuno Mendonça, Olufemi Ezekiel Adeleke, and Gabriela Jorge da Silva. 2017. "Molecular Screening of Antibiotic-Resistant Determinants among Multidrug-Resistant Clinical Isolates of *Proteus Mirabilis* from SouthWest Nigeria." *African Health Sciences* 17 (2): 356–65. <https://doi.org/10.4314/ahs.v17i2.9>.

© 2023. Thermo Fisher Scientific Inc. Sva prava pridržana. ATCC® je zaštitni znak društva ATCC. Svi ostali zaštitni znakovi vlasništvo su društva Thermo Fisher Scientific Inc. i njegovih podružnica. Ove informacije nisu namijenjene za poticanje na uporabu tih proizvoda na način kojim se mogu povrijediti prava intelektualnog vlasništva drugih društava.

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Viitteet

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Acar, Sümeyra, Belkis Levent, and Ekrem Atalan. 2015. "A Molecular Epidemiological Investigation of Multistate Outbreaks of Salmonella Enteritidis from Clinical and Environmental Samples in Turkey, 2000–2010." *Turkish Journal of Medical Sciences* 45 (1): 76–83. <https://doi.org/10.3906/sag-1401-33>.
- Alabi, Olumuyiwa Samuel, Nuno Mendonça, Olufemi Ezekiel Adeleke, and Gabriela Jorge da Silva. 2017. "Molecular Screening of Antibiotic-Resistant Determinants among Multidrug-Resistant Clinical Isolates of Proteus Mirabilis from SouthWest Nigeria." *African Health Sciences* 17 (2): 356–65. <https://doi.org/10.4314/ahs.v17i2.9>.

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Referenzen

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Acar, Sümeýra, Belkis Levent, and Ekrem Atalan. 2015. "A Molecular Epidemiological Investigation of Multistate Outbreaks of Salmonella Enteritidis from Clinical and Environmental Samples in Turkey, 2000–2010." *Turkish Journal of Medical Sciences* 45 (1): 76–83. <https://doi.org/10.3906/sag-1401-33>.
5. Alabi, Olumuyiwa Samuel, Nuno Mendonça, Olufemi Ezekiel Adeleke, and Gabriela Jorge da Silva. 2017. "Molecular Screening of Antibiotic-Resistant Determinants among Multidrug-Resistant Clinical Isolates of *Proteus Mirabilis* from SouthWest Nigeria." *African Health Sciences* 17 (2): 356–65. <https://doi.org/10.4314/ahs.v17i2.9>.

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Atsauces

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Acar, Sümeyra, Belkıs Levent, and Ekrem Atalan. 2015. "A Molecular Epidemiological Investigation of Multistate Outbreaks of Salmonella Enteritidis from Clinical and Environmental Samples in Turkey, 2000–2010." *Turkish Journal of Medical Sciences* 45 (1): 76–83. <https://doi.org/10.3906/sag-1401-33>.
5. Alabi, Olumuyiwa Samuel, Nuno Mendonça, Olufemi Ezekiel Adeleke, and Gabriela Jorge da Silva. 2017. "Molecular Screening of Antibiotic-Resistant Determinants among Multidrug-Resistant Clinical Isolates of *Proteus mirabilis* from SouthWest Nigeria." *African Health Sciences* 17 (2): 356–65. <https://doi.org/10.4314/ahs.v17i2.9>.

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Nuorodos

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Acar, Sümeyra, Belkıs Levent, and Ekrem Atalan. 2015. "A Molecular Epidemiological Investigation of Multistate Outbreaks of Salmonella Enteritidis from Clinical and Environmental Samples in Turkey, 2000–2010." *Turkish Journal of Medical Sciences* 45 (1): 76–83. <https://doi.org/10.3906/sag-1401-33>.
- Alabi, Olumuyiwa Samuel, Nuno Mendonça, Olufemi Ezekiel Adeleke, and Gabriela Jorge da Silva. 2017. "Molecular Screening of Antibiotic-Resistant Determinants among Multidrug-Resistant Clinical Isolates of *Proteus Mirabilis* from SouthWest Nigeria." *African Health Sciences* 17 (2): 356–65. <https://doi.org/10.4314/ahs.v17i2.9>.

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Referencie

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Acar, Sümeýra, Belkis Levent, and Ekrem Atalan. 2015. "A Molecular Epidemiological Investigation of Multistate Outbreaks of *Salmonella* Enteritidis from Clinical and Environmental Samples in Turkey, 2000–2010." *Turkish Journal of Medical Sciences* 45 (1): 76–83. <https://doi.org/10.3906/sag-1401-33>.
- Alabi, Olumuyiwa Samuel, Nuno Mendonça, Olufemi Ezekiel Adeleke, and Gabriela Jorge da Silva. 2017. "Molecular Screening of Antibiotic-Resistant Determinants among Multidrug-Resistant Clinical Isolates of *Proteus* *Mirabilis* from SouthWest Nigeria." *African Health Sciences* 17 (2): 356–65. <https://doi.org/10.4314/ahs.v17i2.9>.

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