	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Piperacillin/Tazobactam Discs 36 µg/110 µg (TZP36/TZP110)

REF CT1616B / CT0725B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the Oxoid AST generic instructions for use supplied with the product and available online.

Piperacillin/Tazobactam Antimicrobial Susceptibility Test (AST) Discs are 6mm discs which contain specific amounts of the antimicrobial agents piperacillin/tazobactam. The discs are labelled on both sides with details of the antimicrobial (TZP) and amount present (µg): TZP36 (36 µg) or TZP110 (110 µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Piperacillin/Tazobactam AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for *in vitro* susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which Piperacillin and Tazobactam have been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling, and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100,^{1ac} or EUCAST² breakpoint tables must be used to interpret zone size results.

Species with published breakpoints according to current FDA literature:

Gram-Negative

- *Enterobacteriaceae*
- *Pseudomonas aeruginosa*
- *Acinetobacter* spp.
- *Haemophilus influenzae* and *parainfluenzae*

Gram-Positive

- N/A

Species with published breakpoints according to current EUCAST literature (all associated with 36 µg):

Gram-Negative

- *Enterobacterales*
- *Pseudomonas* spp.
- *Haemophilus influenzae*
- *Achromobacter xylosoxidans*

Gram-Positive

- N/A

Species with published breakpoints according to current CLSI literature (all associated with 110 µg):

Gram-Negative

- *Enterobacterales*
- *Pseudomonas aeruginosa*
- *Acinetobacter* spp.
- *Haemophilus influenzae* and *parainfluenzae*

Gram-Positive

- N/A

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Oxoid Piperacillin/Tazobactam AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. A suitable therapeutic agent for in vivo use can be determined using filter paper discs impregnated with specified concentrations of antimicrobial agents placed on the surface of a suitable test medium. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognized zone diameter ranges for the specific antimicrobial agent/organism combination under test.

For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology, refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation in the References section below.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Oxoid Piperacillin/Tazobactam AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST.²

Materials Provided

Oxoid Piperacillin/Tazobactam AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Piperacillin/Tazobactam AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with CT1616B & CT0725B	
Component Description	Material Description
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Pale beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its dessiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Piperacillin/Tazobactam AST Disc Reagents	
Reagent	Description of Function
Piperacillin	White to off-white powder. Piperacillin is a penicillin antibiotic that exerts its mechanism of action through inhibiting the synthesis of the bacterial cell wall and septum.
Tazobactam	White to off-white powder. Tazobactam is a beta-lactamase inhibitor antimicrobial. Tazobactam exerts its mechanism of action through inhibiting the bacterial beta-lactamase enzymes, preventing the breakdown of piperacillin by said enzymes.

The concentration of antibiotic on the AST Disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Oxoid Piperacillin/Tazobactam AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Piperacillin/Tazobactam AST Discs have a shelf-life of 3 years if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Raw Data of Zone Sizes in Accordance with EUCAST Methodology Taken from the 3 Most Recent Batches of Piperacillin/Tazobactam (36 µg) AST Discs.

EUCAST Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3278315	<i>Escherichia coli</i>	35218	MHA	21	27	24	22	22	22	2 ± 0 (CV=0)
	<i>Escherichia coli</i>	25922	MHA	21	27	24	23	23	23	-1 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	23	29	26	26	26	25	-0.33 ± 0.58 (CV=-1.7)
	<i>Klebsiella pneumoniae</i>	700603	MHA	14	20	17	17	17	17	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV) (%)
							1	2	3	
3230873	<i>Escherichia coli</i>	35218	MHA	21	27	24	23	23	23	-1 ± 0 (CV=0)
	<i>Escherichia coli</i>	25922	MHA	21	27	24	23	23	23	-1 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	23	29	26	25	25	25	-1 ± 0 (CV=0)
	<i>Klebsiella pneumoniae</i>	700603	MHA	14	20	17	17	17	17	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV) (%)
							1	2	3	
3207306	<i>Escherichia coli</i>	35218	MHA	21	27	24	22	22	22	-2 ± 0 (CV=0)
	<i>Escherichia coli</i>	25922	MHA	21	27	24	23	23	23	-1 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	23	29	26	25	25	25	-1 ± 0 (CV=0)
	<i>Klebsiella pneumoniae</i>	700603	MHA	14	20	17	17	17	17	0 ± 0 (CV=0)

Table 4. Raw Data of Zone Sizes in Accordance with CLSI Methodology Taken from the 3 Most Recent Batches of Piperacillin/Tazobactam (110 µg) AST Discs.

CLSI Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3278316	<i>Staphylococcus aureus</i>	25923	MHA	27	36	31.5	29	29	29	-2.5 ± 0 (CV=0)
	<i>Escherichia coli</i>	25922	MHA	24	30	27	24	24	24	-3 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	25	33	29	29	29	29	0 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	33	38	36	34	34	34	-2 ± 0 (CV=0)
	<i>Escherichia coli</i>	35218	MHA	24	30	27	24	24	24	-3 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV) (%)
							1	2	3	
3245370	<i>Staphylococcus aureus</i>	25923	MHA	27	36	31.5	33	33	33	1.5 ± 0 (CV=0)
	<i>Escherichia coli</i>	25922	MHA	24	30	27	27	27	27	0 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	25	33	29	29	29	30	0.33 ± 0.58 (CV=1.7)
	<i>Haemophilus influenzae</i>	49247	HTM	33	38	36	37	37	37	1 ± 0 (CV=0)
	<i>Escherichia coli</i>	35218	MHA	24	30	27	25	25	26	-1.67 ± 0.58 (CV=-0.35)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV) (%)
							1	2	3	
3226130	<i>Staphylococcus aureus</i>	25923	MHA	27	36	31.5	34	34	36	3.16 ± 1.15 (CV=0.36)
	<i>Escherichia coli</i>	25922	MHA	24	30	27	27	27	27	0 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	25	33	29	28	29	29	-0.33 ± 0.58 (CV=-1.7)
	<i>Haemophilus influenzae</i>	49247	HTM	33	38	36	35	37	37	0.83 ± 1.15 (CV=1.39)
	<i>Escherichia coli</i>	35218	MHA	24	30	27	26	26	26	-1 ± 0 (CV=0)

Clinical Performance Characteristics

The purpose of this study was to determine the causative pathogens and their antimicrobial susceptibility pattern in neonatal sepsis patients.³ The study was done at the Tripoli Medical Centre in Libya from 1997 to 2001. Isolates were identified using routine biochemical techniques and antimicrobial susceptibility testing was performed using the disc diffusion method. A variety

of antibiotics were tested, including the Oxoid Piperacillin/Tazobactam AST Discs, according to the National Committee for Clinical Laboratory Standards (now the CLSI guidelines). The isolates were categorised as resistant or susceptible.

The study collected a total of 1,431 blood cultures of which 801 were positive for bacterial growth. Susceptibility testing results demonstrated 43%, 44%, 26% and 9% of *Serratia marcescens*, *Serratia liquefaciens*, *Klebsiella pneumoniae*, and *Enterobacter agglomerans* isolates to be resistant to piperacillin/tazobactam, respectively. All *E. cloacae* isolates were found to be susceptible to piperacillin/tazobactam. This high susceptibility rate of piperacillin/tazobactam was much higher than most other studied antibiotics (Figure 1).³

The study found piperacillin/tazobactam to be effective against Gram-negative isolates causing neonatal septicaemia especially among the *Enterobacter cloacae* and *E. agglomerans* isolates. Susceptibility rates to piperacillin/tazobactam were much higher than most other studied antibiotics. This study did not report any inconsistent or out of specification results. Therefore, this study has demonstrated that the Oxoid Piperacillin/Tazobactam AST Discs are effective for in vitro analysis and potential clinical treatment of Gram-negative causing neonatal sepsis.

Another study was done to evaluate three combination antimicrobials against extended spectrum beta-lactamase (ESBL) producing Gram negative bacteria in India.⁴

This retrospective study collected a variety of clinical samples from hospitalised patients in New Delhi from September 2003 to August 2004. Gram-negative isolates were identified and their ESBL status was determined. The isolates had their antimicrobial susceptibility tested using a variety of antibiotics, including the Oxoid Piperacillin/Tazobactam (110 µg) TZP110 AST Discs. The susceptibility testing was performed in accordance to NCCLS (now the CLSI guidelines) and isolates were categorised as resistant, intermediate, or sensitive. The following reference strains were used as quality control: *P. aeruginosa* ATCC-27853, *E. coli* ATCC-25922, and *K. pneumoniae* ATCC-700603.

The study collected a total of 80,116 clinical samples of which 13,091 were Gram-negative and 9,004 of these were ESBL positive. The susceptibility testing demonstrated that piperacillin/tazobactam was the most active antimicrobial towards these isolates, with a total of 81.37% of isolates being susceptible (Figure 2).⁴ In comparison, cefoperazone/sulbactam (C/S) and ticarcillin/clavulanic acid (T/C) demonstrated a susceptibility rate of 76.06% and 45.48%, respectively. Piperacillin/tazobactam susceptibility among the *Pseudomonas spp.*, *E. coli* and *Klebsiella spp.* isolates was found to be significantly higher when compared to C/S.

This study has shown that the combination antibiotic piperacillin/tazobactam has a higher in vitro activity towards ESBL producing gram-negative bacteria when compared to other combination antibiotics. This study did not report any out of specification or inconsistent results. Therefore, this study has shown that the Oxoid Piperacillin/Tazobactam (110 µg) TZP110 AST Discs are an effective and consistent method for the in vitro analysis of ESBL positive gram-negative pathogens.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. Disk Diffusion Method Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
 - b. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
3. M. S. Ellabib, *et al.* (2004). Changing Pattern of Neonatal Bacteremia *Saudi Medical Journal* 25(12):1957-1956.
4. S. Mohanty, *et al.* (2005). *The Indian Journal of Medical Research* 122 (5): 425-428.

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of Symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit

	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision Information

Version	Date of issue and modifications introduced
1.0	2023-01-17

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. Disk Diffusion Method Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
 - b. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
3. M. S. Ellabib, et al. (2004). Changing Pattern of Neonatal Bacteremia Saudi Medical Journal 25(12):1957-1956.
4. S. Mohanty, et al. (2005). The Indian Journal of Medical Research 122 (5): 425-428.

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